

USER GUIDE

# Chromium Single Cell DNA Training Kit

FOR USE WITH

Chromium Single Cell DNA Training Kit, 4 rxns PN-1000061

## Notices

### Document Number

CG000171 | Rev B

### Legal Notices

© 2019 10X Genomics, Inc (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/trademarks](http://www.10xgenomics.com/trademarks). 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: [www.10xgenomics.com/patents](http://www.10xgenomics.com/patents). The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at [www.10xgenomics.com/legal-notices](http://www.10xgenomics.com/legal-notices), or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

### Instrument & Licensed Software Updates Warranties

Updates to existing Instruments and Licensed Software may be required to enable customers to use new or existing products. In the event of an Instrument failure resulting from an update, such failed Instrument will be replaced or repaired in accordance with the 10x Limited Warranty, Assurance Plan or service agreement, only if such Instrument is covered by any of the foregoing at the time of such failure. Instruments not covered under a current 10x Limited Warranty, Assurance Plan or service agreement will not be replaced or repaired.

### Chromium Controller Usable Temperature Range

The recommended temperature range for the Chromium Controller or the Chromium Single Cell Controller when running the protocol outlined in this Chromium Single Cell DNA Reagents Kits User Guide is 19-25°C (66-77°F), which is different than the usable temperature range of 18-28°C (64-82°F) stated in the Chromium Controller Specifications (Document CG00020) and the Chromium Single Cell Controller Specifications (Document CG00050). Running the Chromium Controller or the Chromium Single Cell Controller outside the recommended temperature range of 19-25°C (66-77°F) when using the reagent kits and chip kits described herein will invalidate the warranty of these products.

### Support

Email: [support@10xgenomics.com](mailto:support@10xgenomics.com)

10x Genomics

6230 Stoneridge Mall Road

Pleasanton, CA 94588 USA

---

# Table of Contents

Introduction	4
Objective	5
Chromium Single Cell DNA Training Kit	6
Chromium Accessories	7
Recommended Thermal Cyclers	7
Additional Kits, Reagents & Equipment	8
Tips & Best Practices	9
Training Step 1	13
1.0 Cell Bead Generation	14
1.1 Load Chromium Chip C	16
1.2 Run the Chromium Controller & Transfer Cell Beads	18
Training Step 2	19
2.0 Cell Bead Processing	20
2.1 Solvent Exchange	21
2.2 Lysis	22
2.3 Filtration	23
Troubleshooting	24
7.1 Cell Beads	25
7.2 Chromium Controller Errors	27

# Introduction

Objective

Chromium Single Cell DNA Training Kit

Chromium Single Cell DNA Accessories

Chromium Accessories

Recommended Thermal Cyclers

Additional Kits, Reagents & Equipment

## Objective

The purpose of this User Guide is to train new users on:

- How to handle Activation Agent and Cell Matrix
- How to prepare Cell Bead Mix
- How to load a Chromium Chip C with the Cell Bead Mix, CB Polymer, and Partitioning Oil
- How to load a Chromium Chip C into the Chromium Controller and run the Controller
- How to inspect the resulting Cell Beads in the Chromium Chip C
- How to transfer the Cell Beads
- How to process the Cell Beads in preparation for loading Chromium Chip D

Firmware Version 3.0 or higher is required in the Chromium Controller or the Single Cell Chromium Controller used for this protocol.

For guidance on other items, refer to the User Guides cited below:

- For guidance on qualifying the Chromium Controller or Chromium Single Cell Controller, refer to the Chromium Controller Specifications (CG00020) or the Chromium Single Cell Controller Specifications (CG00050), and the Chromium Controller Readiness Test (CG00025) or the Chromium Single Cell Controller (CG00051).
- For guidance on preparing a sample that can be quantified and moved through library generation and sequencing, refer to the applicable Sample Preparation Demonstrated Protocols and to the Chromium Single Cell DNA Reagent Kits User Guide (CG000153).

## Chromium Single Cell DNA Training Kit

### Chromium Single Cell DNA Training Kit, 4 rxns PN-1000061

#### Chromium Single Cell DNA Cell Bead Reagent Kit, 4 rxns PN-1000034 (store at -20°C)

##### Chromium Single Cell DNA Cell Bead Reagents

	#	PN
Activation Agent (4 rxns)	1	2000039
<input type="radio"/> CB Buffer	2	2000040

10xGenomics.com

10x  
GENOMICS

#### Chromium Single Cell DNA Cell Bead Polymer Kit, 4 rxns PN-1000035 (store at 4°C)

##### Chromium Single Cell DNA Cell Bead Polymer

	#	PN
<input checked="" type="radio"/> CB Polymer	1	2000041
Cell Matrix (4 rxns)	1	2000042
<input type="radio"/> Lysis 1	1	2000049
<input checked="" type="radio"/> Lysis 2	1	3000146

10xGenomics.com

10x  
GENOMICS

#### Chromium Chip C Single Cell DNA Kit, 16 rxns PN-1000032 (store at ambient temperature)

##### Chromium Partitioning Oil

	#	PN
Partitioning Oil	2	220088

##### Chromium Recovery Agent

	#	PN
<input type="radio"/> Recovery Agent	2	220016

##### Chromium Chip C & Gaskets

	#	PN
Chip C Single Cell DNA	2	2000030
Gasket	1	3000072

10xGenomics.com

10x  
GENOMICS

## Flowmi™ Filters, 50 rxns, PN-1000055

	#	PN
Flowmi	1	1000055

## Chromium Single Cell DNA Accessory Kit, PN-1000058

Product	#	PN (Orderable)	PN (Item)
10x Magnetic Separator A	1	1000054*	2000067
10x Chromium Chip D Holder	1	1000053*	3000109
CABLE USB A MALE - B MALE 1M BLK	1	-	3000173

\*May be ordered individually or as part of the Chromium Single Cell DNA Accessory Kit. Also available as part of the Chromium Controller Accessory Kit, PN-110204.

## Chromium Accessories

Product	PN (Orderable)	PN (Item)
10x Vortex Adapter	120251	330002
10x Vortex Clip	120253	230002
10x Chip Holder	120252	330019

## Recommended Thermal Cyclers

Thermal cyclers used must support uniform heating of 100 µl emulsion volumes.

Supplier	Description	Part Number
BioRad	C1000 Touch Thermal Cycler with 96-Deep Well Reaction Module	1851197
Eppendorf	MasterCycler Pro	North America 950030010 International 6321 000.01
Thermo Fisher Scientific	Veriti 96-Well Thermal Cycler	4375786

## Additional Kits, Reagents & Equipment

The items in the table below have been validated by 10x Genomics and are highly recommended for the Chromium Single Cell DNA protocol. Substituting materials may adversely affect system performance.

Supplier	Description	Part Number (US)
<b>Plastics</b>		
Eppendorf	PCR Tubes 0.2 ml 8-tube strips	951010022
	DNA LoBind Tubes, 1.5 ml	022431021
	DNA LoBind Tubes, 2.0 ml	022431048
USA Scientific	TempAssure PCR 8-tube strip	1402-4700
Rainin	Tips LTS W-O 200UL Filter RT-L200WFLR	30389241
	Tips LTS 20UL Filter RT-L10FLR	30389226
	Tips LTS 200UL Filter RT-L200FLR	30389240
	Tips LTS 1ML Filter RT-L1000FLR	30389213
<b>Kits &amp; Reagents</b>		
Thermo Fisher Scientific	Nuclease-free Water	AM9937
Corning Cellgro	Phosphate-Buffered Saline (PBS) 1X without calcium and magnesium	21-040-CV
Millipore Sigma	1 N NaOH	S2770-100ML
	Phosphate-Buffered Saline (PBS) with 10% Bovine Albumin	SRE0036
Ricca Chemical Company	Glycerin (glycerol), 50% (v/v) Aqueous Solution	3290-32
Teknova	1 M Tris-HCl, pH 7.2	T1072
<b>Equipment</b>		
VWR	Vortex Mixer	10153-838
	Divided Polystyrene Reservoirs	41428-958
Eppendorf	Eppendorf ThermoMixer C	5382000023
	Eppendorf SmartBlock PCR 96	5306000006
	Eppendorf SmartBlock 1.5 ml, Thermoblock for 24 Reaction Vessels	5360000038
Rainin	Pipet-Lite LTS Pipette L-2XLS	17014393
	Pipet-Lite LTS Pipette L-10XLS	17014388
	Pipet-Lite LTS Pipette L-20XLS	17014392
	Pipet-Lite LTS Pipette L-100XLS	17014384
	Pipet-Lite LTS Pipette L-200XLS	17014391
	Pipet-Lite LTS Pipette L-1000XLS	17014382
	Pipet-Lite Multi Pipette L8-10XLS	17013802
	Pipet-Lite Multi Pipette L8-20XLS	17013803
	Pipet-Lite Multi Pipette L8-50XLS	17013804
	Pipet-Lite Multi Pipette L8-200XLS	17013805



# Tips & Best Practices



## Icons



Tip & Best Practices section includes additional guidance



Signifies critical step requiring accurate execution



Troubleshooting section includes additional guidance

## Images

- The Training Kit is for 4 rxns. However, all images and illustrations shown are for 8 rxns.

Emulsion-safe  
Plastics

- Use 10x Genomics validated emulsion-safe plastic consumables when handling GEMs as some plastics can destabilize GEMs.

General  
Reagent  
Handling

- Fully thaw and thoroughly mix reagents before use.
- Calculate reagent volumes with 10% excess of 1 rxn values.
- Cover Partitioning Oil tubes and reservoirs to minimize evaporation.

50% Glycerol  
Solution

- Purchase 50% glycerol solution from Ricca Chemical Company, Glycerin (glycerol), 50% (v/v) Aqueous Solution, PN-3290-32.
- OR
- Prepare 50% glycerol solution:
    - i. Mix an equal volume of water and 99% Glycerol, Molecular Biology Grade.
    - ii. Filter through a 0.2- $\mu$ m filter.
    - iii. Store at  $-20^{\circ}\text{C}$  in 1-ml LoBind tubes. 50% glycerol solution should be equilibrated to room temperature before use.

## Pipette Calibration


- Follow manufacturer's calibration and maintenance schedules.
- Pipette accuracy is **critical** in **Training Step 2**.

## Chromium Chip Handling

- Minimize exposure of reagents, chips, and gaskets to sources of particles and fibers, laboratory wipes, frequently opened flip-cap tubes, and dusty surfaces.
- Execute steps without pause or delay, unless indicated. When multiple chips are to be used, load, run, and collect the content from one chip before loading the next.
- Fill all unused input wells in rows labeled 1, 2, and 3 on a chip with an appropriate volume of 50% glycerol solution before loading the used wells. **DO NOT** add glycerol to the Recovery Wells.
- Avoid contacting the bottom surface of the chip with gloved hands and other surfaces. Frictional charging can lead to inadequate priming of the channels, potentially leading to either clogs or wetting failures.
- Minimize the distance that a loaded chip is moved to reach the Chromium Controller.
- Keep the chip horizontal to prevent wetting the gasket with oil, which depletes the input volume and may adversely affect the quality of the assay.


## 10x Chip Holders

- 10x Chip Holders encase Chromium Chips.
- The holder lid flips over to become a stand, holding the chip at 45 degrees for optimal Recovery Well content removal.
- Squeeze the black sliders on the back side of the holder together to unlock the lid and return the holder to a flat position.
- Use the 10x Chip Holder to encase Chip C.



## Chromium Chip & Holder Assembly

- Align notch on the chip (upper left corner) and the holder.
- Insert the left-hand side of the chip under the guide. Depress the right-hand side of the chip until the spring-loaded clip engages.
- Close the lid before dispensing reagents into the wells.



### Activation Agent Handling

- Use one tube of Activation Agent per sample. DO NOT puncture foil seal of tubes not used at the time.
- Store unused tubes at  $-20^{\circ}\text{C}$ .



### Cell Matrix Handling

- Use one tube of Cell Matrix per sample. DO NOT puncture foil seal of tubes not used at the time.
- Store unused tubes at  $4^{\circ}\text{C}$ . DO NOT freeze Cell Matrix.



### Chromium Chip C Loading

- Place the assembled chip and holder flat on the bench with the lid closed.
- Dispense at the bottom of the wells without introducing bubbles.
- Wait for the Cell Bead Mix to drain into the bottom of the pipette tips and dispense again to ensure complete volume transfer.
- Refer to [Load Chromium Chip C](#) for specific instructions.



### 10x Magnetic Separator A

- Offers one position of the magnets relative to a tube.
- 10x Magnetic Separator A is different from 10x Magnetic Separator.
- Use for Cell Bead Processing.






# Training Step 1

## Cell Bead Generation

- 1.1 Load Chromium Chip C
- 1.2 Run the Chromium Controller & Transfer Cell Beads

1

## 1.0 Cell Bead Generation

GET STARTED!				
Action	Item	10x PN	Preparation & Handling	Storage
<b>Equilibrate to Room Temperature</b>	1X PBS	-	-	-
<b>Place on Ice</b>	<b>Activation Agent</b> (Cell Bead Reagent Kit)	2000039	Use one tube per sample. Puncture the tube foil seal immediately before use. Return unused tubes to -20°C.	-20°C
				
	<b>Cell Matrix</b> (Cell Bead Polymer Kit)	2000042	Use one tube per sample. Puncture the tube foil seal immediately before use. Return unused tubes to 4°C.	4°C
				
	 <b>CB Polymer</b> (Cell Bead Polymer Kit)	2000041	-	4°C
	<b>Training Sample</b> PBS + 0.04% BSA	-	Training sample (PBS + 0.04% BSA) is used as a substitute for single cell suspension, which is prepared in PBS + 0.04% BSA.	-
<b>Obtain</b>	<b>Partitioning Oil</b>	220088	-	Ambient
	<b>Nuclease-free Water</b>	-	-	-
	<b>Chromium Chip C</b>	2000030	See Tips & Best Practices.	Ambient
	<b>10x Vortex Adapter</b>	330002	See Tips & Best Practices.	Ambient
	<b>10x Gasket</b>	370017	See Tips & Best Practices.	Ambient
	<b>10x Chip Holder</b>	330019	See Tips & Best Practices.	Ambient
	<b>50% glycerol solution</b> If using <8 reactions	-	See Tips & Best Practices.	Ambient



Firmware Version 3.0 or higher is required in the Chromium Controller or the Single Cell Chromium Controller used for the Single Cell DNA protocol.

## 1.1

### Load Chromium Chip C

See Tips & Best Practices for chip handling instructions. When loading the chip, raising and depressing the pipette plunger should each take ~5 sec. When dispensing, raise the pipette tips at the same rate as the liquid is rising, keeping the tips slightly submerged.

#### TIPS

**a. Assemble Chromium Chip C in a 10x Chip Holder** See Tips & Best Practices

**b. Dispense 50% Glycerol Solution into Unused Chip Wells** (if <8 samples per chip)

- i. 75  $\mu$ l into unused wells in row labeled 1.
- ii. 40  $\mu$ l into unused wells in row labeled 2.
- iii. 200  $\mu$ l into unused wells in row labeled 3.

DO NOT add 50% glycerol solution to the top row of Recovery Wells.  
DO NOT use any substitute for 50% glycerol solution.

**c. Prepare Cell Bead Mix**


- i. Remove the Activation Agent and Cell Matrix tube strips from the holders. Centrifuge briefly. Place on ice.
- ii. Using a pipette tip, puncture the foil seal of one Activation Agent and one Cell Matrix tube per sample and enlarge the holes, avoiding contact between the pipette tip and the Activation Agent or the Cell Matrix.
- iii. Add 20  $\mu$ l nuclease-free water to each Activation Agent tube. Pipette mix to dissolve.

DO NOT mix during steps iv-vi.

- iv. Layer 4.7  $\mu$ l Activation Agent solution on top of the Cell Matrix avoiding contact between the Cell Matrix and the pipette tip.  
DO NOT mix.

- v. Layer the 6.5  $\mu$ l PBS on top of the Activation Agent. DO NOT mix.

- vi. Resuspend training sample and immediately layer 4.0  $\mu$ l training sample on top of PBS.  
DO NOT mix.




See Tips & Best Practices for chip handling instructions. When loading the chip, raising and depressing the pipette plunger should each take ~5 sec. When dispensing, raise the pipette tips at the same rate as the liquid is rising, keeping the tips slightly submerged.



#### d. Load Row Labeled 1

Using a **50- $\mu$ l multi-channel pipette** (set to 20  $\mu$ l), gently pipette mix the Cell Bead Mix until homogeneous. Using the same pipette tips, dispense **23  $\mu$ l** Cell Bead Mix into the bottom center of each well in the **row labeled 1** without introducing bubbles.

Wait for the Cell Bead Mix to drain into the bottom of the pipette tips and dispense again to ensure complete transfer.




Dislodge any Cell Bead Mix stuck along the sidewalls of the well using a pipette tip. Accurate input volume is critical for correct Cell Bead volume yield and optimal performance.

#### e. Prepare CB Polymer

Centrifuge the CB Polymer briefly. Vortex **10 sec**. Centrifuge briefly again.

#### f. Load Row Labeled 2

Using a **single-channel pipette**, slowly dispense **40  $\mu$ l** CB Polymer into the bottom center of each well in the **row labeled 2** without introducing bubbles.




Wait **60 sec** for the chip to prime.

#### g. Load Row Labeled 3

Dispense **250  $\mu$ l** Partitioning Oil into each well in the **row labeled 3**.

Failure to add Partitioning Oil can damage the Chromium Controller.



#### h. Attach 10x Gasket

Align the notch with the top left-hand corner. Ensure the gasket holes are aligned with the wells. Avoid touching the smooth gasket surface. **DO NOT** press down on the gasket.



Keep horizontal to avoid wetting the gasket.



## 1.2 Run the Chromium Controller & Transfer Cell Beads

- Press the eject button on the Controller to eject the tray.
- Place the loaded chip with the gasket in the tray. Press the button to retract the tray.
- Confirm program on screen. Press the play button.



- At completion of the run (~11.5 min), the Controller will chime. **Immediately** proceed to the next step. Steps e – i should be executed within **5 min**.

- Press the eject button to remove the chip.
- Discard the gasket. Open the chip holder. Fold the lid back until it clicks to expose the wells at 45 degrees.

- Remove **60 µl** Partitioning Oil from lowest point of the Recovery Wells in the top row using **regular-bore** pipette tips.

- Using **wide-bore** pipette tips (pipette set to 95 µl), over the course of **~30 sec** aspirate remaining Partitioning Oil and Cell Beads from the bottom of each Recovery Well and transfer into a PCR 8-tube strip. Dispense slowly with the pipette tips against the sidewalls without introducing bubbles. Using the same pipette tips, transfer any remaining Partitioning Oil and Cell Beads from the Recovery Wells to the tubes.



DO NOT flick, shake or vortex the tubes to prevent wetting the sidewalls and caps of the tubes.

- Immediately** place the capped tube strip on a thermomixer and shake at **21°C** for **16-24 h** at **1000 rpm**.

Usable Temperature Range  
19-25°C (66-77°F)




Expose Wells



Remove Partitioning Oil



Transfer Cell Beads  
(use wide-bore pipette tips)



Cell Beads (wide-bore pipette tips)



# Training Step 2

## Cell Bead Processing

- 2.1 Solvent Exchange
- 2.2 Lysis
- 2.3 Filtration



## 2.0 Cell Bead Processing

GET STARTED!				
Action	Item	10x PN	Preparation & Handling	Storage
<b>Equilibrate to Room Temperature</b>	1X PBS	-	-	-
	<input type="radio"/> CB Buffer (Cell Bead Reagent Kit)	2000040	-	-20°C
<b>Place on Ice</b>	<input type="radio"/> Lysis 1 (Cell Bead Polymer Kit)	2000049	Mix by inverting tube 5x. DO NOT vortex.	4°C
	<input checked="" type="radio"/> Lysis 2 (Cell Bead Polymer Kit)	3000146	Vortex. Centrifuge briefly.	4°C
<b>Obtain</b>	Recovery Agent	220016	-	Ambient
	0.4 N NaOH Prepare 1 ml for 8 reactions	-	Prepare fresh.	Ambient
	1 M Tris-HCl, pH 7.2	-	-	Ambient
	10x Magnetic Separator A	2000067	Different from the 10x Magnetic Separator.	Ambient
	10x Vortex Clip	230002	-	Ambient
	1.5 ml magnetic tube holder	-	-	Ambient
	Flowmi Filters, 70 µm	1000055	-	Ambient

## 2.1 Solvent Exchange



- a. Remove the tube strip containing Cell Beads from the thermomixer.
- b. Remove **45  $\mu$ l** Partitioning Oil from the bottom of the tube. **5-10  $\mu$ l** Partitioning Oil should remain. **DO NOT** aspirate Cell Beads. If Cell Beads are aspirated, return the solution to the tube, reduce removal volume by **5  $\mu$ l**, and reattempt removal.
- c. Add **95  $\mu$ l** PBS.
- d. Cap the tube strip and place in a 10x Vortex Clip. Vortex **20 sec**. Centrifuge **~1 sec** only. The Cell Bead suspension may not be uniform.
- e. Add **60  $\mu$ l** Recovery Agent. Wait **60 sec**.
- f. Gently invert **~25x**. **DO NOT** vortex.
- g. Centrifuge **10 sec** to separate the phases. The top aqueous layer contains Cell Beads. The bottom pink layer contains Recovery Agent/Partitioning Oil. Place on a 10x Magnetic Separator **A** for **30 sec**.

TIPS

10x Magnetic Separator **A** (magnet **A**) is different from 10x Magnetic Separator.


- h. Insert pipette tips along the sidewalls of the tube, away from the Cell Beads and slowly remove **120  $\mu$ l** Recovery Agent/Partitioning Oil/PBS from the bottom of the tube, retaining the Cell Beads. If the level of samples remaining is not at the etched mark on magnet **A**, adjust by either adding or removing PBS.



A small volume of Recovery Agent/Partitioning Oil may remain in the tube. **DO NOT** aspirate Cell Beads.

- i. Add **120  $\mu$ l** PBS. Remove from magnet **A**.
- j. Gently invert until homogeneous and centrifuge **~1 sec**.
- k. Place on magnet **A** for **30 sec**.
- l. Remove **120  $\mu$ l** from the bottom of the tube, retaining the Cell Beads.
- m. Remove from magnet **A**.


Cell Beads from Thermomixer

Remove Oil (~5-10  $\mu$ l Oil Remains)

Add PBS, Vortex, Centrifuge



Phase Separation



Place on Magnetic Separator A



Retain Cell Beads



## 2.2

### Lysis

- a. Add **120 µl** Lysis 1. Mix by inverting until homogeneous. Centrifuge **~1 sec**.
- b. Place on magnet **A** for **30 sec**.
- c. Remove **120 µl** from the bottom of the tube, retaining the Cell Beads.
- d. Incubate in a thermal cycler using the following protocol.

Lid Temperature	Reaction Volume	Run Time
85°C	75 µl	10 min
Step	Temperature	Time
Incubate	72°C	00:10:00
Hold	4°C	Hold

- e. Add **100 µl** Lysis 1. Mix by inverting until homogeneous. Centrifuge **~1 sec**.
- f. Add **5 µl** Lysis 2. **Immediately** mix by inverting until homogeneous.
- g. Incubate **15 min** at **room temperature**. Gently invert every **5 min** to prevent settling of beads.
- h. Centrifuge **~1 sec**. Place on magnet **A**. Wait **30 sec**.
- i. Remove **105 µl** from the bottom of the tube, retaining the Cell Beads.
- j. Add **120 µl** PBS. Mix by inverting until homogeneous. Centrifuge **~1 sec**.
- k. Place on magnet **A**. Wait **30 sec**.
- l. Remove **120 µl** from the bottom of the tube, retaining the Cell Beads.
- m. **Repeat** steps j – l for a total of 2 washes.



- n. Adjust the liquid level such that the meniscus is at the etched mark on magnet **A** by either adding or removing PBS.



- o. Remove from magnet **A**. Add **75 µl** freshly prepared 0.4 N NaOH. Mix by inverting until homogeneous.
- p. Incubate **5 min** at **room temperature**. Gently invert every **~1 min**. Centrifuge **~1 sec** at the end of incubation.
- q. Add **50 µl** 1 M Tris-HCl, pH 7.2. Mix by inverting until homogeneous. Centrifuge **~1 sec**.
- r. Place on magnet **A**. Wait **30 sec**.



- s. Remove **120 µl** from the bottom of the tube, retaining the Cell Beads.
- t. Add **120 µl** CB Buffer. Mix by inverting until homogeneous. Centrifuge **~1 sec**.
- u. Place on magnet **A**. Wait **30 sec**.
- v. Remove **120 µl** from the bottom of the tube, retaining the Cell Beads.
- w. **Repeat** steps t – v for a total of 2 washes.

## 2.3 Filtration

- a. Add **450  $\mu$ l** CB Buffer to a 1.5-ml microcentrifuge tube for each sample.
- b. Using a **single-channel** pipette (set to 75  $\mu$ l), pipette mix Cell Beads. Using the same pipette tip, transfer **all** Cell Beads to the microcentrifuge tube containing CB Buffer.
- c. Using a **single-channel** pipette (set to 520  $\mu$ l), gently pipette mix the Cell Bead suspension (volume may be less than **520  $\mu$ l**). Using the same pipette tip, aspirate the Cell Bead suspension and attach the pipette tip to a **70  $\mu$ m** Flowmi filter.



- d. **Immediately** pass the Cell Bead suspension slowly through the Flowmi filter to the **first stop**, into a new 1.5-ml microcentrifuge tube. Variation in volume of filtrate is normal.
- e. Place the capped tube containing the Cell Bead suspension on a 1.5-ml magnetic tube holder for **3 min** or until the solution clears.
- f. Remove **260  $\mu$ l** solution, retaining the Cell Beads.
- g. Remove from the magnet. Pipette mix using a **200- $\mu$ l single-channel pipette** (set to 185  $\mu$ l) until homogeneous without introducing bubbles. Using the same pipette tip, transfer exactly **185  $\mu$ l** to a new tube strip. Some Cell Beads may remain in the 1.5-ml microcentrifuge tube. If volume is less than **185  $\mu$ l**, add CB buffer from the CB Buffer bottle to bring volume to **185  $\mu$ l**.
- h. Centrifuge **~1 sec**. Place on magnet **A** for **3 min**.
- i. Using a **200- $\mu$ l single-channel pipette**, slowly remove **159  $\mu$ l** from the bottom of the tube. A variable small volume of Cell Beads may be aspirated.
- j. Using a **single-channel pipette**, add **7  $\mu$ l** CB Buffer from the CB Buffer bottle to the retained Cell Beads. Remove from the magnet and place on ice.

### Attach Flowmi Filter



### Pass through Flowmi Filter



### Remove 159 $\mu$ l







# Troubleshooting








- 7.1** Cell Beads
- 7.2** Controller Errors

## 7.0 Troubleshooting

### 7.1 Cell Beads

STEP	NORMAL	PROBLEM & SOLUTION
1.1 d Mixing the Cell Bead Mix	Cell Bead Mix is uniform with no visible magnetic bead clumps.	Cell Bead Mix has magnetic bead clumps.  <b>Solution:</b> Using a pipette tip, aspirate from the bottom of the tube and dispense at the liquid surface multiple times to break the magnetic bead clumps.
1.1 d Loading Cell Bead Mix into Chip C	Dispense 23 $\mu$ l Cell Bead Mix into the bottom center of each well in row labeled 1 without introducing bubbles.	<23 $\mu$ l Cell Bead Mix aspirated in the pipette tip for dispensing.  <b>Solution:</b> Using a pipette tip, push down Cell Bead Mix stuck to the sidewalls of the tube. Slowly aspirate from the bottom of the tube holding the pipette tips in the solution for ~2 sec after the plunger is fully released.
1.2 i Transferring tubes to the thermomixer	 Cell Beads are uniform and not stuck on the sidewalls or caps of the tubes.	 Cell Beads are stuck on the sidewalls and caps of the tubes.  DO NOT flick, shake or vortex the tubes when transferring to the thermomixer. <b>Solution:</b> Centrifuge briefly to recover beads.
2.1 b After removal of Partitioning Oil	 5-10 $\mu$ l Partitioning Oil remaining in the tubes.	 >10 $\mu$ l Partitioning Oil remaining in tube*.  <b>Solution:</b> Remove oil from the tube* until remaining oil volume is 5-10 $\mu$ l.



STEP	NORMAL	PROBLEM & SOLUTION
<p>2.2 n Removal of PBS from the bottom of the tubes, retaining the Cell Beads</p>	 <p>PBS removed uniformly, retaining Cell Beads.</p>	 <p>Variable removal of PBS.</p> <p><b>Solution:</b> Add or remove PBS until uniform.</p>
<p>2.2 s Cell Beads migrating to the magnets</p>	 <p>Cell Beads migrate to the magnets uniformly.</p>	 <p>Cell Beads migrate to the magnets variably.</p> <p><b>Solution:</b> Use pipette tip to gently mix Cell Beads to enable uniform migration to the magnets.</p>
<p>2.3 d Pass Cell Bead suspension through Flowmi filter</p>	<p>Cell Bead suspension passes through the Flowmi filter into a new 1.5-ml microcentrifuge tube.</p>	 <p>Flowmi filter clogs while passing the Cell Bead suspension.</p> <p><b>Solution:</b> Discard clogged Flowmi filter and proceed with a new Flowmi filter. Some loss in Cell Bead suspension volume is expected.</p>

## 7.2

### Chromium Controller Errors

If the Chromium Controller or the Chromium Single Cell Controller fails to start, an error tone will sound and one of the following error messages will be displayed:

- a. **Chip not read – Try again:** Eject the tray, remove and/or reposition the 10x Chip Holder assembly and try again. If the error message is still received after trying this more than twice, contact [support@10xgenomics.com](mailto:support@10xgenomics.com) for further assistance.
- b. **Check gasket:** Eject the tray by pressing the eject button to check there is a 10x Gasket on the Chromium Chip. In the case when the 10x Gasket installation was forgotten, install and try again. In the case when a 10x Gasket was already installed, remove, reapply, and try again. If the error message is still received after trying either of these more than twice, contact [support@10xgenomics.com](mailto:support@10xgenomics.com) for further assistance.
- c. **Pressure not at Setpoint:**
  - i. If this message is received within a few seconds of starting a run, eject the tray by pressing the eject button and check for dirt or deposits on the 10x Gasket. If dirt is observed, replace with a new 10x Gasket and try again. If the error message is still received after trying this more than twice, contact [support@10xgenomics.com](mailto:support@10xgenomics.com) for further assistance.
  - ii. If this message is received after a few minutes into the run, the Chromium Chip must be discarded. **Do not try running this Chromium Chip again as this may damage the Chromium Controller.**
- d. **CAUTION: Chip Holder not Present:** Eject the tray by pressing the eject button to check there is a 10x Chip Holder encasing the Chromium Chip. In the case when the 10x Chip Holder was forgotten, install with a 10x Gasket in place, and try again. If the error message is still received after a 10x Chip Holder is confirmed as in place, contact [support@10xgenomics.com](mailto:support@10xgenomics.com) for further assistance.
- e. **Invalid Chip CRC Value:** This indicates the Chromium Chip has encountered an error, should not be run, and must be discarded. Contact [support@10xgenomics.com](mailto:support@10xgenomics.com) for further assistance.
- f. **Endpoint Reached Early:**

If this message is received, contact [support@10xgenomics.com](mailto:support@10xgenomics.com) for further assistance.