

## CASE STUDY

# Upper-arm Skin Rejuvenation

Using DermaLab Elasticity and Ultrasound Probes

To compare Rennova Elleva and Sculptra

Objective evaluation of collagen density and skin elasticity following treatment of the upper inner arms


## CASE STUDY



### SOLUTION & METHOD

20 adults underwent a single treatment session in which **one arm** received 18 mL of **Rennova Elleva** injectable, and the **other arm** received 18 mL of **Sculptra (Galderma)**. Skin elasticity and collagen density were evaluated using the DermaLab Elasticity and Ultrasound Probes. Measurements were performed before treatment and at 30, 60, 90, and 120 days post treatment.

### ULTRASOUND PROGRESSION



**Left image:** Ultrasound image of upper arm before injectable. **Right image:** Upper arm 120 days after treatment with Rennova Elleva. The increasing white/yellow coloration reflects an increase in collagen deposition within the dermis **Reference:** Figure lent from Bravo et al.: "Poly-L-lactic acid (PLLA) in arm contouring: a prospective and blinded trial of Rennova® Elleva and Sculptra®", *Frontiers in Surgery*, Vol 13 (2026), Figure 2.



### RESULTS AND CONCLUSIONS

Both products **increased collagen density** in the upper arms. Furthermore, **Skin elasticity improved** significantly in both treatment sides.

No significant differences were observed between the treated sides, suggesting that both Rennova Elleva and Sculptra are effective options for upper arm rejuvenation.

### MEASUREMENT AREA



**Upper inner arms**



### Case information

**Study title:** Poly-L-lactic acid (PLLA) in arm contouring: a prospective and blinded trial of Rennova® Elleva and Sculptra®

**Published in:** *Frontiers in Surgery*, Vol 13 (2026)

**Authors:** Bravo et al.

**Aim:** To evaluate and compare the clinical and aesthetic outcomes of Rennova Elleva and Sculptra in upper arms.

**Please note:** Measurements were performed using the DermaLab Combo; alternatively, collagen density and elasticity measurements may be obtained and visualized using the DermaLab Aesthetic.

