EN Instructions for use

Composition: Cross-linked sodium hyaluronate 25 mg/mL with 0.3% lidocaine hydrochloride. Phosphate buffer pH 6.7–7.3, sodium chloride.

Description: The device is a sterile, biodegradable, viscoelastic, clear, transparent, isotonic and homogenized injectable gel implant in a prefilled syringe. The device consists of cross-linked HA obtained from *Streptococcus equi* bacteria, formalized to a concentration of 25 mg/mL in a physiologic buffer supplemented with 0.3% lidocaine hydrochloride. The auxiliary role of lidocaine hydrochloride is to reduce the pain associated with injections into the skin. Each box contains one prefilled syringe of 1.0 mL gel implant, two 27 G 1/2 disposable sterile needles and an instruction for use. A label including batch number, serial number and UDI is located at the rear side of the box. This label must be detached and affixed to the field "patient label" (patient label) on the "Implant card" (implant card) at the end of the iFU. The implant card must be handed out to the patient to ensure traceability. The prefilled syringe constitutes the sterile barrier of the gel implant. The needle is sterilized in its container constituting the sterile barrier.

Indications: The device is indicated to correct moderate to severe midface volume deficit in the zygomaticomalar region, antero-medial cheek region and submalar region. The device has to be administered by deep subcutaneous and/or supraperiosteal injection.

Intended purpose: The intended purpose of the device is to restore face volume in order to correct moderate to severe midface volume deficit to treat signs of aging.

Contraindications: The device must not be used in: - patients who tend to develop hypertrophic scarring, have pigment disorders or have a susceptibility to keloid formation. | - patients with a history of autoimmune disease or who are receiving immune therapy. | - patients who are known to be hypersensitive to hyaluronic acid, lidocaine hydrochloride or to amide-type local anaesthetics. | patients who have previously received permanent fillers in the area to be treated. | - pregnant or breastfeeding women. | - patients under 18 years of age. Patients on anticoagulant therapy or patients receiving platelet aggregation inhibitors (e.g. acetylsalicylic acid) should not be treated with the device without consulting their doctors. The device must not be used in areas presenting cutaneous, inflammatory and/or infectious processes (e.g. acne, herpes, ...). The device must not be used in association with laser therapy, chemical peeling, dermabrasion or mesotherapy.

Precautions for use: In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection. Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications. Sensitive skin may be pretreated using a local anaesthetic patch or cream. Discard the syringe, remaining product and needle after use in a special container.

Undesired side effects: Physicians must inform the patient that there are potential side effects and/or incompatibilities associated

with implantation of this product, which may occur immediately or may be delayed. Adverse device effects obtained from clinical study and post-marketing surveillance were as follows: facial pain, redness, granulomas, induration, bacterial infection, edema, inflammation at the injected site were observed (discoloration, discoloration (dispigmentation), hematoma, nodules (inflammatory and non-inflammatory), pain, pruritus, swelling. In addition, the following events and reactions have been observed with similar products and were considered as potential risks for the device (adverse events, without at least a suspected causal relationship with the device, are not listed): abscess, angioedema, bleeding/lumpiness, bleeding, bruising, burning, discomfort, erythema, firmness, herpes reactivation, hypersensitivity/allergy, implant-site mass, implant visibility, irritation, itching, lipoed pattern after accidental arterial embolization, malair edema, migration, necrosis due to vascular compromise, numbness, paresthesia, retinal artery occlusion, scar sarcoidosis, scleromyedema (generalized), telangiectasia, tenderness, vasculitis, vasospasm, vasovagal reaction during injection, vision loss due to retinal artery occlusion. Herpes virus reactivations may be provoked in rare cases by direct damage of neuronal axons by the needle or by tissue manipulation and inflammatory reaction after filler injection. The dissection of the subepidermal plane by filler needle use may increase the incidence of local adverse events. Too superficial placement of the filler or an uneven distribution of the injected product may lead to visible, pale nodules in the skin. It is therefore important to take these possible complications into account. Patients must inform their physician as soon as possible about any inflammatory reactions persisting for more than one week or any other occurring secondary effect. The physician should treat these side effects appropriately. Any undesirable side effects associated with injection of the filler must be reported to the distributor and/or to the manufacturer.

Methods of use: The product must be injected into non-inflamed, disinfected, healthy skin. The technique used is essential for the success of the treatment. Use the 27 G 1/2 needle which is provided with the syringe and inject slowly by applying the appropriate injection technique. Inject low volumes in two or more sessions instead of high volumes in one session. The amount injected will depend on the volume deficit to be corrected and the area to be treated. The maximum dosage is 6 mL per treatment session. Do not inject more than 9 mL per patient per year. Inject the product slowly and apply the least amount of pressure necessary. After the injection, doctors may perform a light massage in order to distribute the product uniformly. Let the patient stay in the office for several minutes after the injection to detect blanching caused by arterial occlusions. Based on a clinical study, filling effect of the device gradually decreased over time. Mild to moderate effect lasted up to 24 weeks in 94.9% patients, 36 weeks in 83.8% patients and 52 weeks in 76.3% patients.

Warnings: Verify the integrity of the syringe and the expiry date before use. Do not use a syringe with an open or shifted tip cap or in an opened or damaged blister. As long as the syringe is stored in its originally sealed blister, the gel implant inside the intact syringe is guaranteed to be sterile. Do not use a needle with a damaged container. Do not use any other needle or syringe than provided by the manufacturer. Do not manipulate / bend the

needle. Do not re-use syringe or needle. The re-use of the device creates a potential infection risk for patients or users. If the 27 G 1/2 needle is blocked, do not increase the pressure on the plunger rod, but stop the injection and replace the needle. Introduction of product into the vasculature may lead to embolization, occlusion of vessels, ischemia, or infarction. Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral haemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur. There are no available clinical data (efficacy, tolerance) about injecting the device into an area which has already been treated with any other filling product or botulinum toxin. Do not inject into blood vessels, bones, tendons, ligaments, nerves or muscles. Do not inject the device into naevi. Do not overcorrect. Patients should be advised not to apply any make-up for 12 hours after the injection and to avoid prolonged exposure to sunlight and UV or tanning beds or Turkish baths for one week after the injection. To avoid a possible risk of product mobility, the patient should be advised not to massage the treatment site for a few days following the injection. There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds such as benzalkonium chloride solutions. Therefore the device should never be placed in contact with these substances or with medical-surgical instruments that have been in contact with these substances. The implanted viscoelastic gel is safe in a magnetic resonance environment. The needle and prefilled syringe system cannot be used in a magnetic resonance environment.

Storage: The device should be stored at 2–25°C / 36–77°F in a dry place in the original blister and folding box and protected from light, heat and frost. Handle with care.

ⓐⓑⓒⓓⓔ Instructions for the correct handling of the device: Hold the Luer-Lock adapter as shown in 1. To remove the top cap, twist 2 and pull carefully 3. Following the instruction above will prevent the formation of air bubbles. Hold the syringe as shown in 4. Open the enclosed needle case and insert the needle firmly 5 (do not use any other needle). Tightly secure the needle by twisting it clockwise 6.

ⓕ Explanation of international symbols: 01: Caution 02: Consult instructions for use 03: Keep dry 04: Temperature limitation 05: Keep away from sunlight 06: Do not use if package is damaged 07: Do not re-use 08: Do not restrict 09: Does not contain any latex 10: Use by 11: Batch number 12: Sterile fluid path inside syringe. Sterilized using steam. Syringe outside not sterile. 13: Sterilized using ethylene oxide 14: Single syringe for single use 15: Single needle for single use 16: Manufacturer 17: Date of manufacture 18: Unique device identifier 19: Reference number 20: Serial number

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Prefilled Syringe CROMA-PHARMA GmbH