## Saypha® FILLER Lidocaine

EN

## Instructions for use

**Composition:**

Cross-linked sodium hyaluronate 23 mg/mL with 0.3% lidocaine hydrochloride. Phosphate buffer pH 6.7–7.3 q.s., sodium chloride, BDDE (crosslinking agent). One syringe contains 1.0 mL of Saypha® FILLER Lidocaine. Sterilized using steam.

**Description:**

Saypha® FILLER Lidocaine is a sterile, biodegradable, viscoelastic, clear, transparent, isotonic and homogenized injectable gel implant. Saypha® FILLER Lidocaine consists of cross-linked hyaluronic acid (HA), obtained from *Streptococcus equi* bacteria, formulated to a concentration of 23 mg/mL in a physiologic buffer supplemented with 0.3% lidocaine hydrochloride. Each box contains one syringe of Saypha® FILLER Lidocaine, two 27G ½" disposable sterile needles reserved for injection of Saypha® FILLER Lidocaine and a product leaflet. The ancillary role of lidocaine hydrochloride is to reduce the patient's discomfort during treatment. A set of two labels showing the batch number is situated at the bottom of the box. One of these labels should be attached to the patient's file and the other should be given to the patient to ensure traceability.

**Indications:**

Saypha® FILLER Lidocaine is a viscoelastic solution to correct light to moderate nasolabial folds. It is indicated to be injected into the mid to deep dermis.

**Exclusion criteria:**

Do not inject Saypha® FILLER Lidocaine into the periorbital area (crow's feet, eye circle, eye lids) or glabellar region. Saypha® FILLER Lidocaine must not be used in:
 

- patients who tend to develop hypertrophic scarring. | - patients with a history of autoimmune disease or who are receiving immune therapy. | - patients who are known to be hypersensitive to hyaluronic acid, lidocaine hydrochloride or to amide-type local anaesthetics. | - pregnant or breastfeeding women. | - patients under 18 years of age.

 Patients on anticoagulants or patients receiving platelet aggregation inhibitors (e.g. ASS) should not be treated with Saypha® FILLER Lidocaine without consulting their doctors. Saypha® FILLER Lidocaine must not be used in areas presenting cutaneous, inflammatory and/or infectious processes (e.g. acne, herpes). Saypha® FILLER Lidocaine must not be used in association with laser therapy, chemical peeling or dermabrasion.

**Precautions for use:**

Saypha® FILLER Lidocaine is only indicated for intradermal injections. Sensitive skin may be pretreated using a local anaesthetic patch or cream. Please note that any anaesthesia may cause redness or local hypersensitivity. There are no available clinical data (efficacy, tolerance) about injecting Saypha® FILLER Lidocaine into an area which has already been treated with another filling product. In order to minimise the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection. Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications. Patients should be advised not to apply any make-up for 12 hours after the injection and to avoid prolonged exposure to sunlight and UV or using saunas or Turkish baths for one week after the injection. If the 27G ½" needle is blocked, do not increase the pressure on the plunger rod, but stop the injection and replace the needle. Do not inject into blood vessels. Do not use Saypha® FILLER Lidocaine in bones, tendons, ligaments or muscles. Do not inject Saypha® FILLER Lidocaine into naevi. Do not overcorrect. Discard the syringe and remaining product after use. Do not use if packaging has been damaged.

**Interactions with other agents:**

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds such as benzalkonium chloride solutions. Therefore Saypha® FILLER Lidocaine should never be placed in contact with these substances or with medical-surgical instruments that have been in contact with these substances.

**Undesired side effects:**

Physicians must inform the patient that there are potential side effects and/or incompatibilities associated with the implantation of this device, which may occur immediately or may be delayed. These include (non-exhaustive list):
 

- Inflammatory reactions (e.g. redness, oedema, erythema . . .) which may be associated with itching and pain on pressure after the injection. These reactions may last for a week. |
- Indurations or nodules at the injection site especially in cases of too superficial placement. In a clinical study of similar substances, local mobility of the injected gel and skin lumpiness were observed; this might be due to the use of too much volume and/or sub-optimal injection technique (e.g., large bolus injections). In some cases implant aspiration was used to resolve the event. | - Discolouration of the injection site. | - Poor effect or weak filling effect if Saypha® FILLER Lidocaine is injected incorrectly.

 Cases of skin necrosis, bacterial infection, formation of granulomas, cystic and acneiform lesions and hypersensitivity have all been reported in the literature following hyaluronic acid injections. Cases of visual loss as the most severe uncommon undesired events due to wrong injection technique have also been reported. Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction. Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. It is therefore important to take these possible complications without any delay into account and to communicate them to the patient before treatment. Patients must inform their physician as soon as possible about any inflammatory reactions which persist for more than one week or any other secondary effect which develops. The physician should treat these as appropriate. Any other undesirable side effects associated with the injection of Saypha® FILLER Lidocaine must be reported to the distributor and/or to the manufacturer.

**Methods of use:**

Saypha® FILLER Lidocaine must be injected into non-inflamed, disinfected, healthy skin. Inject the product slowly and apply the least amount of pressure necessary. The technique used is essential for the success of the treatment (e.g. microdroplet technique). Therefore this device must only be used by doctors who have received specific training on the injection technique for filling wrinkles. The area to be treated should be disinfected thoroughly prior to the injection. Use the 27G ½" needle which is provided with the syringe and inject slowly by applying the appropriate injection technique. The amount injected will depend on the wrinkle or depression which is to be corrected. A maximum volume of 2.0 mL and 1.0 mL per nasolabial fold is recommended for initial and touch-up treatments, respectively. After the injection, doctors may apply a light massage in order to distribute the product uniformly. However, patients should be advised not to massage the area by themselves for one week following the injection in order to avoid a possible risk of product mobility. The filling effect lasts up to 6 months based on clinical data gained with the device.

**Warnings:**

Verify the integrity of the syringe and the expiry date before use. Do not use any other needle or syringe than provided by the manufacturer. Do not re-use; quality and sterility can only be guaranteed for an originally closed syringe. The re-use of the product creates a potential infection risk for patients or users. Do not resterilize. After use the needle must be disposed of in a special container. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur. In order to reduce the risk of possible skin necrosis if intravascular injection has happened consider removal of the product by aspiration/ drainage, extrusion or enzymatic degradation (use of Hyaluronidase has been described in scientific publications).

**Storage:**

Store Saypha® FILLER Lidocaine at 2–25 °C / 36–77 °F in a dry place in the original box. Protect from light, heat and frost and handle with care.

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
## Saypha® FILLER Lidocaine

**EN Instructions for the correct handling of the device:**

Hold the Luer-Lock adapter as shown in ❶. To remove the tip cap, twist ❷ and pull carefully ❸. Following the instruction above will prevent the formation of air bubbles. Hold the syringe as shown in ❹. Open the enclosed needle case and insert the needle firmly ❺ (do not use any other needle). Tightly secure the needle by twisting it clockwise ❻.

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**CROMA-PHARMA GmbH**