Disconnecting / Desconectar / Débranchement / Trennen / Disconnessione /
 Loskoppelen / Frakobling / Irrotaminen / Koble fra / Fränkoppling /
 Odlącznie / Odpojení / Αποσύνδεση / Leválasztás / מתנתק / 연결 해제 / 断开
 连接 / الفصل / Desconectar / Odpojenie

Reconnecting / Reconectar / Branchement / Erneutes Verbinden / Riconnessione / Weer aansluiten / Tilkobling /
 Liittäminen uudelleen / Koble til igjen / Tilkoppling / Ponowne podłączenie / Opétovné připojení / Επανασύνδεση /
 Újracsatlakoztatás / מתחבר מחדש / 재연결 / 重新连接 / إعادة التوصيل / Voltar a conectar / Opätovné pripojenie

Contacts / Contactos / Contacts / Kontakt / Contatti /
 Contactinformatie / Kontakter / Yhteystiedot / Kontaktinformasjon /
 Kontakt / Dane kontaktowe / Obsah / Στοιχεία επικοινωνίας /
 Elérhetőség / אמני קשר / 연락처 / 联系人 / جهات الاتصال / Contactos /
 Kontakty

www.tandemdiabetes.com/contact

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
Indications for use

The AutoSoft 90™ infusion set is indicated for the subcutaneous infusion of insulin administered by an external Tandem Diabetes Care™ pump. The infusion set is indicated for single use. The infusion set is designed to be used for patients specified in the instructions for use accompanying the external pump.

INTENDED USE

The infusion set is intended for subcutaneous infusion of insulin in the treatment of diabetes mellitus. The infusion set has been tested to be compatible with insulins approved for subcutaneous infusion.

CLINICAL BENEFITS

- Continuous subcutaneous infusion of insulin from an external pump
- Fewer needle insertions compared with multiple daily injections

TARGET PATIENT POPULATION

The infusion set is designed to be used for patients specified in the instructions for use accompanying the external pump. Please visit the Tandem Diabetes Care website to find the user guide for your Tandem insulin pump.

CONTRAINDICATIONS

AutoSoft 90 is neither intended nor indicated for use in intravenous infusion of blood or blood products.

WARNINGS

- Carefully examine the instructions for use before inserting AutoSoft 90. Failure to follow instructions may result in pain or injury.
- AutoSoft 90 infusion sets with t:lock™ tubing connectors should only be used with Tandem cartridges featuring the t:lock connector.
- AutoSoft 90 is only sterile and nonpyrogenic if the insertion device is undamaged and unopened. Do not use if the insertion device is already opened or has been damaged. Ensure sterility by checking if the sterile paper and tamperproof seal are unbroken.
- When using AutoSoft 90 for the first time, do so in the presence of a healthcare provider. Since different individuals have varying amounts of subcutaneous layer, ensure that an appropriate cannula length is chosen to fit your needs.
- Inaccurate medication delivery, site irritation or infection may result from

improper insertion or maintenance of your infusion site.

- Do not re-insert the introducer needle into the infusion set. This could cause a tearing or damage of the soft cannula and may lead to unpredictable medication flow.
- Reuse of the infusion set may cause site irritation or infection and may lead to inaccurate medication delivery.
- Be sure that the needle guard is removed before insertion is made.
- If the soft cannula bends during insertion, apply a new AutoSoft 90 immediately at a new site.

PRECAUTIONS

- Never try to fill or free a clogged tubing while the infusion set is inserted. This may result in unpredictable medication flow.
- When unwinding, release the tubing with caution as a hard pull of the tubing can result in damage to the infusion device. Ensure that the infusion set is properly in place in the serter when the tubing is fully released (see figure 8 on page 88).
- Do not leave air in the infusion set. Make sure to fill the tubing completely. Consult pump manufacturer's instructions for use.
- AutoSoft 90 is a single-use device and

should be disposed of immediately after use. Do not clean or re-sterilize.

- Do not in any way put disinfectants, perfumes, deodorants, cosmetics or other products containing disinfectants in contact with the connector or the tubing. This may impact the integrity of the infusion set.
- Replace the infusion set if the adhesive tape becomes loose or is displaced from its original site. Since the cannula is soft, it will not cause any pain if it becomes dislodged, and this may take place without notice. The soft cannula must always be completely inserted to receive the full amount of medication.
- Replace the infusion set every two to three days, or per your healthcare provider's instructions.
- If your insertion site becomes irritated or infected, replace infusion set and use new site until the first site has healed.
- Wash your hands before disconnecting and reconnecting the infusion set. Consult your healthcare provider on how to compensate for missed medication when disconnected.
- Never point a loaded insertion device towards any body part, where insertion is

not desired.

PRODUCT DESCRIPTION

AutoSoft 90 is an infusion device which combines an infusion set with an aid for insertion (serter). As the infusion set is an integrated part of the serter, AutoSoft 90 is delivered assembled and ready for use as a single, sterile unit.

INTENDED USERS

Patients, caregivers and healthcare professionals.

USE ENVIRONMENT

Clinical settings (e.g. hospital and clinic) or home environment.

RECOMMENDATIONS

- Check your blood glucose level one to three hours after insertion of the AutoSoft 90. Measure your blood glucose on a regular basis. Discuss this with your healthcare provider.
- Do not change your infusion set just prior to bedtime, unless blood glucose can be checked one to two hours after insertion.
- If your blood glucose level becomes unexplainably high or an occlusion alarm occurs, check for clogs and leaks. If in

doubt, change your infusion set since the soft cannula could be dislodged or partially clogged. Discuss plan to bring down your blood glucose level with your healthcare provider should any of these problems arise. Test your blood glucose level to make sure the problem is corrected.

- Carefully monitor your blood glucose level when disconnected and after reconnecting.
- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

STORAGE & DISPOSAL

- Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight, or in high humidity.
- Store and handle insulin according to the manufacturer's instructions.
- Put the lid back on the insertion device and dispose in a proper sharp container according to local laws, to avoid risk of needle stick.
- Dispose of a used infusion set according to local regulations for biohazardous waste.

Preparation and insertion

Read all instructions carefully before using the AutoSoft 90 infusion set. For information about pump therapy, possible errors, and potential risks and their causes, consult your pump manufacturer’s instructions for use before connecting AutoSoft 90. When using AutoSoft 90, follow proper hygiene procedures. Please see the corresponding illustrations on page 88-90.

CONTENTS

- A. Lid
- B. Needle guard
- C. Soft cannula
- D. Adhesive tape
- E. Cannula housing
- F. Connector needle
- G. Tubing
- H. Introducer needle
- I. Insertion device
- J. Sterile paper
- K. Three raised dots – press here when lifting the lid off (see figure 6)
- L. Lined indentations – press here when loading AutoSoft 90 (see figure 10)

- M. Round indentations – press here simultaneously when inserting AutoSoft 90 (see figure 12)
- N. Disconnect cover (see figure B)

INSTRUCTIONS

Unfold the cover of this booklet to view images that correspond to the steps below.

- 1. Wash your hands before inserting AutoSoft 90.
- 2. Suggested insertion sites. Choose an insertion site as recommended by your healthcare provider. Do not use an area right next to your previous insertion site.
- 3. Clean the area of insertion with a disinfectant as directed by your healthcare provider. Be sure the area is air dry before inserting AutoSoft 90.
- 4. Pull the red tape to remove the seal.
- 5. Remove the sterile paper.
- 6. Press the three raised dots on each side of the lid firmly with one hand and lift up the lid with the other hand.

CAUTION: Be careful not to bend or touch the AutoSoft 90 introducer needle when preparing the serter and infusion set.
IMPORTANT NOTE: If the infusion set is not securely placed in the insertion device with the needle pointing straight ahead

prior to insertion, pain or minor injury may occur. To avoid this, please ensure that the following steps are taken gently so the infusion set is not loosened / removed by accident.

- 7. **Unwinding the tubing:** Carefully twist the beginning of the tubing out of the slot. Pull up gently. Now unwind the tubing from the insertion device by pulling it carefully upwards.
CAUTION: Do not touch the AutoSoft 90 introducer needle when unwinding the tubing. Be careful not to pull too hard when the tubing is fully released. You could accidentally remove the infusion set from the serter.
- 8. Connect AutoSoft 90 to a cartridge and the pump. Prime infusion set. Ensure there are no air bubbles in the tubing. Consult pump manufacturer’s instructions for use.
CAUTION: When priming the infusion set, hold AutoSoft 90 with the needle pointing down to ensure that medication does not come in contact with the adhesive backing paper.
- 9. Pull up gently to remove the adhesive backing paper.

- 10. **Prepare the inserter:** Place fingers on lined indentations. Press lined indentations on each side. Pull spring up until you hear a “click”.
- 11. Carefully remove the needle guard by gently twisting and then pulling it. Check that the soft cannula does not exceed the introducer needle.
CAUTION: Before positioning AutoSoft 90 over a clean insertion site, place the tubing in the slot to ensure that the tubing is not caught underneath the infusion set during insertion.
- 12. **Insert AutoSoft 90:** Position AutoSoft 90 over the insertion site. Simultaneously press the round indentations on both sides of the insertion device to insert AutoSoft 90.
- 13. Push gently on the center of the inserter to secure the adhesive tape on the skin.
- 14. Remove the serter/introducer needle by grasping the center and pulling it gently back. Gently massage the tape onto the skin. Fill the empty space in the cannula with the following amount of insulin:
U-100 INSULIN:
6 mm: 0.3 units (0.003 ml)
9 mm: 0.5 units (0.005 ml)

- 15. Push the lid back in place until you hear a “click.”

DISCONNECT

- A. AutoSoft 90 allows you to temporarily disconnect your pump without having to change the infusion set. Hold the adhesive in place by placing a finger just in front of the cannula housing, and gently squeeze the sides of connector. Pull the connector needle straight out from the cannula housing.
- B. Insert the disconnect cover into the cannula housing until you hear a “click”.

RECONNECT

- C. If necessary, prime the infusion set until medication dispenses from the needle.
CAUTION: When priming the infusion set, hold AutoSoft 90 with the needle pointing down.
- D. Remove the disconnect cover from the cannula housing by pressing the grips. Place a finger on the adhesive tape just in front of the cannula housing.
- E. Push the connector needle straight in until you hear a “click”.

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© 2025 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, Tandem logos, and AutoSoft are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. / Tous droits réservés. Tandem Diabetes Care, logos Tandem et AutoSoft sont des marques déposées ou des marques de commerce de Tandem Diabetes Care, Inc. aux États-Unis et/ou dans d'autres pays. AW-1000184_G