# WEB™

**ANEURYSM EMBOLIZATION SYSTEM**

**INSTRUCTIONS FOR USE**

---

<b>WEB™ ANEURYSM EMBOLIZATION SYSTEM</b>	
Instructions for Use .....	English ..... 1
<b>SYSTÈME D'EMBOUISATION D'ANÉVRISME WEB™</b>	
Mode d'emploi.....	French / Français ..... 5
<b>WEB™ ANEURYSMA-EMBOLISATIONSSYSTEM</b>	
Gebrauchsanweisung.....	German / Deutsch ..... 9
<b>SISTEMA DE EMBOLIZACIÓ DE ANEURISMA WEB™</b>	
Instrucciones de uso .....	Spanish / Español..... 13
<b>SISTEMA DI EMBOLIZZAZIONE DI ANEURISMI WEB™</b>	
Istruzioni per l'uso .....	Italian / Italiano ..... 17
<b>WEB™ SISTEMA DE EMBOLIZAÇÃO DE ANEURISMA</b>	
Instruções de utilização.....	Portuguese / Português..... 21
<b>WEB™-ANEURISMEEMBOLISERINGSSYSTEM</b>	
Brugsanvisning.....	Danish / Dansk ..... 25
<b>WEB™ ANEURYSMA-EMBOLISATIESYSTEEM</b>	
Gebbruksaanwijzing .....	Dutch / Nederlands..... 29
<b>ANEURYSMAN EMBOLISAATIOJÄRJESTELMÄ WEB™</b>	
Käyttöohjeet .....	Finnish / Suomi ..... 33
<b>WEB™ ANEURYSMEMBOLISERINGSSYSTEM</b>	
Bruksanvisning .....	Swedish / Svenska ..... 37
<b>WEB™-ANEURISMEEMBOLISERINGSSYSTEM</b>	
Bruksanvisning .....	Norwegian / Norsk..... 41
<b>ΣΥΣΤΗΜΑ ΕΜΒΟΛΙΣΜΟΥ ΑΝΕΥΡΥΣΜΑΤΟΣ WEB™</b>	
Οδηγίες χρήσης .....	Greek / Ελληνικά ..... 45
<b>WEB™ ANEVRIZMA EMBOLIZASYON SISTEM</b>	
Kullanma Talimatları.....	Turkish / Türkçe ..... 49
<b>WEB™ 動脈瘤栓塞系統</b>	
使用說明 .....	Traditional Chinese / 繁體中文..... 53
<b>WEB™ 동맥류 색전 폐색 시스템</b>	
사용 지침 .....	Korean / 한국어 ..... 56
<b>СИСТЕМА ЗА ЕМБОЛИЗАЦИЯ ПРИ АНЕВРИЗЪМ WEB™</b>	
Инструкции за употреба .....	Bulgarian / Български ..... 60

<b>SUSTAV ZA EMBOLIZACIJU ANEURIZME WEB™</b>	
Upute za uporabu.....	Croatian / Hrvatski..... 64
<b>SYSTÉM PRO EMBOLIZACI ANEURYZMAT WEB™</b>	
Návod k použití.....	Czech / Česky ..... 68
<b>ANEURÛSMI EMBOLISEERIMISSÛSTEEM WEB™</b>	
Kasutusjuhend.....	Estonian / Eesti keel ..... 72
<b>WEB™ ANEURIZMAEMBOLIZÁLÓ RENDSZER</b>	
Használati utasítás.....	Hungarian / Magyar..... 76
<b>WEB™ ANEIRISMAS EMBOLIZĀCIJAS SISTĒMA</b>	
Lietošanas instrukcija.....	Latvian / Latviski..... 80
<b>WEB™ ANEURIZMOS EMBOLIZACIJOS SISTEMA</b>	
Naudojimo instrukcijos.....	Lithuanian / Lietuvių k..... 84
<b>SYSTEM DO EMBOLIZACJI TĘTNIAKÓW WEB™</b>	
Instrukcja obsługi.....	Polish / Polski..... 88
<b>SISTEM PENTRU EMBOLIZAREA ANEVRISMELOR WEB™</b>	
Instrucțiuni de utilizare.....	Romanian / Română..... 92
<b>СИСТЕМА ЭМБОЛИЗАЦИИ АНЕВРИЗМЫ WEB™</b>	
Инструкции по применению.....	Russian / Русский ..... 96
<b>WEB™ SISTEM ZA EMBOLIZACIJU ANEURIZME</b>	
Uputstvo za upotrebu.....	Serbian / Srpski..... 100
<b>SISTEM ZA EMBOLIZACIJO ANEVRIZME WEB™</b>	
Navodila za uporabo.....	Slovenian / Slovenščina..... 104
<b>WEB™ SYSTÉM NA EMBOLIZÁCIU ANEURYZMY</b>	
Návod na použitie.....	Slovak / Slovenčina..... 108
<b>СИСТЕМА ЕМБОЛІЗАЦІЇ АНЕВРИЗМИ WEB™</b>	
Інструкції із застосування.....	Ukrainian / Українська..... 112
<b>نظام الانصمام لتمدد الأوعية الدموية WEB™</b>	
تعليمات الاستخدام.....	Arabic / العربية..... 116
<b>СИСТЕМ ЗА ЕМБОЛИЗАЦИЈА НА АНЕВРИЗМИ WEB™</b>	
Упатство за употреба.....	Macedonian / Македонски..... 120
<b>SISTEMA DE EMBOLIZAÇÃO DE ANEURISMA WEB™</b>	
Instruções de uso.....	Brazilian Portuguese / Português do Brasil..... 124

# English

## WEB™ Aneurysm Embolization System

### Instructions for Use

#### DEVICE DESCRIPTION

The WEB Aneurysm Embolization System consists of an implantable embolization device attached to a delivery device. The delivery device is powered by a hand-held, battery-powered detachment control device designed specifically for the WEB Aneurysm Embolization System. The detachment control device is provided separately.

The WEB embolization device is manufactured from nitinol wires in a braided, self-expanding mesh configuration. The WEB embolization devices are provided in a broad range of sizes (diameters and lengths) to satisfy the needs of the physician. During treatment, the physician selects the appropriate device size based on the size, shape and location of the intracranial aneurysm or other vascular abnormality to be occluded. The WEB embolization device is delivered to the treatment site on the delivery device through standard neuro-interventional wire-reinforced microcatheters with a specified minimum inner diameter (see Table 1 below). An introducer sheath on the outside of the delivery device assists in the placement of the system into the microcatheter.

The WEB device that remains in the patient following detachment consists of:

Table 1 – Quantitative Implant Material Information

Implant Materials		Mass (mg)*
Metallic Components	Nitinol, Platinum, Platinum-Iridium Alloy	≤ 40
Non-metallic Components	Epoxy, PET	≤ 0.1

\*Approximate content

#### INTENDED PURPOSE

The WEB Aneurysm Embolization System is intended for the endovascular embolization of ruptured and unruptured intracranial aneurysms and other neurovascular abnormalities such as arteriovenous fistulae (AVF).

The WEB Aneurysm Embolization System is also intended for vascular occlusion of blood vessels within the neurovascular system to permanently obstruct blood flow to an aneurysm or other vascular malformation.

#### POTENTIAL COMPLICATIONS

Potential complications include but are not limited to the following: hematoma at the site of entry, aneurysm rupture, emboli, vessel perforation, parent artery occlusion, hemorrhage, ischemia, vasospasm, clot formation, device migration or misplacement, premature or difficult device detachment, non-detachment, incomplete aneurysm filling, revascularization, post-embolization syndrome, and neurological deficits including stroke and death.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State or Local Health Authority in which the user and/or patient is established.

#### REQUIRED ADDITIONAL ITEMS

- WEB Aneurysm Embolization System detachment control device
- Wire-reinforced microcatheter with distal tip RO marker (see Table 2)
- Guide catheter compatible with microcatheter
- Steerable guidewire compatible with microcatheter
- Two rotating hemostatic Y valves (RHV)
- One three-way stopcock
- One one-way stopcock
- Sterile saline
- Pressurized sterile saline drip

Table 2 – Microcatheter Sizes

WEB Embolization Device (Diameter) Range	Minimum Microcatheter Inner Diameter (inch)	Recommended VIA™ Microcatheter <sup>1</sup>
W2 – WEB SINGLE 8 – 9 mm	0.027	VIA 27
W2 – WEB SINGLE 10 – 11 mm	0.033	VIA 33
W4 – WEB SINGLE 4 – 7 mm	0.021	VIA 21
W5 – WEB SINGLE 3 – 7 mm	0.017	VIA 17


<sup>1</sup> Use of a different catheter may result in extreme friction and damage to the device

#### WARNINGS AND PRECAUTIONS

- **CAUTION:** This device should be used only by physicians trained in percutaneous, intravascular and neurovascular techniques and procedures at medical facilities with the appropriate fluoroscopic equipment.
- **CAUTION:** The WEB embolization device should be used by physicians who have


received appropriate training for this device.

- **CAUTION:** Using this device in a catheter that is not recommended or required may result in extreme friction and damage to the device.
- The WEB Aneurysm Embolization System (see Figure 1) is provided sterile and non-pyrogenic unless the unit package is opened or damaged. Do not use if the packaging is damaged. Use before expiration date noted on the product packaging.
- The WEB Aneurysm Embolization System is intended for single use only. The detachment control device is intended to be used for one patient. Do not resterilize and/or reuse the device. Reuse and/or resterilization can increase risk of infection, cause a pyrogenic response or other life threatening complications. Reuse and/or resterilization can degrade product performance, leading to device malfunction. Dispose of all devices in accordance with applicable hospital, administrative and/or local government policy.
- The WEB embolization device must be delivered only through a compatible microcatheter with a PTFE inner surface coating. Damage to the embolization and delivery device may occur and necessitate removal of both the device and microcatheter from the patient.
- The operator should be aware that ≥0.021" microcatheters, in distal blood vessels, may increase the risk of thromboembolism.
- Steam shaping 0.021" and greater microcatheters may result in improper WEB embolization device delivery and deployment, depending on the degree of shaping and catheter deflection during WEB embolization device delivery.
- High quality, digital subtraction fluoroscopic road mapping, with orthogonal views, is mandatory to achieve correct placement of the embolization device.
- Advance and retract the device slowly. Do not advance the delivery device with excessive force. Determine the cause of any unusual resistance. Remove the device if excessive friction is noted and check for damage.
- If repositioning is required, take special care to retract or to advance the device under fluoroscopy, including new road map to confirm catheter position.
- Do not rotate the delivery device during or after delivery of the embolization device. Rotating the device may result in damage or premature detachment.
- If an embolization device must be retrieved from the vasculature after detachment, retrieval devices (e.g. alligator and snare) should be used per their manufacturer's instructions.
- The WEB embolization device foreshortens during delivery (~60%) (e.g. see Figure 2a, a properly deployed 11mm wide x 9mm long device will measure ~20 mm long in a 0.032" microcatheter).
- When properly deployed, the two radio-opaque markers should be separated and fluoroscopically visible (e.g. see Figure 2b, depending on working projection and placement in the aneurysm, the distance between the proximal to distal marker should approximate the labeled WEB embolization device length).
- WEB embolization device visibility may vary with diameter; larger sizes may be more visible than smaller sizes. Examples are shown in Figure 2c.
- The pictures in (a) through (c) below illustrate WEB embolization device deployment. Initially, the distal implant marker band exits the microcatheter (a). As the implant is advanced, it begins to expand in diameter (b). When the distance between the catheter marker band and implant tip is about 1/3 of the total implant marker-marker distance, the implant diameter is generally about 1/2 of its fully deployed diameter (b). When the implant distal marker band to catheter distal marker band distance is about 2/3 of the total implant marker-marker distance, the implant has reached about 4/5 of its fully deployed diameter and the distal marker band begins moving into the distal recess (c).


#### Note:

- WEB embolization devices are available in both wide neck and spherical shapes.
- VIA 17 Microcatheters have a proximal marker band not shown in the drawings or photos below. This proximal catheter marker band is not used for WEB embolization device delivery.


- If the radio-opaque markers are clustered (i.e. a shorter distance between markers than expected), retract WEB embolization device into the microcatheter and evaluate the microcatheter/aneurysm position with multiple fluoroscopic angles.
- The embolization device cannot be detached with any other power source other than a MicroVent Inc. detachment control device. Ensure that at least two detachment control devices are available before initiating an embolization procedure.
- Batteries are pre-loaded into the detachment control device. Do not attempt to remove or replace the batteries.
- Do not use in conjunction with radio frequency (RF) devices.
- Patients who are allergic to nickel may have an allergic reaction to this device.

**Figure 1: WEB Aneurysm Embolization System (not to scale)**


**Figure 2: WEB embolization device in catheter partially deployed foreshortened (a) and properly fully deployed (b)**


**(a) WEB embolization device partially in catheter, partially deployed in aneurysm – foreshortened view**

**(b) WEB embolization device fully deployed in aneurysm**

**(c) WEB Embolization Device Visibility**

## PROCEDURE

### Catheterization of the Lesion

1. Using standard interventional procedures, access the vessel with a guide catheter. The guide catheter should have an inner diameter large enough to allow for contrast injection while the microcatheter is in place.
2. Attach a rotating hemostatic valve (RHV) to the hub of the guide catheter. Attach a three-way stopcock to the side arm of the RHV and then connect a line for continuous infusion of flush solution.
3. Select a microcatheter with the appropriate inner diameter (see Table 1).
4. After the microcatheter has been positioned inside the lesion, remove the guidewire.
5. Attach a second RHV to the hub of the microcatheter. Attach a one-way stopcock to the sidearm of the second RHV and connect the flush solution line to the stopcock.
6. Open the stopcock to allow flush through the microcatheter with sterile flush solution.

### Device Size Selection

7. Perform fluoroscopic road mapping.
8. Measure and estimate the size of the lesion to be treated.
9. Select the appropriately sized embolization device. The size should be chosen based upon angiographic assessment of the aneurysm diameter, height, neck and shape.

### Preparation of the Device for Delivery

10. Remove the detachment control device from its protective packaging and place it within the sterile field. Do not use any power source from another medical device manufacturer to detach the embolization device.
11. Remove the delivery device from the packaging hoop by pulling the proximal end until the introducer exits the hoop.
12. Slowly advance the embolization device out of the introducer sheath and inspect for any damage. Do not use the device if there is any damage.
13. Firmly insert the proximal end of the delivery device into the distal end of the detachment control device. Do not press the detachment button. If the light flashes green and the device beeps, the detachment control device is good to use. If not, replace the detachment control device and repeat this step with a new detachment control device. If the light flashes green and beeps with this second controller, proceed to the next step. Otherwise, obtain a new WEB embolization device delivery system.
14. While holding the introducer sheath vertically, gently retract the embolization device back into the introducer sheath.

### Introduction and Deployment of the Device

15. Open the RHV on the microcatheter to accept the introducer sheath.
16. Insert the introducer sheath through the RHV. Seat the distal tip of the introducer sheath at the distal end of the microcatheter hub and close the RHV lightly around the introducer to secure the RHV to the introducer. Do not over-tighten the RHV. Only insert the introducer into the microcatheter hub until a slight resistance is felt. Do not over insert.
17. Push the device into the lumen of the microcatheter. Use caution to avoid catching the embolization device on the junction between the introducer sheath and hub of the microcatheter.
18. Push the device through the microcatheter until the proximal end of the delivery device meets the proximal end of the introducer sheath.
19. Loosen the RHV.
20. Retract the introducer sheath just out of the RHV.
21. Close the RHV around the delivery device.
22. Slide the introducer sheath completely off the delivery device using care not to kink or damage the delivery system.
23. Carefully advance the device until the distal end of the device reaches the last marker on the microcatheter.
24. Reposition the tip of the microcatheter so that it sits just at the neck of the aneurysm. Do not insert the microcatheter completely inside the aneurysm.
25. Under fluoroscopic guidance, slowly advance the embolization device out the tip of the microcatheter. Continue to advance the embolization device into the lesion until optimal deployment is achieved. The following may require repositioning or removal of the embolization device and/or repositioning of the microcatheter:
  - a. If the embolization device size is not appropriate, remove and replace with another device.
  - b. If undesirable movement of the embolization device is noted following placement and prior to detachment, remove the device. Movement of the embolization device may indicate that the device could migrate once it is detached.
  - c. If the embolization device (implant) does not fully open:
    - i. Retrieve the implant, reposition the microcatheter further proximally and re-deploy the implant to allow more room for expansion; **or**
    - ii. Replace the implant with another implant of the same or alternate size.The embolization device should not be retracted and deployed more than twice. After two attempts, remove the embolization device and replace with another device of the same or alternate size. To minimize the potential risk of emboli, **DO NOT** allow an inappropriately sized or non-optimally positioned device to reside in the aneurysm significantly beyond the activated clotting time (ACT). Experience has shown that blood/clot can also prevent the WEB embolization device from full deployment and recapture. To minimize the risks of potential complications, the status of the patient's anti-platelet medication regimen should be considered when making a decision to remove the entire system from the aneurysm for replacement by a new device.
26. Angiographic assessment should always be performed prior to detachment to ensure that the embolization device is not markedly protruding into the parent vessel.
27. The WEB embolization device should be placed with the proximal surface (center of proximal marker) aligned with the aneurysm neck and the proximal marker extending beyond the neck.
28. If redeploying, retract device so that there is no contact with the aneurysm before attempting redeployment.
29. Tighten the RHV to prevent movement of the embolization device.

30. Verify that the distal portion of the delivery device is not under tension or compression prior to detachment. This could cause the microcatheter tip to move resulting in aneurysm or vessel rupture.

### Detachment of the Device

31. The detachment control device is pre-loaded with batteries and will activate when the delivery device is properly connected.
32. Verify that the RHV is firmly locked around the delivery device before attaching the detachment control device to ensure that the embolization device does not move during the connection process.
33. Ensure that the delivery device gold connectors are clean and free from blood or contrast. If necessary, wipe the connectors with sterile water and dry before connecting.
34. Insert the proximal end of the delivery device into the detachment control device. When the delivery device is properly connected, the light will flash green and an intermittent tone will be heard.
35. Verify the embolization device position before pressing the detachment button.
36. Push the detachment button. During firing, the light should be solid green and the beep should be continuous.
37. Verify detachment by first loosening the RHV valve, then pulling back slowly on the delivery device and verifying that there is not embolization device movement. If the embolization device does not detach, push the detachment button again. If the device is still not detached, obtain a new detachment control device and attempt detachment up to two additional times. If it does not detach, remove the delivery device.
38. Verify the position of the embolization device angiographically through the guide catheter.
39. Prior to removing the microcatheter from the treatment site, place an appropriately sized guidewire completely through the microcatheter lumen to ensure that no part of the embolization device remains within the microcatheter.

The physician has the discretion to modify the device deployment technique based on the complexity and variation in embolization procedures. Any modifications must be consistent with the previously described procedures, warnings, precautions and patient safety information in these instructions for use.

### SPECIFICATIONS FOR DETACHMENT CONTROL DEVICE

- Model number: FG00175-001
- Output voltage: 11.2-11.8 VDC
- Input Voltage: 24 VDC
- Batteries: 4 each A23
- Type of applied part: BF
- Equipment not suitable for use in the presence of flammable mixtures
- The detachment control device is a single use device, preloaded with batteries and packaged sterile. No cleaning, inspection or maintenance is required. The detachment control device should not be cleaned, re-sterilized or re-used.
- Batteries are pre-loaded in the detachment control device. Do not attempt to remove or replace the batteries prior to use.
- If the detachment control device does not perform as described in the Detachment section of these instructions for use, discard the detachment control device and replace it with a new unit.

### PACKAGING, STERILIZATION, STORAGE AND OPERATING CONDITIONS

The WEB embolization and delivery device are placed inside a protective dispenser hoop and packaged in a pouch and unit carton. The detachment control device is packaged separately in a protective pouch and unit carton. The devices will remain sterile unless the package is opened, damaged or the expiration date has passed. If the sterile packaging is unintentionally opened or damaged discard the device.

The WEB embolization device and delivery device are sterilized by gamma radiation. A small round indicator label has been affixed to the packaging of the WEB embolization device and delivery device. This indicator turns from yellow to red upon exposure to radiation sterilization and must be red in order to use the device. If the indicator is yellow, do not use the WEB embolization device.

The detachment control device is sterilized by ethylene oxide. A small round indicator label has been affixed to the packaging of the detachment control device. This indicator turns from purple to green upon ethylene oxide sterilization and must be green in order to use the device. If the indicator is purple, do not use the device. Store at a controlled room temperature in a dry place.

The device should be used at a temperature of 20° – 23°C and a relative humidity of 30 – 60%. Atmospheric pressure variations do not impact device functionality.

### MATERIALS

The WEB Aneurysm Embolization System does not contain latex or PVC materials.

### WARRANTY

MicroVention Inc. warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness. Handling, storage, cleaning and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedures and other matters

beyond MicroVention's control directly affect the device and this warranty is limited to the repair and replacement of this device, through its expiration date. MicroVention Inc. shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. MicroVention Inc. neither assumes, nor authorizes any other person to assume or for it, any other or additional liability or responsibility in connection with this device. MicroVention, Inc. assumes no liability with respect to device reused, reprocessed or resterilized and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.

Prices, specification and model availability are subject to change without notice.

## **MRI INFORMATION**

### **MR Conditional**

The WEB Aneurysm Embolization Device was determined to be MR-conditional.

Non-clinical testing demonstrated that the WEB Aneurysm Embolization Device is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

### **Static Magnetic Field**

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 4,000-Gauss/cm (extrapolated) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

### **MRI-Related Heating**

Under the scan conditions defined above, the WEB Aneurysm Embolization Device is expected to produce a maximum temperature rise of +1.4 °C after 15 minutes of continuous scanning.

### **Artifact Information**

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the WEB Aneurysm Embolization Device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. In non-clinical testing, the artifact size (as imaged with a gradient echo pulse sequence and 3 Tesla MRI system) extends approximately 5-mm relative to the size and shape of the WEB Aneurysm Embolization Device.

## **SUMMARY OF SAFETY AND CLINICAL PERFORMANCE**

The Summary of Safety and Clinical Performance (SSCP) for the device will be accessible in the European database on medical devices (Eudamed: <https://ec.europa.eu/tools/eudamed>), when available.

Permanent implant. Follow-up required at the discretion of the physician

eIFU website: [www.microvention.com/eIFU-MicroVention.com](http://www.microvention.com/eIFU-MicroVention.com)

© Copyright 2022 MicroVention, Inc. All rights reserved.

WEB™ and VIA™ are registered trademarks in the United States.

MicroVention™ is a trademark of MicroVention, Inc., registered in the United States and other jurisdictions.

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© Copyright 2022 MicroVention, Inc. Tous droits réservés.
© Copyright 2022 MicroVention, Inc. Alle Rechte vorbehalten.
© Copyright 2022 MicroVention, Inc. Todos los derechos reservados.
© Copyright 2022 MicroVention, Inc. Tutti i diritti riservati.
© Copyright 2022 MicroVention, Inc. Todos os direitos reservados.
© Copyright 2022 MicroVention, Inc. Alle rettigheder forbeholdes.