SOFIA™ Flow 88 Aspiration Catheter

INSTRUCTIONS FOR USE

SOFIA™ Flow 88 Aspiration Catheter	
Instructions for Use	English 1
Cathéter d'aspiration SOFIA™ Flow 88	
Mode d'emploi.....	French / Français 4
SOFIA™ Flow 88 Aspirationskatheter	
Gebrauchsanweisung	German / Deutsch..... 7
Catéter de aspiración SOFIA™ Flow 88	
Instrucciones de uso	Spanish / Español..... 10
Catetere di aspirazione SOFIA™ Flow 88	
Istruzioni per l'uso.....	Italian / Italiano..... 13
Cateter de aspiração SOFIA™ Flow 88	
Instruções de utilização	Portuguese / Português..... 16
SOFIA™ Flow 88-aspirationskateter	
Brugsanvisninger	Danish / Dansk..... 19
SOFIA™ Flow 88-aspiratiekatheter	
Gebruiksaanwijzing.....	Dutch / Nederlands..... 22
SOFIA™ Flow 88 -imukatetri	
Käyttöohjeet.....	Finnish / Suomi 25
SOFIA™ Flow 88 aspirationskateter	
Bruksanvisning	Swedish / Svenska 28
SOFIA™ Flow 88 aspirasjonskateter	
Bruksanvisning	Norwegian / Norsk 31
Καθετήρας αναρρόφησης SOFIA™ Flow 88	
Οδηγίες χρήσης	Greek / Ελληνικά 34
SOFIA™ Flow 88 Aspirasyon Kateterii	
Kullanım Talimatları	Turkish / Türkçe..... 37
SOFIA™ Flow 88 抽吸導管	
使用說明.....	Traditional Chinese / 繁體中文 40
SOFIA™ Flow 88 흡인 카테터	
사용 안내서	Korean / 한국어 43
Аспирационен катетър SOFIA™ Flow 88	
Инструкции за употреба.....	Bulgarian / Български 46
Aspiracijski kateter SOFIA™ Flow 88	
Upute za upotrebu	Croatian / Hrvatski 49
Aspirační katétr SOFIA™ Flow 88	
Návod k použití.....	Czech / Český..... 52
SOFIA™ Flow 88 aspiratsioonikateeter	
Kasutusjuhend	Estonian / Eesti keel..... 55
SOFIA™ Flow 88 aspirációs katéter	
Használati útmutató	Hungarian / Magyar 58
SOFIA™ Flow 88 aspirācijas katetrs	
Lietošanas pamācība.....	Latvian / Latviski 61
Siurbiamasis kateteris „SOFIA™ Flow 88“.	
Naudojimo instrukcija	Lithuanian / Lietuvių k..... 64
Cewnik aspiracyjny SOFIA™ Flow 88	
Instrukcja użycia	Polish / Polski 67
Cateter de aspirație SOFIA™ Flow 88	
Instrucțiuni de utilizare.....	Romanian / Română 70
Аспирационный катетер SOFIA™ Flow 88	
Инструкция по применению ...	Russian / Русский..... 73
SOFIA™ Flow 88 aspiracioni kateter	
Uputstvo za upotrebu	Serbian / Srpski 76
Aspiracijski kateter SOFIA™ Flow 88	
Navodila za uporabo.....	Slovenian / Slovenščina..... 79
Aspiračný katéter SOFIA™ Flow 88	
Návod na použitie.....	Slovak / Slovenčina 82
Аспіраційний катетер SOFIA™ Flow 88	
Інструкція з використання.....	Ukrainian / Українська..... 85
SOFIA™ Flow 88 قسطرة الشفط	
تعليمات الاستخدام	Arabic / العربية 88
Катетер за аспирација SOFIA™ Flow 88	
Упатство за употреба	Macedonian / Македонски..... 91
Cateter de aspiração SOFIA™ Flow 88	
Instruções de uso	Brazilian Portuguese / Português do Brasil . 94

English
SOFIA™ Flow 88 Aspiration Catheter
Instructions for Use

Carefully read all instructions prior to use.

DEVICE DESCRIPTION

The SOFIA Flow 88 Aspiration Catheter is a non-tapered, single-lumen, flexible catheter equipped with coil and braid reinforcement. The distal segment is designed to facilitate vessel selection with 60 cm of distal-shaft hydrophilic coating for navigation through the vasculatures. The radiopaque marker is located at the distal end of the catheter for visualization under fluoroscopy.

CONTENTS

One Catheter
One Introducer Sheath

INTENDED PURPOSE

The SOFIA Flow 88 Aspiration Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA Flow 88 Aspiration Catheter can be used to facilitate introduction of diagnostic or therapeutic agents. The SOFIA Flow 88 Aspiration Catheter is not intended for use in coronary arteries. Moreover, the SOFIA Flow 88 Aspiration Catheter is intended for use in removal/aspiration of emboli and thrombi from selected blood vessels in the arterial system, including the peripheral and neuro vasculatures.

CONTRAINDICATIONS

There are no known contraindications.

CAUTION

Rx Only: Federal (USA) law restricts this device to sale by or on the order of a physician.

Do not use if pouch is opened or damaged.

This device is intended for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose in accordance with hospital, administrative and/or local government policy.

WARNINGS

The SOFIA Flow 88 Aspiration Catheter should only be used by physicians who have received appropriate training in interventional techniques.

The SOFIA Flow 88 Aspiration Catheter is provided sterile and non-pyrogenic. Do not use if the packaging is breached or damaged.

Inspect the SOFIA Flow 88 Aspiration Catheter prior to use. Do not use the device if any damages or irregularities are observed.

Appropriate anti-coagulation and anti-platelet therapy should be administered per standard medical practice.

The SOFIA Flow 88 Aspiration Catheter should be manipulated under fluoroscopic guidance. Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.

Do not use the SOFIA Flow 88 Aspiration Catheter with Ethiodol or Lipiodol contrast media or other such contrast media which includes the components of those agents.

Do not use organic solvents as the device may be damaged.

Do not use automated high-pressure contrast injection equipment with the SOFIA Flow 88 Aspiration Catheter as it may damage the device.

Do not exceed 2070 kPa (300 psi) maximum recommended infusion pressure. Excess pressure may damage the device or injure the patient.

Torquing the SOFIA Flow 88 Aspiration Catheter excessively while kinked may damage the device, resulting in separation of the device. Withdraw the entire device (the device, Compatible Catheter, and guidewire) if the device is severely kinked.

The Introducer Sheath is not intended for use inside the patient body. Ensure that the Introducer Sheath is removed from the SOFIA Flow 88 Aspiration Catheter once the distal shaft of the SOFIA Flow 88 Aspiration Catheter is placed inside the patient body.

Excessive aspiration with the distal tip of the SOFIA Flow 88 Aspiration Catheter covered by the vessel wall may cause vessel injury. Carefully investigate location of the distal tip under fluoroscopy prior to aspiration.

Do not attempt to clear the inner lumen of the SOFIA Flow 88 Aspiration Catheter by infusion while keeping the device in the patient body.

When flow from the lumen stops or becomes stagnant during aspiration, do not attempt to clear the inner lumen of the SOFIA Flow 88 Aspiration Catheter by infusion while keeping the device in the patient body.

PRECAUTIONS

Exercise care in handling the SOFIA Flow 88 Aspiration Catheter to reduce the chance of accidental damage.

Verify compatibility of the SOFIA Flow 88 Aspiration Catheter when using other ancillary devices commonly used in intravascular procedures. The physician must be familiar with percutaneous, intravascular techniques and possible complications associated with the procedure.

Use caution when manipulating the SOFIA Flow 88 Aspiration Catheter in tortuous vasculature to avoid damage. Avoid advancing or withdrawal against resistance until the cause of resistance is determined.

Presence of calcifications, irregularities, or other devices may damage the SOFIA Flow 88 Aspiration Catheter and potentially affect its insertion or removal.

Maintain perfusion of heparinized saline for inner lumen of the SOFIA Flow 88 Aspiration Catheter to prevent thrombus formation.

If removed from the patient, the hydrophilic coating on the SOFIA Flow 88 Aspiration Catheter should be hydrated with heparinized saline. Do not allow the coating to dry.

POTENTIAL COMPLICATIONS

Potential complications include but are not limited to: vessel or aneurysm perforation, vasospasm, hematoma at the site of entry, embolism, ischemia, intracerebral/intracranial hemorrhage, pseudo aneurysm, seizure, stroke, infection, vessel dissection, thrombus formation, and death.


Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State or Local Health Authority in which the user and/or patient is established.

COMPATIBILITY

Refer to product label for device dimensions. Use the information on labeling provided with other devices to determine device compatibility. When using the SOFIA Flow 88 Aspiration Catheter as a single guiding catheter, choose the appropriate size of a femoral sheath, referring to product label.

Guide Sheath Recommendation*	≥ 0.113 in. ID
Catheter/Microcatheter Recommendation	≤ 0.0835 in. OD
Guidewire Recommendation	≤ 0.038 in. OD

* Guide sheath/catheter compatibility recommendation applies only to MicroVention/Terumo brand catheter/sheaths.

PREPARATION FOR USE

- Carefully remove the SOFIA Flow 88 Aspiration Catheter and the Introducer Sheath from the package.
- Inspect the SOFIA Flow 88 Aspiration Catheter for any damage.
WARNING: Do not use the device if any damages or irregularities are observed.
- Flush the lumen of the SOFIA Flow 88 Aspiration Catheter with heparinized saline. Attach a rotating hemostatic valve (RHV) to the proximal hub of the SOFIA Flow 88 Aspiration Catheter. Set up the line for perfusion of heparinized saline through the sidearm of the RHV.
- Hydrate the hydrophilic coating on the SOFIA Flow 88 Aspiration Catheter with heparinized saline before use. Keep the coating hydrated and do not allow the coating to dry.

DELIVERY OF THE SOFIA FLOW 88 ASPIRATION CATHETER

- Go to step 6 or 7, depending on the situation described below and choose appropriate devices for navigation of the SOFIA Flow 88 Aspiration Catheter. Refer to the Compatibility section for Compatible devices.
- Navigation through the vasculature, except for the intracranial vasculature**
 - Prepare Compatible Catheter/Microcatheter and 0.035" or 0.038" Guidewire for navigation of the SOFIA Flow 88 Aspiration Catheter.
 - Insert the Guidewire into the Compatible Catheter. Insert the Compatible Catheter with Guidewire into the SOFIA Flow 88 Aspiration Catheter and advance the Compatible Catheter with Guidewire until the Compatible Catheter and the SOFIA Flow 88 Aspiration Catheter are aligned at the distal end.

- c. Carefully insert the SOFIA Flow 88 Aspiration Catheter and the Compatible Catheter with Guidewire through a hemostatic valve of the femoral sheath. If necessary, use the Introducer Sheath provided in the package to facilitate with the introduction.
 - d. Remove the Introducer Sheath, if used, from the SOFIA Flow 88 Aspiration Catheter once the distal shaft of the SOFIA Flow 88 Aspiration Catheter is placed inside the patient body.
WARNING: Introducer Sheath is not intended for use inside the patient body.
 - e. Under fluoroscopic guidance, advance or withdraw the SOFIA Flow 88 Aspiration Catheter over the Compatible Catheter with Guidewire until desired position is attained or before the intracranial position is achieved.
WARNING: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
WARNING: Torquing the SOFIA Flow 88 Aspiration Catheter excessively while kinked may damage the device resulting in separation of the device. Withdraw the entire device (the device, Compatible Catheter, and Guidewire) if the device is severely kinked.
WARNING: Do not use automated high-pressure contrast injection equipment with the SOFIA Flow 88 Aspiration Catheter as it may damage the device. Do not exceed 2070 kPa (300 psi) maximum recommended infusion pressure. Excess pressure may damage the device or injure the patient.
 - f. Go to step 7 for navigation through the intracranial vasculatures. Otherwise, proceed to step 8.
- 7. Navigation through the intracranial vasculature**
- a. Prepare a Compatible Catheter/Microcatheter and Compatible Guidewire for navigation of the SOFIA Flow 88 Aspiration Catheter.
 - b. Slowly remove, if any, devices previously inserted in the SOFIA Flow 88 Aspiration Catheter. Insert the Compatible Catheter with Guidewire into the SOFIA Flow 88 Aspiration Catheter.
 - c. Under fluoroscopic guidance, advance or withdraw the SOFIA Flow 88 Aspiration Catheter over the Compatible Catheter and Guidewire until desired position is attained.
WARNING: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
WARNING: Torquing the SOFIA Flow 88 Aspiration Catheter excessively while kinked may damage the device resulting in separation of the device. Withdraw the entire device (the device, Compatible Catheter, and Guidewire) if the device is severely kinked.
WARNING: Do not use automated high-pressure contrast injection equipment with the SOFIA Flow 88 Aspiration Catheter as it may damage the device. Do not exceed 2070 kPa (300 psi) maximum recommended infusion pressure. Excess pressure may damage the device or injure the patient.
8. Slowly remove the Compatible Catheter, or Guidewire if necessary. Make sure that continuous perfusion of heparinized saline is maintained through the sidearm of the RHV.
NOTE: The Compatible Catheter used to navigate the SOFIA Flow 88 Aspiration Catheter may be kept for the rest of procedure.

ASPIRATION THROUGH THE SOFIA FLOW 88 ASPIRATION CATHETER WITH SYRINGE

1. Under fluoroscopic guidance, make sure that the distal tip of the SOFIA Flow 88 Aspiration Catheter is located in a straight portion of the vessel. Reposition the distal tip if necessary. Refer to Figure-1 for location of the distal tip.
WARNING: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
WARNING: Torquing the SOFIA Flow 88 Aspiration Catheter excessively while kinked may damage the device resulting in separation of the device. Withdraw the entire device (the device, Compatible Catheter, and Guidewire) if the device is severely kinked.

Figure-1

2. Make sure the distal tip of the SOFIA Flow 88 Aspiration Catheter is engaged with the emboli or thrombi under fluoroscopy.
 3. At the proximal end of the SOFIA Flow 88 Aspiration Catheter, remove the perfusion line and attach a stopcock with 30 cc or 60 cc locking syringe to the sidearm of the RHV.
 4. Tighten the RHV.
 5. With the stopcock closed pull and lock the plunger of the syringe.
 6. To begin aspiration, open the stopcock.
 7. Make sure that the syringe aspirates blood, emboli, or thrombi through the system. Close the stopcock if the syringe does not aspirate any blood, emboli, or thrombi, or slow aspiration is observed. Carefully investigate the cause of the restriction and reposition the distal tip if necessary. If blood, emboli, or thrombi is stuck in the SOFIA Flow 88 Aspiration Catheter remove the entire device (the device, Compatible Catheter, and Guidewire) and clear the inner lumen. If restriction still remains and the distal tip of the SOFIA Flow 88 Aspiration Catheter has been repositioned correctly, close the stopcock, re-attach the syringe, and resume aspiration. Maintain aspiration to make sure the emboli or thrombi remains fully engaged with the distal tip of the SOFIA Flow 88 Aspiration Catheter. With the emboli or thrombi fully engaged pull back the SOFIA Flow 88 Aspiration Catheter out of the patient body.
WARNING: Excessive aspiration with the distal tip of the SOFIA Flow 88 Aspiration Catheter covered by the vessel wall may cause vessel injury.
WARNING: Do not attempt to clear inner lumen of the SOFIA Flow 88 Aspiration Catheter by infusion while keeping the device in the patient body.
 8. After aspiration is completed, remove the SOFIA Flow 88 Aspiration Catheter from the patient body. If re-access to the vasculatures with the same device is desired, flush and clean the inner lumen of the device by infusion. Inspect the device for any damage. Follow steps 5 through 8 in the “Delivery of the SOFIA Flow 88 Aspiration Catheter” section for navigation.
- WARNING:** Do not use the device if any damages or irregularities are observed.

ASPIRATION THROUGH THE SOFIA FLOW 88 ASPIRATION CATHETER WITH PUMP

1. Under fluoroscopic guidance, make sure that the distal tip of the SOFIA Flow 88 Aspiration Catheter is located in a straight portion of the vessel. Reposition the distal tip if necessary. Refer to Figure-1 for location of the distal tip.
WARNING: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
WARNING: Torquing the SOFIA Flow 88 Aspiration Catheter excessively while kinked may damage the device resulting in separation of the device. Withdraw the entire device (the device, Compatible Catheter, and Guidewire) if the device is severely kinked.

Figure-1

2. Make sure the distal tip of the SOFIA Flow 88 Aspiration Catheter is engaged with the emboli or thrombi under fluoroscopy.
3. Attach aspiration tubing to the aspiration pump and turn on the pump (refer to the directions for use of the aspiration tubing and aspiration pump manual). Confirm that the aspiration gauge reads -20 inHg. Ensure the stopcock/valve on the aspiration tubing is in the closed position.
4. Remove the perfusion line from the side arm of the RHV that is connected to the SOFIA Flow 88 Aspiration Catheter and attach the aspiration tubing to the side arm of the RHV.
5. Tighten the RHV.
6. Make sure the distal tip of the SOFIA Flow 88 Aspiration Catheter is engaged with the embolus or thrombus under fluoroscopy.
7. To begin aspiration, turn the aspiration tubing stopcock/valve to the open position and check to see if blood, thrombus, or embolus are aspirated through the system.
8. If after 10 seconds blood is still observed flowing through the system, stop aspiration. To stop aspiration, turn aspiration tubing stopcock/valve to the closed position. Carefully reposition the distal tip of the SOFIA Flow 88 Aspiration Catheter to engage the thrombus and resume aspiration.
9. If flow is restricted or absent, maintain aspiration to make sure any thrombus or embolus is fully engaged with the distal tip of the SOFIA Flow 88 Aspiration Catheter. With the thrombus or embolus fully engaged, slowly pull back the SOFIA Flow 88 Aspiration Catheter and completely withdraw out of the patient.
WARNING: Excessive aspiration with the distal tip of the SOFIA Flow 88 Aspiration Catheter engaged with vessel wall may cause vessel injury.
10. With the SOFIA Flow 88 Aspiration Catheter outside of the patient body, flush the catheter to clear the catheter of any thromboembolic material that may be inside.
WARNING: Do not attempt to clear inner lumen of the SOFIA Flow 88 Aspiration Catheter by infusion while keeping the device in the patient body. Remove the SOFIA Flow 88 Aspiration Catheter from the patient body before attempting to clear the lumen.
11. If repeated access to the vasculature with the same device is desired, flush and clean the inner lumen of the device by infusion. Inspect the device for any damage. Reintroduce the SOFIA Flow 88 Aspiration Catheter into the body and follow steps 5 through 8 in the “Delivery of the SOFIA Flow 88 Aspiration Catheter” section to navigate the catheter to the target site.
WARNING: Do not use the device if any damage or irregularities are observed.

12. Aspiration may be attempted up to 3 times, for a total of 2 minutes of aspiration per attempt with the SOFIA Flow 88 Aspiration Catheter.
13. Once aspiration is completed, remove the SOFIA Flow 88 Aspiration Catheter from the patient body.

The physician has the discretion to modify described manipulations of the SOFIA Flow 88 Aspiration Catheter to accommodate the complexity and variation in procedures. Any technique modification must be consistent with previously described instructions, warnings, precautions, and patient safety information.

STORAGE

Keep dry and away from sunlight. See the product label for the device shelf life. Do not use the device beyond the labeled shelf life.

MATERIALS

The SOFIA Flow 88 Aspiration Catheter is not manufactured with natural rubber latex, polyvinylchloride (PVC), or di-2-ethylhexyl phthalate (DEHP).

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) for the device will be accessible in the European database on medical devices (Eudamed: <https://ec.europa.eu/tools/eudamed>), when available.

WARRANTY

MicroVention, Inc. warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness. Handling, storage, cleaning and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure and other matters beyond MicroVention’s control directly affect the device and the results obtained from its use. MicroVention’s obligation under this warranty is limited to the repair or replacement of this device through its expiration date. MicroVention shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. MicroVention neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. MicroVention assumes no liability with respect to devices reused, reprocessed or resterilized and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.

Prices, specifications, and model availability are subject to change without notice.

© Copyright 2023 MicroVention, Inc. All rights reserved.

MicroVention™ and SOFIA™ are trademarks of MicroVention, Inc., registered in the United States and other jurisdictions.

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© Copyright 2023 MicroVention, Inc. Tous droits réservés.
© Copyright 2023 MicroVention, Inc. Alle Rechte vorbehalten.
© Copyright 2023 MicroVention, Inc. Todos los derechos reservados.
© Copyright 2023 MicroVention, Inc. Tutti i diritti riservati.
© Copyright 2023 MicroVention, Inc. Todos os direitos reservados.
© Copyright 2023 MicroVention, Inc. Alle rettigheder forbeholdes.