LVIS™ EVO™

Intraluminal Support Device

Instructions for Use

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The MicroVention LVIS (Low-Profile Visualized Intraluminal Support) EVO device [Figures 1, 2 and 3a] is a self-expanding Nitinol with Platinum core, single wire braid, compliant, closed-cell design that can be deployed and retrieved by a single operator. The LVIS EVO device is sterile and non-pyrogenic and is packaged as a single unit with an introducer sheath and a detachable push wire.

Figure 1.
LVIS EVO Device - Components



Figure 2.
LVIS EVO Device – Implant Delivery

	Distal Markers	Helical Markers	Proximal Markers
LVIS EVO Device	4	N/A	4

Figure 3a.
LVIS EVO Device Implant Dimensions


Table 1: LVIS EVO Device Product Specification and Undeployed Length, Metal Surface Area %

Table 1: LVIS EVO Device Product Information

LVIS EVO						
Product Code	Undeployed Length † (mm)	Total Length / Working Length* (mm) at device diameter (mm)				
		2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm
LEV2512	20	16 / 15	12 / 11			
LEV2517	29	23 / 22	17 / 16			
LEV2522	38	30 / 29	22 / 21			
LEV2527	47	36 / 35	27 / 26			
LEV3018	34	28 / 27	24 / 23	18 / 17		
LEV3024	44	37 / 36	31 / 30	24 / 23		
LEV3028	54	46 / 45	39 / 38	28 / 27		
LEV3032	60	50 / 49	43 / 42	32 / 31		
LEV3517	32	28 / 27	26 / 25	22 / 21	17 / 16	
LEV3522	44	39 / 38	35 / 34	30 / 29	22 / 21	
LEV3528	56	49 / 48	44 / 43	37 / 36	28 / 27	
LEV3534	67	60 / 59	53 / 52	45 / 44	34 / 33	
LEV4013	22		20 / 19	18 / 17	15 / 14	13 / 12
LEV4018	36		31 / 30	28 / 27	24 / 23	18 / 17
LEV4021	43		36 / 35	33 / 32	28 / 27	21 / 20
LEV4027	56		48 / 47	43 / 42	37 / 36	27 / 26
LEV4031	63		53 / 52	48 / 47	41 / 40	31 / 30
All Sizes Compatible with Headway™ 17 Microcatheter (inner diameter = 0.017" or 0.43 mm), Scepter C™ or Scepter XC™ Occlusion Balloons						
* Total Length (which includes flared ends) = Working Length + 1 mm (0.5 mm each side)						
† Within Headway 17 Microcatheter (inner diameter = 0.017" or 0.43 mm)						

Table 1 (continued): LVIS EVO Device Product Information

LVIS EVO					
Product Code	Device Diameter (mm)				
	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm
	Free Area (%)				
LEV2512	74	74			
LEV2517	74	72			
LEV2522	74	72			
LEV2527	73	72			
LEV3018	76	77	75		
LEV3024	75	77	75		
LEV3028	75	77	74		
LEV3032	76	77	74		
LEV3517	76	78	79	77	
LEV3522	75	78	79	76	
LEV3528	75	78	79	76	
LEV3534	76	78	78	75	
LEV4013		81	83	83	81
LEV4018		79	81	81	79
LEV4021		80	81	81	78
LEV4027		80	81	81	78
LEV4031		81	82	82	79

INDICATIONS FOR USE

The LVIS EVO device is indicated for use with neurovascular embolization coils in patients ≥ 18 years of age for the treatment of wide-neck (neck width ≥ 4 mm or dome to neck ratio < 2) saccular intracranial aneurysms arising from a parent vessel with a diameter ≥ 2.0 mm and ≤ 4.5 mm.

CONTRAINDICATIONS

Use of the LVIS EVO device is contraindicated under these circumstances:

- Patients in whom anticoagulant, anti-platelet therapy or thrombolytic drugs are contraindicated;
- Patients with known hypersensitivity to metal, such as nickel-titanium and metal jewelry;
- Patients with anatomy that does not permit passage or deployment of the LVIS EVO device;
- Patients with an active bacterial infection;
- Patients with a pre-existing stent in place at the target aneurysm.

WARNINGS

Do not use the LVIS EVO device for acutely ruptured intracranial aneurysms within a minimum of 30 days from intracranial aneurysm rupture.

Should unusual resistance be felt at any time during access or removal, the introducer/microcatheter and LVIS EVO device should be removed as a single unit. Applying excessive force during delivery or retrieval of the LVIS EVO device can potentially result in loss or damage to the device and delivery components.

The LVIS EVO device should only be used by physicians trained in endovascular interventional neuroradiology, radiology, neurosurgery or interventional neurology for the treatment of intracranial aneurysms.

Selection of the LVIS EVO device size is important for proper product performance and patient safety and must be based on pre-treatment angiograms for correct and accurate vessel measurements from multiple views.

It is imperative to use the LVIS EVO device with compatible microcatheters. If repeated friction is encountered during LVIS EVO device delivery, verify microcatheter is not kinked or in extremely tortuous anatomy. Confirm that the microcatheter does not ovalize. Confirm that there is adequate sterile flush solution.

Do not reposition the LVIS EVO device in the parent vessel without fully retrieving the device. The LVIS EVO device MUST be retrieved into the microcatheter and re-deployed at the desired target location or removed completely from the patient.

Do not attempt to re-position the LVIS EVO implant after detachment.

Do not shape the tip of the delivery wire.

Do not torque the delivery wire while advancing or retracting the LVIS EVO device. A torque device should not be used.

PRECAUTIONS

The LVIS EVO device is provided sterile for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Carefully inspect the sterile package and the LVIS EVO device prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components, or if the packaging is damaged.

See the product label for the device shelf life. Do not use the device beyond the labeled use by date.

Exercise caution when crossing the deployed/detached LVIS EVO device with adjunctive devices such as guidewires, catheters, microcatheters or balloon catheters to avoid disrupting the device geometry and device placement.

The LVIS EVO device with neurovascular embolization coils may create local field inhomogeneity and susceptibility artifacts during magnetic resonance angiography (MRA), which may degrade the diagnostic quality to assess effective intracranial aneurysm treatment.

The safety and effectiveness of the device has not been established in the treatment of large and giant wide-neck intracranial aneurysms.

The benefits may not outweigh the risks of treatment in patients with wide-neck intracranial aneurysms ≤ 5 mm in size, or reduced life expectancy, in the absence of additional risk factors for intracranial aneurysm rupture.

The safety and effectiveness of the device has not been well established in the posterior circulation.

Ensure that the specific embolization coil models and sizes used are indicated for the embolization of intracranial aneurysms.

Potential Adverse Events

The following potential risks and complications associated with general anesthesia, cerebral angiography, intracranial catheterization, intracranial stent placement or intra-saccular coil deployment have been identified below:

- Allergic reaction, including but not limited to: contrast dye, nitinol metal, and any other medications used during the procedure;

- Aphasia
- Blindness;
- Cardiac Arrhythmia;
- Coil prolapsed or migration into normal vessel adjacent to aneurysm
- Complications of arterial puncture including pain, local bleeding, local infection and injury to the artery, vein or adjacent nerves;
- Cranial neuropathy;
- Death;
- Device fracture, migration or misplacement;
- Dissection or perforation of the parent artery;
- Headache;
- Hemorrhage (i.e., intracerebral hemorrhage (ICH), subarachnoid hemorrhage (SAH), or retroperitoneal (or in other locations));
- Hemiplegia;
- Hydrocephalus;
- Infection;
- Injury to normal vessel or tissue;
- Ischemia;
- Mass effect;
- Myocardial Infarction;
- Neurological deficits;
- Occlusion of non-target side branches;
- Pseudo aneurysm formation;
- Reactions to anti-platelet/anti-coagulant agents;
- Reactions due to radiation exposure;
- Reactions to anesthesia and related procedures;
- Reactions to contrast agents;
- Renal failure;
- Aneurysm rupture;
- Stenosis of stented segment;
- Seizure;
- Stent thrombosis;
- Stroke or TIA (Transient Ischemic Attack);
- Thromboembolic event (T/E);
- Vasospasm;
- Visual impairment.

Potential Risks Associated with X-ray Exposure: The use of the LVIS EVO device requires fluoroscopy, which presents potential risks associated with X-ray exposure. The risks of angiographic and fluoroscopic X-ray radiation doses to the patient include risks such as alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia that increase in probability as procedure time and number of procedures increase. The probability of adverse event occurrence increases as the procedure time and the number of procedures increase. Operators should take all necessary precautions to limit X-ray radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors whenever possible.

MRI Safety Information

The LVIS EVO (Low-profile Visualized Intraluminal Support) device is MR Conditional. A patient with the LVIS EVO (Low-profile Visualized Intraluminal Support) device may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

MR Conditional

Parameter	Condition
Nominal Values of Static Magnetic Field (T)	1.5-T or 3-T, only
Maximum Spatial Field Gradient (T/m and gauss/cm)	40-T/m (4,000-gauss/cm)
Type of RF Excitation	Circularly Polarized (CP) (i.e., quadrature-driven)
Transmit RF Coil Information	There are no transmit RF coil restrictions. Accordingly, the following may be used: body transmit RF coil and all other RF coil combinations (i.e., body RF coil combined with any receive-only RF coil, transmit/receive head RF coil, transmit/receive knee RF coil, etc.)
Operating Mode of MR System	Normal Operating Mode
Maximum Whole Body Averaged SAR	2-W/kg (Normal Operating Mode)
Maximum Head SAR	3.2-W/kg (Normal Operating Mode)
Limits on Scan Duration	Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back to back sequences/series without breaks)
MR Image Artifact	The presence of this implant produces an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant is located in the area of interest.

CLINICIAN USE INFORMATION

Materials

The following parts are required to use the LVIS EVO devices:

- LVIS EVO devices should be introduced only Headway 17 Microcatheter (0.017 inch inner diameter) or Scepter C / Scepter XC Occlusion Balloon (0.0165 inch inner diameter).

Other accessories for performing a procedure and NOT supplied; to be selected based on the physician's experience and preferences:

- Appropriate-sized Guiding catheter for use with selected microcatheter
- Headway 17 microcatheter
- Scepter C / Scepter XC Occlusion Balloon
- Microcatheter-compatible guidewires
- Saline solution/heparin-saline solution continuous flush set
- Contrast solution
- Rotating Hemostatic Valve (RHV)
- Pressurized sterile Infusion solutions – IV stand
- Femoral arterial sheath, compatible with delivery guide catheter
- Femoral artery access device, sterile needle, guidewire

The LVIS EVO device does not contain latex or PVC materials.

PACKAGING AND STORAGE

The LVIS EVO device is placed inside a protective, plastic dispenser coil and packaged in a pouch and unit carton. The LVIS EVO device and dispenser coil will remain sterile unless the package is opened, damaged, or the expiration date has passed. Keep dry and away from sunlight .

SHELF LIFE

See the product label for the device shelf life. Do not use the device beyond the labeled use by date.

PREPARATION FOR USE

Device and Delivery System Selection

Appropriate selection of the LVIS EVO device is important for patient safety. In order to choose the optimal LVIS EVO device model size for any given lesion, examine pre-treatment angiograms for correct and accurate vessel measurements.

HOW SUPPLIED

Sterile: This device is sterilized with E-Beam irradiation. Non-pyrogenic.

Contents: One (1) LVIS EVO device.

Storage: Keep dry and away from sunlight.

DIRECTIONS FOR USE

1. Gain vascular access according to standard practice.
 2. Place guide catheter in the appropriate target vessel.
 3. a. Navigate the microcatheter (.017" MicroVention Headway 17 or Scepter C / Scepter XC Occlusion Balloon) over a guidewire at least 15 mm distal to the aneurysm neck or target location.
b. A second microcatheter can be navigated into the aneurysm sac for future coil deployment steps using the jailing technique (steps 22-24). In this technique, the microcatheter is effectively jailed between the vessel wall and outer surface of the stent and the coils are kept within the aneurysm and outside of the reconstructed vessel lumen.
 4. Remove the guidewire.
 5. Maintain flush through the microcatheter per standard endovascular practice.
 6. Select an appropriately sized LVIS EVO device (Refer to Table 1).
 7. Carefully inspect the LVIS EVO device package for damage to the sterile barrier.
 8. Peel open the pouch using aseptic technique.
 9. Carefully place the dispenser coil into the sterile field.
 10. a. Unclip the molded cap attached to the delivery wire from the dispenser coil. Pull on the proximal end of the delivery wire until the introducer exits the dispenser coil. Hold the delivery wire and introducer together while continuing to remove the entire device. Do not partially deploy the LVIS EVO device from the introducer.
b. After removal from the dispenser coil, carefully push on the delivery wire and in a bowl of saline, partially deploy the LVIS EVO implant up to 5 mm or 50% (whichever occurs first, being careful not to detach the implant) from the distal introducer tip (Refer to Table 1 and Figure 3b). Check for the following:
 - Implant distal marker uniformity
 - Implant distal end shows even displacement with no entanglement
 - Implant tracks smoothly through introducer
- Warning: DO NOT FULLY DEPLOY LVIS EVO device.**
If the device is deployed, DO NOT attempt to reload the device. Use a new device.
- c. With the LVIS EVO implant and introducer sheath positioned and hydrated within the bowl of saline, gently manipulate the LVIS EVO implant within the saline to hydrate the implant and minimize visible air bubbles. Carefully pull back on the delivery wire to fully retrieve the LVIS EVO implant and the delivery wire tip within the introducer.
- Warning: DO NOT CONTINUE** if any defect is observed; return the unit to MicroVention, Inc.
11. Confirm that the tip of the delivery wire is entirely within the introducer.
 12. Confirm that the delivery wire is not kinked and that the introducer tip is not damaged. **DO NOT CONTINUE** if either defect is observed; return the unit to MicroVention, Inc.
Warning: Do not shape the tip of the delivery wire.
 13. Partially insert the distal end of the introducer into the RHV connected to the microcatheter. Tighten the RHV locking ring. Flush the y-connector of the RHV with sterile saline and verify that fluid exits the proximal end of the introducer.
Warning: Purge the LVIS EVO device carefully to avoid the accidental introduction of air into the system. [Figure 4]
 14. Untighten the RHV locking ring and advance the introducer until it is **fully engaged** with the microcatheter hub, then tighten the RHV locking ring.
Warning: Confirm that there are no air bubbles trapped anywhere in the system.
Caution: Confirm that there is no gap between the introducer and the microcatheter hub to enable LVIS EVO device introduction into the microcatheter. [Figure 5]
 15. Advance the delivery wire to transfer the LVIS EVO device from within the introducer into the microcatheter.
Warning: Do not torque the delivery wire while advancing or retracting the LVIS EVO device. A torque device should not be used.
 16. Continue advancing the delivery wire into the microcatheter until the proximal tip of the delivery wire enters the introducer. Loosen the RHV locking ring, remove the introducer, and set it aside.
Note: Fluoroscopy may be used up to this point at the physician's discretion.
Warning: Do not apply undue force. If resistance is encountered at any point during LVIS EVO device delivery or manipulation, withdraw the unit and select a new LVIS EVO device.
 17. Track the LVIS EVO device through the microcatheter to the tip. Carefully advance the LVIS EVO device until the device exit marker on the proximal end of the delivery wire approaches the RHV on the hub of the microcatheter. At this time, fluoroscopic guidance must be initiated.
 18. Position the LVIS EVO device for deployment, ensuring a sufficient length of stent will be deployed on either side of the aneurysm neck, by aligning the LVIS EVO implant distal radiopaque end markers sufficiently past the aneurysm neck. [Figure 6]

Note: A proper push/pull technique, encompassing sufficient delivery wire push force, in addition to an opposing microcatheter withdrawal force, will facilitate properly deploying the LVIS EVO device to achieve full expansion and good vessel apposition.

Note: Slowly advancing the LVIS EVO device while adjusting the microcatheter position will ensure accurate deployment. Maintain simultaneous control of the LVIS EVO device and microcatheter in order to position and expand the device at the proper location.

Caution: Using a rapid microcatheter withdrawal technique to deploy the LVIS EVO device is not recommended and may result in device elongation.

19. If LVIS EVO device positioning is not satisfactory, the LVIS EVO device may be recaptured and repositioned if it is not fully deployed. The LVIS EVO device may be recaptured until the point where the proximal end of the LVIS EVO device markers are aligned 3 mm proximally with the microcatheter distal marker band (approximately 75% deployed). [Figure 7]

Caution: If resistance is felt while recapturing the LVIS EVO device, do not continue to recapture the device. Withdraw the microcatheter slightly to unsheath the LVIS EVO device (without exceeding the recapture limit), and then attempt to recapture the LVIS EVO device.

Caution: The LVIS device must not be re-deployed more than three times.

Note: The LVIS EVO device delivery wire should not be utilized as a guidewire after stent deployment. Do not torque the LVIS EVO device. A torque device should not be used.

20. If LVIS EVO device positioning is satisfactory, carefully retract the microcatheter and advance the delivery wire together, to allow the LVIS EVO device to deploy across the neck of the aneurysm. Ensure the device proximal radiopaque end markers are sufficiently proximal to the aneurysm neck to ensure an adequate landing zone. The LVIS EVO device will expand, and the total length may foreshorten up to 60% from its undeployed length (refer to Table 1) as it exits the microcatheter. Ensure microcatheter is retracted and clear from the proximal flared ends.

Note: Visualize and refer to the implant radiopaque end markers to maintain adequate implant length on each side of the aneurysm neck or target location to ensure appropriate neck coverage. [Figure 8]

Warning: Do not detach the LVIS EVO device if it is not properly positioned in the parent vessel. Observe the delivery wire distal tip to assure it remains within the desired location of the parent vessel.

21. Prior to removing the delivery wire and if necessary, carefully position the microcatheter distal to the LVIS EVO device to maintain access through the LVIS EVO device. Remove and discard the delivery wire.

Warning: The LVIS EVO device delivery wire should not be utilized as a guidewire. Do not torque the LVIS EVO device. A torque device should not be used.

22. a. If applicable, advance a .017" inner diameter (or suitable size) microcatheter over the guidewire.
b. If a second microcatheter has been placed into the aneurysm in step 3b, detachable coils can be delivered into the aneurysm through the second microcatheter (replacing steps 22-24).

Warning: Ensure that the jailed microcatheter does not move while constantly observing LVIS EVO device marker positions during the coiling procedure to ensure that the device does not migrate from its deployed position.

23. Use the guidewire and microcatheter to access the aneurysm through the LVIS EVO device cells.

Warning: Observe LVIS EVO device marker position during placement of the microcatheter into the aneurysm to ensure that the LVIS EVO device does not migrate or dislodge from its deployed position.

Note: Access to the aneurysm may be facilitated by the use of a microcatheter that has been shaped.

24. After the microcatheter is positioned within the aneurysm, detachable coils may be delivered into the aneurysm according to conventional methods.


Warning: Observe LVIS EVO device marker position during the coiling procedure to ensure that the device does not migrate from its deployed position.

25. After placing the last coil, verify that the LVIS EVO device has remained patent and properly positioned. Advance a guidewire, if necessary, to the microcatheter tip and carefully remove the microcatheter.

Note: A microcatheter may be positioned into the aneurysm sac prior to delivery of the LVIS EVO device. The microcatheter will be supported by the LVIS EVO device during delivery of embolic coiling. After completing the coiling, the coiling microcatheter should be carefully removed to avoid dislodging the LVIS EVO device.


26. After completing the procedure, withdraw and discard all applicable accessory devices.

27. **Caution:** Carefully watch the LVIS EVO device distal and proximal markers when passing through the deployed LVIS EVO device with embolic coiling microcatheters to avoid displacing the LVIS EVO device


- Check for the following:
- Implant distal marker uniformity
 - Implant distal end shows even displacement with no entanglement
 - Implant tracks smoothly through introducer
- Warning: DO NOT FULLY DEPLOY LVIS EVO device.

[Figure 3b. Step 10b]

13. Verify that fluid exits the proximal end of the introducer

[Figure 4. Step 13]


[Figure 5. Step 14]

18. Position distal markers sufficiently distal to the aneurysm neck

[Figure 6. Step 18]

19. The LVIS EVO device can be recaptured and repositioned if not yet fully deployed.

[Figure 7. Step 19]

20. Ensure sufficient proximal landing zone from aneurysm neck.

[Figure 8. Step 20]

WARRANTY DISCLAIMER

MicroVention warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for particular purpose. Handling, storage of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure, and other matters beyond MicroVention's control directly affect the device and the results obtained from its use. MicroVention's obligation under this warranty is limited to the replacement of this device through its expiration date. MicroVention shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. MicroVention neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. MicroVention assumes no liability with respect to devices reused, reprocessed or resterilized or used after the expiration date stated on product label and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.

Prices, specifications, and model availability are subject to change without notice.

© Copyright 2025 MicroVention, Inc. All rights reserved.

MicroVention™, LVIS™, EVO™, Headway™, Scepter C™ and Scepter XC™ are trademarks of MicroVention, Inc., registered in the United States and other jurisdictions. All third-party products are trademarks™ or registered © trademarks and remain the property of their respective holders.

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© Copyright 2025 MicroVention, Inc. Tous droits réservés.