# **FRED™ Jr.**

Flow Re-Direction Endoluminal Device

## **INSTRUCTIONS FOR USE**

<b>FRED™ JR. SYSTEM</b>	
Instructions for Use .....	English .....1
<b>SYSTÈME FRED™ JR</b>	
Mode d'emploi.....	Français .....4
<b>FRED™ JR. SYSTEM</b>	
Gebrauchsanweisung.....	Deutsch .....7
<b>SISTEMA FRED™ JR.</b>	
Instrucciones de uso .....	Español .....10
<b>SISTEMA FRED™ JR.</b>	
Istruzioni per l'uso .....	Italiano .....13
<b>SISTEMA FRED™ JR.</b>	
Instruções de Utilização .....	Português .....16
<b>FRED™ JR.-SYSTEM</b>	
Brugsanvisning.....	Dansk .....19
<b>FRED™ JR.-SYSTEEM</b>	
Gebruiksaanwijzing .....	Nederlands .....22
<b>FRED™ JR. -JÄRJESTELMÄ</b>	
Käyttöohjeet .....	Suomi .....25
<b>FRED™ JR.-SYSTEM</b>	
Bruksanvisning .....	Svenska .....28
<b>FRED™ JR. SYSTEM</b>	
Bruksanvisning .....	Norsk .....31
<b>ΣΥΣΤΗΜΑ FRED™ JR.</b>	
Οδηγίες χρήσης .....	Ελληνικά .....34
<b>FRED™ JR. SİSTEMİ</b>	
Kullanım Talimatları.....	Türkçe .....37
<b>СИСТЕМА FRED™ JR.</b>	
Инструкции за употреба .....	български .....40
<b>SUSTAV FRED™ JR.</b>	
Upute za upotrebu.....	Hrvatski.....43
<b>SYSTÉM FRED™ JR.</b>	
Návod k použití.....	Česky.....46

<b>FRED™ JR. SÜSTEEM</b>	
Kasutusjuhend.....	Eesti keel .....49
<b>FRED™ JR. SYSTEM</b>	
Használati útmutató.....	Magyar.....52
<b>FRED™ JR. SISTĒMA</b>	
Lietošanas instrukcija.....	Latviešu valodā.....55
<b>„FRED™ JR.“ SISTEMA</b>	
Naudojimo instrukcija.....	Lietuviškai.....58
<b>SYSTEM FRED™ JR.</b>	
Instrukcja użycia.....	Polski.....61
<b>SISTEMUL FRED™ JR.</b>	
Instrucţiuni de utilizare.....	Română.....64
<b>СИСТЕМА FRED™ JR.</b>	
Инструкция по применению.....	Русский.....67
<b>FRED™ JR. SISTEM</b>	
Uputstvo za upotrebu.....	Srpski.....70
<b>SISTEM FRED™ JR.</b>	
Navodila za uporabo.....	Slovenščina.....73
<b>SYSTÉM FRED™ JR.</b>	
Návod na použitie.....	Slovenčina.....76
<b>СИСТЕМА FRED™ JR.</b>	
Інструкція з використання.....	Українська.....79
<b>ماظن FRED™ JR.</b>	
تعليمات الاستخدام.....	العربية.....82

# English FRED™ Jr. System Instructions for Use

## DEVICE DESCRIPTION

The MicroVention Flow Re-Direction Endoluminal Device (FRED Jr.) system (Figures 1, 2 and 3) is a self-expanding nickel-titanium, single wire braid, compliant closed cell paired-stent design that can be simultaneously deployed and retrieved by a single operator. The FRED Jr. system features integrated dual layer coverage designed to focus mainly at the neck of an aneurysm. The FRED Jr. system has distal and proximal markers on its ends as well as interweaved helical marker strands delineating the inner working length of the stent to provide fluoroscopic visibility. The FRED Jr. system is packaged sterile as a single unit with an introducer sheath and a detachable delivery wire.

Figure 1. FRED Jr. System Setup


Figure 2: FRED Jr. System Marker Nomenclature


Figure 3: Fred Jr. System Implant


Table 1: Qualitative and Quantitative Implant Material Information

Implant Material		Approximate Mass (mg)
Metallic Components	Nitinol, Tantalum	≤ 49
Non-metallic Components	Adhesive	≤ 0.13

The FRED Jr. system may also be used with embolic coils for the treatment of intracranial neurovascular lesions.

## CONTRAINDICATIONS

Use of the FRED Jr. system is contraindicated under these circumstances:

- Patients in whom anticoagulant, antiplatelet therapy or thrombolytic drugs are contraindicated
- Patients with known hypersensitivity to nickel-titanium
- Patients in whom angiography demonstrated inappropriate anatomy that does not permit passage or deployment of the FRED Jr. system

## POTENTIAL COMPLICATIONS

Possible complications include but are not limited to the following:

- Bleeding or Hemorrhage including intracerebral, retroperitoneal or other locations
- Complications of arterial puncture including pain, local bleeding (hematoma) or injury to the artery or adjacent nerves
- Device migration
- Distal Embolization
- Headache
- Incomplete aneurysm occlusion
- Neurologic deficits including stroke and/or death
- Perforation or dissection of the vessel(s)
- Pseudoaneurysm formation
- Rupture or perforation of aneurysm
- Transient ischemic attack (TIA) or ischemic stroke
- Vasospasm
- Vessel occlusion
- Vessel stenosis or thrombosis

## WARNINGS

Should unusual resistance be felt at any time during access or removal, the introducer/guide catheter/microcatheter and FRED Jr. system should be removed as a single unit. Applying excessive force during delivery or retrieval of the FRED Jr. system can potentially result in loss or damage to the device and delivery components.

The FRED Jr. system should only be used by physicians trained in endovascular interventional neuroradiology, radiology, neurosurgery or interventional neurology for the treatment of intracranial aneurysms or other vascular lesions.

It is imperative to use the FRED Jr. system with a Headway™ 21 microcatheter. If repeated friction is encountered during FRED Jr. system delivery, verify microcatheter is not kinked or in extremely tortuous anatomy. Confirm that the microcatheter does not ovalize. Confirm that there is adequate sterile heparinized flush solution.

Do not reposition the FRED Jr. system in the parent vessel without fully retrieving the device. The FRED Jr. system MUST be retrieved/resheathed into the microcatheter and re-deployed at the desired target location or removed completely from the patient.

Do not attempt to re-position the FRED Jr. implant after deployment/detachment.

The safety and performance of the FRED Jr. system for use in conjunction with intravascular medical devices other than neurovascular embolization coils have not been established.

## CAUTIONS

This product should only be used by experienced physicians who have completed endovascular training in the use of the FRED Jr. system. This device is used for percutaneous neurointerventional and peripheral vascular procedures as indicated by a representative from MicroVention-Terumo or a MicroVention-authorized distributor.

The FRED Jr. system does not contain latex or PVC materials.

The FRED Jr. system is provided sterile for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Carefully inspect the sterile package and the FRED Jr. system prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components, or if the package is opened or damaged.

See the product label for shelf life. Do not use the FRED Jr. system beyond the labeled use by date.

Exercise caution when crossing the deployed/detached FRED Jr. system with adjunctive devices such as guidewires, catheters, microcatheters or balloon catheters to avoid disrupting the device geometry and device placement.

## CLINICAL BENEFIT

FRED devices are intended for embolization of intracranial aneurysms as well as treatment of other neurovascular lesions when used with embolic coils. The embolization results in occlusion of intracranial neurovascular aneurysms.

## ADDITIONAL NOTICE TO USER

- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

The FRED Jr. delivery wire contains a hazardous substance, Cobalt (CAS No. 7440-48-4) is classified as a carcinogenic, mutagenic, or toxic for reproduction (CMR) substance of Category 1A or 1B, or as an endocrine disruptor (ED), and is present at a concentration >0.1% weight/weight.

## INTENDED PURPOSE/ INDICATIONS

The FRED Jr. system is intended for endovascular embolization of intracranial neurovascular aneurysms.

- The Summary of Safety and Clinical Performance (SSCP) for the device will be accessible in the European database on medical devices after the launch of the European Database on Medical Devices (EUDAMED) (EUDAMED: <https://ec.europa.eu/tools/eudamed>). The SSCP will be linked to the Basic UDI-DI in the EUDAMED public website.
- Dispose of device in accordance with hospital policy and local regulations for biohazardous waste.
- A patient implant card is included in the package. This card should be completed and provided to the patient.
- Permanent implant. Follow-up required at the discretion of the physician
- The electronic instructions for use (eIFU) is available via MicroVenton website: <https://terumoneuro.com/products/product-use-and-safety>

## MR INFORMATION

Non-clinical testing has demonstrated that the FRED Jr. system is MR Conditional. A patient with this device can be safely scanned in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2,500-gauss/cm (25-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the FRED Jr. system is expected to produce a maximum temperature rise of 2.8 °C at 1.5-Tesla or 3.6 °C at 3-Tesla after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the device extends approximately 4-mm from the FRED Jr. system when imaged with a gradient echo pulse sequence and a 3-Tesla MR system.

MicroVenton, Inc. recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation or equivalent organization. A FRED Jr. system patient implant card is included in the package, which should be completed and provided to the patient.

## CLINICIAN USE INFORMATION

### Materials

The following parts are required to use the FRED Jr. system:

- FRED Jr. system should be introduced only by means of a Headway 21 microcatheter

Other accessories for performing a procedure and NOT supplied; to be selected based on the physician's experience and preferences:

- Appropriate-sized Guiding catheter for use with selected microcatheter
- Headway 21 microcatheter
- Microcatheter-compatible guidewires
- Saline solution/heparin-saline solution continuous flush set
- Contrast solution
- Rotating Hemostatic Valve (RHV)
- Pressurized sterile Infusion solutions – IV stand
- Femoral arterial sheath, compatible with delivery guide catheter
- Femoral artery access device, sterile needle, guidewire

## PACKAGING AND STORAGE

The FRED Jr. system is placed inside a protective, plastic dispenser coil and packaged in a pouch and unit carton. The FRED Jr. system and dispenser coil will remain sterile unless the package is opened, damaged, or the expiration date has passed. Keep Dry and away from Sunlight.

## SHELF LIFE

See the product label for the device shelf life. Do not use the device beyond the labeled use by date.

## PREPARATION FOR USE

### Device and Delivery System Selection

Appropriate selection of the FRED Jr. system is important for patient safety. In order to choose the optimal FRED Jr. system model size for any given lesion, examine pre-treatment angiograms for correct and accurate vessel measurements.

### Directions For Use

1. Gain vascular access according to standard angiographic practice and perform diagnostic angiogram to document target aneurysm and parent vessel to confirm vessel diameter.
2. Place appropriate size guidecatheter according to standard practice.
3. Per physician discretion, if embolization coils are going to be used along with the FRED Jr. system, position a .017 inch (0.43 mm) inner diameter (or suitable size) microcatheter (to be utilized for coil embolization) coaxially into the target aneurysm. If the FRED Jr. system is going to be used without embolization coils, disregard this step.
4. Coaxially navigate a microcatheter (MicroVenton Headway 21 microcatheter) over a guidewire distal to the aneurysm neck or target location. Remove the guidewire.
5. Maintain flush through the microcatheter(s) per standard endovascular practice.
6. Select an appropriate-sized FRED Jr. system (Refer to Figure 3 and FRED JR. System Inservice Guide) according to the size of the parent vessel/aneurysm neck.
 

**Note:** The FRED Jr. system implant foreshortens significantly (up to 50%) as it expands to the diameter of the parent vessel. Take implant foreshortening into account when sizing and deploying the FRED Jr. system.
7. Carefully inspect the package for damage to the sterile barrier. Peel open the pouch using aseptic technique and place the dispenser coil into the sterile field.

- a. Unclip the molded cap attached to the delivery wire from the dispenser coil. Pull on the proximal end of the delivery wire until the introducer exits the dispenser coil. Hold the delivery wire and introducer together while continuing to remove the entire device.
- b. After removal from the dispenser coil, carefully push on the delivery wire and in a bowl of saline, only partially deploy the FRED Jr. implant up to 5 mm or 50% (whatever occurs first, being careful not to detach the implant) from the distal introducer tip (Refer to Figures 1 and 4). Check for the following:
  - Implant distal marker uniformity
  - Implant distal end shows even displacement with no entanglement
  - Implant tracks smoothly through introducer

### Warning: DO NOT FULLY DEPLOY the FRED Jr. system.

c. With the FRED Jr. implant and introducer sheath positioned and hydrated within the bowl of saline, gently manipulate the FRED Jr. implant within the saline to hydrate the implant and minimize visible air bubbles. Carefully pull back on the delivery wire to fully retrieve the FRED Jr. implant within the introducer.

**Warning: DO NOT CONTINUE** if any defect is observed; return the unit to MicroVenton, Inc.

9. Confirm that the device is entirely within the introducer and the introducer tip is not damaged. **DO NOT CONTINUE** if either defect is observed; return the unit to MicroVenton, Inc.
10. Partially insert the distal end of the introducer into the RHV connected to the Headway 21 microcatheter. Tighten the RHV locking ring. Flush the RHV with sterile saline and verify that fluid exits the proximal end of the introducer, hydrating the introducer.

**Warning:** Purge the FRED Jr. system carefully to avoid the accidental introduction of air into the system. [Figure 5]

11. Untighten the RHV locking ring and advance the introducer until it is **fully engaged** with the Headway 21 microcatheter hub, then tighten the RHV locking ring.

**Warning:** Confirm that there are no air bubbles trapped anywhere in the system.

**Caution:** The introducer must be properly engaged with the microcatheter hub to enable FRED Jr. system to be introduced into the microcatheter. [Figure 6]

12. Advance the delivery wire to transfer the FRED Jr. system from within the introducer into the microcatheter.

**Warning: Do not torque the delivery wire while advancing or retracting the FRED Jr. system.**

13. Continue advancing the delivery wire into the microcatheter until the proximal tip of the delivery wire enters the introducer. Loosen the RHV locking ring, remove the introducer, and set it aside.

**Note:** Fluoroscopy may be used up to this point at the physician's discretion.

**Warning:** Do not apply undue force. If resistance is encountered at any point during delivery or manipulation, withdraw the unit and select a new FRED Jr. system.

14. Track the device through the microcatheter to the tip. Carefully advance until the device exit marker on the proximal end of the delivery wire approaches the RHV. At this time, fluoroscopic guidance must be initiated.

15. Position the FRED Jr. system for deployment by aligning the FRED Jr. system implant distal radiopaque end markers approximately 3 mm past the aneurysm neck or approximately half the width of the neck past the neck for larger aneurysms. [Figure 7]

**Note:** A slow, proper push/pull technique, encompassing sufficient delivery wire push force, in addition to an opposing microcatheter withdrawal force, to remove excess microcatheter slack while maintaining the microcatheter tip within the center of the parent vessel, will facilitate properly deploying the FRED Jr. system at the proper location, to achieve full expansion and good vessel apposition.

**Note:** If applicable, verify microcatheter placed into aneurysm in step 3 is still properly positioned for coil delivery.

**Caution:** Using a rapid microcatheter withdrawal technique to deploy the FRED Jr. system is not recommended and may result in device elongation or improper deployment. Be aware of delivery wire tip position during deployment.

16. If FRED Jr. system positioning is not satisfactory, the implant may be recaptured and repositioned if it is not fully deployed. The implant may be recaptured until the point where the distal-most wire marker, collocated distal to the implant proximal markers, is aligned approximately 50% of length proximal to the microcatheter distal marker band. [Figure 8]

**Caution:** If resistance is felt while recapturing the device, do not continue to recapture. Withdraw the microcatheter slightly to unseat the device (without exceeding the recapture limit), and then attempt to recapture again.

**Caution:** The FRED Jr. system must not be re-deployed more than three times.

**Caution: The FRED Jr. system delivery wire should not be utilized as a guidewire. Do not torque the FRED Jr. system. A torque device should not be used.**

17. If FRED Jr. system positioning is satisfactory, carefully advance the delivery wire while retracting the microcatheter as needed to minimize slack, maintaining the microcatheter around the center of the parent vessel, to allow the implant to deploy across the neck of the aneurysm. Ensure the implant proximal radiopaque end markers are approximately 3 mm or approximately half the width of the neck proximal to the aneurysm neck for adequate coverage.

**Note:** The FRED Jr. system will expand and may foreshorten up to 50% from its undeveloped length. Visually verify opening of the proximal end, ensuring that the microcatheter distal tip marker is pulled back, adequately away from the implant proximal end, to allow the proximal end to freely open.

**Note:** Visualize and refer to implant radiopaque end markers to maintain adequate implant length of approximately 3 mm or approximately half the width of the neck on each side of the aneurysm neck/target location to ensure adequate coverage. [Figure 7]

**Warning:** Do not detach the FRED Jr. system if it is not properly positioned in the parent vessel.

**Warning:** If applicable, observe FRED Jr. system marker position during coiling procedure to ensure that the device does not migrate.


18. Should you wish to position the microcatheter distal to the placed FRED Jr. stent; remove the delivery wire and replace with a guidewire that is compatible with Headway 21. Carefully navigate both the guidewire and catheter through the stent so as not to displace the stent from the desired position.  
**Caution:** The FRED Jr. system delivery wire should not be utilized as a guidewire. Do not torque the FRED Jr. system. A torque device should not be used.
19. Carefully inspect the deployed FRED Jr. implant under fluoroscopy to confirm that it is completely apposed to the vessel wall and not kinked. If the implant is not fully apposed or is kinked, consider utilizing a suitable microguidewire and/or occlusion balloon catheter to fully open the implant.  
**Note:** The jailed microcatheter should be carefully removed to avoid dislodging the FRED Jr. implant.
21. After completing the procedure, withdraw and discard all applicable accessory devices.  
**Caution:** Carefully watch the FRED Jr. implant distal and proximal markers when passing through the implanted device with other devices to avoid displacing the implant.

**HOW SUPPLIED**

Sterile: This device is sterilized with E-Beam irradiation. Non-pyrogenic  
 Contents: One (1) FRED Jr. system  
 Storage: Store product in a dry, cool place.


**WARRANTY DISCLAIMER**

MicroVention warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for particular purpose. Handling, storage, cleaning, and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure, and other matters beyond MicroVention's control directly affect the device and the results obtained from its use. MicroVention's obligation under this warranty is limited to the repair or replacement of this device through its expiration date. MicroVention shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. MicroVention neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. MicroVention makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.  
 Prices, specifications, and model availability are subject to change without notice.  
 © Copyright 2025 MicroVention, Inc. All rights reserved.  
 MicroVention™, FRED™ Jr, and Headway™ are trademarks of MicroVention, Inc., registered in the United States and other jurisdictions.  
 All third-party products are trademarks™ or registered® trademarks and remain the property of their respective holders.


[Figure 4. Step 8b. Check for the following:


- Implant distal marker uniformity
  - Implant distal end shows even displacement with no entanglement
  - Implant tracks smoothly through introducer
- Warning: DO NOT FULLY DEPLOY FRED Jr. system]


[Figure 5. Step 10. Verify that fluid exits the proximal end of the introducer]

[Figure 6. Step 11. Seat in microcatheter]

[Figure 7. Steps 15 and 17. Position distal and proximal radiopaque end markers approximately 3 mm distal and proximal, respectively to the aneurysm neck]

[Figure 8. Step 16. The FRED Jr. device can be recaptured and repositioned if not yet fully deployed]

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© Copyright 2025 MicroVention, Inc. Tous droits réservés.
© Copyright 2025 MicroVention, Inc. Alle Rechte vorbehalten.
© Copyright 2025 MicroVention, Inc. Todos los derechos reservados.
© Copyright 2025 MicroVention, Inc. Tutti i diritti riservati.
© Copyright 2025 MicroVention, Inc. Todos os direitos reservados.
© Copyright 2025 MicroVention, Inc. Alle rettigheder forbeholdes.