Precipitating Hydrophobic Injectable Liquid

Liquid Embolic System

INSTRUCTIONS FOR USE

PHIL™ Liquid Embolic System	
Instructions for Use.....	English 1
Système embolique liquide PHIL™	
Mode d'emploi.....	Français 5
PHIL™ Flüssiges Embolisierungssystem	
Gebrauchsanweisung.....	Deutsch..... 9
Sistema embólico líquido PHIL™	
Instrucciones de uso.....	Español 13
Sistema embolico liquido PHIL™	
Istruzioni per l'uso.....	Italiano 17
Sistema embólico líquido PHIL™	
Instruções de utilização.....	Português 21
PHIL™ flydende, embolisk system	
Brugsanvisning.....	Dansk 25
PHIL™ vloeibaar embolisch systeem	
Gebruiksaanwijzing.....	Nederlands..... 29
PHIL™-neste-embolisaatiojärjestelmä	
Käyttöohjeet.....	Suomi..... 33
PHIL™ flytande embolisystem	
Bruksanvisning.....	Svenska 37
PHIL™ flytende embolisk system	
Bruksanvisning.....	Norsk 41
Σύστημα υγρού εμβολικού παράγοντα PHIL™	
Οδηγίες χρήσης.....	Ελληνικά..... 45
PHIL™ Sivi Emboli Sistemi	
Kullanım Talimatları.....	Türkçe 49
Течна емболизираща система PHIL™	
Инструкции за употреба.....	Български 53

Sustav tekućeg sredstva za embolizaciju PHIL™	
Upute za upotrebu	Hrvatski 57
Tekutý embolický systém PHIL™	
Návod k použití	Česky 61
Vedel embooliasüsteem PHIL™	
Kasutusjuhend	Eesti keel..... 65
PHIL™ folyékony embolizálórendszer	
Használati útmutató	Magyar 69
PHIL™ šķidrā embolizācijas līdzekļa sistēma	
Lietošanas instrukcija.....	Latviešu valodā..... 73
PHIL™ skystų embolinių medžiagų sistema	
Naudojimo instrukcija	Lietuvių k..... 77
Płynny system embolizacyjny PHIL™	
Instrukcja użycia	Polski 81
Sistem de embolizare cu lichid PHIL™	
Instrucțiuni de utilizare	Română..... 85
Жидкая эмболическая система PHIL™	
Инструкция по применению.....	Русский 89
PHIL™ sistem tečnog sredstva za embolizaciju	
Uputstvo za upotrebu	Srpski 93
Tekoči sistem za embolizacijo PHIL™	
Navodila za uporabo.....	Slovenščina..... 97
Tekutý embolizačný systém PHIL™	
Návod na použitie	Slovenčina 101
Система рідинної емболізації PHIL™	
Інструкція з використання	Українська..... 105
PHIL™ النظام الصمي السائل	
تعليمات الاستخدام.....	العربية..... 109

English

PHIL™ Liquid Embolic System

Instructions for Use

Carefully read all instructions prior to use.

DEVICE DESCRIPTION

The PHIL device is a non-adhesive liquid embolic agent comprised of a co-polymer dissolved in DMSO (dimethyl sulfoxide). An iodine component is chemically bonded to the co-polymer to provide a radiopaque element during fluoroscopic visualization. A PHIL Liquid Embolic Starter Kit System consists of a sterile, pre-filled, 1 mL syringe of PHIL liquid embolic, a sterile, pre-filled 1 mL syringe of DMSO, and a Universal microcatheter Adaptor. A Re-PHIL Liquid Embolic System consists of two sterile, pre-filled 1 mL syringes of PHIL liquid embolic and a Universal microcatheter Adaptor. A DMSO compatible delivery microcatheter that is indicated for use in the neurovascular or peripheral vasculature is used to access the embolization target site. The PHIL Liquid Embolic System is available in several product formulations: PHIL 35%, PHIL 30%, PHIL 25%, and low viscosity PHIL LV liquid embolics. PHIL LV and PHIL 25% liquid embolic devices travel more distally and penetrate deeper into the nidus due to their lower viscosities compared to PHIL 30% or 35% liquid embolic devices. It is recommended that PHIL LV liquid embolic is used with a flow-arrest technique, such as a balloon catheter, to increase injection control. Final solidification of PHIL material occurs within three minutes for all product formulations.

The PHIL device that remains in the patient following the procedure consists of:

Table 1 – Quantitative Implant Material Information

Implant Materials		Mass (g)*
Metallic Components	None, PHIL does not contain metallic components	0
Non-metallic components	PHIL Polymer	0.3 – 0.5

* Approximate content

PRINCIPLE OF OPERATION

The PHIL device is delivered by slow, controlled injection through a microcatheter into the vascular malformation under fluoroscopic control. The DMSO solvent dissipates into the blood, causing the copolymer to precipitate *in situ* into a coherent embolus. The PHIL device immediately forms a skin as the polymeric embolus solidifies from the outside to the inside, while traveling more distally in the vascular lesion.

CONTENTS

PHIL Liquid Embolic System Starter Kit:

- One 1.0 mL syringe of PHIL liquid embolic
- One 1.0 mL syringe of DMSO
- One Universal Adaptor

Re-PHIL Liquid Embolic System:

- Two 1.0 mL syringe of PHIL liquid embolic
- One Universal Adaptor

INTENDED PURPOSE

The PHIL device is intended for use in the embolization of lesions in the peripheral and neurovasculature, including arteriovenous malformations and hypervascular tumors.

CONTRAINDICATIONS

The use of the PHIL device is contraindicated when any of the following conditions exist:

- When patient has severe iodine allergy.
- When optimal microcatheter placement is not possible.
- When provocative testing indicates intolerance to the occlusion procedure.
- When vasospasm stops blood flow.
- Not for use with premature infants (<1,500 g) or individuals with significant liver and kidney function impairment.

CAUTIONS

- Do not use if packaging is opened or damaged.
- This device is intended for single use only. Do not reuse, reprocess, or re-sterilize. Reuse, reprocessing, or re-sterilization may compromise the integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- After use, dispose in accordance with hospital, administrative and/or local government policy.

WARNINGS

- Performing embolization to occlude blood vessels is a high-risk procedure. This device should be used only by physicians with neuro or peripheral interventional training and a thorough knowledge of the vascular pathology to be treated, vascular architecture, angiographic techniques, and super-selective embolization techniques.
- Vascular malformation embolization may influence or change blood flow patterns, thereby subjecting arteries supplying the vascular malformation or the normal surrounding tissues and perivascular space around the malformation to increased pressures, or by causing an increase in the intra-nidal pressure. These conditions could result in hemorrhagic complications.
- Care must be taken to avoid venous outflow occlusion prior to nidal or arterial feeders as this may also produce malformation hemorrhage. If the PHIL device extravasates outside the vascular space, secondary to vessel wall compromise, a sub-acute inflammatory response to the material may occur in the surrounding vascular space which may lead to tissue damage.
- There may be some topical hypersensitivity and/or release of histamines from Dimethyl Sulfoxide (DMSO).
- Therapeutic embolization should not be performed when high blood flow precludes safe delivery of the embolic agent to prevent non-target embolization.
- **The microcatheter tip should be placed as distal as possible, and as close to the target vascular lesion as possible to prevent any non-target embolization of normal surrounding tissue or cranial nerves.**
- Premature solidification of the PHIL device may occur if the microcatheter or luer hub comes in contact with any solution such as saline, blood or contrast.

- Use only DMSO compatible microcatheters indicated for use in the neuro or peripheral vasculature. Other microcatheters or syringes may not be compatible with DMSO and their use can result in thromboembolic events due to microcatheter degradation.
- Use only the MicroVenton pre-filled syringes to inject DMSO and the PHIL device. Other syringes may not be compatible with DMSO.
- Rapid injection of DMSO into the vasculature space may lead to vasospasm and/or angioneurosis.
- In the event of microcatheter occlusion, excessive force to advance the plunger may result in microcatheter rupture due to over pressurization.
- Do not allow more than **1 cm of the PHIL device** to reflux back over microcatheter tip. Excessive reflux may result in difficult microcatheter removal.
- After using a microcatheter with the PHIL device, do not attempt to clear or inject any material through it. Such attempts may lead to embolus or non-target embolization.
- STOP injection if the PHIL device is not visualized exiting the microcatheter tip. If the microcatheter becomes occluded, over-pressurization and vessel rupture can occur. During injection, continuously verify under fluoroscopy that the PHIL device is exiting the microcatheter tip.
- STOP injection if increased resistance is observed. If increased resistance occurs, determine the cause (e.g. occlusion in microcatheter lumen) and replace the microcatheter if needed. Do not attempt to clear or overcome resistance by applying increased injection pressure, as use of excessive pressure may result in microcatheter or vessel rupture and embolization of non-target areas.
- Only use thumb pressure to inject the PHIL device. Using the palm of hand to advance plunger may result in microcatheter or vessel rupture due to over pressurization in the event of microcatheter occlusion.
- DO NOT interrupt the PHIL device injection for longer than three minutes prior to re-establishing injection. This may cause solidification of the PHIL device within the microcatheter tip resulting in microcatheter occlusion. Use of excessive pressure to clear the microcatheter may result in microcatheter or vessel rupture and embolization of non-target areas.
- The PHIL Liquid Embolic System is provided sterile and non-pyrogenic unless the unit package is opened or damaged. Do not use if the packaging is damaged. Use before expiration date noted on the product packaging.
- The PHIL Liquid Embolic System is intended for single use only. The kit is intended to be used for one patient. Do not resterilize and/or reuse the device. Reuse and/or resterilization can increase risk of infection, cause a pyrogenic response or other life threatening complications.
- Caution is to be used for nickel-sensitive and/or pregnant female patients, as the device packaged in glass syringes may contain Nickel.
- Inspect product packaging prior to use. Do not use if sterile barrier is open or damaged.
- Use prior to expiration date.
- Verify that the microcatheters and accessories (see directions for use) used in direct contact with the PHIL polymer are clean and compatible with the material and do not trigger polymerization or degrade with contact. Use only DMSO compatible microcatheters indicated for use in the neuro or peripheral vasculature, and MicroVenton sterile, pre-filled syringes. Other microcatheters or syringes may not be compatible with DMSO and their use can result in thromboembolic events due to microcatheter degradation. Refer to the Warnings and Directions for Use sections.
- Upon completion of the PHIL device injection, wait three minutes, slightly aspirate the syringe, and then gently pull the microcatheter to separate it from the PHIL cast. Failure to wait this recommended time to retrieve the microcatheter after injection may result in fragmentation of the PHIL device into non-target vessels.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to:

- Hematoma at the puncture site
- Non-target arterial thrombosis
- Ischemic events due to embolic migration, vasospasm, thrombosis
- Hemorrhagic accidents: vascular rupture – perforation
- Hemodynamic changes induced by the embolization may result in hemorrhagic complications
- These ischemic or hemorrhagic complications may result in various functional neurological deficits, stroke, and possibly death.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State or Local Health Authority in which the user and/or patient is established.

TRAINING

Serious, including fatal, consequences could result with the use of the PHIL Liquid Embolic System without adequate training. Contact your MicroVenton representative for information on training courses.

HOW SUPPLIED

The PHIL device is provided sterile and non-pyrogenic. The PHIL Liquid Embolic System is available in several product formulations: PHIL 35%, PHIL 30%, PHIL 25%, and PHIL LV. With all formulations, care must be taken prior to injection to identify any normal non-target vessels so that the PHIL device does not occlude normal vessels or vascular territories.

- PHIL 35% liquid embolic is recommended for higher flow scenarios and larger fistulous components.
- PHIL 30% liquid embolic is recommended for moderate to high flow scenarios and when feeding pedicle injections are conducted close to the nidus.
- PHIL 25% liquid embolic is recommended for lower flow scenarios and when more distal travel is needed.
- PHIL LV liquid embolic is recommended for reaching distal vessels and to penetrate deeper into the nidus or malformation. The formulation is less viscous, and the use of flow arrest technique such as balloon microcatheter is recommended to increase delivery control.

Final solidification occurs within three minutes for all product formulations.

PRECAUTIONS

- The safety and effectiveness have not been studied in the following patient populations:
 - Pregnant and nursing women
 - Individuals less than 18 years old
 - Individuals with feeding pedicle aneurysms not associated with the nidus, or distal feeders to the nidus or vascular malformation.
- Some data indicate that dimethyl sulfoxide potentiates other concomitantly administered medications.
- A garlic-like taste may be noted by the patient with use of the PHIL device due to the DMSO component. This taste may last several hours. An odor on the breath and skin may be present.

PACKAGING, STERILIZATION, STORAGE AND OPERATING CONDITIONS

The PHIL liquid embolic syringe, DMSO syringe and universal adaptor are placed inside separate polycarbonate trays and packaged in a unit carton. The devices will remain sterile unless the package is opened, damaged or the expiration date has passed. If the sterile packaging is unintentionally opened or damaged discard the device.

The PHIL liquid embolic syringe, DMSO syringe and universal adaptor are sterilized by steam sterilization. A small round indicator label has been affixed to the packaging of the PHIL liquid embolic syringe tray, DMSO syringe tray and universal adaptor tray. This indicator turns from blue to pink upon exposure to steam sterilization and must be pink in order to use the device. If the indicator is blue, do not use the PHIL embolic device.

Store the PHIL Liquid Embolic System at room temperature in a dry place. Keep away from direct sunlight. If product freezes due to exposure to colder temperatures, thaw at room temperature before use.

DIRECTIONS FOR USE

Verify that adequate sedation is used throughout the embolization procedure. Insufficient sedation may result in patient discomfort or movement. Patient movement during embolic agent injection may result in embolization of a non-target vessel.

1. Confirm microcatheter placement with injection of contrast agent per institutional procedure.

Warning: Use only DMSO compatible microcatheters indicated for use in the neuro or peripheral vasculature, other microcatheters may not be compatible with DMSO and their use can result in thromboembolic events due to microcatheter degradation.

Warning: Use PHIL LV liquid embolic with a flow arrest technique, such as a balloon catheter (ie. Scepter™ Occlusion Balloon catheters) to increase injection control. Use of PHIL LV liquid embolic without a flow arrest technique may result in non-target embolization.

2. Completely clear the microcatheter with approximately 10 mL of saline
3. Holding the syringe upwards, inject DMSO using the pre-filled DMSO syringe into the delivery microcatheter in sufficient volume to fill the microcatheter dead space. Refer to delivery microcatheter labeling for dead space volume.

Warning: Use only the MicroVention pre-filled syringes to inject DMSO and the PHIL device. Other syringes may not be compatible with DMSO.

4. The plunger of the PHIL syringe may stick slightly initially. Pull back slightly on the plunger of the capped PHIL syringe to free up its movement.
5. Remove the cap from the PHIL syringe. Attach the Universal Adaptor to the PHIL syringe and purge air.
6. Remove the DMSO syringe from the microcatheter and overfill and wash the luer hub with the balance of the DMSO.
7. Connect the PHIL syringe and adaptor assembly to the microcatheter hub, making sure there is no air in the hub during the connection. Ensure that the hypotube is seated at the bottom most point of the microcatheter hub, and tighten all connections before use.

Warning: Premature solidification of the PHIL device may occur if microcatheter luer contacts any solution such as saline, blood or contrast.

8. Begin injecting the PHIL device to displace the DMSO.

Warning: Only use thumb pressure to inject the PHIL device. Using the palm of hand to advance plunger may result in microcatheter or vessel rupture due to over pressurization in the event of microcatheter occlusion.

9. Monitor the volume injected to anticipate the first infusion of liquid embolic into the vasculature.

Warning: Do not allow more than 1 cm of the PHIL device to reflux back over microcatheter tip. Excessive reflux may result in difficult microcatheter removal.

10. After using a microcatheter with the PHIL device, do not attempt to clear or inject any material through it. Such attempts may lead to embolus or embolization of a non-target area.

Warning: STOP injection if increased resistance is observed. If increased resistance occurs, determine the cause (e.g. occlusion in microcatheter lumen) and replace the microcatheter if needed. Do not attempt to clear or overcome resistance by applying increased injection pressure, as use of excessive pressure may result in microcatheter or vessel rupture and embolization of non-target areas.

Warning: DO NOT interrupt the PHIL device injection for longer than three minutes prior to re-injection. Solidification of the PHIL device may occur at the microcatheter tip resulting in microcatheter occlusion and use of excessive pressure to clear the microcatheter may result in microcatheter rupture.

11. Upon completion of the PHIL device injection, wait three minutes, slightly aspirate the syringe, and then gently pull the microcatheter to separate it from the PHIL cast.

Difficult microcatheter removal or microcatheter entrapment may be caused by any of the following:

- Angioarchitecture of the vascular lesion, very distal afferent, lengthened, and tortuous pedicle.
- Vasospasm of the feeding pedicle.
- Reflux of the PHIL device over the distal microcatheter tip/shaft.

Should microcatheter removal become difficult, the following will assist in microcatheter retrieval:

- Carefully pull the microcatheter to assess any resistance to removal.
- If resistance is felt, remove any slack or redundancy in the microcatheter.
- Gently apply traction to the microcatheter (approximately 3-4 cm of stretch to the microcatheter).
- Hold this traction for a few seconds and release. Assess traction on vasculature to minimize risk of vessel rupture or hemorrhage. This should be performed under fluoroscopic control.
- This process can be repeated intermittently until microcatheter is retrieved.

For entrapped microcatheters:

- Under some difficult clinical situations, rather than risk rupturing the malformation and consequent hemorrhagic complications by applying too much traction on an entrapped microcatheter, it may be safer to leave the microcatheter in the vascular system.
- This is accomplished by stretching the microcatheter and cutting the shaft near the entry point of vascular access allowing the microcatheter to remain in the artery.
- If the microcatheter breaks during removal, distal migration or coiling of the microcatheter may occur. Same day surgical resection should be considered to minimize the risk of thrombosis.

MR INFORMATION

The PHIL device has been determined to be MR Safe as defined in ASTM F2503-08.

Image Artifact Information

MR testing has demonstrated that the PHIL device is MR Safe and shows no artifacts on MR images in relation to the size and shape of this device.

MATERIALS

The PHIL device is not manufactured with natural rubber latex, polyvinylchloride (PVC), or di 2 ethylhexyl phthalate (DEHP).

WARRANTY

MicroVention warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for particular purpose. Handling, storage, cleaning and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure and other matters beyond MicroVention's control directly affect the device and the results obtained from its use. MicroVention's obligation under this warranty is limited to the repair or replacement of this device through its expiration date and MicroVention shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. MicroVention neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. MicroVention assumes no liability with respect to devices reused, reprocessed or re-sterilized and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.

Prices, specifications and model availability are subject to change without notice.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) for the device will be accessible in the European database on medical devices (EUDAMED: <https://ec.europa.eu/tools/eudamed>), when available.

Permanent implant. Follow-up required at the discretion of the physician.

eIFU website: <https://www.terumoneuro.com/products/product-use-and-safety>

Copyright 2025 MicroVention, Inc. All rights reserved.

MicroVention™, PHIL™, and Scepter™ are trademarks of MicroVention, Inc., registered in the United States and other jurisdictions.

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