# **FRED<sup>TM</sup>** <sup>TM</sup>

Flow Re-direction Endoluminal Device

## **INSTRUCTIONS FOR USE**

<b>FRED™ X™ SYSTEM</b>	
Instructions for Use .....	English .....1
<b>SYSTÈME FRED™ X™</b>	
Mode d'emploi.....	Français .....4
<b>FRED™ X™ SYSTEM</b>	
Gebrauchsanweisung.....	Deutsch .....7
<b>SISTEMA FRED™ X™</b>	
Instrucciones de uso .....	Español.....10
<b>SISTEMA FRED™ X™</b>	
Istruzioni per l'uso .....	Italiano .....13
<b>SISTEMA FRED™ X™</b>	
Instruções de Utilização .....	Português .....16
<b>FRED™ X™-SYSTEM</b>	
Brugsanvisning.....	Dansk .....19
<b>FRED™ X™-SYSTEEM</b>	
Gebruiksaanwijzing .....	Nederlands .....22
<b>FRED™ X™ -JÄRJESTELMÄ</b>	
Käyttöohjeet .....	Suomi .....25
<b>FRED™ X™-SYSTEMET</b>	
Bruksanvisning .....	Svenska .....28
<b>FRED™ X™ SYSTEM</b>	
Bruksanvisning .....	Norsk .....31
<b>ΣΥΣΤΗΜΑ FRED™ X™</b>	
Οδηγίες χρήσης.....	Ελληνικά.....34
<b>FRED™ X™ sistemi</b>	
Kullanım Talimatları.....	Türkçe .....37
<b>СИСТЕМА FRED™ X™</b>	
Инструкции за употреба .....	български .....40
<b>SUSTAV FRED™ X™</b>	
Upute za upotrebu.....	Hrvatski.....43
<b>SYSTÉM FRED™ X™</b>	
Návod k použití.....	Česky .....46

<b>FRED™ X™ SÜSTEEM</b>	
Kasutusjuhend.....	Eesti keel .....49
<b>FRED™ X™ SYSTEM</b>	
Használati útmutató.....	Magyar.....52
<b>FRED™ X™ SISTĒMA</b>	
Lietošanas instrukcija.....	Latviešu valodā.....55
<b>„FRED™ X™“ SISTEMA</b>	
Naudojimo instrukcija.....	Lietuviškai.....58
<b>SYSTEM FRED™ X™</b>	
Instrukcja użycia.....	Polski.....61
<b>SISTEMUL FRED™ X™</b>	
Instrucţiuni de utilizare.....	Română.....64
<b>СИСТЕМА FRED™ X™</b>	
Инструкция по применению.....	Русский.....67
<b>SISTEM FRED™ X™</b>	
Uputstvo za upotrebu.....	Srpski.....70
<b>SISTEM FRED™ X™</b>	
Navodila za uporabo.....	Slovenščina.....73
<b>SYSTÉM FRED™ X™</b>	
Návod na použitie.....	Slovenčina.....76
<b>СИСТЕМА FRED™ X™</b>	
Інструкція з використання.....	Українська.....79
<b>نظام FRED™ X™</b>	
تعليمات الاستخدام.....	العربية.....82

**English**  
**FRED™ System**  
**Instructions for Use**

**DEVICE DESCRIPTION**

The FRED X System includes:

- A self-expanding, single wire braid, closed-cell, surface treated nitinol stent.
- An integrated dual layer construct designed to focus mainly at the neck of an aneurysm.
- Radiopaque markers on its proximal and distal ends as well as interwoven helical marker strands delineating the inner working length of the stent to provide fluoroscopic visibility.
- 7 different implant diameters ranging from 2.5 to 5.5 mm and in various implant lengths ranging from 13 to 45 mm.
- The FRED X 2.5 - 3.0 mm diameter implants are compatible with the Headway™ 21 Microcatheter. The FRED X System 3.5 - 5.5 mm diameter implants are compatible with the Headway 27 Microcatheter.

**CONTENTS**

- The FRED X System is packaged sterile as a single unit with an implant, introducer sheath and a detachable delivery wire.

**INTENDED PURPOSE/ INDICATIONS**

- The FRED X System is intended for endovascular embolization of intracranial neurovascular aneurysms.
- The FRED X System may also be used with embolic coils for the treatment of intracranial neurovascular lesions.


**CONTRAINDICATIONS**

- Patients with known hypersensitivity to nickel-titanium
- Patients in whom angiography demonstrated inappropriate anatomy that does not permit passage or deployment of the FRED X System
- Patient in whom anticoagulant, antiplatelet therapy or thrombolytic drugs are contraindicated


**POTENTIAL COMPLICATIONS**

Possible complications include but are not limited to the following:


- Bleeding or Hemorrhage including intracerebral, retroperitoneal or other locations
- Complications of arterial puncture including pain, local bleeding (hematoma) or injury to the artery or adjacent nerves
- Device migration
- Distal Embolization
- Headache
- Incomplete aneurysm occlusion
- Neurologic deficits including stroke and/or death
- Perforation or dissection of the vessel(s)
- Pseudoaneurysm formation
- Rupture or perforation of aneurysm
- Transient ischemic attack (TIA) or ischemic stroke
- Vasospasm
- Vessel occlusion
- Vessel stenosis or thrombosis


**Figure 2: FRED X Implant Nomenclature**

**Figure 3: FRED X Marker Nomenclature 2.5 - 3.0 mm Outer Diameter Devices**

**Figure 4: FRED X Marker Nomenclature 3.5 - 5.5 mm Outer Diameter Devices**

**Figure 1. FRED X System Setup**

**Table 1: FRED X System Models and Dimensions**

Device	Outer Diameter (mm)	Total Length (mm)	Working Length (mm)
FRED X System	2.5	13 to 30	8 to 26
	3.0	13 to 32	9 to 27
	3.5	13 to 40	7 to 36
	4.0	13 to 44	7 to 38
	4.5	15 to 45	8 to 39
	5.0	15 to 36	9 to 29
	5.5	22 to 32	14 to 26

**Table 2: Qualitative and Quantitative Implant Material Information**

Implant Material		Approximate Mass (mg)
Metallic Components	Nitinol, Tantalum	≤ 49
Non-metallic Components	Adhesive, Polymer	≤ 0.14

The FRED X delivery wire contains a hazardous substance, Cobalt (CAS No. 7440-48-4) is classified as a carcinogenic, mutagenic, or toxic for reproduction (CMR) substance of Category 1A or 1B, or as an endocrine disruptor (ED), and is present at a concentration >0.1% weight/weight.

## WARNINGS

- Should unusual resistance be felt at any time during access or removal, the introducer/guide catheter/microcatheter and FRED X System should be removed as a single unit. Applying excessive force during delivery or retrieval of the FRED X System can potentially result in loss or damage to the device and delivery components.
- The FRED X System should only be used by physicians trained in endovascular interventional neuroangiography, radiology, neurosurgery or interventional neurology for the treatment of intracranial aneurysms or other vascular lesions.
- The FRED X System should only be delivered through the appropriate Headway microcatheter, Headway 27 for FRED X System sizes 3.5 - 5.5 mm and Headway 21 for FRED X System sizes 2.5 - 3.0 mm.
- If repeated friction is encountered during FRED X System delivery, verify microcatheter is not kinked or in extremely tortuous anatomy. Confirm that the microcatheter does not ovalize. Confirm that there is adequate sterile heparinized flush solution.
- Do not reposition the FRED X System in the parent vessel without fully retrieving the device. The FRED X System MUST be retrieved/resheathed into the microcatheter and re-deployed at the desired target location or removed completely from the patient.
- Do not attempt to re-position the FRED X implant after deployment/detachment.
- The safety and performance of the FRED X System for use in conjunction with intravascular medical devices other than neurovascular embolization coils have not been established.

## CAUTIONS

- This product should only be used by experienced physicians who have completed endovascular training in the use of the FRED X System. This device is used for percutaneous neurointerventional and peripheral vascular procedures as indicated by a representative from MicroVenton-Terumo or a MicroVenton-authorized distributor.
- The FRED X System does not contain latex or PVC materials.
- The FRED X System is provided sterile for single use only. Do not reuse, reprocess or sterilize. Reuse, reprocessing or sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Carefully inspect the sterile package and the FRED X System prior to use to verify that neither has been damaged during shipment. Do not use kinked or damaged components, or if the package is opened or damaged.
- See the product label for shelf life. Do not use the FRED X System beyond the labeled use by date.
- Exercise caution when crossing the deployed/detached FRED X System with adjunctive devices such as guidewires, catheters, microcatheters or balloon catheters to avoid disrupting the device geometry and device placement.

## MR INFORMATION

Non-clinical testing has demonstrated that the FRED X System is MR Conditional. A patient with this device can be safely scanned in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2,500-gauss/cm (25-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the FRED X System is expected to produce a maximum temperature rise of 2.8 °C at 1.5-Tesla and 3.6 °C at 3-Tesla after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the device extends approximately 4-mm from the FRED X System when imaged with a gradient echo pulse sequence and a 3-Tesla MR system.

MicroVenton, Inc. recommends that the patient register the MR conditions disclosed in this IFU with the MediciAlert Foundation or equivalent organization. A FRED X System patient implant card is included in the package, which should be completed and provided to the patient.

## CLINICIAN USE INFORMATION

The following accessories are required to use the FRED X System:

- Note: Other accessories for performing a procedure are NOT supplied; should be selected based on the physician's experience and preferences

- Appropriate-sized Guiding catheter for use with selected microcatheter
- Headway 27 microcatheter (3.5 - 5.5 mm diameter devices)
- Headway 21 microcatheter (2.5 - 3.0 mm diameter devices)
- Microcatheter-compatible guidewires
- Saline solution/heparin-saline solution continuous flush set Contrast solution
- Contrast solution
- Rotating Hemostatic Valve (RHV)
- Pressurized sterile Infusion solutions – IV stand
- Femoral arterial sheath, compatible with delivery guide catheter
- Femoral artery access device, sterile needle, guidewire

## PACKAGING AND STORAGE

The FRED X System is placed inside a protective, plastic dispenser coil and packaged in a pouch and unit carton. The FRED X System and dispenser coil will remain sterile unless the package is opened, damaged, or the expiration date has passed. Keep Dry and away from Sunlight.

## SHELF LIFE

See the product label for the device shelf life. Do not use the device beyond the labeled use by date.

## PREPARATION FOR USE

Device and Delivery System Selection

Appropriate selection of the FRED X System is important for patient safety. In order to choose the optimal FRED X System model size for any given lesion, examine pre-treatment angiograms for correct and accurate vessel measurements.

## DIRECTIONS FOR USE

1. Gain vascular access according to standard angiographic practice and perform diagnostic angiogram to document target aneurysm and parent vessel to confirm vessel diameter.
2. Place appropriate size guide catheter according to standard practice. Per physician discretion, if embolization coils are going to be used along with the FRED X System, position a suitable size microcatheter to be utilized for coil embolization coaxially into the target aneurysm. If the FRED X System is going to be used without embolization coils, disregard this step.
3. Coaxially navigate the appropriate Headway microcatheter (Headway 21 for 2.5 - 3.0 mm diameter devices and Headway 27 for 3.5 - 5.5 mm devices) over a guidewire at least 15 mm distal to the aneurysm neck or target location. Remove the guidewire.
4. Maintain flush through the microcatheter(s) per standard endovascular practice.
5. Select an appropriate-sized FRED X device according to the size of the parent vessel/aneurysm neck. Note: The FRED X System implant foreshortens (up to 60%) as it expands to the diameter of the parent vessel. Take implant foreshortening into account when sizing and deploying the FRED X System.
6. Carefully inspect the package for damage to the sterile barrier. Peel open the pouch using aseptic technique and place the dispenser coil into the sterile field.
7. a. Unclip the molded cap attached to the delivery wire from the dispenser coil. Pull on the proximal end of the delivery wire until the introducer exits the dispenser coil. Hold the delivery wire and introducer together while continuing to remove the entire device.  
b. After removal from the dispenser coil, carefully push on the delivery wire and in a bowl of saline, only partially deploy the FRED X implant up to 5 mm or 50% (whatever occurs first, being careful not to detach the implant) from the distal introducer tip.

### Check for the following:

- Implant distal marker uniformity
- Implant distal end shows even displacement with no entanglement
- Implant tracks smoothly through introducer

**Warning: DO NOT FULLY DEPLOY** the FRED X system.

c. With the FRED X implant and introducer sheath positioned and hydrated within the bowl of saline, gently manipulate the FRED X implant within the saline to hydrate the implant and minimize visible air bubbles. Carefully pull back on the delivery wire to fully retrieve the FRED X implant and delivery wire tip within the introducer.

**Warning: DO NOT CONTINUE** if any defect is observed; return the unit to MicroVenton, Inc.

8. Confirm that the device is entirely within the introducer, the tip of the delivery wire is not kinked, and the introducer tip is not damaged. DO NOT CONTINUE if either defect is observed; return the unit to MicroVenton, Inc.
9. Partially insert the distal end of the introducer into the RHV connected to the compatible Headway microcatheter. Tighten the RHV locking ring. Flush the RHV with sterile saline and verify that fluid exits the proximal end of the introducer, hydrating the introducer.  
**Warning: Purge** the FRED X System carefully to avoid the accidental introduction of air into the system
10. Tighten the RHV locking ring and advance the introducer until it is fully engaged with the Headway microcatheter hub, then tighten the RHV locking ring.

**Warning:** Confirm that there are no air bubbles trapped anywhere in the system.  
**Caution:** The introducer must be properly engaged with the microcatheter hub to enable FRED X System to be introduced into the microcatheter.

11. Advance the delivery wire to transfer the FRED X System from within the introducer into the microcatheter.

**Warning:** Do not torque the delivery wire while advancing or retracting the FRED X System.

12. Continue advancing the delivery wire into the microcatheter until the proximal tip of the delivery wire enters the introducer. Loosen the RHV locking ring, remove the introducer, and set it aside.

**Note:** Fluoroscopy may be used up to this point at the physician's discretion.

**Warning:** Do not apply undue force. If resistance is encountered at any point during delivery or manipulation, withdraw the unit and select a new FRED X System.

13. Track the device through the microcatheter to the tip. Carefully advance until the device exit marker on the proximal end of the delivery wire approaches the RHV. At this time, fluoroscopic guidance must be initiated.

14. Position the FRED X System for deployment by aligning the FRED X System implant distal radiopaque end markers past the aneurysm neck allowing for adequate distal and proximal device landing zones.

**Note:** A slow, proper push/pull technique, encompassing sufficient delivery wire push force, in addition to an opposing microcatheter withdrawal force, to remove excess microcatheter slack while maintaining the microcatheter tip within the center of the parent vessel, will facilitate properly deploying the FRED X System at the proper location, to achieve full expansion and good vessel apposition.

**Note:** If applicable, verify microcatheter placed into aneurysm in step 3 is still properly positioned for coil delivery.

**Caution:** Using a rapid microcatheter withdrawal technique to deploy the FRED X System is not recommended and may result in device elongation or improper deployment. When deploying 3.5-5.5 mm diameter devices be aware of the delivery wire tip position.

15. If FRED X system positioning is not satisfactory, the implant may be recaptured and repositioned if it is not fully deployed. The FRED X device may be recaptured up to approximately 75% of its deployed length.

**Caution:** If resistance is felt while recapturing the device, do not continue to recapture. Withdraw the microcatheter slightly to unsheath the device (without exceeding the recapture limit), and then attempt to recapture again.

**Caution:** The FRED X System must not be re-deployed more than three times.

16. If FRED X System positioning is satisfactory, carefully advance the delivery wire while retracting the microcatheter as needed to minimize slack, maintaining the microcatheter around the center of the parent vessel, to allow the implant to deploy across the neck of the aneurysm. Ensure the implant proximal radiopaque end markers are in the advised position (see step 15) proximal to the aneurysm neck for adequate coverage. **Note:** The FRED X System will expand and may foreshorten up to 60% from its undeployed length. Visually verify opening of the proximal end, ensuring that the microcatheter distal tip marker is pulled back, adequately away from the implant proximal end, to allow the proximal end to freely open. Push forward on the delivery wire to assist in maintaining access within the implant as needed.

**Note:** Visualize and refer to implant radiopaque end markers to maintain adequate implant length on each side of the aneurysm neck/target location to ensure appropriate coverage.

**Warning:** Do not detach the FRED X System if it is not properly positioned in the parent vessel.

**Warning:** If applicable, observe the FRED X System marker position during coiling procedure to ensure that the device does not migrate.

17. Prior to removing the delivery wire and if necessary, position the microcatheter distal to the implanted device to maintain access through the implanted device. Remove and discard the delivery wire.

**Caution:** The FRED X System delivery wire should not be utilized as a guidewire. Do not torque the FRED X System. A torque device should not be used.

18. Carefully inspect the deployed FRED X implant under fluoroscopy to confirm that it is completely apposed to the vessel wall and not kinked. If the FRED X implant is not fully apposed or is kinked, consider utilizing a suitable micro guidewire and/or occlusion balloon catheter to fully open the implant.

19. If applicable, detachable coils may be delivered into the aneurysm sac following conventional methods, utilizing the jailed microcatheter from step 3. Verify that the FRED X implant remains patent and properly positioned.

**Note:** The jailed microcatheter should be carefully removed to avoid dislodging the FRED X implant.

20. After completing the procedure, withdraw and discard all applicable accessory devices.

**Caution:** Carefully watch the FRED X implant distal and proximal markers when passing through the implanted device with other devices to avoid displacing the implant.

## WARRANTY DISCLAIMER

MicroVention warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or

otherwise, including, but not limited to, any implied warranties of merchantability or fitness for particular purpose. Handling, storage, cleaning, and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure, and other matters beyond MicroVention's control directly affect the device and the results obtained from its use. MicroVention's obligation under this warranty is limited to the repair or replacement of this device through its expiration date. MicroVention shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. MicroVention neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. MicroVention assumes no liability with respect to devices reused, reprocessed or resterilized and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device. Prices, specifications, and model availability are subject to change without notice.

© Copyright 2025 MicroVention, Inc. All rights reserved.

✕ is a trademark of MicroVention, Inc.

MicroVention™, FRED™, and Headway™ are trademarks of MicroVention, Inc., registered in the United States and other jurisdictions.

All third-party products are trademarks™ or registered® trademarks and remain the property of their respective holders.

## CLINICAL BENEFIT

FRED devices are intended for embolization of intracranial aneurysms as well as treatment of other neurovascular lesions when used with embolic coils. The embolization results in occlusion of intracranial neurovascular aneurysms.

## ADDITIONAL NOTICE TO USER

- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established
- The Summary of Safety and Clinical Performance (SSCP) for the device will be accessible in the European database on medical devices after the launch of the European Database on Medical Devices (EUDAMED) (EUDAMED: <https://ec.europa.eu/tools/eudamed>). The SSCP will be linked to the Basic UDI-DI in the EUDAMED public website.
- Dispose of device in accordance with hospital policy and local regulations for biohazardous waste.
- A patient implant card is included in the package. This card should be completed and provided to the patient.
- Permanent implant. Follow-up required at the discretion of the physician
- The electronic instructions for use (eIFU) is available via MicroVention website: <https://terumoneuro.com/products/product-use-and-safety>

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© Copyright 2025 MicroVention, Inc. Tous droits réservés.
© Copyright 2025 MicroVention, Inc. Alle Rechte vorbehalten.
© Copyright 2025 MicroVentión, Inc. Todos los derechos reservados.
© Copyright 2025 MicroVention, Inc. Tutti i diritti riservati.
© Copyright 2025 MicroVention, Inc. Todos os direitos reservados.
© Copyright 2015 MicroVention, Inc. Alle rettigheder forbeholdes.