SOFIA™ Catheter

Instructions for Use

Carefully read all instructions prior to use.

DEVICE DESCRIPTION

The SOFIA Catheter is a non-tapered, single-lumen, flexible catheter equipped with the coil and the braid reinforcement. The distal segment is steam-shapeable to facilitate vessel selection and also has a hydrophilic coating for navigation through the vasculatures. The radiopaque marker is located at the distal end of the catheter for visualization under fluoroscopy.

CONTENTS

One Catheter
One Introducer Sheath
One Shaping Mandrel

INDICATIONS FOR USE/INTENDED PURPOSE

The SOFIA Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA Catheter can be used to facilitate introduction of diagnostic or therapeutic agents. The SOFIA Catheter is not intended for use in coronary arteries. Moreover, the SOFIA Catheter is intended for use in removal/aspiration of emboli and thrombi from selected blood vessels in the arterial system, including the peripheral and neuro vasculatures.

CONTRAINDICATIONS

There are no known contraindications.

CAUTION

Rx Only: Federal (USA) law restricts this device to sale by or on the order of a physician.

Do not use if pouch is opened or damaged.

This device is intended for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose in accordance with hospital, administrative and/or local government policy.

WARNINGS

The SOFIA Catheter should only be used by physicians who have received appropriate training in interventional techniques.

The SOFIA Catheter is provided sterile and non-pyrogenic. Do not use if the packaging is breached or damaged.

Inspect the SOFIA Catheter prior to use. Do not use the device if any damages or irregularities are observed.

Appropriate anti-coagulation and anti-platelet therapy should be administered per standard medical practice.

The SOFIA Catheter should be manipulated under fluoroscopic guidance. Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.

Do not use the SOFIA Catheter with Ethiodol or Lipiodol contrast media or other such contrast media which includes the components of those agents.

Do not use organic solvents as the device may be damaged.

Do not exceed 2070 kPa (300 psi) maximum recommended infusion pressure. Excess pressure may damage the device or injure the patient. Carefully monitor placement of the distal tip when using a power injector to infuse.

Do not make more than 90 degree angle on the Shaping Mandrel. Steaming of the distal tip with more than 90 degree angle may result in damage to the device.

Do not steam the same device more than once, which may result in damage to the device.

Torquing the SOFIA Catheter excessively while kinked may damage the device resulting in separation of the device. Withdraw the entire device (the device, microcatheter, and guidewire) if the device is severely kinked.

The Introducer Sheath is not intended for use inside the patient body. Ensure that the Introducer Sheath is removed from the SOFIA Catheter once the distal shaft of the SOFIA Catheter is placed inside the patient body.

Excessive aspiration with the distal tip of the SOFIA Catheter covered by the vessel wall may cause vessel injury. Carefully investigate location of the distal tip under fluoroscopy prior to aspiration.

Do not attempt to clear the inner lumen of the SOFIA Catheter by infusion while keeping the device in the patient body.

When flow from the lumen stops or becomes stagnant during aspiration, do not attempt to clear the inner lumen of the SOFIA Catheter by infusion while keeping the device in the patient body.

PRECAUTIONS

Exercise care in handling the SOFIA Catheter to reduce the chance of accidental damage.

Verify compatibility of the SOFIA Catheter when using other ancillary devices commonly used in intravascular procedures. The physician must be familiar with percutaneous, intravascular techniques and possible complications associated with the procedure.

Use caution when manipulating the SOFIA Catheter in tortuous vasculature to avoid damage. Avoid advancing or withdrawal against resistance until the cause of resistance is determined.

Presence of calcifications, irregularities, or other devices may damage the SOFIA Catheter and potentially affect its insertion or removal.

Maintain perfusion of heparinized saline for inner lumen of the SOFIA Catheter to prevent thrombus formation.

If removed from the patient, the hydrophilic coating on the SOFIA Catheter should be hydrated with heparinized saline. Do not allow the coating to dry.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to: vessel or aneurysm perforation, vasospasm, hematoma at the site of entry, embolism, ischemia, intracerebral/intracranial hemorrhage, pseudo aneurysm, seizure, stroke, infection, vessel dissection, thrombus formation, and death.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State or Local Health Authority in which the user and/or patient is established.

COMPATIBILITY

Refer to product label for device dimensions. Use the information on labeling provided with other devices to determine device compatibility. When using the SOFIA Catheter as a single guiding catheter, choose the appropriate size of a femoral sheath, referring to product label for catheter dimensions.

PREPARATION FOR USE

1. Carefully remove the SOFIA Catheter and the Introducer Sheath from the package.
2. Inspect the SOFIA Catheter for any damage.

WARNING: Do not use the device if any damages or irregularities are observed.

If steam shaping is desired, use the technique outlined in step 3. Otherwise proceed to step 4.

3. STEAM SHAPING

- a. Bend the Shaping Mandrel for desired shape.
WARNING: Do not make more than 90 degree angle on the Shaping Mandrel. Steaming of the distal tip with more than 90 degree angle may result in damage to the device.
- b. Carefully insert the Shaping Mandrel into the distal tip of the SOFIA Catheter.
- c. Hold the distal segment together with the Shaping Mandrel and steam it for 30 seconds.
- d. Immediately place the shaped distal segment into heparinized saline to set the shape.
- e. Inspect the distal shaft for any damage.
WARNING: Do not use the device if any damages or irregularities are observed.
- f. Remove the Shaping Mandrel from the SOFIA Catheter.
Do not use the device if any damages or irregularities are observed.
WARNING: Do not steam the same device more than once, which may result in damage to the device.

4. Flush the lumen of the SOFIA Catheter with heparinized saline. Attach a rotating hemostatic valve (RHV) to the proximal hub of the SOFIA Catheter. Set up the line for perfusion of heparinized saline through the sidearm of the RHV.
5. Hydrate the hydrophilic coating on the SOFIA Catheter with heparinized saline before use. Keep the coating hydrated and do not allow the coating to dry.

DELIVERY OF THE SOFIA CATHETER

6. Go to step 7 or 8, depending on the situation described below and choose appropriate devices for navigation of the SOFIA Catheter.
7. **Navigation through the vasculature, except for the intracranial vasculature**
 - a. Prepare 0.035" or 0.038" Guidewire for navigation of the SOFIA Catheter.
 - b. Insert the guidewire into the SOFIA Catheter and advance the Guidewire until the Guidewire and the SOFIA Catheter are aligned at the distal end.
 - c. Using the Introducer Sheath provided in the package, carefully insert the SOFIA Catheter and the Guidewire through a hemostatic valve of the femoral sheath.
 - d. Remove the Introducer Sheath from the SOFIA Catheter once the distal shaft of the SOFIA Catheter is placed inside the patient body.
WARNING: Introducer Sheath is not intended for use inside the patient body.
 - e. Under fluoroscopic guidance, advance or withdraw the SOFIA Catheter over the guidewire until desired position is attained or before the intracranial position is achieved. Select vessels by slowly torquing the SOFIA Catheter if necessary.
WARNING: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
WARNING: Torquing the SOFIA Catheter excessively while kinked may damage the device resulting in separation of the device. Withdraw the entire device (the device, Microcatheter, and Guidewire) if the device is severely kinked.

WARNING: Do not exceed 2070 kPa (300 psi) maximum recommended infusion pressure. Excess pressure may damage the device or injure the patient. Carefully monitor placement of the distal tip when using a power injector to infuse.

Catheter	Approximate Nominal Flow Rates at 100 and 300 psi Infusion Pressure					
	Saline		60% Contrast		76% Contrast	
SOFIA 5F 115cm	100 psi	300 psi	100 psi	300 psi	100 psi	300 psi
	11.5 mL/sec	20.6 mL/sec	9.9 mL/sec	17.9 mL/sec	7.7 mL/sec	16.2 mL/sec
SOFIA 5F 125cm	100 psi	300 psi	100 psi	300 psi	100 psi	300 psi
	11.9 mL/sec	20.7 mL/sec	9.5 mL/sec	17.0 mL/sec	7.7 mL/sec	16.4 mL/sec
SOFIA 6F 115cm	100 psi	300 psi	100 psi	300 psi	100 psi	300 psi
	19.0 mL/sec	37.0 mL/sec	18.0 mL/sec	32.0 mL/sec	16.0 mL/sec	30.0 mL/sec
SOFIA 6F 125cm	100 psi	300 psi	100 psi	300 psi	100 psi	300 psi
	19.0 mL/sec	36.0 mL/sec	18.0 mL/sec	32.0 mL/sec	15.0 mL/sec	30.0 mL/sec
SOFIA 6F 131cm	100 psi	300 psi	100 psi	300 psi	100 psi	300 psi
	19.0 mL/sec	36.0 mL/sec	18.0 mL/sec	32.0 mL/sec	15.0 mL/sec	30.0 mL/sec
SOFIA 6F 135cm	100 psi	300 psi	100 psi	300 psi	100 psi	300 psi
	19.0 mL/sec	36.0 mL/sec	17.0 mL/sec	31.0 mL/sec	15.0 mL/sec	28.0 mL/sec

f. Go to step 8 for navigation through the intracranial vasculatures. Otherwise proceed to step 9.

8. Navigation through the intracranial vasculature

- Prepare Microcatheter and compatible Guidewire for navigation of the SOFIA Catheter.
- Slowly remove, if any, devices previously inserted in the SOFIA Catheter. Insert the Microcatheter with the Guidewire into the SOFIA Catheter.
- Under fluoroscopic guidance, advance or withdraw the SOFIA Catheter over the Microcatheter and the Guidewire until desired position is attained. Select vessels by slowly torquing the SOFIA Catheter if necessary.

WARNING: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.

WARNING: Torquing the SOFIA Catheter excessively while kinked may damage the device resulting in separation of the device. Withdraw the entire device (the device, Microcatheter, and Guidewire) if the device is severely kinked.

WARNING: Do not exceed 2070 kPa (300 psi) maximum recommended infusion pressure. Excess pressure may damage the device or injure the patient. Carefully monitor placement of the distal tip when using a power injector to infuse.


- Slowly remove the Guidewire or the Microcatheter if necessary. Make sure that continuous perfusion of heparinized saline is maintained through the sidearm of the RHV.

NOTE: The Microcatheter used to navigate the SOFIA Catheter may be kept for the rest of procedure.

ASPIRATION THROUGH THE SOFIA CATHETER WITH SYRINGE

- Under fluoroscopic guidance, make sure that the distal tip of the SOFIA Catheter is located in a straight portion of the vessel. Reposition the distal tip if necessary. Refer to **Figure-1** for location of the distal tip.
WARNING: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
WARNING: Torquing the SOFIA Catheter excessively while kinked may damage the device resulting in separation of the device. Withdraw the entire device (the device, Microcatheter, and Guidewire) if the device is severely kinked.

Figure-1


- Make sure the distal tip of the SOFIA Catheter is engaged with the emboli or thrombi under fluoroscopy.
- At the proximal end of the SOFIA Catheter, remove the perfusion line and attach a stopcock with 30 cc or 60 cc locking syringe to the sidearm of the RHV.
- Tighten the RHV.
- With the stopcock closed pull and lock the plunger of the syringe.
- To begin aspiration, open the stopcock.
- Make sure that the syringe aspirates blood, emboli, or thrombi through the system. Close the stopcock if the syringe does not aspirate any blood, emboli, or thrombi, or slow aspiration is observed. Carefully investigate the cause of the restriction and reposition the distal tip if necessary. If blood, emboli, or thrombi is stuck in the SOFIA Catheter remove the entire device (the device, Microcatheter, and Guidewire) and clear the inner lumen. If restriction still remains and the distal tip of the SOFIA Catheter has been repositioned correctly, close the stopcock, re-attach the syringe, and resume aspiration. Maintain aspiration to make sure the emboli or thrombi remains fully engaged with the distal tip of the SOFIA Catheter. With the emboli or thrombi fully engaged pull back the SOFIA Catheter out of the patient body.
WARNING: Excessive aspiration with the distal tip of the SOFIA Catheter covered by the vessel wall may cause vessel injury.
WARNING: Do not attempt to clear inner lumen of the SOFIA Catheter by infusion while keeping the device in the patient body.
- After aspiration is completed, remove the SOFIA Catheter from the patient body. If re-access to the vasculatures with the same device is desired, flush and clean the inner lumen of the device by infusion. Inspect the device for any damage. Follow steps 6 through 9 in the delivery of the SOFIA Catheter section for navigation.
WARNING: Do not use the device if any damages or irregularities are observed.

ASPIRATION THROUGH THE SOFIA CATHETER WITH PUMP

- Under fluoroscopic guidance, make sure that the distal tip of the SOFIA Catheter is located in a straight portion of the vessel. Reposition the distal tip if necessary. Refer to **Figure-1** for location of the distal tip.
WARNING: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
WARNING: Torquing the SOFIA Catheter excessively while kinked may damage the device resulting in

separation of the device. Withdraw the entire device (the device, Microcatheter, and Guidewire) if the device is severely kinked.

Figure-1

2. Make sure the distal tip of the SOFIA Catheter is engaged with the emboli or thrombi under fluoroscopy.
3. Attach aspiration tubing to the aspiration pump and turn on the pump (refer to the directions for use of the aspiration tubing and aspiration pump manual). Confirm that the aspiration gauge reads -20 inHg. Ensure the stopcock on the aspiration tubing is in the closed position.
4. Remove the perfusion line from the side arm of the RHV that is connected to the SOFIA Catheter and attach the aspiration tubing to the side arm of the RHV.
5. Tighten the RHV.
6. Make sure the distal tip of the SOFIA Catheter is engaged with the embolus or thrombus under fluoroscopy.
7. To begin aspiration, turn the aspiration tubing stopcock to the open position and check to see if blood, thrombus, or embolus are aspirated through the system.
8. If after 10 seconds blood is still observed flowing through the system, stop aspiration. To stop aspiration, turn aspiration tubing stopcock to the closed position. Carefully reposition the distal tip of the SOFIA Catheter to engage the thrombus and resume aspiration.
9. If flow is restricted or absent, maintain aspiration to make sure any thrombus or embolus is fully engaged with the distal tip of the SOFIA Catheter. With the thrombus or embolus fully engaged, slowly pull back the SOFIA Catheter and completely withdraw out of the patient.
WARNING: Excessive aspiration with the distal tip of the SOFIA Catheter engaged with vessel wall may cause vessel injury.
10. With the SOFIA Catheter outside of the patient body, flush the catheter to clear the catheter of any thromboembolic material that may be inside.
WARNING: Do not attempt to clear inner lumen of the SOFIA Catheter by infusion while keeping the device in the patient body. Remove the SOFIA Catheter from the patient body before attempting to clear the lumen.
11. If repeated access to the vasculature with the same device is desired, flush and clean the inner lumen of the device by infusion. Inspect the device for any damage. Reintroduce the SOFIA Catheter into the body and follow steps 6 through 9 in the "Delivery of the SOFIA Catheter" section to navigate the catheter to the target site.
WARNING: Do not use the device if any damage or irregularities are observed.
12. Aspiration may be attempted up to 3 times, for a total of 2 minutes of aspiration per attempt with the SOFIA Catheter.
13. Once aspiration is completed, remove the SOFIA Catheter from the patient body.

The physician has the discretion to modify described manipulations of the SOFIA Catheter to accommodate the complexity and variation in procedures. Any technique modification must be consistent with previously described instructions, warnings, precautions and patient safety information.

STORAGE

Avoid exposure to water, sunlight, extreme temperatures and high humidity during storage. Store the SOFIA Catheter under controlled room temperature. See the product label for the device shelf life. Do not use the device beyond the labeled shelf life.

MATERIALS

The SOFIA Catheter is not manufactured with natural rubber latex, polyvinylchloride (PVC), or di-2-ethylhexyl phthalate (DEHP).

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) for the device will be accessible in the European database on medical devices (Eudamed: <https://ec.europa.eu/tools/eudamed>), when available.

SYMBOLS

	Lot Number		Do Not Reuse
	Catalog Number		Medical Device
	Contents		Use-by Date
	Sterilized Using Ethylene Oxide		Date of Manufacture
	CE Mark		Manufacturer
	Consult Instructions for use		Non-pyrogenic
	For Prescription Use Only		Do not use if package is damaged
	Single Sterile barrier system with protective packaging inside		Country of manufacture

WARRANTY

MicroVention warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness. Handling, storage, cleaning and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure and other matters beyond MicroVention's control directly affect the device and the results obtained from its use. MicroVention's obligation under this warranty is limited to the repair or replacement of this device and MicroVention shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. MicroVention neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. MicroVention assumes no liability with respect to devices reused, reprocessed or resterilized and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.

Prices, specifications and model availability are subject to change without notice.

© Copyright 2021 MicroVention, Inc. All rights reserved.

SOFIA™ is a registered trademark of MicroVention, Inc. in the United States and other jurisdictions.

eIFU website: www.microvention/eIFU-MicroVention.com

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© Copyright 2021 MicroVention, Inc. Tous droits réservés.
© Copyright 2021 MicroVention, Inc. Alle Rechte vorbehalten.
© Copyright 2021 MicroVention, Inc. Todos los derechos reservados.
© Copyright 2021 MicroVention, Inc. Todos os direitos reservados.
© Copyright 2021 MicroVention, Inc. Alle rettigheder forbeholdes.
© Copyright 2021 MicroVention, Inc. Alle rechten voorbehouden.