**Summary of Safety and Clinical Performance
for
ERIC™ Retrieval Device**

SSCPPT1101204

MicroVention, Inc.
35 Enterprise
Aliso Viejo, California 92656, USA

DOCUMENT CHANGE HISTORY

SSCP Revision	Change Description	NB approved/verified
A	Initial Release	☐ Yes, Validation language: ☒ No*

*Annual entries must be included. If a revision is not required, an entry stating such must be added.

Sign Page

Document Author:

Signature:

Date:

RA Approver:

Signature:

Date:

Legal Approver:

Signature

Date:

TABLE OF CONTENTS

1	SUMMARY OF SAFETY AND CLINICAL PERFORMANCE [PATIENT VERSION].....	5
1.1	Device Identification and General Information	5
1.2	Intended Use of the Device	6
1.3	Device Description.....	6
1.4	Risks and Warnings	6
1.5	Summary of Clinical Evaluation and Post-Market Clinical Follow-up.....	8
1.6	Possible Diagnostic or Therapeutic Alternatives	8
1.7	Suggested Training for Patients.....	8

LIST OF TABLES

Table 1.1: Device Identification and General Information.....	5
Table 1.2: Intended Use	6
Table 1.3: Device Description	6

1 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE [PATIENT VERSION]

Document Revision: A
Date Issued: 05/23/2025

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of its safety and clinical performance prepared for healthcare professionals is found in the first part of this document.

The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your healthcare professional in case you have questions about your medical condition or about the use of the device in your situation. This SSCP is not intended to replace an Implant card or the Instructions For Use to provide information on the safe use of the device.

1.1 Device Identification and General Information

Table 1.1: Device Identification and General Information

Device Names	
Device Trade Name	ERIC Retrieval Device
Device Class	III
Basic UDI-DI	37015174ERICHX
Year when first certificate (CE) was issued	2013
Legal Manufacturer	
Name & Address	MicroVention, Inc. 35 Enterprise Aliso Viejo, California, 92656 USA
Authorized Representative	
Name & Address	MicroVention Europe SARL 30 bis, rue du Vieil Abreuvoir 78100 Saint-Germain-en-Laye, France
Notified Body	
Name & Address	DQS Medizinprodukte GmbH August-Schanz-Straße 21 D-60433 Frankfurt am Main Germany

1.2 Intended Use of the Device

Table 1.2: Intended Use

Intended Use	
Intended Purpose	This device is used to restore blood flow in the blood vessels that supply to brain.
Indications for Use	This device is used to restore blood flow in the blood vessels that supply to brain.
Intended Patient Group(s)	These implants are used in patients who suffered from sudden blood cut off to a part of brain.
Contraindications and/or Limitations	• Patients with known hypersensitivity to nickel-titanium. • Patients with abnormal narrowing of blood vessels close to the blood clot site that may prevent safe removal of the device. • Patients with image evidence of tear in the carotid artery in the neck.

1.3 Device Description

Table 1.3: Device Description

Device Description	
Description of the Device	This device is used to restore the blood supply by removing the solid mass that forms in the blood vessel in brain. This device consists of retrieval segments attached to a pusher wire and is delivered through a small catheter to the block site.
Materials or substances in contact with the patient's tissues	The device was made of Nickle and Titanium. These materials are safe to use and have been used on humans for a long time.
Information about medicinal substances in the device	The device does not contain any medicinal substances.
Description of how device achieves its intended mode of action	This device helps remove the solid blood mass and restore the blood flow.
Description of Accessories	The device is introduced into the patient body using a small catheter.
Description of other Devices or Products intended to be used in combination	None

1.4 Risks and Warnings

Contact your healthcare professional if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

- **How potential risks have been controlled or managed**

The company uses a standard process to predict device risk. The process gives a thorough estimate of what might happen when the device is used. A list of possible harms and what might cause them is compiled. The Instructions of Use also describes any warning or precautions. These may be associated with remaining risk. Reports of any harms are tracked in a complaint database. Reports can come from users or other health care professionals or from publications. The database is assessed on a regular basis. Harms are investigated if rates are increasing or at a certain level. Actions such as changes in labelling or recall can be taken if needed.

- **Remaining risks and undesirable effects**

There is always a risk of undesirable side effects when you have any type of surgery. It can be difficult to know the exact source of some side effects.

The device Instructions of Use identified the following known potential side effects:

- Blood vessel rupture or Blood vessel bulge rupture
- Narrowing of the blood vessels
- Bad Bruise at the entry site of the device
- Blockage in the blood vessel caused by an air bubble, piece of tissue or blood clot
- Lack of blood and oxygen supply to part of the brain
- Bleeding in the brain
- False bulge in the brain vessels
- Seizure
- Stroke
- Infection
- Tear on the inside of the blood vessel
- Reduced blood flow through a vessel caused by clot
- Death

Undesirable side effects have been reported after the use of the device. The events could be related to the device. They could also be related to the surgery or to the patient. The rate of the events we know about are listed below. Some events are published in medical journals. Some events are reported to the company that makes the device. The benefits of the device are greater than the identified risks. You should discuss these events with your doctor if you have questions.

Reported Undesirable Effects	Typical Reported Rates or Range of Rates (%)
Narrowing of the Blood vessels	1.1
Blood vessel rupture or Blood vessel bulge rupture	0.0 - 4.4
Blockage in the blood vessel caused by an air bubble, piece of tissue or blood clot	6.5
Bleeding in the brain	0.0 - 5.6
Stroke	2.1
Lack of blood and oxygen supply to part of the brain	0.01
Death	2.3 - 30.0

- **Warnings and precautions**

There are some warnings and precautions related to the use of ERIC Retrieval device. All of these relate to the surgical steps performed by the surgeon. Your surgeon will let you know of any actions you must take before surgery.

- **Summary of any field safety corrective action, (FSCA including FSN) if applicable**

FSCA is an action taken by manufacturer to reduce the risk of death or serious health issues associated with the use of device. There was one FSCA action for incorrect expiration date. It is in the process of closing. Two Corrective and Preventive Actions (CAPA) were identified. One is about Expired Materials were not identified and removed from the place of manufacturing and another one is about wrong shelf lifetime. One CAPA is closed, and the other one is in the process of closing.

1.5 Summary of Clinical Evaluation and Post-Market Clinical Follow-up

- **Clinical background of the device**

These devices were first placed in the market in 2013. There is sufficient data to confirm the safety and performance of the device.

- **The clinical evidence for the CE-marking**

A thorough search of journals was done to find studies where the ERIC Retrieval device was used. The date range of the search was 01 April 2013 to 17 February 2025. The search found 15 publications using the device in 672 patients. Several outcomes were studied to evaluate how well the device work. Outcomes measured how well the device cleared the blood clot and restored the blood flow to brain.

- **Safety**

When compared ERIC Retrieval device was safe and successful as other similar devices. Rates of undesirable side effects were low and similar. No new or unknown risks were found, Benefit outweigh the risk for the use of implants when used as intended.

1.6 Possible Diagnostic or Therapeutic Alternatives

Talk to your doctor about alternative to this device. They will know what is best for your overall health and condition.

- **General description of therapeutic alternatives**

Treatment other than surgical interventions may be best for some patients. You and your doctor will decide what is right for you. When symptoms are more severe and non-surgical interventions does not work, a mechanical embolectomy may be an option. There are advantages and disadvantages to each treatment option.

1.7 Suggested Training for Patients

This device is not handled directly by the patient. No training is required.

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