Summary of Safety and Clinical Performance
for
Precipitating Hydrophobic Injectable Liquid
PHIL™ System
SSCPPT23-0016

MicroVention, Inc.
35 Enterprise
Aliso Viejo, California 92656, USA

MicroVention Europe SARL
30 bis, rue du Vieil Abrevoir
78100 Saint-Germain-en-Laye, France

DOCUMENT CHANGE HISTORY

SSCP Revision	Change Description	NB approved/verified
A	Initial release	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No* Validation language:

*Annual entries must be included. An entry stating such must be added if a revision is not required.

Sign Page

Document Author:

Signature:

Date:

RA Approver:

Signature:

Date:

Legal Approver:

Signature

Date:

TABLE OF CONTENTS

1	SUMMARY OF SAFETY AND CLINICAL PERFORMANCE [PATIENT VERSION].....	5
1.1	Device Identification and General Information	5
1.2	Intended Use of the Device	6
1.3	Device Description.....	6
1.4	Risks and Warnings	8
1.5	Summary of Clinical Evaluation and Post-Market Clinical Follow-up.....	9
1.6	Possible Diagnostic or Therapeutic Alternatives	10
1.7	Suggested Training for Users.....	10

LIST OF TABLES

Table 2.1	Device Identification and General Information	5
Table 2.2	Intended Use	6
Table 2.3	Device Description	6

1 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE [PATIENT VERSION]

Document Revision:

Date Issued:

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of its safety and clinical performance prepared for healthcare professionals is found in the first part of this document.

The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your healthcare professional in case you have questions about your medical condition or about the use of the device in your situation. This SSCP is not intended to replace an Implant card or the Instructions For Use to provide information on the safe use of the device.

1.1 Device Identification and General Information

Table 1.1 Device Identification and General Information

Device Names	
Device Trade Name	Precipitating Hydrophobic Injectable Liquid (PHIL) Liquid Embolic System
EMDN Code	C010402020302 Embolization Fluids
Device Class	III
Basic UDI-DI	37015174PHILKC
Year when first certificate (CE) was issued for the device	2014
Legal Manufacturer	
Name & Address	MicroVention Europe SARL (MVE) 30 bis, rue du Vieil Abreuvoir 78100 Saint Germain-en-Laye, France
Manufacturer SRN	MVE: FR-MF-000004449
Notified Body	
Name & Address	DQS Medizinprodukte GmbH August-Schanz-Straße 21 D-60433 Frankfurt am Main Germany

1.2 Intended Use of the Device

Table 1.2 Intended Use

Intended Use	
Intended Purpose	The PHIL device is used to block off problem areas in blood vessels in the body and brain, like when there are too many blood vessels or a tumor that has too much blood.
Indications for Use	The PHIL device is used to block off problem areas in blood vessels in the body and brain, like when there are too many blood vessels or a tumor that has too much blood.
Intended Patient Group(s)	Doctors can use the PHIL device to help treat many problems with blood vessels in the brain and body, like tangled blood vessels, small holes between vessels, weak spots in vessels, and tumors that have too much blood.
Contraindications and/or Limitations	<p>Doctors should not use the PHIL device if:</p> <ul style="list-style-type: none"> • The person is very allergic to iodine (a type of medicine). • The tiny tube can't be put in the right spot. • A test shows the person's body can't handle blocking a blood vessel. • Blood can't flow because the vessels are too tight. • The person is a very small baby (under 1,500 grams) or has serious liver or kidney problems.

1.3 Device Description

Table 1.3 Device Description

Device Description					
Description of the Device	<p>PHIL is a special liquid used by doctors to block off problem blood vessels. It is made from soft plastic and has iodine in it so doctors can see it on X-ray pictures.</p> <p>Doctors use a tiny tube (called a microcatheter) to guide the PHIL liquid into the right spot inside the body while watching on a screen. PHIL stays liquid in its special medicine (called DMSO), but when it touches blood, it turns solid. This helps block the blood vessel.</p> <p>How well PHIL works depends on how fast the doctor pushes it in, how fast the blood is moving, and how big the blood vessels are. When PHIL is used, it quickly makes a soft shell on the outside while the inside stays liquid for a little longer. This helps it go deeper into the vessel before it fully hardens.</p> <p>PHIL stays inside the blood vessel for a long time to keep it closed and help the patient.</p> <table border="1" data-bbox="537 1602 1370 1820"> <thead> <tr> <th>Key Functional Element</th> <th>Materials</th> </tr> </thead> <tbody> <tr> <td>Flushing Solution</td> <td>Dimethyl sulfoxide (DMSO)</td> </tr> </tbody> </table>	Key Functional Element	Materials	Flushing Solution	Dimethyl sulfoxide (DMSO)
Key Functional Element	Materials				
Flushing Solution	Dimethyl sulfoxide (DMSO)				

	<p>PHIL Embolic Liquid</p>	<p>PHIL is made from special soft plastics mixed with a liquid called DMSO to help it flow. It also has iodine so doctors can see it on X-rays.</p>
	<p>Syringe</p>	<p>The barrel (the tube part) is made of strong glass called borosilicate glass, which doesn't break easily</p> <p>The Luer (the small metal tip where the needle or tube connects) is made from strong stainless steel that doesn't rust</p> <p>The stopper is made from a soft rubber-like material that helps seal the liquid inside and keeps it from leaking</p> <p>The plunger is made from stainless steel, a strong metal that doesn't rust and helps push the liquid out</p> <p>The cap is made from strong plastic that helps cover and protect the device</p> <p>The finger grip is made from strong plastic (yellow or white) that helps you hold and control the device easily</p>
	<p>Universal Adaptor</p>	<p>The hub is made from strong plastic and stainless steel, which helps connect parts and keeps them sturdy and rust-free</p> <p>The Tuohy Borst body is also made from polypropylene, a strong plastic that helps hold the parts in place</p> <p>The Tuohy Borst cap is made from polypropylene, a strong plastic that helps close and hold parts together</p> <p>The Tuohy gasket is made from silicone, a soft, rubbery material that helps keep things sealed tight</p>
	<p>Tray Packaging</p>	<p>The tray is made from polycarbonate, a tough plastic that doesn't break easily</p> <p>The tray lid is made from Tyvek, a strong, paper-like material that keeps things clean and safe.</p> <p>The insert (the paper piece inside the box) is made from the same strong, white cardboard as the box.</p> <p>The box (carton) is made from strong, white cardboard called bleached paperboard</p>
<p>Materials or substances in contact with the patient's tissues</p>	<p>Two special chemicals are mixed in a liquid called DMSO. One help make things stick, and the other helps hold water</p>	

Information about medicinal substances in the device	The PHIL Embolic System does not have any medicine, animal parts, or blood in it.
Description of how device achieves its intended mode of action	Doctors use a tiny tube to send the PHIL liquid to the problem spot in the blood vessel. They watch where it goes using a special X-ray. The liquid works because it can mix with one kind of liquid (called DMSO), but not with water. When the PHIL liquid touches blood, it starts to turn solid. This helps block the blood vessel in just the right place. How well it works depends on how fast the liquid is given, how fast the blood is moving, and how big the blood vessel is. The outside of the liquid gets hard first, while the inside stays soft to help it go farther if needed.
Description of Accessories	<ol style="list-style-type: none"> 1. DMSO is a special liquid used to mix with PHIL. It helps the PHIL stay in liquid form and has iodine in it so doctors can see it on an X-ray. Doctors use a tiny tube filled with DMSO to guide the PHIL liquid to the problem spot while watching it on a screen. 2. The universal adaptor is a small part that connects the syringe to the tiny tube (microcatheter). It helps the PHIL liquid go in smoothly without getting mixed too much with DMSO, so it works better when it gets to the right spot.
Description of other Devices or Products intended to be used in combination	<ol style="list-style-type: none"> 1. Tiny metal coils can be put in a blood vessel to slow down the blood. This helps doctors add a special liquid to block the blood vessel. 2. Doctors use tiny tubes called microcatheters. Some new ones have tips that can come off, which helps doctors take the tube out more easily after the liquid gets hard. 3. Some tiny tubes have little balloons on them. These helps keep healthy blood vessels safe and let the special liquid go deeper into the problem spot. Special balloon tubes, like the Scepter Mini, can send the liquid into small, far-away blood vessels and help keep it from going the wrong way.

1.4 Risks and Warnings

Contact your healthcare professional if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

- How potential risks have been controlled or managed
The company uses a standard process to predict device risk. The process gives a thorough estimate of what might happen when the device is used. A list of possible harms and what might cause them is compiled. The Instructions of Use also describes any warning or precautions. These may be associated with remaining risk. Reports of any harms are tracked in a complaint database. Reports can come from users or other health care professionals or from publications. The database is assessed on a regular basis. Harms are investigated if rates are increasing or at a certain level. Actions such as changes in labelling or recall can be taken if needed.

- Remaining risks and undesirable effects
Sometimes, problems can happen like bruising, blood clots, or bleeding if a blood vessel breaks. These problems might lead to serious issues like blood moving, a stroke, or even death.
- Warnings and precautions

Warning:

Only trained doctors should use this device because it can be risky and might cause bleeding or other problems if not used carefully. Doctors must use the right tools, go slowly, and watch closely to make sure the liquid goes to the right place and doesn't hurt healthy parts.

Cautions:

This device should only be used one time—don't reuse or try to clean it. Using it again could break it or make someone sick

Precautions:

This device hasn't been tested on kids, pregnant women, or people with some blood vessel problems, and it might cause a garlic taste or smell. Always check the package, use the right tools, and gently remove the tube after waiting 3 minutes to keep everything safe.

- Summary of any field safety corrective action, (FSCA including FSN) if applicable

During this time, 3 reports about PHIL devices were looked at and all were fixed successfully. There were also 12 improvement actions started- 9 are finished, 2 are being worked on, and 1 is still being checked.

1.5 Summary of Clinical Evaluation and Post-Market Clinical Follow-up

- Clinical background of the device
These devices were first placed in the market in 2014.
- The clinical evidence for the CE-marking
In 22 studies with 406 patients, doctors used the PHIL device to treat brain and body blood vessel problems, and almost all treatments worked very well. There were a few problems and 2 patient deaths, but none were caused by the PHIL device, showing that it works safely and effectively.
- Safety
PHIL is as safe as other similar treatment. No new or unknown risks were found, benefits will be more than the risks when used as they should be.

1.6 Possible Diagnostic or Therapeutic Alternatives

When considering alternative treatments, it is recommended to contact your healthcare professional who can take into account your individual situation.

- General description of therapeutic alternatives

There are other liquid materials like PHIL that doctors use to block off blood vessels. These include:

- n-BCA (a kind of medical glue), like TruFill
- Onyx, made by a company called Medtronic

These are used to treat brain blood vessel problems (like bAVMs) and other body blood vessels with too much blood, like AVMs or tumors.

- n-BCA is sticky like glue.
- Onyx and PHIL are not sticky — they're more like rubber.

Doctors today mostly use non-sticky ones like PHIL or Onyx because they are safer. These types are less likely to trap the tube inside and give doctors more time to control how the liquid spreads.

n-BCA gets hard really fast when it touches blood, but this can sometimes cause it to move to the wrong place or form clots.

Onyx can take longer to show up on the X-ray and may cause too much radiation or even cause the tube to get stuck. It also uses DMSO, which can make blood vessels squeeze too much (vasospasm).

All these liquids have been used for many years, and doctors have learned a lot about how well they work and the risks they bring. That's how we know PHIL and other similar agents can help with bAVM treatment.

1.7 Suggested Training for Users

Only specially trained doctors who know a lot about blood vessel problems and how to treat them should use the PHIL liquid. They need to understand how blood vessels work, how to take special X-ray pictures, and how to carefully block tiny blood vessels.

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