SOFIA™ 88 Catheter

INSTRUCTIONS FOR USE

SOFIA™ 88 Catheter	
Instructions for Use	English..... 1
 Cathéter SOFIA™ 88	
Mode d’emploi.....	Français..... 3

English
SOFIA™ 88 Catheter
Instructions for Use

Carefully read all instructions prior to use.

DEVICE DESCRIPTION

The SOFIA 88 Catheter is a non-tapered, single-lumen, flexible catheter equipped with coil and braid reinforcement. The distal segment is designed to facilitate vessel selection with 60 cm of distal-shaft hydrophilic coating for navigation through the vasculature. The radiopaque marker is located at the distal end of the catheter for visualization under fluoroscopy.

CONTENTS

- One Catheter
- One Introducer Sheath

INDICATIONS FOR USE

The SOFIA 88 Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA 88 Catheter can be used to facilitate introduction of diagnostic agents or therapeutic devices. The SOFIA 88 Catheter is not intended for use in coronary arteries.

CONTRAINDICATIONS

There are no known contraindications.

CAUTION

Rx Only: Federal (USA) law restricts this device to sale by or on the order of a physician.

Do not use if pouch is opened or damaged.

This device is intended for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose in accordance with hospital, administrative and/or local government policy.

WARNINGS

The SOFIA 88 Catheter should only be used by physicians who have received appropriate training in interventional techniques.

The SOFIA 88 Catheter is provided sterile and non-pyrogenic. Do not use if the packaging is breached or damaged.

Inspect the SOFIA 88 Catheter prior to use. Do not use the device if any damages or irregularities are observed.

Appropriate anti-coagulation and anti-platelet therapy should be administered per standard medical practice.

The SOFIA 88 Catheter should be manipulated under fluoroscopic guidance. Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.

Limit the exposure to X-ray radiation doses to patients and physicians by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors when possible.

Do not use the SOFIA 88 Catheter with Ethiodol or Lipiodol contrast media or other such contrast media which includes the components of those agents.

Do not use organic solvents as the device may be damaged.

Do not use the SOFIA 88 Catheter for delivery of liquid embolic agents, including those containing dimethyl sulfoxide (DMSO) or n-butyl cyanoacrylate (n-BCA).

Do not use automated high-pressure contrast injection equipment with the SOFIA 88 Catheter as it may damage the device.

Do not exceed 2070 kPa (300 psi) maximum recommended infusion pressure. Excess pressure may damage the device or injure the patient.

Torquing the SOFIA 88 Catheter excessively while kinked may damage the device, resulting in separation of the device. Withdraw the entire device (the device, compatible catheter, and guidewire) if the device is severely kinked.

The Introducer Sheath is not intended for use inside the patient body. Ensure that the Introducer Sheath is removed from the SOFIA 88 Catheter once the distal shaft of the SOFIA 88 Catheter is placed inside the patient body.

Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

PRECAUTIONS

Exercise care in handling the SOFIA 88 Catheter to reduce the chance of accidental damage.

Verify compatibility of the SOFIA 88 Catheter when using other ancillary devices commonly used in intravascular procedures. The physician must be familiar with percutaneous, intravascular techniques and possible complications associated with the procedure.

Use caution when manipulating the SOFIA 88 Catheter in tortuous vasculature to avoid damage. Avoid advancing or withdrawal against resistance until the cause of resistance is determined.

Presence of calcifications, irregularities, or other devices may damage the SOFIA 88 Catheter and potentially affect its insertion or removal.

Maintain perfusion of heparinized saline for inner lumen of the SOFIA 88 Catheter to prevent thrombus formation.

If removed from the patient, the hydrophilic coating on the SOFIA 88 Catheter should be hydrated with heparinized saline. Do not allow the coating to dry as this may impact the coating safety and performance.

Avoid pre-soaking devices for long durations when the device is not in use as this may impact the coating safety and performance.

Avoid wiping the device with dry gauze as this may damage the device coating. Avoid excessive wiping of the coated device.

POTENTIAL COMPLICATIONS

Potential complications include but are not limited to: vessel or aneurysm perforation, vasospasm, hematoma at the site of entry, embolism, ischemia, intracerebral/intracranial hemorrhage, pseudo aneurysm, seizure, stroke, infection, vessel dissection, thrombus formation, and death.

Potential complications related to angiographic and fluoroscopic X-ray radiation doses include, but are not limited to, alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia. The probability of occurrence of complications may increase as procedure time and number of procedures increase.

COMPATIBILITY

Refer to product label or below diagram for device dimensions. Use the information on labeling provided with other devices to determine device compatibility. When using the SOFIA 88 Catheter as a single guiding catheter, choose the appropriate size of a femoral sheath, referring to product label.

Guide Sheath Recommendation*	≥ 0.113 in. ID
Catheter/Microcatheter Recommendation	≤ 0.0835 in. OD
Guidewire Recommendation	≤ 0.038 in. OD

* Guide sheath/catheter compatibility recommendation applies only to MicroVention/Terumo brand catheter/sheaths.


DIAGRAM OF SOFIA 88 CATHETER

* 120cm length is not cleared or marketed in the United States.

PREPARATION FOR USE

- Carefully remove the SOFIA 88 Catheter and the Introducer Sheath from the package.
- Inspect the SOFIA 88 Catheter for any damage.
WARNING: Do not use the device if any damages or irregularities are observed.
- Flush the lumen of the SOFIA 88 Catheter with heparinized saline. Attach a rotating hemostatic valve (RHV) to the proximal hub of the SOFIA 88 Catheter. Set up the line for perfusion of heparinized saline through the sidearm of the RHV.
- Hydrate the hydrophilic coating on the SOFIA 88 Catheter with heparinized saline before use. Allow at least 30 seconds for hydration. Keep the coating hydrated and do not allow the coating to dry as this may impact the coating safety and performance.

DELIVERY OF THE SOFIA 88 CATHETER

5. Go to step 6 or 7, depending on the situation described below and choose appropriate devices for navigation of the SOFIA 88 Catheter. Refer to the Compatibility section for compatible devices.
6. **Navigation through the vasculature, except for the intracranial vasculature**
- a. Prepare Compatible Catheter/Microcatheter and Guidewire for navigation of the SOFIA 88 Catheter.
 - b. Insert the guidewire into the Compatible Catheter. Insert the Compatible Catheter with Guidewire into the SOFIA 88 Catheter and advance the Compatible Catheter with Guidewire until the Compatible Catheter and the SOFIA 88 Catheter are aligned at the distal end.
 - c. Carefully insert the SOFIA 88 Catheter and the Compatible Catheter with Guidewire through a hemostatic valve of the femoral sheath. If necessary, use the Introducer Sheath provided in the package to facilitate with the introduction.
 - d. Remove the Introducer Sheath, if used, from the SOFIA 88 Catheter once the distal shaft of the SOFIA 88 Catheter is placed inside the patient body.
WARNING: Introducer Sheath is not intended for use inside the patient body.
 - e. Under fluoroscopic guidance, advance or withdraw the SOFIA 88 Catheter over the Compatible Catheter with Guidewire until desired position is attained or before the intracranial position is achieved. Select vessels by slowly torquing the SOFIA 88 Catheter if necessary.
WARNING: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
WARNING: Torquing the SOFIA 88 Catheter excessively while kinked may damage the device resulting in separation of the device. Withdraw the entire device (the device, Compatible Catheter, and Guidewire) if the device is severely kinked.
WARNING: Do not use automated high-pressure contrast injection equipment with the SOFIA 88 Catheter as it may damage the device. Do not exceed 2070 kPa (300 psi) maximum recommended infusion pressure. Excess pressure may damage the device or injure the patient.
 - f. Go to step 7 for navigation through the intracranial vasculatures. Otherwise, proceed to step 8.
7. **Navigation through the intracranial vasculature**
- a. Prepare Compatible Catheter/Microcatheter and Guidewire for navigation of the SOFIA 88 Catheter.
 - b. Slowly remove, if any, devices previously inserted in the SOFIA 88 Catheter. Insert the Compatible Catheter with the Guidewire into the SOFIA 88 Catheter.
 - c. Under fluoroscopic guidance, advance or withdraw the SOFIA 88 Catheter over the Compatible Catheter and the Guidewire until desired position is attained. Select vessels by slowly torquing the SOFIA 88 Catheter if necessary.
WARNING: Do not advance or withdraw the device when excessive resistance is met until the cause of resistance is determined.
WARNING: Torquing the SOFIA 88 Catheter excessively while kinked may damage the device resulting in separation of the device. Withdraw the entire device (the device, Compatible Catheter, and Guidewire) if the device is severely kinked.
WARNING: Do not use automated high-pressure contrast injection equipment with the SOFIA 88 Catheter as it may damage the device. Do not exceed 2070 kPa (300 psi) maximum recommended infusion pressure. Excess pressure may damage the device or injure the patient.
8. Slowly remove the Compatible Catheter or Guidewire if necessary. Make sure that continuous perfusion of heparinized saline is maintained through the sidearm of the RHV.
NOTE: The Compatible Catheter used to navigate the SOFIA 88 Catheter may be kept for the rest of procedure.

The physician has the discretion to modify described manipulations of the SOFIA 88 Catheter to accommodate the complexity and variation in procedures. Any technique modification must be consistent with previously described instructions, warnings, precautions and patient safety information.

STORAGE

Avoid exposure to water, sunlight, extreme temperatures and high humidity during storage. Store the SOFIA 88 Catheter under controlled room temperature. See the product label for the device shelf life. Do not use the device beyond the labeled shelf life.

MATERIALS

The SOFIA 88 Catheter is not manufactured with natural rubber latex, polyvinylchloride (PVC), or di-2-ethylhexyl phthalate (DEHP).

WARRANTY

MicroVention, Inc. warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness. Handling, storage, cleaning and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure and other matters beyond MicroVention's control directly affect the device and the results obtained from its use. MicroVention's obligation under this warranty is limited to the repair or replacement of this device through its expiration date. MicroVention shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. MicroVention neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. MicroVention assumes no liability with respect to devices reused, reprocessed or resterilized and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.

Prices, specifications and model availability are subject to change without notice.

© Copyright 2024 MicroVention, Inc. All rights reserved.

MicroVention™ and SOFIA™ are trademarks of MicroVention, Inc., registered in the United States and other jurisdictions.

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© Copyright 2024 MicroVention, Inc. Tous droits réservés.
MicroVention™ et SOFIA™ sont des marques de MicroVention, Inc., déposées aux États-Unis et dans d'autres juridictions.