ERICTM

Retrieval Device

INSTRUCTIONS FOR USE

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English ERIC™ Retrieval Device Instructions for Use

Carefully read all instructions prior to use.

DEVICE DESCRIPTION

The ERIC Retrieval Device is a resheathable mechanical thrombectomy device for restoring blood flow by removing clots from occluded vasculature. The ERIC Retrieval Device consists of retrieval segments attached to a pusher wire and is delivered through a Microcatheter to the occlusion site. Radiopaque markers provide visualization under fluoroscopic guidance to aid navigation. The device exit markers on the proximal end of the pusher wire indicate the initiation of fluoroscopic guidance.

Reference Appendix for Design Diagrams

CONTENTS

One Retrieval Device

INTENDED PURPOSE

The ERIC Retrieval Device is intended for use in the revascularization of acute ischemic stroke caused by intracranial occlusive vessels in patients.

CONTRAINDICATIONS

- Patients with known hypersensitivity to nickel-titanium.
- Patients with stenosis proximal to the thrombus site that may prevent safe recovery of the ERIC Retrieval Device.
- Patients with angiographic evidence of carotid dissection.

CAUTION

- Do not use if pouch is opened or damaged. Visually inspect all sterile barrier systems that are labeled as sterile, immediately prior to use in order to determine if breaches in sterile barrier system integrity are evident.
- This device is intended for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- After use, dispose in accordance with hospital, administrative and/or local government policy.

WARNINGS

- The ERIC Retrieval Device should only be used by physicians who have received appropriate training in interventional techniques.
- The ERIC Retrieval Device is provided sterile and non-pyrogenic. Do not use if the packaging is breached or damaged.
- Appropriate anti-coagulation and anti-platelet therapy should be administered per standard medical practice.
- Do not advance or withdraw the ERIC Retrieval Device when excessive resistance is observed. Assess the source of resistance using fluoroscopic means. If needed, reseat the ERIC Retrieval Device into Microcatheter and remove the entire system under aspiration. If resistance is encountered during reseat, stop reseat and remove the entire system under aspiration.
- Position the distal tip marker of Microcatheter just proximal to the retrieval segments after deployment of the ERIC Retrieval Device and

maintain the distal tip marker in the same position during withdrawal to reduce risk of device fracture.

- Do not perform more than three (3) retrieval attempts per device.
- Do not apply excessive force to ERIC Retrieval Device when cleaning the device for additional retrieval attempts.
- Ensure during cleaning of the ERIC Retrieval Device that all foreign components (clot, fibers, etc) are fully removed before reinsertion.
- Users and /or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State or Local Health Authority in which the user and /or patient is established.

PRECAUTIONS

- Exercise care in handling the ERIC Retrieval Device to reduce the chance of accidental damage.
- Use of organic solvents may damage the ERIC Retrieval Device.
- Verify that the inner lumen of Microcatheter is compatible with the ERIC Retrieval Device prior to use.
- Verify ERIC Retrieval Device compatibility when using other ancillary devices commonly used in intravascular procedures. Physician must be familiar with percutaneous, intravascular techniques and possible complications associated with the procedure.
- Use caution when manipulating the ERIC Retrieval Device in tortuous vasculature to avoid damage to the vasculature or the device. Avoid advancing or withdrawing against resistance until the cause of resistance is determined.
- Presence of calcifications, irregularities, or other devices may damage the ERIC Retrieval Device and potentially affect its insertion or removal.
- Maintain saline perfusion between the ERIC Retrieval Device Microcatheter and Guiding Catheter to prevent thrombus formation.

POTENTIAL COMPLICATIONS

Potential complications include but are not limited to vessel or aneurysm perforation, vasospasm, hematoma at the site of entry, embolism, ischemia, intracerebral/intracranial hemorrhage, pseudo aneurysm, seizure, stroke, infection, vessel dissection, thrombus formation, and death.

COMPATIBILITY

Refer to Table 2 for product compatibility with Microcatheters.

All ERIC Retrieval Devices are compatible with a Microcatheter with an inner diameter of 0.017 inch or larger.

Table 2: Recommended Microcatheter Sizes

ERIC Retrieval Device Diameter	Minimum Microcatheter Inner Diameter (inch)
3mm	0.017
4mm	0.017
6mm	0.017

PREPARATION FOR USE

1. Administer anti-coagulation and anti-platelet medications per standard medical practice.
2. By performing angiography, determine the location and size of the target anatomy to be revascularized.
3. Select suitable Guiding Catheter for thrombus removal and navigate it to an appropriate place. Attach the rotating hemostatic valve (RHV) to proximal end of the Guiding Catheter and connect

Table 1 – ERIC Retrieval Device: Product Specifications and Sizing Recommendations

Size (mm)	Recommended Microcatheter ID (inch)	Retrieval Segment Diameter (mm)	Retrieval Length (mm)	Additional Radiopaque Markers		Pusher Wire Length (cm)	Length from Distal Tip to Fluoro safe Marker (cm)
	Min.			Distal	Prox.		Max.
3x15	.017	3.0	15	1	1	203	145
3x20	.017	3.0	20				
4x24	.017	4.0	24				
4x30	.017	4.0	30				
6x35	.017	6.0	35				
6x44	.017	6.0	44				

- saline perfusion line to the side port of the RHV. Keep maintaining saline perfusion for the rest of procedure.
- Choose the suitable ERIC Retrieval Device and remove the device by pulling it from the dispenser tube. Select appropriate Microcatheter.
 - Connect the RHV to proximal end of the Microcatheter. Attach saline perfusion line to the side port of the RHV. Keep maintaining saline perfusion for the rest of procedure.
 - Navigate the Microcatheter using a suitable Guidewire. Under Fluoroscopic guidance, advance the Microcatheter through the thrombus and tighten the RHV to fix the Microcatheter to the Guiding Catheter.

DELIVERY OF THE ERIC RETRIEVAL DEVICE

- Insert the distal end of the Introducer Sheath into the RHV attached to the Microcatheter. Tighten the RHV around the Introducer Sheath and verify saline flushes back from the proximal end of the Introducer Sheath.
- Loosen the RHV and advance the distal end of the Introducer Sheath until the distal end reaches the end of the hub of the Microcatheter. Tighten the RHV around the Introducer Sheath and make sure there is no bubble observed within the entire system.
- Gently push the pusher wire of the ERIC Retrieval Device into the Microcatheter until the distal end of the pusher wire approaches to the proximal end of the Introducer Sheath. Loosen the RHV and remove the Introducer Sheath from the pusher wire and tighten the RHV around the pusher wire.
- Keep pushing the pusher wire until the device exit marker approaches to the RHV. Fluoroscopic guidance must be initiated at this moment.
- Carefully advance the ERIC Retrieval Device until the distal marker of the Microcatheter and the distal radiopaque tip of the ERIC Retrieval Device are aligned.

Note: The distal radiopaque tip of the ERIC Retrieval Device is approximately 1mm.

WARNING: Do not advance the ERIC Retrieval Device when excessive resistance is observed. Assess the source of resistance using fluoroscopic means. The devices may be damaged and cause injury if the ERIC Retrieval Device is advanced.

- Align the proximal marker and the proximal interface of the thrombus to allow the device to engage the thrombus in the proximal segment and extend beyond the clot. Loosen the RHV around the Microcatheter. Withdraw the Microcatheter to deploy the ERIC Retrieval Device while maintaining the retrieval segments in the same location. Refer to Figures 2 & 3 for device deployment.
- Once the ERIC Retrieval Device is fully deployed, place the distal tip of the Microcatheter close to the proximal end of the retrieval segments. Tighten the RHV around the Microcatheter to fix the system. Refer to Figure-2 for device deployment.

RETRACTION OF THE ERIC RETRIEVAL DEVICE

- Make sure that the distal tip of the Microcatheter is positioned just proximal to the retrieval segments.

WARNING: Maintain the distal tip marker in the same position during withdrawal to reduce risk of device fracture.
- Loosen the RHV around the Microcatheter just enough to allow retraction of the Microcatheter while maintaining seal for aspiration.
- Replace saline perfusion line with 60 cc syringe to the side port of the RHV at the proximal end of the Guiding Catheter. Make sure that the RHV maintains seal within the Guiding Catheter by applying aspiration. Adjust the RHV if necessary.
- Start aspiration to the Guiding Catheter using a 60 cc syringe and slowly retract the ERIC Retrieval Device together with the Microcatheter to retrieve thrombus. While maintaining aspiration to the Guiding Catheter, continue retracting the ERIC Retrieval Device and the Microcatheter until the retrieval segments arrive in the proximal end of the Guiding Catheter.

WARNING: Do not withdraw the ERIC Retrieval Device when excessive resistance is observed. Assess the source of resistance using fluoroscopic means. If needed resheath the ERIC Retrieval Device with Microcatheter and remove the entire system under aspiration. If resistance is encountered during resheathing, stop resheathing and remove the entire system under aspiration.

WARNING: Do not perform more than three (3) retrieval attempts using the same ERIC Retrieval Device.
- Open the RHV and remove the ERIC Retrieval Device and the Microcatheter from the Guiding Catheter.
- Aspirate the Guiding Catheter to make sure inner lumen of the Guiding Catheter is clear from any thrombus.
- Replace the 60 cc syringe with saline perfusion line to the side port of the RHV at the proximal end of the Guiding Catheter.

- If additional retrieval attempts are desired with the same device, inspect the ERIC Retrieval Device for any damage. Do not use the same device if any damage is observed and use a new device. Clean the ERIC Retrieval Device using saline. Ensure during cleaning of the ERIC Retrieval Device that all foreign components (clot, fibers, etc) are fully removed before reinsertion. Use extra caution when cleaning the ERIC Retrieval Device. Use a Microcatheter to navigate the ERIC Retrieval Device for subsequent retrieval attempts following the same steps from 1 through 21 above.

WARNING: Do not perform more than three (3) retrieval attempts using the same ERIC Retrieval Device.

WARNING: Do not apply excessive force to the ERIC Retrieval Device when cleaning the device for additional retrieval attempts.

STORAGE

Keep dry and away from sunlight. See the product label for the device shelf life. Do not use the device beyond the labeled shelf life.

MATERIALS

The ERIC Retrieval Device is not manufactured with natural rubber latex, synthetic latex, or polyvinylchloride (PVC).

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) for the device will be accessible in the European database on medical devices (EUDAMED: <https://ec.europa.eu/tools/eudamed>), when available.

WARRANTY

MicroVention, Inc. warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness. Handling, storage, cleaning and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure and other matters beyond MicroVention's control directly affect the device and the results obtained from its use. MicroVention's obligation under this warranty is limited to the repair or replacement of this device through its expiration date. MicroVention shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. MicroVention neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. MicroVention assumes no liability with respect to devices reused, reprocessed or resterilized and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.

Prices, specifications and model availability are subject to change without notice.

© Copyright 2025 MicroVention, Inc. All rights reserved.

MicroVention™ and ERIC™ are trademarks of MicroVention, Inc., registered in the United States and other jurisdictions.

All third-party products are trademarks™ or registered® trademarks and remain the property of their respective holders.

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© Copyright 2025 MicroVention, Inc. Tous droits réservés.
© Copyright 2025 MicroVention, Inc. Alle Rechte vorbehalten.
© Copyright 2025 MicroVention, Inc. Todos los derechos reservados. MicroVention™ y ERIC™ son marcas comerciales de MicroVention, Inc. registradas en los Estados Unidos y en otras jurisdicciones.
© Copyright 2025 MicroVenton, Inc. Tutti i diritti riservati.
© Copyright 2025 MicroVention, Inc. Todos os direitos reservados.
© Copyright 2025 MicroVention, Inc. Alle rettigheder forbeholdes.