VIA™

Microcatheter

INSTRUCTIONS FOR USE

VIA™ Microcatheter		
Instructions for Use	English	1
Microcathéter VIA™		
Mode d'emploi.....	French / Français	3
VIA™ Microcatheter		
Gebrauchsanweisung.....	German / Deutsch	5
Microcatéter VIA™		
Instrucciones de uso	Spanish-Europe / Español.....	7
Microcatetere VIA™		
Istruzioni per l'uso	Italian / Italiano	9
Microcateter VIA™		
Instruções de utilização.....	Portuguese-Europe / Português.....	11
VIA™-mikrokateter		
Brugsanvisning.....	Danish / Dansk	13
VIA™-mikrokatheter		
Gebbruksaanwijzing	Dutch / Nederlands.....	15
VIA™ Microcatheter		
Käyttöohjeet	Finnish / Suomi	17
VIA™ Microcatheter		
Bruksanvisning	Swedish / Svenska	19
VIA™ Mikrokateter		
Bruksanvisning	Norwegian / Norsk.....	21
Μικροκαθετήρας VIA™		
Οδηγίες χρήσης	Greek / Ελληνικά	23
VIA™ Mikrokateter		
Kullanma Talimatları.....	Turkish / Türkçe	25
VIA™ 微导管		
使用说明	Chinese / 简体中文	27
VIA™ 마이크로카테터		
사용 지침	Korean / 한국어	29
VIA™ Microcatheter		
Инструкции за употреба	Bulgarian / Български	31
VIA™ Microcatheter		
Upute za uporabu.....	Croatian / Hrvatski.....	33

VIA™ Microcatheter		
Návod k použití.....	Czech / Čeština	35
VIA™ mikrokateeter		
Kasutusjuhend.....	Estonian / Eesti keel	37
VIA™ mikrokatéter		
Használati utasítás.....	Hungarian / Magyar	39
VIA™ mikrokatetrs		
Lietošanas instrukcija	Latvian / Latviski	41
VIA™ mikrokateris		
Naudojimo instrukcijos	Lithuanian / Lietuvių k.....	43
Mikrocewnik VIA™		
Instrukcja obsługi	Polish / Polski	45
VIA™ Microcatheter		
Instrucțiuni de utilizare.....	Romanian / Română.....	47
Микрокатетер VIA™		
Инструкции по применению	Russian / Русский	49
VIA™ mikrokateeter		
Uputstvo za upotrebu	Serbian / Srpski	51
Mikrokatéter VIA™		
Návod na použitie.....	Slovak / Slovensky	53

English

VIA Microcatheter

Instructions for Use

DEVICE DESCRIPTION

The VIA™ Microcatheter is a single lumen catheter designed to be introduced over a steerable guidewire into the vasculature. The physician inserts the catheter into the vein or artery through the skin (percutaneous) using a sheath or guidewire. The device can then be navigated to the treatment site. The distal 100 cm of the catheter contains a hydrophilic coating on the surface which aids navigation and assists with manipulation while in the vasculature. Throughout the procedure the physician can obtain the position of the catheter by the radiopaque marker bands using fluoroscopic techniques (VIA17 Microcatheter has 2 radiopaque marker bands). Diagnostic and interventional devices can be delivered through the lumen of the catheter to the treatment site. The proximal end of the catheter incorporates a standard luer adapter to facilitate attachment of accessories.

INTENDED USE / INDICATIONS FOR USE

VIA 17, 17 Preshaped 45°, 17 Preshaped 90° - The VIA Microcatheter is intended for the introduction of interventional devices (such as WEB device / coils / stents) and infusion of diagnostic agents (such as contrast media) into the neuro, peripheral, and coronary vasculature.

CONTRAINDICATIONS

The VIA Microcatheter is contraindicated for use with liquid embolic materials, such as n-butyl 2-cyanoacrylate or ethylene vinyl alcohol & DMSO (dimethyl sulfoxide).

POTENTIAL COMPLICATIONS

Potential complications include but are not limited to the following: hematoma at the site of entry, aneurysm rupture, emboli, vessel perforation, parent artery occlusion, hemorrhage, ischemia, vasospasm, vascular thrombosis and neurological deficits including stroke and death.


Potential X-ray radiation exposure related adverse events include but are not limited to: alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, delayed neoplasia, tissue necrosis, and risks associated with contrast dye.

COMPATIBILITY

Microcatheter	Recommended Guide Catheters	Compatible Guidewires
VIA 17 PRESHAPED 45°	5.0F/0.055"(1.40MM) ID OR LARGER	0.014"(0.36MM) OD OR SMALLER
VIA 17 PRESHAPED 90°	5.0F/0.055"(1.40MM) ID OR LARGER	0.014"(0.36MM) OD OR SMALLER
VIA 17	5.0F/0.055"(1.40MM) ID OR LARGER	0.014"(0.36MM) OD OR SMALLER

WARNINGS

CAUTION: This device should be used only by physicians trained in percutaneous, intravascular and neurovascular techniques and procedures at medical facilities with the appropriate fluoroscopic equipment.

-  The VIA Microcatheter is provided sterile and non-pyrogenic unless the unit package is opened or damaged.
-  Do not use if the packaging is damaged. Use before expiration date noted on the product packaging.
-  The VIA Microcatheter is intended for single use only.
-  Do not resterilize and/or reuse the device. Reuse and/or resterilization can increase risk of infection, cause a pyrogenic response or other life-threatening complications. Reuse and/or resterilization can degrade product performance, leading to device malfunction. Dispose of all devices in accordance with applicable hospital, administrative and/or local government policy.
- Never advance or withdraw a device against resistance until the cause of the resistance is determined by fluoroscopy. Excessive force against resistance may result in damage to the device or vessel perforation.
- Always monitor infusion rates when using the microcatheter.
- When injecting contrast for angiography, ensure the catheter is not kinked or occluded.
- Do not exceed 300 psi maximum recommended infusion pressure. Excess pressure may result in catheter damage or patient injury.
- Steam shaping may result in improper device delivery and deployment, depending on the degree of shaping and catheter deflection during device delivery.
- Using the VIA Microcatheter with guide catheters smaller than what is recommended (see compatibility table above) may result in damaging the hydrophilic coating.
- The VIA Microcatheter has not been evaluated for use in the pediatric population (<22 yrs of age).

PRECAUTIONS

- Prior to use, examine the VIA Microcatheter to verify that it has not been damaged during shipment.
- The VIA Microcatheter has a lubricious hydrophilic coating on the outside of the catheter. It must be kept hydrated in order to be lubricious. This can be accomplished by attaching a Y-connector to a continuous saline drip.
- High quality, digital subtraction fluoroscopic road mapping, with orthogonal views, is mandatory to achieve correct placement of the microcatheter and embolization device.
- If repositioning is required, take special care to retract or to advance the device under fluoroscopy.
- The operator should be aware that microcatheters, in distal blood vessels, may increase the risk of thromboemboli.
- Limit the exposure to X-ray radiation doses to patients and physicians by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors, when possible. The risk of X-ray radiation exposure complications may increase as procedure time and number of procedures increase.
- If removed from the patient, the hydrophilic coating on the VIA Microcatheter should be hydrated with heparinized saline. Do not allow the coating to dry as this may impact the coating safety and performance.
- Avoid pre-soaking devices for long durations when the device is not in use as this may impact the coating safety and performance.
- Avoid wiping the device with dry gauze as this may damage the device coating. Avoid excessive wiping of the coated device.

ACCESSORIES

- 1 Shaping Mandrel
- 1 Introducer Sheath

PROCEDURE

Catheterization of the Lesion

1. Using standard interventional procedures, access the vessel with a guide catheter. The guide catheter should have an inner diameter large enough to allow for contrast injection while the microcatheter is in place.
2. Attach a rotating hemostatic valve (RHV) to the hub of the guide catheter. Attach a three-way stopcock to the side arm of the RHV and then connect a line for continuous infusion of flush solution.
3. Select the appropriate VIA Microcatheter size depending on the size of the device that will be delivered.
4. Gently remove the VIA Microcatheter dispenser coil or tray from the pouch.
 - a. **For straight tip Microcatheter configuration only (Dispenser Coil):** Flush the dispenser coil with sterile flush solution through the female luer attached to the coil. Once hydrated, do not allow the catheter to dry as this may impact coating safety and performance; place in basin of sterile saline solution if needed.
 - b. **For Preshaped tip Microcatheter configuration only (Tray):** Carefully remove the lid. Carefully remove the Microcatheter tip from the slot of the tip retainer by pinching the two opposite sides and discard the tip retainer. Lift the Microcatheter hub from the tray and remove the Microcatheter. Fully hydrate the Microcatheter with sterile saline solution in basin. Once hydrated, do not allow the catheter to dry as this may impact coating safety and performance; place in basin of sterile saline solution if needed.
5. If desired, carefully introduce the Microcatheter through the introducer sheath.
6. Prepare an appropriately sized guidewire and insert into the microcatheter per the manufacturer's instructions.
7. Introduce the guidewire and microcatheter as a unit into the guiding catheter until the distal tip of the guide catheter has been reached. Alternatively advance the guidewire and microcatheter until the desired site has been reached. Verify catheter position using fluoroscopy. Gently tighten RHV, as needed, without crushing the microcatheter.
8. Peel away introducer sheath by pulling on the tab, if used.
9. After the microcatheter has been positioned inside the lesion, remove the guidewire.

- Flush the ID of the microcatheter with sterile flush solution by attaching a syringe to the catheter hub
- Attach a second RHV to the hub of the microcatheter. Attach a one-way stopcock to the sidearm of the second RHV and connect the flush solution line to the stopcock.
- Open the stopcock to allow flush through the microcatheter with sterile flush solution.

MicroVention, Inc. is a wholly-owned subsidiary of Terumo Americas Holding, Inc. which is a wholly-owned subsidiary of Terumo Corporation.

VIA™ is a registered trademark of Sequent Medical, Inc. in the United States.

Removal of the VIA Microcatheter

- Under fluoroscopic guidance, withdraw the VIA Microcatheter until the entire device has been removed from the patient.

SHAPING MANDREL

WARNING: The steam shaping mandrel is not intended for use in the human body.

Use only a steam source to shape the catheter tip. Do not use other heat sources.

Prior to use, inspect the tip of the catheter for any damage that may have resulted from steam shaping. Do not use a catheter that has been damaged in any way.

VIA 17 Microcatheters are supplied with one straight shaping mandrel.

Microcatheter models with a preshaped tip may be steam shaped if desired.

Follow steps below for using straight shaping mandrel:

- Remove shaping mandrel from card and insert into distal tip of catheter.
- If desired, remove introducer sheath from the card and carefully introduce the Microcatheter through the introducer sheath.
- Carefully bend catheter tip and shaping mandrel into desired shape. A slight over exaggeration may be required to account for catheter relaxation.
- Shape the catheter by holding the shaped portion approximately 1 inch (2.5 cm) from the steam source for approximately 30 seconds. DO NOT EXCEED 30 SECONDS.
- Allow catheter tip to cool in air or saline before removing the mandrel. Remove the mandrel and discard. Multiple shaping is not recommended.
- Inspect the tip for any damage that may have resulted from steam shaping. If any damage is found, do not use the catheter.

FLOW RATES

Flowrates at 100 and 300 psi	Saline		50% Saline/50% Omnipaque 350	
	100 psi (690 kPa)	300 psi (2070 kPa)	100 psi (690 kPa)	300 psi (2070 kPa)
VIA 17 & PRESHAPED	0.5 mL/sec	0.9 mL/sec	0.2 mL/sec	0.5 mL/sec

PACKAGING, STERILIZATION AND STORAGE

The VIA Microcatheter is placed inside a protective dispenser coil (straight tip) or tray (preshaped tip) and packaged in a pouch and unit carton. The device will remain sterile unless the package is opened, damaged or the expiration date has passed.

STERILE EO The VIA Microcatheter is sterilized by ethylene oxide. A small round indicator label has been affixed to the packaging of the VIA Microcatheter. This indicator turns from purple to green upon ethylene oxide sterilization and must be green in order to use the device. If the indicator is purple, do not use the device.

Store at a controlled room temperature in a dry place.

MATERIALS

The VIA Microcatheter is not made with natural rubber, latex, or PVC.

WARRANTY

MicroVention, Inc. warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness. Handling, storage, cleaning and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure and other matters beyond MicroVention's control directly affect the device and the results obtained from its use. MicroVention's obligation under this warranty is limited to the repair or replacement of this device and MicroVention shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. MicroVention neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. MicroVention assumes no liability with respect to devices reused, reprocessed or resterilized and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Prices, specification and model availability are subject to change without notice.

© Copyright 2023 MicroVention, Inc. All rights reserved.

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© Copyright 2023 MicroVention, Inc. Tous droits réservés.
© Copyright 2023 MicroVention, Inc. Alle Rechte vorbehalten.
© Copyright 2023 MicroVention, Inc. Todos los derechos reservados.
© Copyright 2023 MicroVention, Inc. Tutti i diritti riservati.
© Copyright 2023 MicroVention, Inc. Todos os direitos reservados.
© Copyright 2023 MicroVention, Inc. Alle rettigheder forbeholdes.