AZURTM
Vascular Plug

Instructions for Use

AZUR™ Vascular Plug Instructions for Use.....	English.....	1
Bouchon vasculaire AZUR™ Mode d'emploi	Français.....	5
AZUR™ Gefäßverschlusskörper Gebrauchsanweisung	Deutsch	9
Tapón vascular AZUR™ Instrucciones de uso	Español	13
Tappo vascolare AZUR™ Istruzioni per l'uso	Italiano	17
Tampão vascular AZUR™ Instruções de utilização	Português	21
AZUR™ Vascular Plug Brugsanvisning.....	Dansk	25
AZUR™ Vasculaire plug Gebruiksaanwijzing	Nederlands	29
AZUR™ -verisuonitulppa Käyttöohjeet	Suomi	33
AZUR™ -kärlplugg Bruksanvisning	Svenska	37
AZUR™ vaskulær plugg Bruksanvisning	norsk	41
Αγγειακό βύσμα AZUR™ Οδηγίες χρήσης	Ελληνικά	45
AZUR™ Vasküler Tıkaç Kullanım Talimatları	Türkçe.....	49
Васкуларна тапа AZUR™ Инструкции за употреба	Български.....	53
Vaskularni čep AZUR™ Upute za uporabu	Hrvatski	57
Cévní zátka AZUR™ Návod k použití	Česky	61

Vaskulaarne kork AZUR™ Kasutusjuhend	Eesti keel	65
AZUR™ érdugó Használati útmutató	Magyar	69
AZUR™ vaskulārais aizbāznis Lietošanas instrukcija.....	Latviešu valodā	73
„AZUR™“ kraujagyslių kištukas Naudojimo instrukcija.....	Lietuvių k.	77
Zatyczka naczyniowa AZUR™ Instrukcja użycia.....	Polski	81
Dispozitiv vascular AZUR™ Instrucţiuni de utilizare	Română	85
Сосудистая заглушка AZUR™ Инструкция по применению	Русский	89
Vaskularni čep AZUR™ Uputstvo za upotrebu	Srpski	93
Vaskularni čep AZUR™ Navodila za uporabo	Slovenščina	97
Cievna zátka AZUR™ Návod na použitie	Slovenčina	101
Судинна заглушка AZUR™ Інструкція з використання	Українська	105
AZUR™ سدادة الأوعية الدموية تعليمات الاستخدام.....	العربية	109

English
AZUR™ Vascular Plug
Instructions for Use

DEVICE DESCRIPTION

The AZUR Vascular Plug is an embolization device consisting of a conformable, self-expanding nitinol braided wire frame surrounding a flexible, occlusive membrane. The AZUR Vascular Plug is deployed in an appropriately sized vessel to reduce or block the flow of blood. The implant will have radiopaque markers to provide visual confirmation of deployment location during the interventional treatment. The implant is deliverable through a catheter on a detachable delivery system. The delivery pusher is powered by an AZUR Detachment Controller to detach the implant. The catheter and AZUR Detachment Controller are provided separately.

The implantable plug component that remains in the patient following detachment consists of:

Table 1

Implant Material		AZUR Vascular Plug*
Metallic Substances	<ul style="list-style-type: none"> • Nitinol Wire • Platinum, Iridium 	≤ 0.016 g
Non-metallic sub-stances	<ul style="list-style-type: none"> • ePTFE, PET • polyolefin, adhesive 	≤ 0.003 g
* Approximate content		

INDICATIONS FOR USE

The AZUR Vascular Plug is intended to reduce or block the rate of blood flow in vessels of the peripheral vasculature.

INTENDED USER

This device should only be used by physicians who have undergone training in the use of the AZUR system for embolization procedures as prescribed by a representative from Terumo or a Terumo-authorized distributor.

CONTRAINDICATIONS

Use of the AZUR Vascular Plug is contraindicated in any of the following circumstances:

- When patient has known hypersensitivity to nickel-titanium.
- When end arteries lead directly to nerves.
- When arteries supplying the lesion to be treated are not large enough to accept emboli.
- In the presence of severe atheromatous disease.
- In the presence of vasospasm (or likely onset of vasospasm).

POTENTIAL COMPLICATIONS

Potential complications include but are not limited to: hematoma at the site of entry, vessel perforation, unintended parent artery occlusion, incomplete filling, vascular thrombosis, hemorrhage, ischemia, vasospasm, edema, implant migration or misplacement, premature or difficult implant detachment, clot formation, revascularization, post-embolization syndrome, and neurological deficits including stroke and possibly death.

The physician should be aware of these complications and instruct patients when indicated. Appropriate patient management should be considered.

REQUIRED ADDITIONAL ITEMS

- AZUR Detachment Controller
- 0.027" – 0.047" ID reinforced catheter for delivery of the AZUR system
- Guide Catheter compatible with delivery catheter, if applicable
- Steerable guidewires compatible with catheter
- 2 rotating hemostatic Y valves (RHV)
- 1 three-way stopcock
- Sterile Saline and/or lactated Ringer's Injection
- Pressurized sterile saline drip
- 1 one-way stopcock

Note : The PG Pro and Headway 27 microcatheters manufactured by MicroVenture were used for 0.027"ID catheter compatibility assessment.

WARNINGS AND PRECAUTIONS

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State or Local Health Authority in which the user and/or patient is established.

- Refer to instructions supplied with all interventional devices to be used with the AZUR Vascular Plug for their intended uses, contraindications, and potential complications.
- This device is intended for single use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Angiography is required for pre-embolization evaluation, operative control, and post-embolization follow up.
- Do not advance the delivery pusher with excessive force. Determine the cause of any unusual resistance, remove the AZUR system, and check for damage.
- Advance and retract the AZUR system slowly and smoothly. Remove the entire AZUR system if excessive friction is noted. If excessive friction is noted with a second AZUR system, check the catheter for damage or kinking.
- The implant must be properly positioned within a maximum of 3 positioning attempts. If the implant cannot be properly positioned after 3 attempts, simultaneously remove the device and the catheter.
- If repositioning is necessary, take special care to retract the implant under fluoroscopy in a one-to-one motion with the delivery pusher. If the implant does not move in a one-to-one motion with the delivery pusher, or if repositioning is difficult, gently remove and discard the entire device.
- Tortuosity or complex vessel anatomy may affect accurate placement of the implant.
- The long-term effect of this product on extravascular tissues has not been established so care should be taken to retain this device in the intravascular space.
- Always ensure that at least two AZUR Detachment Controllers are available before starting an AZUR system procedure.
- The implant cannot be detached with any power source other than an AZUR Detachment Controller.
- Do NOT place the delivery pusher on a bare metallic surface.
- Always handle the delivery pusher with surgical gloves.
- Do NOT use in conjunction with radio frequency (RF) devices.


Figure 1 – Diagram of AZUR System Setup


Figure 2 – Diagram of AZUR Vascular Plug Packaging

System Size	Model No.	Treatable Vessel Range
Small	45-180500	2.5 mm – 4.5 mm
Medium	45-180800	4.5 mm – 6.5 mm
Large	45-181000	6.5 mm – 8.0 mm

Figure 3 – AZUR Vascular Plug Dimensions

PREPARATION FOR USE

1. Refer to Figure 1 for the set-up diagram.
2. Attach a rotating hemostatic valve (RHV) to the hub of the guiding catheter. Attach a 3-way stopcock to the side arm of the RHV and then connect a line for continuous infusion of flush solution.
3. Attach a second RHV to the hub of the catheter. Attach a 1-way stopcock to the sidearm of the second RHV and connect the flush solution line to the stopcock.
4. Open the stopcock and flush the catheter with sterile flush solution and then close the stopcock. To minimize the risk of thromboembolic complications, it is critical that a continuous infusion of appropriate sterile flush solution be maintained into the guide catheter, the femoral sheath, and the catheter.

CATHETERIZATION OF THE LESION

5. Access the parent vessel or vascular lesion using standard interventional procedures.
6. After the catheter has been positioned at the target site, remove the guidewire.

IMPLANT SIZE SELECTION

7. Perform fluoroscopic road mapping.
8. Measure and estimate the size of the lesion to be treated.
9. Select an implant size based on Figure 3.
10. Correct implant size selection increases effectiveness and patient safety. In order to choose the optimum implant size for any given lesion, examine the pre-treatment angiograms. The appropriate implant size should be chosen based upon angiographic assessment of the diameter and length of the target vessel.

PREPARATION OF THE AZUR SYSTEM FOR DELIVERY

11. Remove the AZUR Detachment Controller from its protective packaging. Pull the white pull-tab from the side of the detachment controller. Discard the pull-tab and place the detachment controller in the sterile field. The AZUR Detachment Controller is packaged separately as a sterile device. **Do not use any power source other than the AZUR Detachment Controller to detach the implant. The AZUR Detachment Controller is intended to be used on one patient. Do not attempt to re-sterilize or otherwise re-use the AZUR Detachment Controller.**
12. Peel open the AZUR Vascular Plug pouch and remove the dispenser hoop from the pouch. Prior to using the device, remove the wedge from the dispenser hoop. While keeping the AZUR Vascular Plug implant and introducer in place, pull the pusher out of the dispenser hoop such that the proximal end of the delivery pusher is exposed. Use care to avoid contaminating this end of the delivery pusher with foreign substances such as blood or contrast. Firmly insert the proximal end of the delivery pusher into the funnel section of the AZUR Detachment Controller. See Figure 4. **Do not push the detachment button at this time.**


Figure 4 – Testing the AZUR Detachment Controller

13. Wait three seconds and observe the indicator light on the detachment controller.
 - If the green light does not appear or if a red light appears, replace the device.
 - If the light turns green, then turns off at any time during the three-second observation, replace the device.
 - If the green light remains solid green for the entire three-second observation, continue using the device.
14. While the implant is still in the loading cover, inspect the implant for any irregularities or damage. **If any damage to the implant or delivery pusher is observed, DO NOT use the device.**
15. Flush the implant by filling the loading cover with at least 5-ml of saline with a syringe or by submerging in a saline filled bowl. **DO NOT REMOVE THE LOADING COVER DURING THIS STEP.**

16. With one hand securing the green introducer, grab the delivery pusher and gently pull proximally to begin retracting the implant. See Figure 5.


Figure 5 – Pull Delivery Pusher in Proximal Direction to Retract Implant

17. Retract the implant completely into the introducer such that the distal tip is approximately 1 cm inside the introducer. **IF THE IMPLANT IS PULLED TOO FAR INTO THE INTRODUCER, IT MAY INCREASE RESISTANCE DURING INTRODUCTION INTO THE DELIVERY CATHETER.**
18. Once the implant is retracted approximately 1 cm into the introducer, pull the introducer until it disengages from the loading cover. Continue pulling proximally until it is completely removed from the dispenser hoop. **The implant should remain in its position within the introducer during this step.**

INTRODUCTION AND DEPLOYMENT OF THE AZUR SYSTEM

19. Open the RHV on the catheter just enough to accept the introducer sheath of the AZUR system.
20. Insert the distal end of the introducer sheath of the AZUR system through the RHV and into the hub of the catheter until the sheath is firmly seated. Tighten the RHV **lightly** around the introducer sheath to secure the RHV to the introducer. **Do not over-tighten the RHV around the introducer sheath. Excessive tightening could damage the device.**
21. Visually verify the flushing solution is infusing normally. Saline should be seen exiting the proximal end of the introducer sheath.
22. Advance the delivery pusher and implant into the lumen of the catheter until the proximal end of the delivery pusher extends approximately 12" past the proximal end of the introducer sheath. Loosen the RHV. Retract the introducer sheath proximally just out of the RHV. Close the RHV around the delivery pusher. Slide the introducer sheath proximally until it is completely off of the delivery pusher. Use care not to kink the delivery system.
23. Discard the introducer sheath. The AZUR system cannot be re-sheathed after introduction into the catheter.
24. Under fluoroscopic guidance, slowly advance the implant out the tip of the catheter. Continue to advance the implant into the desired site until the proximal radiopaque marker on the implant is aligned or slightly distal of catheter distal tip RO marker, positioning the detachment zone just outside the catheter tip. See Figure 6. Reposition if necessary. If the implant size is not suitable, remove and replace with another device. If undesirable movement of the implant is observed under fluoroscopy following placement and prior to detachment, remove the implant and replace with another more appropriately sized implant. Movement of the implant may indicate that the implant could migrate once it is detached. **DO NOT** rotate the delivery pusher during or after delivery of the implant into the vasculature. Rotating the delivery pusher may result in damage to the implant or premature detachment of the implant from the delivery pusher, which could result in implant migration. Angiographic assessment should also be performed prior to detachment to ensure that the implant is not extending into undesired vasculature.
25. Complete the deployment and any repositioning. **If the implant cannot be properly positioned within 3 attempts, simultaneously remove the device and the catheter.**
26. Tighten the RHV to prevent movement of the implant.
27. Verify that the distal shaft of the delivery pusher is not under stress before implant detachment. Axial compression or tension could cause the tip of the catheter to move during implant delivery. Catheter tip movement could cause the vessel to perforate.


Figure 6 – Position of Marker Bands for Detachment

DETACHMENT OF THE IMPLANT

28. The AZUR Detachment Controller is pre-loaded with battery power and will activate when a delivery pusher is properly connected. It is in a "power off"

mode when no delivery pusher is attached. It is not necessary to push the button on the side of the AZUR Detachment Controller to activate it.

29. Verify that the RHV is firmly locked around the delivery pusher before attaching the AZUR Detachment Controller to ensure that the implant does not move during the connection process.
30. Although the delivery pusher's gold connectors are designed to be compatible with blood and contrast, every effort should be made to keep the connectors free of these items. If there appears to be blood or contrast on the connectors, wipe the connectors with sterile water or saline solution before connecting to the AZUR Detachment Controller.
31. Connect the proximal end of the delivery pusher to the AZUR Detachment Controller by firmly inserting the proximal end of the delivery pusher into the funnel section of the AZUR Detachment Controller. See Figure 4.
32. When the AZUR Detachment Controller is properly connected to the delivery pusher, a single audible tone will sound and the light will turn green to signal that it is ready to detach the implant. If the detachment button is not pushed within 30 seconds, the solid green light will slowly flash green. Both flashing green and solid green lights indicate that the device is ready to detach. If the green light does not appear, check to ensure that the connection has been made. If the connection is correct and no green light appears, replace the AZUR Detachment Controller.
33. Verify the implant position before pushing the detachment button.
34. Push the detachment button. When the button is pushed, an audible tone will sound and the light will flash green.
35. At the end of the detachment cycle, three audible tones will sound and the light will flash yellow three times. This indicates that the detachment cycle is complete. If the implant does not detach during the detachment cycle, leave the AZUR Detachment Controller attached to the delivery pusher and attempt another detachment cycle when the light turns green.
36. The light will turn red after the number of detachment cycles specified on the AZUR Detachment Controller labeling. DO NOT use the AZUR Detachment Controller if the light is red. Discard the AZUR Detachment Controller and replace it with a new one when the light is red.
37. Verify detachment of the implant by first loosening the RHV valve, then pulling back slowly on the delivery system and verifying that there is no implant movement. If the implant did not detach, do not attempt to detach it more than two additional times. If it does not detach after the third attempt, remove the delivery system.
38. After detachment has been confirmed, slowly retract and remove the delivery pusher. **Advancing the delivery pusher once the implant has been detached involves the risk of vessel perforation. Do NOT advance the delivery pusher once the implant has been detached.**
39. Verify the position of the implant angiographically.

The physician has the discretion to modify the implant deployment technique to accommodate the complexity and variation in embolization procedures. Any technique modifications must be consistent with the previously described procedures, warnings, precautions and patient safety information.

SPECIFICATIONS FOR AZUR DETACHMENT CONTROLLER

- Output voltage: 8 ± 1 VDC
- Cleaning, preventative inspection, and maintenance: The AZUR Detachment Controller is a single use device, preloaded with battery power, and packaged sterile. No cleaning, inspection, or maintenance is required. If the device does not perform as described in the Detachment section of these Instructions, discard the AZUR Detachment Controller and replace it with a new unit.
- The AZUR Detachment Controller is a single use device. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Batteries are pre-loaded into the AZUR Detachment Controllers. Do not attempt to remove or replace the batteries prior to use.
- After use, dispose of the AZUR Detachment Controller in a manner consistent with local regulations.

PACKAGING AND STORAGE

The AZUR system is placed inside a protective, plastic dispenser hoop and packaged in a pouch and unit carton. The AZUR system and dispenser hoop will remain sterile unless the package is opened, damaged, or the expiration date has passed. Store at a controlled room temperature in a dry place.

The AZUR Detachment Controller is packaged separately in a protective pouch and carton. The AZUR Detachment Controller has been sterilized; it will remain sterile unless the pouch is opened, damaged, or the expiration date has passed. Store at a controlled room temperature in a dry place.

After use, dispose of the delivery system and detachment controller in accordance with hospital, administrative and/or local government policy.

SHELF LIFE

See the product label for the device shelf life. Do not use the device beyond the labeled shelf life.

MR INFORMATION

Non-clinical testing demonstrated that the implant is **MR conditional**. A patient can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only.
- Maximum spatial gradient magnetic field of 2,500-gauss/cm (25-T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the First Level Controlled Operation Mode.

MRI-Related Heating

Under the scan conditions defined, the implant is expected to produce a maximum temperature rise of 3.2°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

Image Artifact Information

In non-clinical testing, the image artifact caused by the implant extends approximately 2-mm from the device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system.

Terumo Corporation recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation or equivalent organization.

MATERIALS

The AZUR system does not contain latex or PVC materials.

DISPOSAL

Dispose of device/syringe in accordance with hospital policy and local regulations for biohazardous waste.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) for the device will be accessible in the European database on medical devices after the launch of the European Database on Medical Devices (Eudamed) (Eudamed: <https://ec.europa.eu/tools/eudamed>). The SSCP will be linked to the Basic UDI-DI in the Eudamed public website. A patient implant card is included in the package. This card should be completed and provided to the patient.

WARRANTY

MicroVention and Terumo warrant that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for particular purpose. Handling, storage, cleaning and sterilization of the device as well as factors relating to the patient, diagnosis, treatment, surgical procedure and other matters beyond MicroVention or Terumo's control directly affect the device and the results obtained from its use. MicroVention or Terumo's sole obligation under this warranty is limited to the repair or replacement of this device through its expiration date, and MicroVention and Terumo shall not be liable for any incidental, indirect, special or consequential loss, damage or expense directly or indirectly arising from the use of this device. Neither MicroVention nor Terumo assumes, or authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. Neither MicroVention nor Terumo assumes any liability with respect to devices reused, reprocessed or resterilized and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use, with respect to such device.

Prices, specifications, and model availability are subject to change without notice.

MicroVention™ and Headway™ are trademarks of MicroVention, Inc., registered in the United States and other jurisdictions.

© Copyright 2025 Terumo Corporation. All rights reserved.

All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties.

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© Copyright 2025 Terumo Corporation. Tous droits réservés.
© Copyright 2025 Terumo Corporation. Alle Rechte vorbehalten.
© Copyright 2025 Terumo Corporation. Todos los derechos reservados.
© Copyright 2025 Terumo Corporation. Tutti i diritti riservati.
© Copyright 2025 Terumo Corporation. Todos os direitos reservados.
© Copyright 2025 Terumo Corporation. Alle rettigheder forbeholdes.