Clinell Universal Wipes



Clinell Universal Wipes can be used for surface disinfection and cleaning of non invasive medical devices

Why?

The most common pathogens may well survive or persist on surfaces for days, weeks and months and can thereby be a continuous source of transmission if no regular preventive surface cleaning and disinfection is performed.

What to clean?

Recommended equipment to be cleaned between every patient:

- Dental chairs
- Local work surfaces
- Curing lamps
- ✓ Inspection lights & handles
- Spittoons
- Aspirators
- ✓ Trolleys/delivery units
- X-ray units

Supports compliance with HTM 01-05¹, WHTM 01-05², NIPCM³ and UK IPC COVID-19 guidance⁴

Please ensure that you follow the manufacturer's guidelines.

GAMA Healthcare Ltd..

The Maylands Building, Maylands Avenue, Hemel Hempstead, Hertfordshire, HP2 7TG, UK. T: +44 (0)20 7993 0030 E: info@gamahealthcare.com www.gamahealthcare.com

 Department of Health (2013) Health Technical Memorandum 01-05: Decontamination in primary care dental practices.
DH: London

How?

After the completion of the patient treatment, you should discard and replace any disposable materials. Observe any post AGP downtime.

Perform hand hygiene, risk assess and don appropriate PPE.



Take one wipe from the pack. Working from clean to dirty, wipe in an "S" shaped pattern, taking care to not go over the same area twice.



Change wipe if it becomes dry or soiled and dispose of in the appropriate waste bin.



Let the surface air dry. Dispose of PPE in the appropriate waste bin and perform hand hygiene. Discard and replace all disposable fittings.

Ready for next patient.

The frequency of cleaning should be risk assessed with particular attention given to frequently touched surfaces.

Product range

Clinell Universal Wipes replace the traditional two-step system of a detergent clean followed by disinfection with a combined detergent/disinfectant wipe.

- Removal of surface contamination and bio-matter
- Prevention of microbial resistance
- Savings in staff time and patient care time
- Effectively kills 99.99% of norovirus and MRSA in 10 seconds and kills the COVID-19 virus in 30 secs
- Alcohol-free, preventing equipment damage and supports compliance with HTM 01-05¹, WHTM 01-05², NIPCM³ and UK IPC COVID-19 guidance⁴
- Provides a pre-measured dose of disinfectant, preventing dilution errors

For product information and codes:

www.gamahealthcare.com/products/universal-range





Health Protection Scotland Updated July 2018

- Decontamination of reusable non-invasive care equipment.

4. Public Health England (2020) COVID-19: infection prevention and control dental appendix Updated 10 August 2021

