


Glassware Etching in Laboratory Glassware Washers: Causes, Impacts, and Mitigation Strategies



Introduction

Glassware washers, essential components in modern scientific laboratories, simplify the cleaning process of various scientific apparatuses due to their convenience, efficiency, and effectiveness. However, labs occasionally encounter the issue of glassware etching, which can compromise the integrity and functionality of the glassware. Etching is a form of glassware corrosion that manifests as a dull and frosted appearance on the glass surface. The physical and chemical alteration of the glass often results in irreversible damage. This degradation not only impairs the aesthetics, but also jeopardizes the accuracy and precision of experimental procedures involving etched glassware. This article provides a comprehensive examination of the causes, consequences, and possible mitigation strategies for glassware etching within laboratory washers.

Causes of Etching:

The process of etching in glassware can be attributed to multiple factors.

Chemical Interaction: The leading cause involves the interaction between glass and substances such as acetone, especially when coupled with high heat. For example, high pH detergents can extract silica from the glass surface, causing it to become rough and etched¹. Additionally, if residual chemicals are not adequately rinsed before washing, they can react with certain materials on the glassware, leading to etching.

Water Quality: Contaminants in the water can have a profound effect on the glassware²:

- **Particulates:** Elements like iron, copper, and manganese, have the potential to form deposits on rinsed items, increasing the risk of discoloration or corrosion.
- **pH:** Fluctuations in pH levels due to various water contaminants can affect detergent performance,

particularly with enzymatic detergents.

- **Organics:** Organic compounds might hinder detergent effectiveness. High levels of Total Organic Carbon (TOC) for rinsing can leave water residues, causing discoloration.
- **Hardness ions:** Calcium and magnesium compounds can reduce detergent efficiency by interacting with detergent surfactants, hindering dirt dispersion. These compounds are less soluble in warm water compared to cold water, potentially leading to scaling. Insufficient removal of these ions may result in visible water marks on the rinsed items.
- **Silica:** Silica buildup on the glassware or tools being washed, washer interior, and heating components can cause etching.

Temperature Fluctuations: Extreme temperature variations during the wash and rinse cycles can cause thermal shock, leading to microcracks and surface corrosion on glassware.

By addressing these factors, the risk of etching can be reduced, ensuring the preservation of the glassware's integrity and the accuracy of scientific results.

Impacts of Glassware Etching:

Glassware etching has wide-ranging impacts for laboratory glassware. Etched, rough surfaces can adsorb chemicals unpredictably, leading to inaccuracies in reagent concentrations and skewing quantitative data. The obscured, frosted glass impairs the observation of reactions inside the glassware, limiting its utility in procedures that rely on visual monitoring. Furthermore, etched glassware becomes more fragile and susceptible to breakage, leading to reduced utility and higher replacement costs.

Mitigation Strategies:

Efficient mitigation strategies can minimize the adverse impacts of etching.

Detergent Selection: Utilizing a final working pH that is mildly alkaline or neutral can decrease the corrosive effects on glassware³. Liquid detergents are often less likely to etch glassware than powdered alternatives⁴. Avoid using excessive detergent concentrations, and only use the least amount of detergent necessary for thorough cleaning as high concentrations of detergent can contribute to glassware etching.

Water Treatment: Implementing water treatment equipment to reduce water contaminants can reduce glassware etching. For example, the Labconco WaterPro[®] PS Polishing System or WaterPro[®] RO System can provide high purity type I or III water. An inline water softener or conditioning unit may be necessary.

Temperature Management: Control temperature fluctuations during the wash and drying cycles to prevent thermal shock. A high water temperature makes glassware more susceptible to etching than a cooler water temperature. For glassware longevity, wash glassware with the lowest temperature required for thorough cleaning.

Routine Maintenance: Regular inspection, cleaning, and maintenance of the washer, including checking rinse aid, are crucial to prevent etching.

Manufacturer's Approach: Labconco offers comprehensive solutions for glassware etching, incorporating features like optimized wash cycles, automated detergent dispensing, precise temperature control, and pure water

Visit [labconco.com](https://www.labconco.com) to request information for all Labconco products.

LABCONCO

inlets. Innovations like Labconco's CleanPoint™ conductivity monitoring, help reduce unnecessary water rinses, further aiding in the mitigation process.

Conclusion:

Glassware etching in laboratory washers poses significant challenges to the integrity and lifespan of laboratory glassware. By carefully managing elements such as detergents, water quality, temperature, and routine maintenance, laboratories can improve the life and reliability of their glassware. Through the insights and strategies provided in this paper, a comprehensive framework can be formed to tackle glassware etching, preserving the integrity of experimental results, and optimizing laboratory processes.

References

1. Fuchs, R.J. *Formulation of Household Automatic Dishwasher Detergents*. Journal of the American Oil Chemists' Society 58.4 (1981): 366A-370A.
2. ANSI/AAMI ST108:2023. *Water for the processing of medical devices*. Association of Medical Instrumentation (2023).
3. Rutkowski, B.J. *The Cause and Prevention of Glassware Corrosion in Domestic Mechanical Dishwashers*. Journal of the American Oil Chemists' Society 48.4 (1971): 166-169.
4. Schwartz, Charles, and Bernard H. Gilmore. *Sodium metaphosphate in mechanical dishwashing*. Industrial & Engineering Chemistry 26.9 (1934): 998-1001.

Visit [labconco.com](https://www.labconco.com) to request information for all Labconco products.

LABCONCO