


Step-by-step guide: iTero™ scanner procedure for Denture & Removable

01 New Scan

- Tap New Scan.



02 Patient details

- Fill in the new patient's name or search for an existing patient.

03 Order section

- Tap Procedure and select Denture/Removable.
- Tap Type and select from the list.

Settings for each Denture Type vary. See pages 2-3 for instructions on each type.

04 Scan Options

- Toggle the Denture Copy Scan.
- Toggle NIRI (Near Infra-Red Imaging) Capture, if necessary.
- Confirm attaching a new sleeve to the wand.

05 Notes

- Enter any free text regarding this Rx.

Your Denture & Removable is now ready for scanning.



Step-by-step guide: iTero™ scanner procedure for Denture & Removable

Denture & Removable types

Partial Denture/Framework

- When Partial Denture/Framework Type is selected, an optional Denture Details window will be open.
- Select a Mould and a Shade System.
- Tap each relevant tooth on the tooth diagram and select if the tooth is Missing or will be Clasped. To undo, select Regular.



The Upper and Lower Denture options will toggle automatically based on the Tooth Diagram.



Immediate Denture

- When Immediate Denture Type is selected, an optional Denture Details window will be open.
- Select a Mould and a Shade System.
- Tap each relevant tooth on the tooth diagram and select if the tooth is Missing or will be Clasped. To undo, select Regular.

The Upper and Lower Denture options will toggle automatically based on the indications defined in the Tooth Diagram.




Step-by-step guide: iTero™ scanner procedure for Denture & Removable

Denture & Removable types - continue



Full Denture Tissue Based

- When Full Denture Tissue Based Type is selected, an optional Denture Details window will be open.
- Select a Stage, Mould, and Shade System.
- Toggle the Upper and Lower Denture options to indicate the relevant arch.



Full Denture Implant Based

- When Full Denture Implant Based Type is selected, an optional Denture Details window will be open.
- Select a Stage, Mould, and Shade System.
- Tap each relevant tooth on the Tooth Diagram and select Missing or Implant Based.
- An Implant Based window will open for each Implant on the Tooth Diagram.
 - Select a Scan Body Manufacturer, Implant Manufacturer, and Implant Type/Scan Body.
 - Save settings by selecting the back arrow.
 - Repeat for all Implant Based teeth.



Upper and Lower Denture options will toggle automatically based on the indications defined in the Tooth Diagram.

To save settings, tap the back arrow.

Continue with step 4 as described on page 1.

Denture/Removable procedure offers a straightforward workflow that enables precise planning and fabrication of partial and full dentures.

01 Clinical case review

Clinical history

- Six months ago, four implants were placed in the maxillary jaw.
- After ensuring a successful osteointegration, the patient can proceed with a fixed denture fabrication.


Details to complete the Rx

- Denture type - implant-based.
- Treatment stage – try in.
- Implants - 13,11,21,23.
- Mould - Square Ovoid.
- Shade body - A3.
- Gingival colour - Light pink
- Material - Ceramic: Zirconia.



02 Complete the Rx

+ Optional step



+ Toggle denture copy scan if scanning of the pre-existing denture is required.

8 Indicate implant position.

9 Add implant characteristics: scan body manufacturer, implant manufacturer implant type/scan body.

+ Add additional information for the lab.

The screenshot shows the software interface with the following sections:

- Scan Options:** Includes 'NIRI Capture' (disabled), 'Denture Copy Scan' (disabled), and 'New Sleeve Attached' (checked).
- Tooth Diagram:** A grid of tooth icons with numbers 18-28 and 48-38. A 'Missing' box is placed over tooth 15, and an 'Implant Based' box is placed over tooth 14.
- Implant Based Detail View:** Shows a dropdown for 'Scan Body Manufacturer' (Straumann), 'Implant Manufacturer' (Straumann), and 'Implant Type/Scan Body' (TissueLevel RN). A 'Delete' button is visible.
- Treatm Table:**

| Tooth No. | Scan Body Manufacturer | Implant Manufacturer | Implant Type/Scan Body | Actions |
|-----------|------------------------|----------------------|------------------------|--------------|
| 18 | | | | |
| 13 | | | | |
| 21 | Implant Based | | | Show Details |
| 23 | Implant Based | | | Show Details |
- Notes:** A text area labeled 'Add Note'.

03 Scan mode

+ Optional step

+ Denture copy scan.

4 Bite scan.

The diagram shows a top-down view of a dental arch with a denture. Three numbered steps are indicated:

- Step 1: A box around a single implant area.
- Step 2: A box around the entire upper arch.
- Step 3: A box around the entire opposing arch.

1 Scan body scan (repeat for every implant).

2 Upper arch scan.

3 Opposing arch scan.

*Note that the A.I. Cleanup is disabled by default in this procedure, and soft tissue will not be automatically deleted. To enable A.I. Cleanup, press the screen and hold; then tap Enable A.I. Cleanup.

04 View mode



Check if all anatomical parts of the soft tissue are captured. Ensure the scan bodies are accurately captured, including its location and angulation are captured accurately.



If preparation requires adjustments use the Eraser tool to erase, and to further rescan the modified areas.

05 Send mode



Send the completed scan to the lab.

