

3.0 mm Cannulated Screw System

Surgical Technique



3.0 mm Cannulated Screw System

The 3.0 mm Cannulated Screw System is part of a series of cannulated screws ranging from 3.0 to 8.0 mm. Matching a combination of screw options with necessary instrumentation, the systems provide convenience and flexibility for the orthopaedic surgeon and the OR staff. Consistent with all Zimmer Biomet implants, the screws offer TiMAX[®] surface treatment which has been shown in mechanical testing, to have increased fatigue strength compared to 316L electropolished stainless steel, type I anodized titanium, and machined titanium.¹ These systems are truly designed with the surgeon in mind.




Figure 1




Figure 2

The following surgical technique describes the use of a 3.0 mm cannulated screw for a distal osteotomy of the first metatarsal.

Guide Wire Insertion

Position the 1.1 mm end of the 2.0 mm Drill Guide (Cat. No. 110008413) at the entry point and insert a 1.1 mm Guide Wire (Cat. No. 14-450350) by drilling to the desired depth (Figure 1), and verify with fluoroscopy.

Countersinking and Measurement

Insert the 3.0 mm Depth Gauge with Countersink (Cat. No. 110008459) over the inserted guide wire. Using the incorporated handle, rotate the depth gauge back and forth until sufficient countersinking has been achieved, and note the screw length measurement on the calibrated gauge (Figure 2).

ⓘ **Note:** To prevent excessive bone removal, countersinking under power is not recommended.




Figure 3




Figure 4

Pre-Drilling

Connect the 2.0 mm Cannulated Drill (Cat. No. 110008402) to an AO adapter and slide the drill over the guide wire. Drill to appropriate depth and remove drill, leaving the guide wire in place (Figure 3), and verify with fluoroscopy.

Screw Insertion

Advance the screw over the guide wire using the BT10 Cannulated Hexalobular Screwdriver (Cat. No. 110008449) attached to the Screwdriver Handle (Cat. No. 110017406) until the head of the screw is completely seated in the bone (Figure 4), and then remove k-wire.




Figure 5

Screw Insertion (cont.)

A 3.0 mm Flat Washer (Cat. No. 110008345) may be used in conjunction with the screw for osteoporotic bone or where the cortex is thin, increasing the surface area of the screw head (Figure 5).

ⓘ **Note:** Washers should not be used if countersinking has already occurred.

ⓘ **Note:** Screw removal or additional adjustment when not using a guide wire in the cannulation of the screw can be achieved by utilizing the BT10 SOLID Hexalobular Screwdriver (Cat. No. 110009961), which is identified by a laser etched 5 mm black band.

In addition, a sterile 3 mm Screw Extractor (Cat. No. 110008497) is available for screw removal if a screw head is broken.

Equipment Maintenance

Special care should be exercised during the cleaning of cannulated instrumentation to ensure that all foreign material is removed from the inside of the instrument. Regular inspection and careful maintenance will increase the useful life of all instrumentation.

Caution

Guide wires should be replaced after each surgery. Use of scored or bent guide wires increases the risk of unwanted guide wire advancement or breakage during screw insertion.

During cannulated screw insertion the surgeon must periodically check for undesired advancement of the guide wire.

The text and schematic drawings included in this protocol are general instructions for the use of Biomet instrumentation only. Placement of screws and plates for the stabilization of specific fractures must be determined using the best judgment of the surgeon.

Implants

3.0 mm Partially-Threaded Cannulated Screws

Non-Sterile	Sterile No.	Description	Quantity
110007492	110009774	3.0 mm x 10 mm Partially-Threaded Cannulated Screw	2
110007494	110009775	3.0 mm x 12 mm Partially-Threaded Cannulated Screw	2
110007496	110009776	3.0 mm x 14 mm Partially-Threaded Cannulated Screw	2
110007498	110009777	3.0 mm x 16 mm Partially-Threaded Cannulated Screw	2
110007500	110009778	3.0 mm x 18 mm Partially-Threaded Cannulated Screw	2
110007502	110009779	3.0 mm x 20 mm Partially-Threaded Cannulated Screw	2
110007504	110009780	3.0 mm x 22 mm Partially-Threaded Cannulated Screw	2
110007506	110009781	3.0 mm x 24 mm Partially-Threaded Cannulated Screw	2
110007508	110009782	3.0 mm x 26 mm Partially-Threaded Cannulated Screw	2
110007510	110009783	3.0 mm x 28 mm Partially-Threaded Cannulated Screw	2
110007512	110009784	3.0 mm x 30 mm Partially-Threaded Cannulated Screw	2
110007514	110009785	3.0 mm x 32 mm Partially-Threaded Cannulated Screw	2
110007516	110009786	3.0 mm x 34 mm Partially-Threaded Cannulated Screw	2
110007518	110009787	3.0 mm x 36 mm Partially-Threaded Cannulated Screw	2
110007520	110009788	3.0 mm x 38 mm Partially-Threaded Cannulated Screw	2
110007522	110009789	3.0 mm x 40 mm Partially-Threaded Cannulated Screw	2
–	110009790	3.0 mm x 45 mm Partially-Threaded Cannulated Screw	0
–	110009792	3.0 mm x 50 mm Partially-Threaded Cannulated Screw	0
–	110009794	3.0 mm x 55 mm Partially-Threaded Cannulated Screw	0
–	110009796	3.0 mm x 60 mm Partially-Threaded Cannulated Screw	0
–	110009798	3.0 mm x 65 mm Partially-Threaded Cannulated Screw	0
–	110009800	3.0 mm x 70 mm Partially-Threaded Cannulated Screw	0

Washers

Non-Sterile	Sterile No.	Description	Quantity
110008345	110008488	Flat Washer 3.0 mm	3

Disposables

Non-Sterile	Sterile No.	Description	Quantity
14-450350	–	1.1 x 150 mm K-Wire Trc Tip (Pk/5)	1
14-450375	–	1.1 x 150 mm K-Wire Thd Tip (Pk/5)	1
–	110008402	2.0 mm Cann Drill w/ AO	2

Instruments

Non-Sterile	Sterile No.	Description	Quantity
110008413	–	1.1 / 2.0 mm Drill Guide	1
110008449	–	BT10 Cann Hexalobular Driver AO	2
110009961	–	BT10 Solid Hexalobular Driver AO	1
110008459	–	3.0 mm Depth Gauge w/ Countersink	1
110008500	–	Guide Wire Pusher 1.1 mm	1
13571	–	Screw Forceps	1
–	110008497	Screw Extractor 3 mm	1
110017406	–	Mini Ratchet Handle AO	1

Case/Tray

Non-Sterile	Sterile No.	Description	Quantity
110008469	–	Cann Screw Case 3.0 mm	1

*Reference IFU 01-50-4202 for sterilization parameters of the Cannulated Screw instrument cases and contents.

Indications & Contraindications

INDICATIONS

Small Cannulated Screws (3.0 mm) are intended for use in:

1. Fixation of small bones, including those in the foot, patella, ankle, wrist and elbow.
2. Arthrodesis of the foot, wrist and elbow.
3. Small and long bone osteotomies.
4. Fracture fixation of small bones, small bone fragments and long bones.

CONTRAINDICATIONS

1. Infections.
2. Patient conditions including blood supply limitations, insufficient quantity or quality of bone.
3. Patients with mental or neurologic conditions who are unwilling or incapable of following postoperative care instructions.
4. Foreign body sensitivity where material sensitivity is suspected or unknown, testing is to be completed prior to implantation of the device.

References

1. Data on file at Biomet. Test #DVA-107504-DVER. Mechanical testing not necessarily indicative of clinical performance.

This material is intended for health care professionals and the Zimmer Biomet sales force only. Distribution to any other recipient is prohibited. All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated. This material must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

Check for country product clearances and reference product specific instructions for use. For complete product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the package insert and www.zimmerbiomet.com, or contact your local Zimmer Biomet representative.

This technique was prepared in conjunction with a licensed health care professional. Zimmer Biomet does not practice medicine. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. The surgeon is responsible for determining the appropriate device(s) and technique(s) for each individual patient and this document does not replace the need for examination and/or advice in whole or in part by a surgeon.

©2016 Zimmer Biomet



Authorized Representative

Biomet UK Ltd.
Waterton Industrial Estate
Bridgend, South Wales
CF31 3XA
UK



ZIMMER BIOMET

Your progress. Our promise.™



Legal Manufacturer

Biomet Orthopedics
P.O. Box 587
56 E. Bell Drive
Warsaw, Indiana 46581-0587
USA

www.zimmerbiomet.com

CE 0086

CE mark on a surgical technique is not valid unless there is a CE mark on the product label.

BMET1048.1-GBL • REV0616