

Performance Characteristics of the Curative HT SARS-CoV-2 Test

Note: This assay has not been FDA approved or cleared but has been submitted to FDA for Emergency Use Authorization (EUA) and is currently under review.

- **Limit of Detection (LoD) - Analytical Sensitivity:**

Table 1. LoD Data

CONCENTRATION	PASSING RATE (%)	
	Manual Extraction	KingFisher Extraction
1,000 copies/mL	6/6 (100%)	6/6 (100%)
500 copies/mL	6/6 (100%)	6/6 (100%)
200 copies/mL	6/6 (100%)	6/6 (100%)
100 copies/mL	6/6 (100%)	5/6 (83%)
50 copies/mL	5/6 (83%)	1/6 (17%)
20 copies/mL	4/6 (67%)	2/6 (33%)

Based on the above results, the Curative HT SARS-CoV2 Test was determined to have an LoD of 200 copies/mL.

For the LoD Confirmation, 100% of samples processed with manual extraction passed, and 100% of samples processed with KingFisher extraction passed.

Table 2. LoD Passing Rate

Concentration	Passing Rate (%)	
	Manual Extraction	KingFisher Extraction
200 copies/mL	20/20 (100%)	20/20 (100%)

- **Inclusivity:**

Based on a January 10, 2022, GSAID database analysis, 6,900,922 viral target genomes were in the GSAID database. Of these genomes 4,914 have 2 or more mismatch bases at any point within the sequences targeted by the Curative HT SARS-CoV-2 Test forward and reverse primers and probe. Treating these numbers as unlinked, randomly assorted, and assuming the 2 mismatches are adjacent and at either primer 3' end would be a worst-case scenario where

4,914/6,900,922 or 0.07% of SARS-CoV-2 genomes in the database are at risk of not being detected by the assay.

In addition, the ability of the Curative HT SARS-CoV-2 Test to detect the SARS-CoV-2 omicron variant, B.1.1.529, was evaluated by testing 20 omicron-positive patient samples and 20 negative simulated patient samples as shown in Table 3 below.

Table 3. Inclusivity of the HT SARS-CoV-2 Test for the Detection of Omicron

		Reference Standard (Sequenced and Contrived Samples)	
		Positive	Negative
The Curative HT SARS-CoV-2 Test	Positive	20	0
	Negative	0	20
	PPA	100.0% (95% C.I. 83.2%-100.0%)*	
	NPA	100.0% (95% C.I. 83.2%-100.0%)*	

*Clopper-Pearson confidence intervals

Patient samples that tested positive for the omicron variant were identified through Curative’s internal sequencing project which provides genomic surveillance and variant tracking for SARS-CoV-2 using de-identified residual positive patient samples. The data collected in this study indicates that the Curative HT SARS-CoV-2 Test detects the omicron variant with 100.0% PPA and 100.0% NPA.

- **Cross-reactivity (Analytical Specificity):**

To assess microbial interference, microorganisms were spiked into contrived SARS-CoV-2-positive samples with a concentration of SARS-CoV-2 that was 2x LoD. No interference was observed.

Table 4. Microorganisms Tested for Cross-Reactivity and Microbial Interference

Microorganism	Test Concentration	Cross-Reactivity	Interference
Adenovirus 7A	10 ⁵ TCID/mL	NONE	NONE
Human Metapneumovirus	10 ⁵ TCID/mL	NONE	NONE
Parainfluenza virus 1	10 ⁵ TCID/mL	NONE	NONE
Parainfluenza virus 2	10 ⁵ TCID/mL	NONE	NONE
Parainfluenza virus 3	10 ⁵ TCID/mL	NONE	NONE
Parainfluenza virus 4	10 ⁴ TCID/mL	NONE	NONE
Influenza A	10 ⁵ TCID/mL	NONE	NONE
Influenza B	10 ⁵ TCID/mL	NONE	NONE
Enterovirus EV68	10 ⁵ TCID/mL	NONE	NONE
Respiratory Syncytial Virus	10 ⁵ TCID/mL	NONE	NONE
Rhinovirus	10 ⁵ TCID/mL	NONE	NONE
<i>Chlamydia pneumoniae</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Haemophilus influenzae</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Legionella pneumophila</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Mycobacterium tuberculosis</i>	10 ⁶ CFU/mL	NONE	NONE

Microorganism	Test Concentration	Cross-Reactivity	Interference
<i>Streptococcus pneumoniae</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Streptococcus pyogenes</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Bordetella pertussis</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Mycoplasma pneumoniae</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Pneumocystis jirovecii</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Candida albicans</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Pseudomonas aeruginosa</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Staphylococcus epidermis</i>	10 ⁶ CFU/mL	NONE	NONE
<i>Streptococcus salivarius</i>	10 ⁶ CFU/mL	NONE	NONE
Human coronavirus 229E	10 ⁴ TCID ₅₀ /mL	NONE	NONE
Human coronavirus OC43	10 ⁵ TCID ₅₀ /mL	NONE	NONE
Human coronavirus HKU1	10 ⁵ copies/mL	NONE	NONE
Human coronavirus NL63	10 ⁵ TCID ₅₀ /mL	NONE	NONE
MERS-coronavirus	10 ⁵ TCID ₅₀ /mL	NONE	NONE
Pooled human nasal wash	N/A	NONE	NONE

- **Interfering Substances:**

Table 5. Potentially Interfering Endogenous and Exogenous Substances

Substance	Active Ingredient	Test Concentration	Interference
Human DNA	Human DNA	1,666.67 ng/L	NONE
Afrin® Original nasal spray	Oxymetazoline hydrochloride (0.05%)	10% v/v	NONE
		15% v/v	YES
Mucin	Purified mucin protein	2.5 mg/mL	NONE
Cepacol® lozenges	Benzocaine, menthol	3 mg/mL	NONE
Chloraseptic® sore throat spray	Phenol, glycerin	5% v/v	NONE
Robitussin® cough syrup	Dextromethorphan, guaifenesin	5% v/v	NONE
Nicotine	Nicotine	0.03 mg/mL	NONE
Human whole blood	Blood	1% v/v	NONE
Inhaled bronchodilators	Albuterol sulphate	0.83 mg/mL	NONE
Antibiotic nasal ointment	Mupirocin	10 mg/mL	NONE
Antibiotic ointment	Tobramycin	4 µg/mL	NONE
Saline nasal spray	Sodium chloride (0.65%)	15% v/v	NONE
Neosynephrine® Original nasal drops	Phenylephrine hydrochloride (0.5%)	15% v/v	NONE
Zycam® nasal gel	Oxymetazoline hydrochloride	15% v/v	NONE
Nasal corticosteroid	Fluticasone propionate	5 µg/mL	NONE

- **Sample Stability:**

Patient samples may be stored in Specimen Transport Buffer for 144 hours (six days) at ambient temperature prior to being tested with the Curative HT SARS-CoV-2 Test.

- **Clinical Evaluation:**

When comparing symptomatic participants' results on the Curative HT SARS-CoV-2 Test and on a highly sensitive FDA-authorized (EUA) molecular comparator test, the Positive Percent Agreement (PPA) was 95.0% (95% confidence interval [95% CI]: 83.1% to 99.4%) and the Negative Percent Agreement (NPA) was 98.2% (95% CI: 93.5% to 99.8%). In total, 27.5% (11/40) of the positive samples had a low viral load as defined by FDA and as measured by the comparator test.

Table 6. Performance of the Curative HT SARS-CoV-2 Test in Symptomatic Individuals

	Comparator SARS-CoV-2		
Curative HT SARS-CoV-2	Positive	Negative	Total
Positive	38	2	40
Negative	2	107	109
Total	40	110	149

Positive Percent Agreement (PPA): 95.0% (95% CI: 83.1%, 99.4%)
 Negative Percent Agreement (NPA): 98.2% (95% CI: 93.5%, 99.8%)

When comparing the Curative HT SARS-CoV-2 Test to the comparator test among asymptomatic participants, the PPA was 95.0% (95% confidence interval [95% CI]: 75.1% to 99.9%) and the NPA was 99.1% (95% CI: 95.0% to 100.0%). In total 30% (6/20) of the positive samples had a low viral load as defined by FDA and as measured by the comparator test.

Table 7. Performance of the Curative HT SARS-CoV-2 Test in Asymptomatic Individuals

	Comparator SARS-CoV-2		
Curative HT SARS-CoV-2	Positive	Negative	Total
Positive	19	1	20
Negative	1	109	110
Total	20	110	130

Positive Percent Agreement (PPA): 95.0% (95% CI: 75.1%, 99.9%)
 Negative Percent Agreement (NPA): 99.1% (95% CI: 95.0%, 100.0%)