

# Guide to FoundationOne®CDx and FoundationOne®Liquid CDx Reports

## FDA-Approved Claims

Any FDA-approved claims will always appear at the beginning of the report, starting on page one.

## Professional Services

The Professional Services section provides information for all reported biomarker and genomic findings. This section is not reviewed or approved by the FDA.

**FOUNDATIONONE®CDx**

PATIENT: DISEASE: Breast invasive ductal carcinoma (IDC)

PHYSICIAN: ORDERING PHYSICIAN: NAME, MEDICAL FACILITY, ADDITIONAL RECIPIENT, MEDICAL FACILITY ID, PHYSICIAN

SPECIMEN: SPECIMEN SITE, SPECIMEN ID, SPECIMEN TYPE, DATE OF COLLECTION, SPECIMEN RECEIVED

**Companion Diagnostic (CDx) Associated Findings**

**1 GENOMIC FINDINGS DETECTED**

PIK3CA E545K

**FDA-APPROVED THERAPEUTIC OPTIONS**

Piqray® (Alpelisib)

**2 OTHER ALTERATIONS & BIOMARKERS IDENTIFIED**

Results reported in this section are not prescriptive or conclusive for labeled use of any specific therapeutic product. See professional services section for additional information.

Microsatellite status: MS-Stable<sup>§</sup>

Tumor Mutational Burden: 4 Muts/Mb<sup>§</sup>

BRCA2: P628fs\*16

FANCG: FANCG (NM\_004629) rearrangement intron 11<sup>§</sup>

MAP2K4: loss<sup>§</sup>

STAG2: splice site 2674-1G>C

STK11: F231L

TP53: C275F

§§ Refer to appendix for limitation statements related to detection of any copy number alterations, gene rearrangements, BRCA1/2 alterations, LOH, MSI, or TMB results in this section.

Please refer to appendix for Explanation of Clinical Significance Classification and for variants of unknown significance (VUS).

© 2021 Foundation Medicine, Inc. All rights reserved.

ABOUT THE TEST: FoundationOne®CDx is the first FDA-approved broad companion diagnostic for solid tumors.

Electronically signed by Julia A. Davis, M.D., Ph.D. | 21 June 2021

Sample Preparation: 150 Second SL, 1st Floor, Cambridge, MA 02141 • CLIA: 2200027301

Sample Analysis: 150 Second SL, 1st Floor, Cambridge, MA 02141 • CLIA: 2200027301

Post-Sequencing Analysis: 150 Second SL, 1st Floor, Cambridge, MA 02141 • CLIA: 2200027301

Foundation Medicine, Inc. | 1.888.999.3639

FDA APPROVED CLAIMS — PAGE 1 of 2

**FOUNDATIONONE®CDx**

PATIENT: DISEASE: Breast invasive ductal carcinoma (IDC)

PHYSICIAN: ORDERING PHYSICIAN: NAME, MEDICAL FACILITY, ADDITIONAL RECIPIENT, MEDICAL FACILITY ID, PHYSICIAN

SPECIMEN: SPECIMEN SITE, SPECIMEN ID, SPECIMEN TYPE, DATE OF COLLECTION, SPECIMEN RECEIVED

**ABOUT THE TEST:** FoundationOne®CDx is a next-generation sequencing (NGS) based assay that identifies genomic findings within hundreds of cancer-related genes. Interpretive content on this page and subsequent pages is provided as a professional service, and is not reviewed or approved by the FDA.

**3 Report Highlights**

- Targeted therapies with NCCN categories of evidence in this tumor type: Alpelisib + Fulvestrant (p. 10), Olaparib (p. 11), Talazoparib (p. 12), Everolimus (p. 13)
- Variants that may inform nontargeted treatment approaches (e.g., chemotherapy) in this tumor type: BRCA2 P628fs\*16 (p. 5)
- Evidence-matched clinical trial options based on this patient's genomic findings: (p. 15)
- Variants in select cancer susceptibility genes to consider for possible follow-up germline testing in the appropriate clinical context: BRCA2 P628fs\*16 (p. 5)

**4 2 Disease relevant genes with no reportable alterations: BRCA1, ERBB2**

1 See About the Test in appendix for details.

BIOMARKER FINDINGS	THERAPY AND CLINICAL TRIAL IMPLICATIONS
Microsatellite status - MS-Stable	No therapies or clinical trials, see Biomarker Findings section
Tumor Mutational Burden - 4 Muts/Mb	No therapies or clinical trials, see Biomarker Findings section

GENOMIC FINDINGS	THERAPIES WITH CLINICAL RELEVANCE (IN PATIENT'S TUMOR TYPE)	THERAPIES WITH CLINICAL RELEVANCE (IN OTHER TUMOR TYPE)
BRCA2 - P628fs*16	Olaparib	Niraparib
10 Trials see p. 15	Talazoparib	Rucaparib
PIK3CA - E545K - subclonal	Alpelisib + Fulvestrant	Everolimus
10 Trials see p. 17	none	Temsirolimus
STK11 - F231L - subclonal	none	none
8 Trials see p. 19		

© 2021 Foundation Medicine, Inc. All rights reserved.

ABOUT THE TEST: FoundationOne®Liquid CDx is the first FDA-approved broad companion diagnostic for solid tumors.

Electronically signed by Julia A. Davis, M.D., Ph.D. | 21 June 2021

Sample Preparation: 150 Second SL, 1st Floor, Cambridge, MA 02141 • CLIA: 2200027301

Sample Analysis: 150 Second SL, 1st Floor, Cambridge, MA 02141 • CLIA: 2200027301

Post-Sequencing Analysis: 150 Second SL, 1st Floor, Cambridge, MA 02141 • CLIA: 2200027301

Foundation Medicine, Inc. | 1.888.999.3639

FDA APPROVED CLAIMS — PAGE 1 of 2

### 1 FDA-Approved CDx Claims

List of FDA-approved companion diagnostics associated with your patient's findings.

### 2 Other Alterations and Biomarkers Identified

For FoundationOne CDx reports, all other genomic and biomarker findings without companion diagnostic claims will appear here. For FoundationOne Liquid CDx reports, alterations shown here are limited to short variants and select rearrangements and copy number alterations. The complete list of genomic and biomarker findings can be found in the Professional Services section.

### 3 Report Highlights

This feature distills important genomic insights in one easy-to-find place, helping you focus on the key actionable results to inform your patient's treatment plan.

### 4 Pertinent Negatives

Identifies important negative results that can be used for patient management. Pertinent negatives do not appear for FoundationOne Liquid CDx.

### 5 Therapies with Clinical Benefit

Therapies for each associated genomic finding are listed in the therapy table. On the left are therapies within your patient's tumor type, and on the right are those with proven clinical benefit in other tumor types. Therapy resistance based on your patient's genomic profile will also be indicated.

### 6 National Comprehensive Cancer Network® (NCCN®) Categories of Evidence and Consensus<sup>1</sup>

Associated NCCN Category that has been assigned to the therapy listed within your patient's tumor type.

### 7 Clinical Trials

Identifies number of trials based on your patient's unique genomic profile with page number for quick reference.

Note: The images shown on this piece are of a sample report and do not represent actual test results. This information is intended to educate healthcare providers on the FoundationOneCDx and FoundationOne Liquid CDx reports and should not be used for patient diagnosis or treatment decisions.

Sample report images last updated December 2021.

## Professional Services Continued

FOUNDATIONONE <sup>®</sup> CDx		PATIENT	TUMOR TYPE Breast invasive ductal carcinoma (IDC)	REPORT DATE
ORDERED TEST #		GENOMIC FINDINGS		
1		<p><b>GENE</b> <b>BRCA2</b></p> <p><b>ALTERATION</b> PGRS1*36</p> <p><b>TRANSCRIPT ID</b> NM_000059</p> <p><b>CODING SEQUENCE EFFECT</b> BRD34C</p> <p><b>VARIANT ALLELE FREQUENCY (%) VAR</b> 60.6%</p> <p><b>POTENTIAL TREATMENT STRATEGIES</b> — Targeted Therapies — Alterations that inactivate BRCA1 or BRCA2 may confer sensitivity to PARP inhibitors<sup>64,65</sup> or ATR inhibitors<sup>67,68</sup>. Clinical responses to PARP inhibitors have been reported for patients with reduced cell viability and increased DNA damage during ATR treatment further support the sensitivity of BRCA2-deficient cells to ATR inhibitors. The WEE1 inhibitor adavosertib alone (ORR: 23%; DCR: 63%); however, in the BRCA mutated cohort, no significant difference in clinical benefit was observed between the combination (ORR: 19%) and monotherapy (ORR: 20%) treatments<sup>69</sup>. In a Phase 1 monotherapy trial of adavosertib that included 9 patients with BRCA2 mutated solid tumors, a patients with BRCA1-mutated cancers (1 with ovarian serous carcinoma and 1 with oral squamous cell carcinoma) achieved PRs, and a third patient with ovarian serous carcinoma harboring mutations in</p> <p>v.4.02.01. For BRCA1 and BRCA2 mutation carriers, the risk of developing breast cancer by age 70 has been found to be approximately 57-65% and 39-49%, respectively, and a lifetime risk of up to 90% has also been reported<sup>67,69</sup>.</p> <p><b>FINDING SUMMARY</b> The BRCA2 tumor suppressor gene encodes a protein that regulates the response to DNA damage<sup>66</sup>. Inactivating mutations in BRCA2 can lead to the inability to repair DNA damage and loss of cell cycle checkpoints, which can lead to tumorigenesis<sup>66</sup>. Alterations such as seen here may disrupt BRCA2 function or expression<sup>66,69,70</sup>.</p> <p><b>POTENTIAL GERMLINE IMPLICATIONS</b> One or more of the BRCA2 variants observed here has been described in the ClinVar database as a likely pathogenic or pathogenic/germline mutation (by an expert panel or multiple submitters)</p>		
FOUNDATIONONE <sup>®</sup> CDx		PATIENT	TUMOR TYPE Breast invasive ductal carcinoma (IDC)	REPORT DATE
2		<p><b>NOTE</b> Clinical trials are ordered by gene and prioritized by age range inclusion criteria (or pediatric patients, depending on ordering medical facility, later trial phase, and verification of trial information within the last two months. While every effort is made to ensure the accuracy of the information contained below, the information available in the public domain is continually updated and should be investigated by the physician or research staff. This is not a comprehensive list of all available clinical trials. Foundation Medicine displays a subset of trial options and make them in order of descending priority: Qualification (or pediatric trial) + Geographical proximity + Later trial phase. Clinical trials listed here may have additional enrollment criteria that may require medical screening to determine final eligibility. For additional information about clinical trials or to conduct a search for additional trials, please see <a href="https://www.foundationmedicine.com/genomic-testing#support-services">clinicaltrials.gov</a>. Or, visit <a href="https://www.foundationmedicine.com/genomic-testing#support-services">https://www.foundationmedicine.com/genomic-testing#support-services</a></p> <p><b>GENE</b> <b>BRCA2</b></p> <p><b>ALTERATION</b> P628fs*16</p> <p><b>NCT02849496</b></p> <p>Olaparib With or Without Atezolizumab in Treating Patients With Locally Advanced or Metastatic Non-HR23-Positive Breast Cancer</p> <p><b>LOCATIONS:</b> Missouri, Illinois, Tennessee, Texas, Ohio, Michigan</p> <p><b>NCT02595931</b></p> <p>ATR Kinase Inhibitor VX-970 and Irinotecan Hydrochloride in Treating Patients With Solid Tumors That Are Metastatic or Cannot Be Removed by Surgery</p> <p><b>LOCATIONS:</b> Missouri, Tennessee, Pennsylvania, Florida, Connecticut, Massachusetts, California</p> <p><b>NCT03221400</b></p> <p>PEN-866 in Patients With Advanced Solid Malignancies</p> <p><b>LOCATIONS:</b> Arkansas, Tennessee, Oklahoma, Michigan, South Carolina, Colorado, Virginia, Maryland, Florida</p> <p><b>NCT04041218</b></p> <p>PARP Inhibition During Pre-surgical Window in Breast/Ovary Cancer</p> <p><b>LOCATIONS:</b> Texas</p> <p><b>NCT04191135</b></p> <p>Study of Olaparib Plus Pembrolizumab Versus Chemotherapy Plus Pembrolizumab After Induction With First-Line Chemotherapy Plus Pembrolizumab in Triple Negative Breast Cancer (TNBC) (MKI-7339-009/KEYLYNK-009)</p> <p><b>LOCATIONS:</b> Oklahoma, Illinois, Texas, Minnesota, Michigan, Georgia, Toronto (Canada)</p> <p><b>PHASE 2</b> PD-1, PARP</p> <p><b>PHASE 1</b> <b>TARGETS</b> ATR</p> <p><b>PHASE 1/2</b> <b>TARGETS</b> PARP, HSP90</p> <p><b>PHASE NULL</b> <b>TARGETS</b> PARP</p> <p><b>PHASE 2/3</b> <b>TARGETS</b> PD-1, PARP</p>		
FOUNDATIONONE <sup>®</sup> CDx		CLINICAL TRIALS		

This content provided as a professional service by Foundation Medicine, Inc. has not been reviewed or approved by the FDA.

© 2022 Foundation Medicine, Inc. All rights reserved.

Financially supported by Douglas Inc. MD, I

John H. Kim, M.D., Ph.D., Laboratory Director (318-220-0275)

Quinn Henderson, M.D., Ph.D., M.Sc., Laboratory Director (318-240-0439)

Foundation Medicine, Inc. | 1.888.688.3639

Sample Preparation: 7010 K2 Cook Road, Morrisville, NC 27560-4109

Sample Analysis: 7010 K2 Cook Road, Morrisville, NC 27560-4109

Post-Sequencing Analysis: 85 Second St., 1st Floor, Cambridge, MA 02140-1200

PROFESSIONAL SERVICES - PAGE 15 OF 20

**FOUNDATIONONE<sup>®</sup> LIQUID CDx**

ORDERED TEST #

PATIENT

TUMOR TYPE

**Ovary epithelial carcinoma (NOS)**

REPORT DATE

---

3

Variant Allele Frequency  
Percentage (VAF%)

HISTORIC PATIENT FINDINGS	VAF%
<b>Blood Tumor Mutational Burden</b>	1 Muts/Mb
<b>Microsatellite status</b>	Cannot Be Determined
<b>Tumor Fraction</b>	Cannot Be Determined
<b>ATM</b>	● Q1537* 0.15%
<b>BRCA2</b>	● N570fs*19 51.9%
<b>TP53</b>	● L194R 0.18%
	● R280K 0.28%
<b>DNMT3A</b>	● splice site 1015-2A>G 1.2%
	● V622fs*7 1.0%

- 1 Biomarker and Genomic Findings**  
Following the initial pages of the report, the professional services section goes into more detail about your patient's findings.
- 2 Clinical Trial Information**  
Detailed information about the clinical trials your patient has been matched to, ranked for the patient based on location and trial phase.
- 3 FoundationOne Liquid CDx Variant Allele Frequency Percentage (VAF%) Graph and Table**  
Shows the detected VAF% and where applicable in the patient's biomarkers and/or genomic signatures. Up to 5 previous tests may be shown. For FoundationOne CDx reports, VAF values are displayed in the Genomic Findings section of Professional Services, alongside other variant information.

## Medical Case Consulting

For additional help with report interpretation, select the “Ask An Expert” feature on your provider portal or contact client services at (888) 988-3639

To learn more about our FDA-approved portfolio, go to [foundationmedicine.com/portfolio](https://foundationmedicine.com/portfolio)

1. Referenced with permission from the National Comprehensive Cancer Network, Inc. © National Comprehensive Cancer Network, Inc. 2021. All rights reserved. To view the most recent and complete version of the recommendations, go online to NCCN.org. NCCN makes no warranties of any kind whatsoever regarding their content, use or application and disclaims any responsibility for their application or use in any way.

FoundationOne®CDx and FoundationOne®Liquid CDx are qualitative next-generation sequencing based *in vitro* diagnostic tests for advanced cancer patients with solid tumors and are for prescription use only. FoundationOne CDx utilizes FFPE tissue and analyzes 324 genes as well as genomic signatures. FoundationOne Liquid CDx analyzes 324 genes utilizing circulating cell-free DNA and is FDA-approved to report short variant in 311 genes. The tests are companion diagnostics to identify patients who may benefit from treatment with specific therapies in accordance with the therapeutic product labeling. Additional genomic findings may be reported and are not prescriptive or conclusive for labeled use of any specific therapeutic product. Use of the tests does not guarantee a patient will be matched to a treatment. A negative result does not rule out the presence of an alteration. Some patients may require a biopsy for testing with FoundationOne CDx when archival tissue is not available which may pose a risk. Patients who are tested with FoundationOne Liquid CDx and are negative for companion diagnostic mutations should be reflexed to tumor tissue testing and mutation status confirmed using an FDA-approved tumor tissue test, if feasible.

For the complete label, including companion diagnostic indications and important risk information, please visit [www.F1CDxLabel.com](http://www.F1CDxLabel.com) and [www.F1LCDxLabel.com](http://www.F1LCDxLabel.com).