Ultra-Blend™ Plus

Calcium Hydroxide Dentin Liner and Protective Base

Description: Ultra-Blend™ plus is a light-cured, radiopaque, dentin-adhering liner and base material containing calcium hydroxide and calcium hydroxyapatite in a urethane dimethacrylate (UDMA) base. Ultra-Blend plus has a higher filler load to minimize shrinkage and a high-viscosity to facilitate ease of placement. It can be easily contoured with a higher- or low-speed handpiece and is available in dentin and opaque white shades.

Indications: Ultra-Blend plus is self-adhering and used as a cavity liner and/or base material prior to placement of restorative materials.

Directions for Use:

A. General

1. Attach Black Micro™ tip securely onto syringe. For larger preparations, use the Black Mini™ tip.
2. Prepare area with a quality bonding system according to manufacturer's instruction.

3. Apply a 0.5-2mm layer of Ultra-Blend plus to the desired area.
NOTE: When applying to an area with difficult access, express a small amount to form a ball at the tip. Use the tip as a ball applicator to prevent over-dispensing.

4. Light cure with VALO™ LED Curing Light on Standard Power for 10 seconds (for lights with output 600 mW/cm^2 cure 20 seconds or for lights with output 600 mW/cm^2 cure for 10 seconds).

NOTE: For base build up, start with a 0.5mm layer and light cure. Follow with 2mm increments, curing between each increment.

5. Ultra-Blend plus will bond to surfaces prepared with dentin bonding agents and to dry dentin.

B. Dentin Liner

1. Rinse exposed dentin with water and air dry.
2. Apply a quality bonding system according to manufacturer's instruction.
3. Place Ultra-Blend plus precisely and in exact volume desired with a Black Micro tip.

4. Light cure with VALO™ LED Curing Light on Standard Power for 10 seconds (for lights with output 600 mW/cm^2 cure 20 seconds or for lights with output 600 mW/cm^2 cure for 10 seconds).

C. Protective Base

1. Rinse exposed dentin with water and air dry.

2. Apply a quality bonding system according to manufacturer's instruction. Mix with a Black Micro™ tip to syringe if base area is small; use a Black Mini™ tip for larger areas.

4. Apply Ultra-Blend plus and light cure with VALO™ LED Curing Light on Standard Power for 10 seconds (for lights with output 600 mW/cm^2 cure 20 seconds or for lights with output 600 mW/cm^2 cure for 10 seconds).

5. If final base is to be thicker than 2mm, apply 0.5mm (or less) as the first layer and light cure 10 seconds. Add additional layers of 2mm thicknesses as needed.

Note: Ultra-Blend plus light cures to a depth of 3mm.

D. Pulp Capping

NOTE: Use Ultra-Blend plus for near pulp (pink) and small mechanical non-hyperemic exposure. If larger exposure and/or if hyperemic pulp, endodontic therapy should be initiated.

1. Apply Consepis™ antibacterial solution with plastic Mini Dento-Infusor™ tip for 60 seconds using a gentle scrubbing action. Do not rinse. Gently dry air.

2. Apply Ultra-Blend plus with the Black Micro tip to the dry dentin at the border of the exposure or then over the exposure or over the thin pink dentin of the near exposure.

NOTE: Minimize covering dentin to maximize area available for dentine bonding techniques.

3. Light cure with VALO™ LED Curing Light on Standard Power for 10 seconds (for lights with output 600 mW/cm^2 cure 20 seconds or for lights with output 600 mW/cm^2 cure for 10 seconds).

4. Apply a quality bonding system according to manufacturer's instruction and restore with product of choice.

E. Opaquer

1. Roughen metal with bur or microabrasion.

2. Rinse exposed metal to be opaqued with water and air dry.
3. Apply a quality bonding system according to manufacturer's instruction and light cure.

4. Apply Ultra-Blend plus over metal incrementally with no layer thicker than 2mm. Light cure with VALO™ LED Curing Light on Standard Power for 10 seconds (for lights with output 600 mW/cm^2 cure 20 seconds or for lights with output 600 mW/cm^2 cure for 10 seconds).

5. Restore with composite or product of choice.

Precautions:

WARNING: Chemical may set in a few to a few hours. To avoid clogging, remove tip and re-cap with Luer Lock caps between each use. Always verify material flow prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.

1. Carefully read and understand all instructions and SDS before using Ultra-Blend plus. Use only as directed.

2. Clinician and patient should wear UV protective eyewear when curing Ultra-Blend plus.

3. For daily use, store at room temperature. Refrigerate bulk storage.

4. Redirect overhead light to prevent premature polymerization.

5. Resins can be sensitizing. Avoid repeated contact of uncured dental resin with skin. Do not use on patients with a known sensitivity to acrylates or other resins. If allergic reaction, dermatitis or rash develops, consult a physician.

6. Dispose of used tips and empty syringes properly.

7. Keep out of reach of children.

8. To avoid cross contamination, do not re-use tips.

9. Do not use after expiration date noted on containers.

10. To avoid cross contamination, disinfect and wipe syringe with an intermediate level disinfectant between uses and/or use a disposable syringe cover. If syringe cover is used, remove top, re-cap, and discard tip and syringe cover.

11. Keep Ultra-Blend plus out of heat/direct sunlight.

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