EN Chromacone® PVS Impression Material

Description/Indications: Chromacone® PVS Impression Material is a superhydrophilic vinyl polysiloxane impression material and is an exciting new advancement in prosthodontics. Chromacone PVS can be utilized for all final impressions for fixed or removable prosthodontics including crowns, bridges, dentures, partial dentures, and implants or whenever impression accuracy is of the prime concern. Chromacone PVS putty or Clone Bite can be used for all bite registrations.

Procedure:

CARTRIDGES:

- Load the cartridge in the dispensing gun.
- Immediately prior to taking the impression remove the cartridge cap and discard.
- Gently squeeze the dispenser to bleed enough material from the cartridge to ensure even flow of the base and catalyst. (About 1/4 inch of material).
- Attach the mixing tip.
- Squeeze the dispenser in order to extrude the product.
- After use, store cartridges on end with the used mixing tips attached to protect the remaining materials from contamination.
- Always bleed the cartridge before each use (see number 3 above).

TIPS:

Use Ultraident® Impression Mixing Tips attached to impression cartridge with intraoral tip to mix material and inject around preparation or use Ultraident® Impression FX Tip and impression syringe to inject and burnish material around preparation.

PLEASE NOTE:

- The temperature of Chromacone PVS Impression Material affects the working and setting times. Ideally it should be kept at room temperature (5-27° C) or (41-80° F).
- Avoid contaminating Chromacone PVS Impression Material with chemicals that could inhibit the set (e.g. acrylic and methacrylate residues, latex rubber and sulphur compounds).

ULTRADENT UNIVERSAL TRAY ADHESIVE:

Use of the tray adhesive is very important. For an accurate impression, coat tray with a thin layer including the edge. Let adhesive dry for 2 minutes before making impression.

TRAYS:

We recommend applying a thin layer of tray adhesive to all trays.

DISPENSING AND GUN ASSEMBLY:

The Ultraident Impression Gun can be sterilized using heat sterilizers, including STATIM Autoclaves, conventional autoclaves, chemical vapor, or dry heat sterilizers.

DISINFECTION:

The impressions should be rinsed under cold running water after removal from the mouth. Soaking for 30 minutes in liquid disinfectants such as 2% gluteraldehyde (acidic), 0.5% - 1% sodium hypochlorite, 0.1% - 1.5% providone-iodone (halogenated phenol) or 0.13% glutaraldehyde (neutral) will not affect the surface or dimensions.

CASTING:

The impression may be poured immediately after removal from the mouth following disinfection. Surface tension will be reduced if the impression is briefly washed out with a surfactant. Rinse thoroughly afterwards with water. Ultraident Impression Materials are compatible with all dental stones, epoxy resins, and polyurethanes.

CONTRA-INDICATIONS

There are no known contra-indications provided when the product is used as directed. There are no known

harmful reactions or secondary effects on patients and/or dental personnel using vinyl polysiloxane impression materials.

IMPORTANT NOTES AND PRECAUTIONS

The temperature of Chromacone PVS Impression Material affects working time. Ideally it should be kept at room temperature between 5-27° C (41-80° F). Even within those boundaries, changes in temperature will affect the working time. A 2° C temperature increase will shorten working time 10-15 seconds; conversely, a 2° C temperature decrease will lengthen working time 10-15 seconds. You may refrigerate to increase working and setting times. Always "bleed" the cartridge before each use by extruding a quarter inch of material before applying a fresh mixing tip. This will ensure even dispensing of fresh base and catalyst pastes. For even more precise application, Ultraident's Intraoral tips are available.

1. Certain chemical agents like superoxol or hydrogen peroxide may interfere with the set of Chromacone PVS, care should be taken to remove these residual materials from impression site. Latex gloves (or anything that contains sulfur) and hand lotions can inhibit or contaminate setting of all Chromacone PVS impression materials. Powder free or polyethylene gloves are recommended or wash hands thoroughly and mix putty bare handed.

2. Proper care should always be taken to ensure good preparation, cleanliness, and moisture control including hemorrhage and saliva in order to achieve an adequate impression.

3. There is no bubbling or gasing with Chromacone PVS. 4. Disinfection can be accomplished with any appropriate disinfectant following the original manufacturer's procedures.

5. Multiple casts can be poured from an impression and Chromacone PVS is compatible with all gypsum products, electroplating and refractory casts for up to 15 days.

6. Temperature affects all set times. All testing for working and setting times was done at 72° F/22° C. Shelf life is at least 3 years.

7. Do not wear latex gloves or hand lotions while mixing the putty. Wash hands and mix equal amounts of catalyst and base.

COMPATIBILITY

Chromacone PVS Impression Material is addition-reaction type silicone impression material. It may not be combined with standard or condensation type silicone materials or poly ether materials.

GENERAL PRECAUTIONS:

- For Professional use only.
- Review instructions, precautions, and SDS before beginning treatment. Use only as directed.
- Keep products out of heat/sunlight.
- Avoid skin exposure to resins.
- Isolate strong chemicals to area of treatment.
- Confirm that patient has no known allergies to treatment materials.
- Test flow of materials from syringe and tip before using intraorally.
- Never force syringe plungers.
- Use Ultraident® Syringe Covers and/or clean and disinfect syringes between patients.

	Technical Data:	
Heavy, Medium, Light/Wash, Super Light Regular Set	Working Time	2:30 min
	Total Setting Time	5:00 min
Heavy, Medium, Light/Wash, Super Light Fast Set Total	Working Time	1:15 min
	Total Setting Time	2:45 min

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