opal bond flow Adhesive

EN

Description/Indications:

Opal Bond Flow adhesive is a low viscosity light cure bonding adhesive used for bonding to etched enamel. It is recommended for permanent retainers, appliances, bondable bite turbos/ramps or similar applications.

Opal Bond Flow adhesive is designed with ideal flowable properties and may be used for direct and indirect bonding. The wear and strength of Opal Bond Flow adhesive is suitable for short term or permanent bonding.

Prepare Enamel for bonding:

1. Clean enamel surfaces to be bonded with pumice or a similar water-based, mild abrasive.
2. Etch, rinse and dry all tooth surfaces prior to bonding.
3. Apply a thin layer of orthodontic primer to surface of etched tooth.

Bonding lingual retainers:

1. Affix the retainer to the prepared enamel surface and express Opal Bond Flow adhesive to desired locations and desired thickness.

2. Ensure a smooth finish to the material before curing for patient comfort.

Note: Take care not to express Opal Bond Flow adhesive interproximally.

3. Cure following light manufacturer's recommended curing methods.

Bondable bite turbos/ramps/buttons:
1. Express material into molds or directly onto etched and primed enamel.

Note: When using molds, do not overfill molds as it may cause flash.

2. Cure following manufacturer's recommended curing methods.

3. To remove, gently squeeze until the composite comes loose. Clean remainder with a fluted bur.

Bonding brackets or appliances:

1. Apply a thin layer of Opal Bond Flow adhesive onto the bracket pad ensuring complete coverage and mesh penetration.

Note: Opal Bond Flow adhesive is designed with a flowable viscosity. For tackier adhesive try Opal Bond MV adhesive or Opal Bond adhesive.

2. Place bracket or appliance coated with Opal Bond Flow adhesive onto etched/primed tooth.

3. Remove any excess composite from circumference of bracket base.

4. Cure following light manufacturer's recommended curing methods

Note: Recommended curing time for halogen lights is 15 seconds on the mesial and 15 seconds on the distal of each tooth. Recommended curing time for LED lights can vary. Refer to curing light manufacturer's recommended curing times.

Precautions:

1. For Professional use only.
2. Review instructions, precautions, and SDS before beginning treatment. Use only as directed.
3. Keep products out of heat/sunlight. Refrigerate.
4. Use quality curing lights.
5. Do not use on patients with contraindications to orthodontic treatment.
6. Cure bonding resin completely, as uncured resin materials will be deficient in physical properties, stress and wear resistance.

7. Any resin can be sensitizing to skin or can cause allergic reactions if individual has been previously sensitized to the allergen. If seen, wash area of tissue exposed to resin completely. If dermatitis, rash or any other allergic reaction is observed, consult a physician.

8. When not in use, ensure cap is secured over the opening of the Opal Bond Flow adhesive syringe.
9. Spritznanzette sind nur für den Einmalgebrauch bestimmt.
9. Tips are single use only
10. Clean and disinfect syringes in between patients.

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© Copyright 2022 Ultradent Products, Inc. All Rights Reserved.