

Opalescence<sup>®</sup>  
whitening toothpaste



Now available direct from  
the whitening experts!

- Actively removes surface stains for whiter teeth
- Formulated to strengthen enamel and help prevent cavities
- Developed by a dentist
- Gentle enough for everyday use
- Unique cool mint flavor for fresh breath that lasts
- Available in a Sensitivity Relief formula to treat sensitive teeth

Ask your dental professional about  
Opalescence Whitening Toothpaste or  
**Buy online at [opalescence.com](http://opalescence.com)**



▶ Enter code  
**GoToothpaste**  
for \$2 off your first order.

For an SDS of this product, please see [www.ultradent.com](http://www.ultradent.com)  
© 2019 Ultradent Products, Inc. All rights reserved.

Manufactured by Ultradent Products, Inc. 505 West Ultradent Drive,  
South Jordan, Utah 84095 [www.opalescence.com](http://www.opalescence.com) Made in USA 1007486.1 021519



Enjoy your  
free sample!

Opalescence  
 go<sup>®</sup>

Professional whitening to go

**Get more confidence in your smile**

Your smile says more about you than any other physical feature.  
It's the first thing people notice. Wouldn't it be great to have  
a simple, affordable method for whitening your smile that's  
convenient enough to take with you anywhere?

Introducing Opalescence Go whitening, the easiest, fastest, and  
most convenient way to brighten your smile. Opalescence Go  
whitening gives you professional results in the convenience of  
pre-loaded disposable trays, so it is ready to go when you are.

Instructions

## Instructions



1. Product must be refrigerated. Before using, allow product to warm to room temperature. This will make it easier for the tray to adhere to the teeth. Remove from package.



3. Gently suck down on the tray. Do not press tray against teeth as it pushes the whitening gel against the gums, causing tissue irritation.



5. Suck down.



2. Center tray on arch.



4. Remove outer tray.



6. Wear Opalescence Go whitening 10% for 30–60 minutes, and Opalescence Go whitening 15% for 15–20 minutes per day. When finished, remove whitening tray and brush teeth.

10% | 30–60

15% | 15–20

For professional use only

## Precautions

1. If you have any questions regarding the appropriate use of this product, including how long it will take to whiten your teeth, please consult your dentist.
2. Regular dental checkups and cleanings are important before and after whitening to maintain a healthy smile.
3. Some patients experience sensitivity with whitening. Should this occur, stop treatment and contact your dentist.
4. Some patients have noticed temporary discomfort of the gums, lips, throat, or tongue. Should any of these symptoms persist more than two days or progressively worsen, call your dentist. These side effects will usually subside within 1–3 days after treatment is discontinued.
5. Pregnant or breastfeeding women, or patients with serious health concerns, should not whiten. Consult your primary care provider prior to treatment.
6. If you have a known allergy or chemical sensitivity to any of the ingredients in Opalescence whitening, do not use. Mint or Melon flavor can be irritating to some patients. If you have a known sensitivity to flavorings, please use the non-flavored Opalescence® PF whitening.
7. Do not use household bleach to whiten teeth.
8. Do not use tobacco or eat while whitening trays are in place.
9. Teeth are naturally darker along the gum line. These areas may require more time to lighten than the rest of the tooth surface and usually remain slightly darker.
10. White spots that appear on the enamel during the whitening process may blend with continued whitening.
11. Crowns, bridges, partial dentures, veneers, and composite fillings will not whiten.
12. Coffee, tobacco, and other products can re-stain your teeth over time. Should this occur, teeth can be re-whitened.
13. Some amalgam fillings may leave a dark purple color in your whitening tray; this is normal.
14. Check expiration date and store according to instructions.
15. Avoid swallowing large amounts of gel or water used to rinse gel residue from the teeth.
16. Avoid contact with eyes.
17. Do not reuse to avoid cross contamination.



Health Hazard



Do not re-use to avoid cross contamination



Keep away from heat/sunlight



Refrigerate



See instructions



Keep out of reach of children