

Opalescence®

tooth whitening systems



Manufactured by

ULTRADENT
PRODUCTS, INC.

505 West Ultradent Drive (10200 South)
South Jordan, UT 84095, USA
www.ultradent.com



Follow your dentist's instructions on how to load gel into your custom bleaching tray. Use no more than 1/3 to 1/2 of syringe.




Brush teeth before inserting tray. Lightly tap tray to adapt tray sides to teeth. Wear Opalescence 10% for 8-10 hours or overnight, Opalescence 15% for 4-6 hours, Opalescence 20% for 2-4 hours, and Opalescence 35% for 30-60 minutes, and Opalescence Quick 45% for 15-30 minutes.



Remove excess gel with clean finger or soft toothbrush. **Gently rinse twice; do not swallow rinsed gel.**



If significant sensitivity occurs, stop treatment and consult dentist.



Clean tray with soft brush and cool tap water. Store tray in case provided.

Bleaching Precautions:

1. Do not use household bleach to whiten teeth.
2. Not intended for use by pregnant women. Consult a health-care practitioner if you are breastfeeding or pregnant.
3. If you have any questions regarding the appropriate use of this product, including how long it will take to bleach your teeth, please consult your dentist as soon as possible.
4. Do not use tobacco or eat while bleaching trays are in place.
5. White spots on the enamel or that appear during bleaching may blend during continued bleaching process.
6. Teeth are naturally darker along the gumline. These areas may require more time to lighten than the rest of the tooth surface and usually remain slightly darker.

General Precautions:

1. A small percentage of patients experience sensitivity with bleaching. Should this occur, remove the tray and contact your dentist.
2. Foods and juices high in citric acid can cause sensitivity to the teeth while bleaching.
3. Some patients have noticed temporary discomfort of the gums, lips, throat, or tongue. Should any of these symptoms persist more than two days or progressively worsen, call your dentist. These side effects will usually subside within 1-3 days after treatment is discontinued.
4. Coffee, tobacco, and other products can re-stain your teeth over time. Should this occur, the teeth can be re-whitened within a few nights using Opalescence.
5. Regular dental checkups and cleanings are important before and after bleaching to maintain a healthy smile.
6. Some old amalgam or "silver" fillings may leave a dark purple color in your bleaching tray; this is normal.
7. Crowns, bridges, partial dentures, veneers, and composite fillings will not bleach.
8. Store bleach out of the sun and heat. Do not freeze.
9. Keep Opalescence out of reach of small children.
10. Do not swallow gel or rinsed gel. Product contains peroxide and may contain fluoride and potassium nitrate; swallowing large amounts can be harmful.
11. Do not use Opalescence if patient has a known allergy or chemical sensitivity to peroxides, Carbopol, glycerin, etc., mint or melon flavor can be irritating to some patients. If a patient has a known sensitivity to mint or melon, another flavor or non-flavored variant should be selected.

Opalescence® whitening toothpaste



Cool mint with fluoride

Indications:

Aids in the prevention of dental decay.

Directions:

Adults and children over 6 years of age: Brush teeth thoroughly, preferably after each meal or at least twice a day, or as directed by a dentist or physician. Children 6 years and younger: Use a pea-sized amount for supervised brushing to minimize swallowing. In case of intake of fluoride from other sources, consult a dentist or physician. Children under 2 years of age: Consult a dentist or physician.

Warning:

As with all fluoride toothpastes, keep out of the reach of children under 6 years of age. If you accidentally swallow more than used for brushing, seek professional assistance or contact a Poison Control Center immediately. Regulations on administering fluoride containing products may differ from country to country.

Toothpaste Ingredients: Contains Sodium Fluoride 0.25% w/w, Glycerin, Water (Aqua), Silica, Sorbitol, Xylitol, Flavor (Aroma), Poloxamer, Sodium Lauryl Sulfate, Carbomer, FD&C Blue#1 (CI 42090), FD&C Yellow#5 (CI 19140), Sodium Benzoate, Sodium Hydroxide, Sparkle (CI 77019, CI 77891), Sucralose and Xanthan Gum.



1oz (20ml) Opalescence
Whitening Toothpaste



Ultradent® Pocket Tray Case



Lot Number



Shade Guide



1.2ml Opalescence



Catalog Number



Recycle



Health hazard



Kosher Certified



See Instructions



Keep away from heat/sunlight

Opalescence
10%, 35% PF, 20% PF, 45%



Refrigerate



Use by date



Keep out of reach of children

Opalescence
10% PF, 15% PF, 20% PF



Store at room temperature

For professional use only.
All Opalescence Products are Gluten Free