

 **ultrapro™ Tx**  
cordless



# Instructions for Use





# Contents

---





<b>1</b>	Preface .....	5
<b>2</b>	Product Description .....	5
<b>3</b>	Indications for Use/Intended Purpose .....	5
<b>4</b>	Contraindications.....	6
<b>5</b>	Warnings and Precautions .....	6
<b>6</b>	Stepwise Instructions.....	7
<b>7</b>	Maintenance .....	9
<b>8</b>	Processing .....	10
<b>9</b>	Storage and Disposal .....	13
<b>10</b>	Technical Considerations .....	14
<b>11</b>	Miscellaneous Information .....	20
<b>12</b>	Symbol Instruction .....	21



## 1. PREFACE

---

Thanks for purchasing the Ultradent Ultrapro™ Tx Cordless hygiene handpiece. For the discerning dental professional who values an untethered alternative to cleaning and polishing, the Ultrapro Tx Cordless handpiece cuts the cord and allows the user to perform prophylaxis flexibly while providing the patient a comfortable and high-quality cleaning experience. It is recommended that you read the accompanying documentation to take full advantage of the features and benefits of the Ultrapro Tx Cordless handpiece.

## 2. PRODUCT DESCRIPTION

---

The Ultrapro Tx Cordless handpiece has an ergonomic and lightweight design for maximum clinician comfort and use. It is equipped with an on/off button and a wireless foot pedal used for dental prophylaxis procedures and is compatible with nearly all disposable prophylaxis angles on the market.

### Product Components:

- 1 – Motor component
- 3 – Autoclavable outer sheaths
- 1 – Charging station
- 2 – Power adapters
- 1 – Foot pedal
- 100 – Disposable barrier sleeves
- 4 – Disposable prophylaxis angles
- 6 – Black O-rings
- 2 – Green O-rings

### Overview of Controls:

The manufacturer accepts no liability for any damage resulting from the improper use of this unit and/or any purpose other than those covered by these instructions.

For all products described, carefully read and understand all instructions in this manual and battery SDS information prior to use.

## 3. INDICATIONS FOR USE/INTENDED PURPOSE

---

The Ultrapro Tx Cordless handpiece is a cordless prophylaxis handpiece equipped with a control button and a wireless foot pedal. It is intended for use with disposable prophylaxis angles in the hygiene operatory to perform cleaning and polishing procedures on tooth surfaces and fillings.

## 4. CONTRAINDICATIONS

---

- Doctors or patients with a pacemaker or hemophilia are prohibited from using this device.
- For patients or users with allergy concerns, refer to product allergen document available at [www.ultradent.com](http://www.ultradent.com). If allergic reaction is observed, rinse exposed area thoroughly with water and have the patient consult their physician.

## 5. WARNINGS AND PRECAUTIONS

---

- Do not place device near a heat source.
- Operate and store device in a safe dry place.
- Do not use lubrication on the motor component.
- Use caution for patients with heart conditions, pregnant women, and children.
- This device requires special precautions regarding electromagnetic compatibility (EMC) and must be in strict accordance with the EMC information for installation and use. Do not use this equipment near fluorescent lamps, radio transmitting devices, remote control devices, or handheld or mobile high-frequency communication devices.
- Long term use of equipment may result in overheating of the motor handpiece and should be left to cool until next use. Please contact Ultradent if any issues arise.
- To avoid damage or malfunction of device, use the original charging base.
- Any changes to the device may violate safety regulations, causing harm to patients. The manufacturer will not accept any liability for a modified device.
- To avoid damage to lithium battery and control circuit, use original power adapter.
- If you think a battery replacement is necessary, please contact Ultradent. Do not attempt to replace the battery.
- To avoid damage, ensure that the device is in the off position before removing or installing the outer sheath or disposable prophylaxis angle.
- To avoid damage to the handpiece and disposable prophylaxis angle, only use disposable prophylaxis angles compatible with this device and be sure that they are correctly installed before starting the handpiece.
- Wireless charging will generate heat. The surface temperature of the charging station and motor handpiece will rise.
- According to IEC 60601-1/UL60601-1, this device must not be used in the presence of a flammable anesthetic gas mixed with air, oxygen, or nitrous oxide. (NOTE: Nitrous oxide by itself is not a flammable anesthetic gas.)
- Do not use handpiece while charging.

## 6. STEPWISE INSTRUCTIONS

---

### Wireless Connection of Foot Pedal

Ensure that the foot pedal is fully charged before first use. If the handpiece is turned on, the Bluetooth will be paired automatically. The Bluetooth indicator will flash while the handpiece and the foot pedal are connecting. If the connection has failed, the Bluetooth indicator will keep flashing. There may be a 3-second delay to the handpiece when the foot pedal is activated. If there is no response to the handpiece after depressing the foot pedal, check the battery and recharge if necessary. After ensuring full connection, step on the foot pedal to control the speed of the handpiece.

### Charging the Handpiece and Foot Pedal

1. Plug the power adapter into an appropriate electrical outlet.
2. Place the motor component on the charging station connected to the power adapter.
3. When the device is charging, the indicator light on the charging station and foot pedal will flash green. Once fully charged, the battery indicator light will illuminate green.

#### Note:

- Handpiece and foot pedal are provided partially charged (charge level 50%–70%). Fully charge both before first use.


Do not depress the foot pedal while charging.


### Preparation

1. Install disposable barrier sleeve onto the motor component.
2. Install the sterilized outer sheath onto the motor component covered with the disposable barrier sleeve.  
**Note:** Before installation, check whether the disposable barrier sleeve is intact. Do not use if damaged.
3. Attach the disposable prophylaxis angle onto the cordless handpiece.

### Use

1. Press the button on the handpiece once to turn device on. Depress wireless foot pedal to control speed of motor handpiece.  
**Note:** There may be a 3-second delay during the Bluetooth connection of the foot pedal and handpiece.
2. Use according to clinical need.
3. After use, discard the disposable prophylaxis angle and disposable barrier sleeve. Disinfect and sterilize outer sheath between each use. Do not immerse sheath in disinfection solution or ultrasonic bath. Only disinfect the motor component, do not sterilize.

INDICATOR PICTURE	CHARGE LEVEL	LED COLOR
	Battery fully charged	3 green
	Charged 2/3	2 bottom indicators green and 1 top indicator off
	Charged 1/3	1 bottom indicator green and 2 top indicators off
	Battery low	1 bottom indicator orange and 2 top indicators off. Unit will soon shut off automatically, charge immediately.



BATTERY INDICATOR'S COLOR	BATTERY OF FOOT PEDAL STATUS
Green light flashes	Charging
Green light illuminates	Fully charged or the battery is sufficient for operation
Orange	Low battery, unit will soon shut off automatically, charge immediately.

### Power Off

The handpiece will automatically shut down once placed on the charging station or after 5 minutes of device not being in use.



## 7. MAINTENANCE

---

### Repair

1. Do not use if product or accessories are damaged, smashed, detached, corroded, or bent.
2. This device does not include accessories for repair usage. Please contact Ultradent for repairs.
3. Repairs are only to be performed by authorized service personnel.



#### Warning:

- Keep the equipment in a dry storage area.
- Do not throw, beat, or shock the equipment.
- Contact Ultradent if the battery needs repair or replacement.

### Replacing the O-ring



#### Warnings

The O-ring should be replaced regularly every 6 months or sooner if needed.

- Remove the disposable prophyl angle and outer sheath from the handpiece.
- Hold the O-ring between your thumb and index finger and apply firm pressure to disengage the ring from the metal channel.
- Pull the old O-ring off and replace with a new O-ring.

### Warranty

Ultradent hereby warrants that this instrument shall, for a period of 2 years\*, conform in all material respects to the specifications therefore as set forth in Ultradent's documentation accompanying the product and be free from any defects in materials/or workmanship. Batteries are covered for a period of 1 year. This warranty applies solely to the original purchaser and is not transferable. All defective products are to be returned to Ultradent. There are no user service components of the Ultrapro Tx Cordless handpiece. Tampering with Ultrapro Tx Cordless handpiece will void its warranty. The Ultrapro Tx Cordless handpiece warranty does not cover customer damage. For example, if an Ultrapro Tx Cordless handpiece is misused or dropped and breaks, the customer would be responsible for paying for any necessary repairs.

Consumable items contained herein are not covered by warranty.

### Shipping container information:

The shipping container you received with your Ultrapro Tx Cordless handpiece was specially designed to safely transport the device. In the unlikely event that you need to return the handpiece for service or repair, please retain the original shipping container.

\*With sales receipt indicating the date of sale to the dentist.

## 8. PROCESSING

The objective of the information provided in this section is to reduce the potential for patient cross-contamination when using the device. Please follow all cleaning, disinfection, and sterilization procedures as recommended.

- Sterilize the outer sheath prior to each use. Additional outer sheaths are available for purchase to accommodate uninterrupted patient procedures.
- Do not place the outer sheath in a disinfectant solution or in an ultrasonic bath.
- Do not use bleach or chloride disinfectant materials.
- For your own safety, please wear personal protective equipment (gloves, glasses, mask).
- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterility.

PART NAME	INFECTION PREVENTION
<b>Disposable Prophy Angle/DPA</b> Prophy angle that helps clinician easily access hard-to-reach areas of the mouth quickly for optimal dental hygiene.	Single use Discard after use
<b>Outer Sheath</b> The outer sheath is a sterilizable protective cover that acts as a sleeve to the internal motor component. It also houses the bearings and acts as an interface between DPA and the motor component. The outer sheath features an on-board user centralized control button for power of the device.	Manual cleaning followed by steam autoclave
<b>Disposable barrier sleeve</b> The disposable barrier sleeve is designed to protect the motor component from debris.	Single use Discard after use
<b>Motor Component</b> The motor component house, the motor, power supply (lithium-ion battery), and 3 indicator lights. The motor component is not sterilizable and is to be used to provide speed and torque.	Manual cleaning and disinfection

PART NAME	INFECTION PREVENTION
<p><b>Power Adapter</b> Converts AC current into DC current required for charging the handpiece.</p>	Manual cleaning and disinfection
<p><b>Foot Pedal</b> Provides a wireless control for rotation speed control of handpiece. It has two indicator lights.</p>	Manual cleaning and disinfection
<p><b>Charging Station</b> The charging station uses an induction charge mechanism, allowing the handpiece to be placed into the charger without regard to any particular alignment.</p>	Manual cleaning and disinfection

### Instruction for Outer Sheath Cleaning and Sterilizing

- The outer sheath can be cleaned manually or by an automated dental instrument washer prior to sterilization.
  - Note: DO NOT put the outer sheath into an ultrasonic bath or immerse in any type of fluid or chemical.
  - The outer sheath is not sterile upon receipt and must be steam autoclaved prior to use in accordance with the following instructions.
  - Wear appropriate PPE (e.g., mask, protective eyewear, and gown).
1. Disassemble by removing the disposable prophylaxis angle (DPA) from the outer sheath and remove the outer sheath along with barrier sleeve from the motor component.

#### 2.A Steps for Manual Cleaning of Outer Sheath

- Rinse without submerging the outer sheath using tap water, purified, distilled, or deionized water at room temperature to remove any visible surface contaminants.
- Prepare enzymatic solution per manufacturer's recommendation.
- Remove visible debris by thoroughly wiping outer sheath with enzymatic solution. If visible soil is seen, repeat the cleaning procedure.

## 2.B Steps for Automated Dental Instrument Washing of Outer Sheath

- Use only enzymatic solution that meets the FDA requirements for cleaning process and follow the manufacturer's instructions for correct use.
3. Rinse the outer sheath to remove enzymatic solution residue.
  4. After manual or automated cleaning, use a dry cloth or paper towel to wipe the surface to remove excess water from the outer sheath.
  5. Package the cleaned and dried outer sheath in a sterilization pouch to maintain sterility after sterilization.

### Sterilization of Outer Sheath

1. Place bagged outer sheath into a steam autoclave per the autoclave manufacturer's instruction.
2. It is recommended to use these validated sterilization parameters: A minimum sterilization time of 4 minutes at a temperature of 132°C/134°C and a pressure of 185kPa to 190kPa with a drying time of 20 minutes.
3. The highest sterilization temperature is 138°C.
4. The outer sheath should remain bagged until ready for use.

Note: The outer sheath has been designed for a large number of sterilization cycles. The materials used in manufacturing were selected accordingly. However, with every renewed preparation for use, thermal and chemical stresses will result in aging of the products. Discoloration and breakdown may occur around 500 uses. End of life is typically determined by wear and damage due to use.

### Motor Component and Accessory Cleaning and Disinfection

NOTE: The motor component, foot pedal, power adapter, and charging station are not waterproof. To prevent damage to the equipment, contamination, or bodily injury, do not immerse or soak any of these components in water or a chemical solution. Do not spray disinfectant or other fluids directly onto the motor component, foot pedal, power adapter, or charging station to avoid liquid from pooling on these handpiece components. The user should spray solution onto a cloth or use a wipe to disinfect these items per the instructions.

### Manual Cleaning Steps

1. Place the motor component, foot pedal, power adapter, and charging station onto a clean surface. Wet a soft cloth with distilled water or deionized water. Thoroughly wipe all the surfaces of the motor component, foot pedal, power adapter, and charging station, etc. until the surface is clean. **DO NOT OVERSATURATE THE CLOTH.**
2. Wipe the surface of components with a dry, soft cloth.
3. Perform the above steps until all visible soil is removed.

## Manual Disinfection Steps

1. Dampen a soft cloth with 70% isopropyl alcohol.
2. Wipe all surfaces of the motor component, foot pedal, power adapter, and charging station with a wet soft cloth for at least 3 minutes or recommended disinfection time of the cleaner.
3. Wipe the surface of the motor component with a dry, clean cloth.
4. Visually inspect to ensure all contaminants have been removed and inspect power supply cord for damage.
5. After cleaning and disinfecting the motor component, install a disposable barrier sleeve before use.

## 9. STORAGE AND DISPOSAL

---

Store in a clean, dry, ventilated, non-corrosive atmosphere with a humidity of 10% to 93%, an atmospheric pressure of 70 KPa to 106 KPa, and a temperature of -20 °C to +55 °C .

After sterilization, the product should remain bagged until ready for use. Product storage must be batched, marked, and recorded.

Dispose of waste according to local rules, guidelines, and regulations.

### Disposal of components

- Dispose of the charging station, foot pedal, and power adapter as industrial waste according to local laws and regulations.
- Dispose of the outer sheath, disposable sleeve, and disposable prophylaxis angle as medical waste according to local laws and regulations.
- The motor component and foot pedal include a battery pack (lithium-ion battery). Do not dispose of in a fire or source of heat. Doing so may cause rupture of the battery pack, scattering of the battery fluid, overheating, smoking, or explosion.

## 10. TECHNICAL CONSIDERATIONS

### Technical Information/Data




ITEMS	SPECIFICATION
Power source	Handpiece: rechargeable lithium battery 3.6V, 750mAh; Foot pedal: rechargeable lithium battery 3.6V, 750mAh
Rated input	100V-240V~, 50Hz/60Hz, 0.4A Max
Rated output	DC5V, 1A
Weight	Handpiece: 121.6 g Foot pedal: 196.5 g
Dimension	Handpiece: 27.6 mm Dia×192 mm Foot pedal: 128 mm×128 mm×34 mm
Type of protection against electric shock	Class II
Degree of protection against electric shock	Type B applied part
Degree of protection against harmful ingress of liquid	<b>Handpiece: IPX3; Foot pedal: IPX1; Charging station: IPX0; Power adapter: IPX0</b>
Operating environment	Ambient temperature: +5 °C ~ +40 °C Relative humidity: 30% ~ 75% Atmospheric pressure: 70 kPa ~106 kPa
Transport and storage condition	Ambient temperature: -20 °C ~ +55 °C Relative humidity: 10% ~ 93% Atmospheric pressure: 70 kPa ~106 kPa
Handpiece speeds (+/-10%)	Speed range: 500 rpm–4000 rpm
Foot pedal speeds (+/-10%)	Speed is controlled and adjusted through varying pressure on the foot pedal. The adjustable range is 500 rpm to 4000 rpm
Maximum torque (+/-10%)	1.2 Ncm
Service life	5 years
Charge time	Handpiece: Approximately 2.5 hours Foot pedal: Approximately 2 hours
Auto-off time	The foot pedal and handpiece will automatically shut down if the standby time exceeds 5 minutes

## Troubleshooting

FAILURE	SOLUTION
Disposable prophy angle is not rotating	<ol style="list-style-type: none"> <li>1. Ensure the outer sheath and the disposable prophy angle are attached securely.</li> <li>2. Ensure the motor component and the outer sheath are attached securely.</li> <li>3. Ensure the disposable prophy angle is not damaged.</li> <li>4. Ensure the outer sheath is not damaged.</li> </ol>
Excessive noise or vibration during operation	<ol style="list-style-type: none"> <li>1. Ensure the outer sheath and the disposable prophy angle are attached securely.</li> <li>2. Ensure the motor component and the outer sheath are attached securely.</li> <li>3. Ensure the disposable prophy angle is not damaged.</li> <li>4. Ensure the outer sheath is not damaged.</li> <li>5. Ensure the O-ring on the motor component or outer sheath is not damaged.</li> </ol>
Difficulty removing or installing the outer sheath	<ol style="list-style-type: none"> <li>1. Ensure the outer sheath is not damaged.</li> <li>2. Check and clean the outer sheath.</li> </ol>
Motor component cannot be charged	<ol style="list-style-type: none"> <li>1. Check for foreign matter between charging station and motor component.</li> <li>2. Check whether the original power adapter is used.</li> </ol>
The foot pedal is unable to connect to the handpiece	<ol style="list-style-type: none"> <li>1. Place the handpiece into the wireless charging base until the indicator light of the charging base flashes. Long press the on/off button on the handpiece while in the charging state. At this time, you will hear a beep; continue pressing the handpiece button.</li> <li>2. Using a thin paperclip-like wire, long press into the hole marked "📶" at the bottom of foot pedal for about 3–5 seconds until the Bluetooth indicator turns from a flashing to steady light. Immediately release the pressed button in the "📶" hole of the foot pedal and the pressed handpiece button to ensure a Bluetooth connection. In the case that this does not work on the first attempt, repeat the steps until there is a Bluetooth connection.</li> </ol>

## EMC-Declaration of Conformity

The device has been tested and homologated in accordance with EN 60601-1-2 for EMC. This does not guarantee in any way that this device will not be affected by electromagnetic interference. Avoid using the device in a high electromagnetic environment.

	Appliance compliance WEEE directive
<b>IPX0</b>	Protection Class IPX0 - IPX0 Classification of ingress of water for charging station – not protected
<b>IPX1</b>	Protection Class IPX1 - IPX1 Classification of ingress of water for foot pedal
<b>IPX3</b>	Protection Class IPX3 - IPX3 Classification of ingress of water for handpiece - Protected against falling spray
	Temperature limitation
	Atmospheric pressure for storage



## Technical Description Concerning Electromagnetic Emission

Table 1: Declaration - Electromagnetic Emissions

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS		
The Ultrapro Tx Cordless handpiece is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultrapro Tx Cordless handpiece should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Ultrapro Tx Cordless handpiece uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	The Ultrapro Tx Cordless handpiece is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

## Technical Description Concerning Electromagnetic Immunity

Table 2: Guidance & Declaration - Electromagnetic Immunity

GUIDANCE & DECLARATION - ELECTROMAGNETIC IMMUNITY			
The Ultrapro Tx Cordless handpiece is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultrapro Tx Cordless handpiece should ensure it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4- 2	±8kV contact ±2, ±4, ±8, ±15kV air	±8kV contact ±2, ±4, ±8, ±15kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4- 4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4- 5	±0.5, ±1kV line to line ±0.5, ±1, ±2kV line to earth	±0.5, ±1kV line to line ±0.5, ±1, ±2kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4- 11	<5 % UT (>95% dip in UT.) for 0.5 cycle  <5 % UT (>95% dip in UT.) for 1 cycle  70% UT (30% dip in UT) for 25 cycles  <5% UT (>95 % dip in UT) for 250 cycles	<5 % UT (>95% dip in UT.) for 0.5 cycle  <5 % UT (>95% dip in UT.) for 1 cycle  70% UT (30% dip in UT) for 25 cycles  <5% UT (>95 % dip in UT) for 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ultrapro Tx Cordless handpiece requires continued operation during power mains interruptions, it is recommended that the Ultrapro Tx Cordless handpiece be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4- 8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

*NOTE: UT is the a.c. mains voltage prior to application of the test level.*

Table 3: Guidance & Declaration - Electromagnetic Immunity Concerning Conducted RF and Radiated RF


GUIDANCE & DECLARATION - ELECTROMAGNETIC IMMUNITY			
The Ultrapro Tx Cordless handpiece is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultrapro Tx Cordless handpiece should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC61000-4- 6 Conducted RF IEC61000-4- 6 Radiated RF IEC61000-4- 3	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM frequency band 3 V/m 80 MHz to 2.7 GHz	3V 6V 3V/m	Portable and mobile RF communications equipment should be used no closer to any part of the Ultrapro Tx Cordless handpiece, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2 \times P^{1/2}$ 80 MHz to 800MHz $d=2.3 \times P^{1/2}$ 800 MHz to 2.7GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a) should be less than the compliance level in each frequency range. b) Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			
a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Ultrapro Tx Cordless handpiece is used exceeds the applicable RF compliance level above, the Ultrapro Tx Cordless handpiece should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Ultrapro Tx Cordless handpiece.			
b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

Table 4: Recommended separation distances between portable and mobile RF communications equipment and the Ultrapro Tx Cordless handpiece.

Recommended separation distances between portable and mobile RF communications equipment and the Ultrapro Tx Cordless handpiece.			
The Ultrapro Tx Cordless handpiece is intended for use in electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Ultrapro Tx Cordless handpiece can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Ultrapro Tx Cordless handpiece as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d=1.2 \times P^{1/2}$	80MHz to 800MHz $d=1.2 \times P^{1/2}$	800MHz to 2,7GHz $d=2.3 \times P^{1/2}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			





















## 11. MISCELLANEOUS INFORMATION

### Transportation

- Prevent excessive shock and vibration during transportation and handle with care.
- Device should not be mixed with dangerous goods during transportation.
- Avoid exposure to direct sunlight or water during transportation to avoid damage.

Report any serious incident to the manufacturer and the competent authority.  
www.ultradent.com / 1.800.552.5512 / 801.572.4200

## 12. SYMBOL INSTRUCTION

	Follow Instructions for Use		Serial number
	Date of manufacturing		Manufacturer
	Type B applied part		Class II equipment
	Handle with care		Keep dry
	Humidity limitation: indicates the acceptable upper and lower limits of relative humidity for transport and storage		Atmospheric pressure limitation: indicates the acceptable upper and lower limits of atmospheric pressure for transport and storage
	Temperature limitation: indicates the maximum and minimum temperature limits at which the item shall be stored and transported		Caution: Federal law restricts this device to sale by or on the order of a licensed dental professional
	Degree of protection against harmful ingress of liquid		Electromagnetic interference from the device is under limits approved by the Federal Communications Commission
	Dispose of in accordance with the Waste Electrical and Electronic Equipment Directive 2002/96/EC for product and accessory disposal		Sterilize up to temperature specified
	Warning		For indoor use only
	Do not reuse: for disposable prophylactic angles and disposable barrier sleeves		Serial number barcode





 ultrapro™ Tx  
cordless



 **ultrapro™ Tx**  
cordless