## ULTRADENT PRODUCTS, INC. FORMA™

Light-cured Composite (2101330)

**ENI Product Description:** FORMA™ restorative composite contains bis-GMA, TEGDMA, Bis-EMA, and UDMA. FORMA™ composite is designed for restorations in anterior and posterior teeth and is suitable for occlusal surfaces. Its formula contains a combination of inorganic filler based on zirconia/silica and barium glass and radiopaque ytterbium trifluoride filler. All shades of FORMA™ composite are radiopaque. FORMA™ composite has a radiopacity value greater than 2 mm of aluminum. Aluminum has a radiopacity equivalent to that of dentine. Thus 1 mm of material having a radiopacity equivalent to 1 mm of aluminum has a radiopacity equivalent to that of dentine and 2 mm of aluminum is equivalent to enamel.

Filler load is of 66% - 72% by volume for the body, enamel and dentin shades. FORMA™ composite is a nano-hybrid composite with nanoparticles in the 5–50 nm range. For all products described, carefully read and understand all instructions and SDS information prior to use.

### Presentation Form(s):

FORMA™ Syringe

**Indications for Use/Purpose of Use:** FORMA™ is an aesthetic composite resin for use in restorations of posterior and anterior teeth.

### Contraindications:

Do not use on patients with a known sensitivity to acrylates or other resins.

### Precautions and warnings:

1. Read and understand all instructions before using product.
2. Use a quality curing light to ensure maximum physical properties of the composite. Uncured material will result in an inferior restoration and may also cause pulpal irritation.
3. Clinician and patient should wear UV and blue light protective eyewear when curing composite.
4. To lessen pre-polymerization of the composite, minimize exposure to bright lights during the procedure.

5. To optimize bond to each succeeding layer, do not remove oxygen inhibition layer between incremental layers of composite.
6. Ensure composite is totally polymerized between layers and interproximally.
7. Syringes are not for intraoral use.
8. Resins can be sensitizing. Avoid repeated contact of uncured dental resin with skin. Do not use on patients with a known sensitivity to acrylates or other resins. If allergic reaction, dermatitis or rash develops, consult a physician.
9. To avoid pre-polymerization of the composite by ambient light, make sure to always cover the tip of the syringe with its cap after composite has been removed.
10. To avoid cross-contamination, use a clean instrument every time you remove the composite from the syringe.
11. Disinfect contaminated syringes with an intermediate level disinfectant between uses.
12. Keep away from extreme heat and direct sunlight.
13. Do not use after expiration date noted on containers.
14. Keep out of reach of children.

### Step by Step Instruction:

#### Preparation:

1. Apply dentin bonding agent to preparation according to manufacturer's instructions.

#### Direct Restoration:

1. Open cap and twist handle (stem) clockwise to express desired amount of composite. Resistance will lessen with each twist.
2. Place FORMA™ composite in incremental layers applying the first layer approximately 0.5 – 1.5 mm thick and cure. Follow with up to 2 mm increments; curing after each increment. For curing lights with an intensity of  $\geq 800$  mW/cm<sup>2</sup> in the wavelength region of 385 – 515 nm, cure each layer for 10 seconds. A final cure of 20 seconds is recommended.
3. After final increment is placed, sculpt and light cure.
4. Shape, finish, and adjust occlusion using preferred method.
5. Immediately replace cap after dispensing composite to prevent polymerization from ambient light.

#### Indirect Restorations:

1. Make impression of prepared tooth and opposing dentition.
2. Prepare master cast and articulate opposing cast using preferred technique.
3. Place FORMA™ composite in incremental layers applying the first layer approx-

imately 0.5 – 1.5 mm thick and cure. Follow with up to 2 mm increments; curing after each increment. For curing lights with an intensity of  $\geq 800$  mW/cm<sup>2</sup> in the wavelength region of 385 – 515 nm, cure each layer for 10 seconds.

4. Use final increments to create proximal contours and ensure marginal adaptation, a final cure of 20 seconds is recommended.
5. Remove restoration. Clean off residual die stone and trim marginal flash. Additions can be made by roughening restoration surface, applying adhesive, and placing FORMA™ composite as needed.
6. Sandblast internal aspect of restoration with 50  $\mu$ m alumina abrasive.
7. Seat restoration onto master cast, evaluate fit, occlusion, and contours. Shape, finish, and adjust occlusion using preferred method.
8. Bond restoration intraorally per adhesive manufacturer's instructions.
9. Following cementation, make necessary occlusal adjustment and polish adjusted surfaces.

### Storage, Transportation, Handling and Disposal:

Store at temperature within 2–27°C.

Properly dispose of syringes, according to current legislation.

**NOTE:** for product SDS and additional information about FORMA composite, please see website: [www.ultradent.com.br](http://www.ultradent.com.br).

**FOR PROFESSIONAL USE ONLY**

DO NOT RE-USE TO AVOID CROSS CONTAMINATION

KEEP AWAY FROM HEAT/SUNLIGHT

FOR PROFESSIONAL USE ONLY

RECOMMENDED STORAGE TEMPERATURE (2-27°C)

REFER TO INSTRUCTIONS FOR INFORMATION ON HOW TO USE

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Ultradent Australia Pty Ltd.