

VISCOSTAT™ HEMOSTATIC



3093 - ViscoStat Dento-Infusor Syringe Kit
4 x 1.2 ml syringes
20 x Metal Dento-Infusor tips



647 - ViscoStat Dento-Infusor IndiSpense Syringe Kit
1 x 30 ml IndiSpense syringe
20 x Metal Dento-Infusor tips
20 x 1.2 ml empty syringes



3070 - ViscoStat Syringe 20pk
1.2 ml syringes



645 - ViscoStat IndiSpense Syringe 1pk
30 ml syringe



1278 - ViscoStat Empty Syringe 20pk
3096 - ViscoStat Empty Syringe 100pk
1.2 ml empty syringes

ASTRINGEDENT™ HEMOSTATIC



111 - Astringedent Bottle 1pk
686 - Astringedent IndiSpense Syringe 1pk
30 ml



2160 - Astringedent Spot Remover 1pk
30 ml bottle

ASTRINGEDENT™ X HEMOSTATIC



112 - Astringedent X Bottle 1pk
690 - Astringedent X IndiSpense Syringe 1pk
30 ml

METAL DENTO-INFUSOR™ TIPS

LOK-TITE™	100pk	500pk
19 ga Metal Dento-Infusor	2559	2560



LOK-TITE™	mm	100pk	500pk
Blue Mini Dento-Infusor	1.20	1086	1440



! Note: Do not use epinephrine preparations with ferric sulfate products (ViscoStat, Astringedent), as blue/black precipitate will occur.

For Profound Hemostasis

Tissue Management

Unparalleled tissue management starts with rapid, profound hemostasis. For more than 40 years, dentists have trusted the immediate hemostatic power, detailed margins, and elimination of surface bleeding and sulcular fluid provided by Ultradent's tissue management products.



Tissue Management

Technique Guide



1 Bleeding subgingival tissues adjacent to preparation. Bleeding must be controlled before starting any procedure.



2 The Metal Dento-Infusor™ tip has a padded brush end ideal for infusing bleeding capillaries with hemostatic agents. Burnish hemostatic agents firmly against sulcus until the bleeding stops.



3 Apply a firm air/water spray to remove residual coagulum.



4 Check tissue for quality, profound hemostasis. Repeat steps 2-3 if bleeding reoccurs.



5 Soak Ultrapak™ knitted retraction cord in hemostatic agent. Pack soaked cord into sulcus and leave in place for 1-3 minutes. This assists with tissue displacement.



6 If a double cord technique is used, remove the upper cord, apply a firm air/water spray, and dry.



VISCOSTAT™ HEMOSTATIC

20% Ferric Sulfate

Viscous formula allows for superior control

- Provides unsurpassed kindness to tissues
- Creates rapid, profound hemostasis
- Viscous solution won't run, drip, or migrate to other teeth



ASTRINGEDENT™ HEMOSTATIC

15.5% Ferric Sulfate

Known as the "classic" hemostatic agent—Ultradent's first product

- Stops moderate bleeding in seconds
- Eliminates sulcular fluid contamination for optimal bonding



ASTRINGEDENT™ X HEMOSTATIC

12.7% Iron Solution

Use for difficult surface capillary bleeding


- Effectively treats difficult-to-stop bleeding
- Is less acidic than competitive iron subsulfate solutions

Use Astringedent X hemostatic solution when a stronger, more potent hemostatic is required and when the attainment of quality hemostasis may be more challenging (e.g., in cases of difficult-to-stop, problem bleeding).

TEST YOUR BURNISHING SKILLS

Are you using adequate pressure to achieve rapid, profound hemostasis?

Directions: Use your Dento-Infusor™ tip to scrub the scratch-off area.



METAL DENTO-INFUSOR™ TIP

- Effectively places hemostatic agents and removes superficial coagulum
- Has a blunt, bent cannula with a padded "scrub brush" for infusing and cleaning the cut sulcus



BLUE MINI™ DENTO-INFUSOR™ TIP

- Offers the same tissue management benefits as the Metal Dento-Infusor tip
- Allows controlled flow of drop-sized quantities