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
Opalescence PF 10% Mint:
 Ingredients: Glycerin, Water (Aqua), Xylitol, Urea (Carbamide) Peroxide, Carbomer, PEG-6, Sodium Hydroxide, EDTA, Peppermint Oil (Aroma), Potassium Nitrate, Sodium Fluoride

Opalescence PF 10% Melon:
 Ingredients: Glycerin, Water (Aqua), Xylitol, Urea (Carbamide) Peroxide, Carbomer, PEG-6, Sodium Hydroxide, Watermelon Flavor (Aroma), EDTA, Potassium Nitrate, Sodium Fluoride

Opalescence F 15% Melon:
 Ingredients: Glycerin, Water (Aqua), Urea (Carbamide) Peroxide, Xylitol, PEG-6, Carbomer, Sodium Hydroxide, Watermelon Flavor (Aroma), EDTA, Sodium Fluoride

Opalescence PF 15%, 16% and 20% Regular:
 Ingredients: Glycerin, Water (Aqua), Urea (Carbamide) Peroxide, Xylitol, Carbomer, PEG-6, Sodium Hydroxide, EDTA, Potassium Nitrate, Sodium Fluoride

Opalescence PF 15%, 16% and 20% Mint:
 Ingredients: Glycerin, Water (Aqua), Urea (Carbamide) Peroxide, Xylitol, Carbomer, PEG-6, Sodium Hydroxide, EDTA, Peppermint Oil (Aroma), Potassium Nitrate, Sodium Fluoride

Opalescence PF 15%, 16% and 20% Melon:
 Ingredients: Glycerin, Water (Aqua), Urea (Carbamide) Peroxide, Xylitol, Carbomer, PEG-6, Sodium Hydroxide, Watermelon Flavor (Aroma), EDTA, Potassium Nitrate, Sodium Fluoride

Opalescence 35% Regular:
 Ingredients: Glycerin, Water (Aqua), Urea (Carbamide) Peroxide, Carbomer, PEG-6, Hydrogen Peroxide, Xylitol, Sodium Hydroxide, EDTA, Potassium Nitrate, Sodium Fluoride, Sucralose

Opalescence 35% Mint:
 Ingredients: Glycerin, Water (Aqua), Urea (Carbamide) Peroxide, Carbomer, PEG-6, Hydrogen Peroxide, Xylitol, Sodium Hydroxide, EDTA, Peppermint Oil (Aroma), Potassium Nitrate, Sodium Fluoride, Sucralose

Opalescence 35% Melon:
 Ingredients: Glycerin, Water (Aqua), Urea (Carbamide) Peroxide, Carbomer, PEG-6, Hydrogen Peroxide, Xylitol, Sodium Hydroxide, Watermelon Flavor (Aroma), EDTA, Potassium Nitrate, Sodium Fluoride, Sucralose

Opalescence Quick 45%:
 Ingredients: Water, Glycerin, Carbamide Peroxide, Hydrogen Peroxide, Carbomer, Polyethylene Glycol, Xylitol, Sodium Hydroxide, Potassium Nitrate, Sodium Fluoride, EDTA, Peppermint Oil (Aroma), Sucralose.

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

EN - PATIENT INSTRUCTIONS

1. If you have any questions regarding the appropriate use of this product, including how long it will take to whiten your teeth, please consult your dentist.

2. Regular dental checkups and cleanings are important before and after whitening to maintain a healthy smile.

3. Some patients experience sensitivity with whitening; should this occur, stop treatment and contact your dentist.

4. Some patients have noticed temporary discomfort of the gums, lips, throat, or tongue. Should any of these symptoms persist more than two days or progressively worsen, call your dentist. These side effects will usually subside within 1-3 days after treatment is discontinued.

5. Pregnant or breastfeeding women, or patients with serious health concerns should not whiten. Consult their primary care provider prior to treatment.

6. If patient has a known allergy or chemical sensitivity to any of the ingredients in Opalescence whitening, do not use. Mint or melon flavor can be irritating to some patients. If a patient has a known sensitivity to flavorings, the non-flavored variant should be selected.

7. Do not use household bleach to whiten teeth.

8. Do not use tobacco or eat while whitening trays are in place.

9. Teeth are naturally darker along the gum line. These areas may require more time to lighten than the rest of the tooth surface and usually remain slightly darker.

10. White spots that appear on the enamel during the whitening process may blend with continued whitening.

11. Crowns, bridges, partial dentures, veneers, and composite fillings will not whiten.

12. Coffee, tobacco, and other products can re-stain your teeth over time. Should this occur, teeth can be re-whitened.

13. Some amalgam fillings may leave a dark purple color in your whitening tray; this is normal.

14. Check expiration date and store according to instructions.

15. Avoid swallowing large amounts of gel or water used to rinse gel residue from the teeth.

16. Avoid contact with eyes.