**EDTA 18% Solution**  
**Filing and Chelating Agent**

Description & Indications:  
Ultradent EDTA is an 18% colorless, low-viscosity, aqueous, dentin-chelating solution. It assists in root canal preparation by decalcifying dentin and removing the smear layer. It is recommended to use EDTA as a final rinse before obturating canal with EndoREZ™ or any other sealer. EDTA is used as a final irrigant for 60 seconds prior to obturation.

Preliminary preparation steps one minute or best results. Quick and effective.  
Prepare IndSpense™ System delivery:  
1. An empty 1.2ml or 5ml syringe is attached to the end of the large, no-waste IndSpense syringe by turning the Luer Lock of the unit dose syringe snugly onto the male thread of the IndSpense as shown in Fig. 1.  
2. With palm grip (Fig. 2), the plunger of the IndSpense syringe is depressed while the unit dose syringe is stabilized at the plunger with the non-dominant hand. Slightly pull out on plunger of unit dose syringe to facilitate loading. Do not express contents of unit dose syringe back into IndSpense syringe. To prevent cross-contamination, a used syringe should not be reattached to the IndSpense syringe.  
3. Remove Luer Lock cap from syringe. Twist a Navifx™ Tip or Navitip™ firmly onto syringe. Gauge Navifx Tip with the rubber stop at least 2mm short of the root terminus as show in Fig. 3.  
Procedure:  
1. Make sure Navifx Tip is loose in canal and at least 2mm short of terminus. Gently deliver the EDTA with Navifx Tip. Move brush stop up and down to agitate and brush solution, cleaning walls with EDTA as shown in Fig. 4, alternating with sodium hypochlorite after each file during instrumentation. Do not express EDTA when canula tip is at full length or if it fits tight against the canal walls.  
2. Perform a final rinse with clean, bacteria-free water, irrigant anesthetic, or Plain EDTA for separation steps one minute or best results. Quick and effective removal of most of the EDTA solution and other irrigants is facilitated by using a Luer Vacuum Adapter with a Capillary Tip or 10mm Micro Capillary. One or two paper points are usually all that are required for final drying of canals, if canals are first vacuumed out as per step 2 above.  
4. Fit master cone, then obturate using EndoREZ™ hydrophilic sealer and gutta percha, or your preferred obturation method.  
Precautions:  
Always use a rubber dam for endodontic therapy. Prevent ingestion of EDTA solution.  
Prevent solution from getting in the eyes or leaving it on the soft tissues for an extended period of time.  
Do not extend tip past root apex or express EDTA past apex.  
Do not express EDTA from canula tip at the root terminus or if the tip fits tight against canal walls. This may force solution past the apex.  
Do not use this product on patients with a known sensitivity to EDTA. The tips and delivery syringes are disposable and intended for single use.

General Precautions:  
For professional use only.  
Review instructions, precautions, and MSDS before beginning treatment.  
Use only as directed.  
Keep products out of heat/sunlight.  
Avoid skin exposure to resins.  
Isolate straw chemicals to area of treatment.  
Confirm that patient has no known allergies to treatment materials.  
Test food of materials from syringe and tip before using intraorally.  
Never force syringe plungers.  
Use Ultradent™ syringe covers and/or clean and disinfect syringes between patients.

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