

Opalescence™
whitening toothpaste



Keep Your Smile Healthy and Bright

- Actively removes surface stains for whiter teeth
- Formulated to strengthen enamel and help prevent cavities
- Developed by a dentist
- Gentle enough for everyday use
- Unique Cool Mint flavor for fresh breath that lasts
- Available in a Sensitivity Relief formula to treat sensitive teeth

Ask your dental professional about
Opalescence Whitening Toothpaste
Learn more at opalescence.com



For an SDS of this product, please see www.ultradent.com
© 2025 Ultradent Products, Inc. All rights reserved.
Manufactured by Ultradent Products, Inc. 505 West Ultradent Drive,
South Jordan, Utah 84095 www.opalescence.com
Made in the USA from globally sourced materials.



1007486AR08 092525

Opalescence



Professional whitening to go

Get more confidence in your smile

Your smile says more about you than any other physical feature. It's the first thing people notice. Wouldn't it be great to have a simple, affordable method for whitening your smile that's convenient enough to take with you anywhere?

Introducing Opalescence Go whitening, the easiest, fastest, and most convenient way to brighten your smile. Opalescence Go whitening gives you professional results in the convenience of pre-loaded disposable trays, so it is ready to go when you are. Contact your dental professional to purchase Opalescence Go whitening.

Instructions

Instructions



1. Remove product from packaging.
 'U' • Upper whitening tray
 'L' • Lower whitening tray



2. Position upper tray on teeth.



3. Bite firmly, then suck on tray for 2 seconds.



4. Remove colored outer tray, leaving the white inner tray on teeth. Repeat the process for the lower tray.



5. After indicated wear time, remove whitening trays and brush teeth.

 10% | 30–60 min.

 15% | 15–20 min.

Precautions

1. Limit contact with soft tissue.
2. Avoid direct contact of the active surface of the tooth whitening product with the gums and/or salivary flow.
3. Some patients may experience temporary discomfort of the gums, lips, throat, tongue, or teeth. Should any of these symptoms persist or worsen after discontinuing use, or in cases of severe discomfort, call your dental professional.
4. Avoid swallowing large amounts of gel. Ingestion of small amounts of gel during use is expected.
5. If the product is ingested, drink large amounts of water and contact a physician. If the product is aspirated, seek immediate medical attention.
6. Restorative materials will not whiten. We recommend whitening dark teeth before restorative placement.
7. Wait at least 7–10 days after whitening procedures before placing adhesive restorations.
8. Those with serious health concerns should consult their primary care provider prior to treatment.
9. Avoid contact with eyes. Rinse immediately if product comes into contact with them.
10. To prevent additional irritation, gingival and general oral health should be evaluated prior to initiating treatment.
11. Monitor patient's whitening progress to prevent whitening beyond the desired shade or translucency.
12. Do not eat, drink, or use tobacco products while whitening trays are in place.
13. Do not wear whitening tray while sleeping.
14. Bring treatment tray to room temperature before use.
15. Pregnant or nursing women should not use teeth whitening products.
16. For patients or users with allergy concerns, refer to product allergen document available at www.ultradent.com. If allergic reaction is observed, rinse exposed area thoroughly with water and have the patient consult their physician.
17. Contains sodium fluoride.

Storage and Disposal

- Store according to outer package label.
- Dispose of waste according to local rules, guidelines, and regulations.

Miscellaneous Information

- Hypocalcified areas (white spot lesions), which may not be visible to the naked eye, will whiten faster (thereby becoming more obvious during whitening). Continue whitening treatment until the unaffected tooth surface blends. Re-evaluate two weeks after whitening treatment when tooth color has stabilized.
- For a small percentage of the population, their teeth will not appreciably whiten
- Some amalgam fillings may leave a dark purple color in your whitening tray. This is normal
- Teeth can re-stain over time. So, do touch up teeth whitening treatments as necessary.

Report any serious incident to the manufacturer and the competent authority.
ultradent.com / 1.800.552.5512 / 801.572.4200



Do not re-use to avoid cross contamination



Keep away from heat/sunlight



Refrigerate



See instructions



Keep out of reach of children