

Opalescence[™]
tooth whitening systems



Ask your dentist today about your
teeth whitening options!



Get a brighter, more confident smile.



Health Hazard



Do not re-use to avoid
cross contamination



5% 40%



Keep away from
heat/sunlight



Recycle



Keep out of reach
of children



Kosher



See instructions

For an SDS of this product, please see www.ultradent.com

© 2022 Ultradent Australia Pty Ltd. P.O. Box 111 Coogee, NSW 2034 All Rights Reserved. Made in the USA from U.S. and globally sourced materials

Manufactured by Ultradent Products, Inc.

505 West Ultradent Drive (10200 South), South Jordan, Utah 84095

1007778AR03 120822



Opalescence[™]
● ● ● go[™]

Professional whitening to go

Get more confidence in your smile

Your smile says more about you than any other physical feature. It's the first thing people notice. Wouldn't it be great to have a simple, affordable method for whitening your smile that's convenient enough to take with you anywhere?

Introducing Opalescence Go whitening, the easiest, fastest, and most convenient way to brighten your smile. Opalescence Go whitening gives you professional results in the convenience of pre-loaded disposable trays, so it is ready to go when you are.

Instructions

Instructions



1. Before using, allow product to warm to room temperature. This will make it easier for the tray to adhere to the teeth. Remove from package.



2. Center tray on arch.



3. Gently suck down on the tray. Do not press tray against teeth as it pushes the whitening gel against the gums, causing tissue irritation.



4. Remove outer tray.



5. Suck down.



6. Wear Opalescence Go whitening 6% for 60–90 minutes, 10% for 30–60 minutes, and 15% for 15–20 minutes per day. When finished, remove whitening tray and brush teeth.

For professional use only

6% | 60–90

10% | 30–60

15% | 15–20

Precautions

1. Limit contact with soft tissue.
2. Avoid direct contact of the active surface of the tooth whitening product with the gums and/or salivary flow.
3. Some patients may experience temporary discomfort of the gums, lips, throat, tongue, or teeth. Should any of these symptoms persist or worsen after discontinuing use, or in cases of severe discomfort, call your dental professional.
4. Avoid swallowing large amounts of gel. Ingestion of small amounts of gel during use is expected.
5. If the product is ingested, drink large amounts of water and contact a physician. If the product is aspirated, seek immediate medical attention.
6. Restorative materials will not whiten. We recommend whitening dark teeth before restorative placement.
7. Wait at least 7–10 days after whitening procedures before placing adhesive restorations.
8. Those with serious health concerns should consult their primary care provider prior to treatment.
9. Avoid contact with eyes. Rinse immediately if product comes into contact with them.
10. To prevent additional irritation, gingival and general oral health should be evaluated prior to initiating treatment.
11. Monitor patient's whitening progress to prevent whitening beyond the desired shade or translucency.
12. Do not eat, drink, or use tobacco products while whitening trays are in place.
13. Do not wear whitening tray while sleeping.
14. Bring treatment tray to room temperature before use.
15. Pregnant or nursing women should not use teeth whitening products.
16. For patients or users with allergy concerns, refer to product allergen document available at www.ultradent.com. If allergic reaction is observed, rinse exposed area thoroughly with water and have the patient consult their physician.
17. Contains sodium fluoride.

Mint Ingredients: Glycerin, Water (Aqua), Hydrogen Peroxide, PVP, Carbomer, Potassium Hydroxide, Silica, Dipotassium Phosphate, Xylitol, Potassium Nitrate, Flavor (Aroma), Citric Acid, Sucralose, Sodium Fluoride, EDTA, Sodium Lauryl Sulfate

Melon Ingredients: Glycerin, Water (Aqua), Hydrogen Peroxide, PVP, Carbomer, Potassium Hydroxide, Flavor (Aroma), Silica, Dipotassium Phosphate, Xylitol, Potassium Nitrate, Citric Acid, Sucralose, Sodium Fluoride, EDTA, Sodium Lauryl Sulfate