opal® bond EN

Description/Indications:

Opal® Bond and Opal® Bond MV are light cure orthodontic adhesives designed for bonding brackets and other bondable appliances to etched enamel.

Bonding Procedure:

1. Clean enamel surfaces to be bonded with pumice or a similar water-based, mild abrasive.
2. Etch, rinse and dry all tooth surfaces prior to bonding.
3. Apply a thin layer of orthodontic primer to surface of etched tooth.
4. Apply a thin layer of light cure Opal Bond onto bonding surface of bracket or other bondable appliance.
5. Place bracket or appliance coated with Opal Bond onto etched tooth.
6. Remove any excess bonding material from circumference of bracket base.
7. Recommended curing time for halogen lights is 15 seconds on the mesial and 15 seconds on the distal of each tooth. Recommended curing time for LED lights can vary. Refer to curing light manufacturer's recommended curing times.

Precautions:

1. For Professional use only.
2. Review instructions, precautions and SDS before beginning treatment. Use only as directed.
3. Keep products out of heat/sunlight.
4. Use quality curing lights.
5. Do not use on patients with contraindications to orthodontic treatment.
6. Cure bonding resins completely, as uncured resin materials will be deficient in physical properties, stress and wear resistance.
7. Any resin can be sensitizing to skin or can cause allergic reactions if individual is already sensitized. If seen, wash area of tissue exposed to resin completely. If dermatitis, rash or any other allergic reaction is observed, consult a physician.
8. When not in use, ensure cap is secured over the opening of the Opal Bond syringe or compule.
9. Clean and disinfect syringes and compules between patients.

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