

PEAK™

UNIVERSAL ADHESIVE SYSTEM

THE GOLD STANDARD IN ADHESIVE DENTISTRY



4554 - Peak Universal Bond Self-Etch Syringe Intro Kit
1 x 1.2 ml Peak Universal Bond syringe
1 x 1.0 ml Peak SE Primer syringe
20 x Black Mini Brush tips
20 x Inspiral Brush tips



4551 - Peak Universal Bond Total-Etch Syringe Intro Kit
1 x 1.2 ml Peak Universal Bond syringe
1 x 1.2 ml Ultra-Etch syringe
20 x Blue Micro tips
20 x Inspiral Brush tips



4541 - Peak Universal Bond Self-Etch Bottle Kit
1 x 4 ml Peak Universal Bond bottle
4 x 1.0 ml Peak SE Primer syringes
40 x Black Mini Brush tips
50 x Mixing wells
50 x Micro Applicator brushes



4542 - Peak Universal Bond Total-Etch Bottle Kit
1 x 4 ml Peak Universal Bond bottle
4 x 1.2 ml Ultra-Etch syringes
40 x Blue Micro tips
50 x Mixing wells
50 x Micro Applicator brushes

REFRIGERATE

4543 - Peak Universal Bond Bottle 1pk
4 ml bottle



REFRIGERATE

4553 - Peak Universal Bond Syringe 4pk
4552 - Peak Universal Bond Syringe 20pk
1.2 ml syringes



REFRIGERATE

5006 - Peak Universal Bond Unit Dose 50pk
0.2 ml unit dose



REFRIGERATE

5135 - Peak SE Primer Syringe 4pk
1.0 ml syringes



4548 - Micro Applicator Brush 400pk
200 x Each color



4545 - Mixing Wells 100pk



Peak Universal Bond adhesive ensures a strong, lasting bond whether you choose a self-etch or a total-etch technique.

- Highest bond strength on both dentin and enamel¹
- Bonds to dentin, enamel, porcelain, metal, composite, and zirconia
- Ideal for indirect and direct bonding, as well as post and core procedures
- Works with self-etch and total-etch techniques
- Contains 0.2% chlorhexidine to help ensure long-term bond strength^{2,3}




SCAN TO VIEW PRICING

1. Data on file. 2. Breschi L, Maravic T, Comba A, et al. Chlorhexidine preserves the hybrid layer in vitro after 10-years aging. *Dent Mater.* 2020;36(5):672-680. doi:10.1016/j.dental.2020.03.009 3. Yaghmoor RB, Jamal H, Abed H, Allan E, Ashley P, Young A. Incorporation of MMP inhibitors into dental adhesive systems and bond strength of coronal composite restorations: A systematic review and meta-analysis of in vitro studies. *Jpn Dent Sci Rev.* 2022;58:298-315. doi:10.1016/j.jdsr.2022.09.004 4. Carrilho MR, Geraldini S, Tay F, et al. In vivo preservation of the hybrid layer by chlorhexidine. *J Dent Res.* 2007;86(6):529-533. 5. Hebling J, Pashley DH, Tjaderhane L, Tay FR. Chlorhexidine arrests subclinical degradation of dentin hybrid layers in vivo. *J Dent Res.* 2005;84(8):741-746.




HIGHEST BOND STRENGTH

Weak bonds can lead to failed restorations, sensitivity, unhappy patients, and damaged reputations for clinicians. Peak Universal Bond adhesive's unique Dymetech monomer blend has been shown to help create long-lasting, strong bonds that measure at or near the original strength of dentin.¹



DYMETECH™ MONOMER BLEND

Dymetech™ phosphate monomer blend is a proprietary blend of three phosphate monomers and twelve cross-linking methacrylate groups that provides time-tested strength and versatility across virtually all substrates.¹ It is the backbone of the strongest bonding agent we have ever tested.¹



The 10-MDP monomer (Image 1) is a long-chained monomer that consists of a methacrylate group on one end and a phosphate group on the other. This can limit the number of bond interactions. The Dymetech monomer blend (Image 2) has more functional groups (phosphates and monomers) than 10-MDP which creates a greater amount of bond interactions.

MULTIPLE DELIVERY OPTIONS

Prepare your tray with any of these methods for strong bonding.

Syringe: easily switch tips for multiple uses, syringe also limits the air that comes into contact with the adhesive, reducing solvent evaporation

Bottle: economic, and drop dispensing eliminates waste

Unit dose: convenience of single-use with just the right amount, reduces risk of cross-contamination

CHLORHEXIDINE FACILITATES BETTER BONDS

Chlorhexidine is present in both Consepsis™ antibacterial solution and Peak Universal Bond adhesive due to its important benefits. In vivo studies have shown that restorations not treated with chlorhexidine (CHX) exhibited a significant decrease in the structural integrity of the collagen network and in bond strength (38% bond strength degradation vs. no degradation in CHX-treated teeth).^{4,5} Add Consepsis solution into your bonding routine with either a self-etch or total-etch method.

NO RINSE	ETCH AND RINSE	WHEN DO YOU USE CONSEPSIS ANTIBACTERIAL SOLUTION?	
		SELF-ETCH - Etch and no rinse Peak SE + Peak Universal Bond	TOTAL-ETCH - Etch and rinse Ultra-Etch + Peak Universal Bond
		1. Peak™ SE primer	1. Ultra-Etch™ etchant
		Consepsis™ solution - RECOMMENDED	Consepsis™ solution - RECOMMENDED
		2. Peak™ Universal Bond adhesive	2. Peak™ Universal Bond adhesive