



Our NaviTip™ tips are treated to create flexible tips with a smooth, round end—ensuring they will easily snake through canals to deliver irrigants precisely where they are needed. This treatment can create a golden hue you may see on some of our NaviTip tips.

### NaviTip™ Tips 29 ga with Single Sideport

- Designed to direct the flow of chemistry through the sideport of the tip before flowing down into the area of the apex, thus reducing the risk of product extrusion
- Flexible, stainless steel cannula easily navigates curved canals
- Bendable tip 0.33 mm

29 ga delivers paste materials: MTAFlow™, EndoREZ™, and UltraCal™ XS.

Note: UltraCal™ XS calcium hydroxide paste should only be used with NaviTip 29 ga Single Sideport tips.

	Tip length	20pk
29 ga - 0.0129" NaviTip	27 mm	4989
29 ga - 0.0129" NaviTip	25 mm	4990
29 ga - 0.0129" NaviTip	21 mm	4991
29 ga - 0.0129" NaviTip	17 mm	4992
29 ga-29 ga NaviTips	27-17 mm	5143

### NaviTip™ Tips 31 ga with Double Sideport Irrigator

- Double sideports deliver irrigants safely, minimizing the possibility of chemicals being expressed past the apex
- The world's smallest cannula navigates the most intricate canal spaces

Designed for: ChlorCid™, ChlorCid™ V, ChlorCid™ Surf, Ultradent™ EDTA 18% Solution, Consepis™, and Citric Acid 20%.

5121 - 21 mm 20pk	5123 - 27 mm 20pk
5122 - 21 mm 50pk	5124 - 27 mm 50pk

### NaviTip™ Tips 29 ga and 30 ga

- Provide controlled delivery to the apex
- Flexible, stainless steel cannula easily navigates curved canals
- Available in 17, 21, 25, and 27 mm lengths and 29 or 30 ga cannulae


29 ga delivers paste materials.

30 ga delivers solutions and gels.

	29 ga - 20pk	30 ga - 20pk	29 ga - 50pk	30 ga - 50pk
17 mm - White	5112	1249	1378	1421
21 mm - Yellow	5113	1349	1374	1422
25 mm - Blue	5114	1250	1376	1423
27 mm - Green	5115	1354	1377	1424
Assorted 17-27 mm	5116	1351	1379	3319

LOK-TITE™

new



LOK-TITE™

### NaviTip™ FX™ Tips 30 ga


- One-of-a-kind brush cleans, scrubs, and irrigates simultaneously
- Rigid 30 ga cannula

Designed for: Ultradent™ Citric Acid 20% Solution.

Listed as an "EXCELLENT" product by a prominent independent research institute.<sup>1</sup>

1452 - 17 mm 20pk
1454 - 25 mm 20pk

1. Clinical Research Associates Newsletter. Volume 29, Issue 1, January 2005.



### Endo-Eze™ Irrigator Tip 27 ga

- Provides ideal reach reducing risk of expressing chemicals past the apex
- Comes with a flexible, blunt cannula with a unique, anti-obturator end
- Non-sterile

Designed for: Ultradent™ 5 ml syringe.

207 - 25 mm 20pk



### Endo-Eze™ Tips

- Work for a variety of endodontic procedures
- Come with flexible, strong cannulae in a variety of sizes
- Bend easily

Designed for: luting materials and air/water delivery.

Use with: TriAway™ Adapter, PermaFlo™ DC (20 ga), and other Ultradent syringes.

Gauge	Diameter	20pk	100pk
18 ga	0.049"; 1.25 mm	345	1428
19 ga	0.042"; 1.06 mm	346	1429
20 ga	0.035"; 0.90 mm	347	1430
22 ga	0.028"; 0.70 mm	348	1431



### Capillary Tips

- Evacuates canals and substantially minimizes use of paper points
- Narrow, flexible taper accesses curved canals
- Do not use to deliver irrigating materials or endodontic chemistries

Designed for: Ultradent™ Luer Vacuum Adapter for moisture removal from endodontic canals.

Diameter	20pk	50pk
0.014"; 0.36 mm	341	3099
0.019"; 0.48 mm	186	1425

LOK-TITE™



Ultradent Tips

# DESIGNED TO DELIVER



The chemistries you use are different. Some are chemically activated—needing to be mixed immediately before delivery. Others have varying viscosities. Some work in pits and fissures, some inside canals, and some on enamel surfaces. Each chemistry you use is designed for a specific purpose. Shouldn't the same be true for your tips?

Ultradent makes tips **designed to deliver** each chemistry we create. Whether you're delivering a flowable composite, a viscous gel, or thick impression material we make the perfect tip for the job. And since our tips are engineered on-site, we test each design to ensure it works perfectly with the chemistry it's intended for.

Check out our tips with **LOK-TITE™** and **COMFORT HUB™**

Luer Lock tips with Lok-Tite feature double threads that lock the tip into place for increased security and wings for easy attachment and removal.  
Comfort Hub tips include larger ergonomic wings that provide a secure, comfortable grasp of the tip.

Black Mini™ Tip

- Dispenses large volumes
- Opaque plastic preserves the flow of light-cured materials

Designed for: Ultra-Blend™ plus, Ultradent™ LC Block-Out Resin, PermaFlo™, TriAway Adapter, UltraTemp™, Opalescence™ Boost™, Ultradent™ Diamond Polish Mint, OpalDam™, Opalescence™ Endo, and OraSeal™ Caulking.

514 - 100pk
1433 - 500pk

LOK-TITE™



Black Mini™ Brush Tip

- Precise, controlled delivery of aqueous materials
- Tight, adjustable brush fibers minimize bubbles
- Unique to Ultradent

Designed for: Consepsis™, Peak™ SE Primer, Peak™ -ZM, Seek™/Sable™ Seek™, Ultradent™ Silane, and Ultradent™ Universal Dentin Sealant.

1169 - 100pk
1432 - 500pk

LOK-TITE™



Black Micro™ FX™ Tip

- Accommodates various viscosities
- Flocked tip fans out to spread materials in a thin, uniform layer
- 22 ga tip

Designed for: PrimaDry™ and PermaSeal™.

1357 - 100pk
1434 - 500pk

LOK-TITE™



Black Micro™ Tip · Blue Micro™ Tip

- Pinpoint precision
- Narrow cannulae accurately deliver materials
- Black Micro tip: 22 ga, Blue Micro tip: 25 ga

Black Micro tip designed for: Ultra-Blend™ plus and EtchArrest™.  
Blue Micro tip designed for: Ultra-Etch™, PrimaDry™, and Opal-Etch™.

1085 - Black 100pk	127 - Blue 100pk
1435 - Black 500pk	1436 - Blue 500pk



Intraoral Impression Tip

- Allows precise placement
- Attaches to Impression Mixing tip

Designed for: Thermo Clone™ VPS and Chromaclone™ PVS.

2906 - 50pk
-------------

LOK-TITE™



Impression Mixing Tips

- Enables direct delivery of impression materials
- Automixing, disposable, and color coded

Designed for: Thermo Clone VPS and Chromaclone PVS.

2902 - Yellow 50pk
2903 - Pink 50pk
2904 - Green 50pk



SoftEZ™ Tip

- Tip fibers provide visible, controlled delivery
- Brush fibers facilitate smooth application

Designed for: Enamelast™.

4712 - 50pk
4711 - 100pk

LOK-TITE™



White Mac™ Tip · White Mini™ Tip

- Dispense large volumes
- All-plastic delivery tips

Designed for: OraSeal™ Caulking, OraSeal™ Putty, Consepsis™ Scrub, and Opalustre™.

1361 - White Mac 100pk
1480 - White Mac 500pk
1247 - White Mini 20pk



Intraoral Tip

- Allows precise placement
- Attaches to dual-barrel mixing tips

Designed for: PermaFlo™ DC.

5922 - 20pk
-------------



LOK-TITE™

Micro 20 ga Tip

- Large-gauge cannula enables consistent flow
- Standard flowable composite delivery tip

Designed for: MTAFlo™, PermaFlo, PermaFlo™ Purple, OpalDam, and OpalDam™ Green.

1252 - 100pk
1437 - 500pk



COMFORT HUB™ LOK-TITE™

Metal Dento-Infusor™ Tip with Comfort Hub™

- Comfort Hub™ design for secure, comfortable control while attaching or removing
- Places hemostatic agents precisely and effectively removes superficial coagulum
- Blunt, bent 19 ga cannula with padded brush enables gentle pressure on the sulcus
- Ultradent's first tip, the "MDI" remains paramount for successful tissue management

Designed for: Astringedent™, Astringedent™ X, ViscoStat™, ViscoStat™ Clear, and Peak™ Universal Bond.

4955 - 100pk
4956 - 500pk




Blue Mini™ Dento-Infusor™ Tip

- Offers the same tissue management benefits as the Metal Dento-Infusor tip
- Allows controlled flow of drop-sized quantities
- Outer diameter 1.2 mm

Designed for: Hemostatics and Consepsis.

1086 - 100pk
1440 - 500pk

LOK-TITE™



Metal Opal™ Seal Tip

- Soft fibers are ideal for applying a consistent, even, thin layer to etched enamel
- Metal cannula provides ideal flow and ease of expression

Designed for: Opal™ Seal.

500-094 - 20pk
500-095 - 100pk
500-097 - 500pk

LOK-TITE™



Ultradent™ Mixing Tip

- Mixes and delivers in one action

Designed for: UltraTemp™, PermaFlo™ DC, and EndoREZ™.

5920 - 20pk
-------------



SST™ (Surgical Suction Tip)

- Ideal for delicate surgeries
- Large-diameter tip opening

Designed for: Ultradent™ Luer Vacuum Adapter for small periodontic or endodontic procedures and controlled suction of Opalescence Boost.

1248 - 20pk
-------------

LOK-TITE™



Inspiralm™ Brush Tip

- Delivers viscous or filled materials smoothly via an internal helical channel and ridge
- Tight, adjustable brush fibers minimize bubbles

Designed for: Composite Wetting Resin, Peak Universal Bond, Ultradent™ Porcelain Etch, Ultra-Etch, UltraSeal™ XT plus, and UltraSeal™ XT hydro.

123 - 100pk
1033-I - 500pk

LOK-TITE™



Micro Capillary™ Tips

- Bright colors are easily identified against soft tissues
- The world's smallest molded tips (inside diameter: 0.008")

Designed for: Perio materials, endodontics, and Ultradent™ Luer Vacuum Adapter.

1120 - 5 mm 20pk	1121 - 10 mm 20pk
1802 - 5 mm 100pk	-

LOK-TITE™



Opal™ Seal Tip

- Soft fibers are ideal for applying a consistent, even, thin layer to etched enamel
- Specially designed to deliver filled resins with internal helical channel and ridge

Designed for: Opal™ Seal.

500-064 - 20pk
500-034 - 100pk

LOK-TITE™

