Pernia™

Fluoride Composite

Product Description

Pernia™ composite is a light-cure, radiopaque, methacrylate-based flowable composite. Pernia® composite contains ~1µm average particle size with narrow upper limit particle distribution. Pernia® composite is available in multiple shades matching the VITA shade guide. It provides superior polishability (particularly for a hybrid) and affords excellent translucency for matching, blending, and depth of cure. It is isotropic; 68% filled mass to mass, and has gum modulus of 1.2 MPa.

Pernia composite is intended to be used by dental professionals on patients of all ages and conditions requiring replacement for hybrid and tooth restorations. Any restrictions regarding patient populations can be found in the Contraindications section.

Pernia composite is a bonding agent and understand all instructions and SOS information prior to use.

Delivery Form(s)

- 1.2 ml delivery syringe
- Delivery tip:
 - Micro 20 ga tip

Indications for Use/Intended Purpose

Pernia composite is approved for anterior and posterior restorations. It is highly recommended for Class III and V restorations, or for the initial lower layers when using indirect building of composites.

Contraindications

- Do not use on patients with known sensitivity to methacrylates or other resins. If dermatitis, eczema, or other allergic reactions are present, expose area thoroughly and have patient consult their physician.
- If a patient has with allergy concerns, refer to product allergen document available at: www.ultradent.com/resources/product-allergens.
- If allergic reaction is observed, discontinue use thoroughly with water and have the patient consult their physician.

Warnings and Precautions

- Resins can be sensitizing. Do not use on patients with a known sensitivity to acrylates or other resins. If allergic reaction, dermatitis, or rash develops, have patient consult their physician.
- To avoid cross-contamination, use clean gloves and a clean instrument
- Avoid curing light to prevent premature polymerization

- To optimize bond to each succeeding layer, do not remove oxygen inhibition layer between incremental layers of composite
- Ensure composite is totally polymerized between layers

Stepwise Instructions

Preparation

1. Remove lube lock cap and securely attach a recommended delivery device. The delivery device can be used when intimate adaptation is required in narrow angles (e.g., Class I, V, and proximal angles Class III).
2. Verify flow and consistency prior to applying intraorally. If resistance is met, replace tip and try again.
3. Refer to bonding agent preparation according to manufacturer's instructions.

Use

1. Apply first layer of Pernia® composite 0.5-1.5 mm (layer 0.5 mm or less when using dentin opaque) to ensure complete adaptation of composite and light cure. If not following with a packable composite, continue with 2 mm increments (0.5 mm or less when using dentin opaque), curing after each increment. For curing lights with an intensity of ≥ 800 mW/cm², the wavelength range of 385-515 nm, cure each layer for 10 seconds. A final cure of 20 seconds is recommended.
2. Shape, finish, and adjust occlusion using preferred method.

Processing

To avoid cross-contamination, wipe syringe with an alcohol wipe. Disinfect device with a disposable barrier sleeve. If a barrier sleeve is not used, syringe should be disposed as single-use. Delivery tips may be reused.

Storage and Storage

Store according to outer package label
 Dispose of waste according to local rules, guidelines, and regulations.

Additional Information

Report any serious incident to the manufacturer and the competent authority.
www.ultradent.com / 800.552.5512 / 801.572.4400
 The Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/en/medicines/screen/search-device>

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