User Manual

Handpiece

J12109

ultrawave[™]
PIEZO ULTRASONIC 

U2 **ULTRADENT**
PRODUCTS, INC.

Contents

1 Documentation

1.1 Associated documentation

1.2 Electronic documentation

2 Required information

2.1 Indications for use

2.2 Operating principle

2.3 Date of inclusion of EC marking

2.4 Latest document update

2.5 Repairing or modifying the device

2.6 Accessory usage conditions

3 Removal from packaging, installation, connections

3.1 Unpacking the device

4 Description of the medical device

4.1 Attach a tip or a file

4.2 Connecting and disconnecting accessories during use

4.3 Connecting the handpiece

4.4 Handpiece cord

5 Cleaning, disinfecting and sterilizing

6 Monitoring and maintenance of the medical device

7 Maintenance

7.1 Identifying incorrect operation

7.1.1 No spray

7.1.2 Ultrasounds not working

7.1.3 Water leakage

7.1.4 Replacing the handpiece seal

8 Technical specifications for the medical device

8.1 Identification

8.2 Handpiece

8.3 Environmental characteristics

8.4 Environmental restrictions

8.5 Main performance characteristics

9 Regulations and standards

9.1 Official Texts

9.2 Medical class of the device

9.3 Standardized Symbols

9.4 Manufacturer identification

9.5 Distributor

9.6 Disposal and recycling

9.7 Disposing of accessories

10 Index

Foreword

The Ultradent Products Inc. medical device that you are about to install and use in your practice is a medical device designed for professional use. It comprises the chosen tool with which you will provide treatment within the context of your work.

To ensure optimum safety for yourself and your patients, comfort in your daily practice and to benefit fully from the technology of your medical device, please read the documentation provided carefully.

If you have received this medical device by mistake, please contact the supplier to arrange for it to be collected.

Please refer to the accessory cleaning, disinfection and sterilization protocols and the handpiece cleaning, disinfection and sterilization protocols for information about the following:

- preparation of parts for sterilization;
- detailed manual and automatic protocols;
- information concerning the sterilization;
- recommendations for the inspection of parts.

1 Documentation

This document contains the following information:

- indications for use;
- description of the medical device;
- installation of the medical device;
- use of the medical device;
- preparation for cleaning and disinfection of the medical device;
- monitoring and general maintenance of the medical device;
- maintenance to be performed by the user.

1.1 Associated documentation

This document must be used in association with the following documents:

Document title	References
Cleaning, disinfection and sterilization protocols for Wrenches	J81009

Document title	References
Cleaning, disinfection and sterilization protocols for Tips	J02009
Cleaning, disinfection and sterilization protocols for Ultradent Products Inc. Handpieces	J12949
General instructions relating to the complete range of Ultradent Products Inc. dental ultrasonic generators	J00039
Ultrawave™ Handpiece User Manual	J12109

1.2 Electronic documentation

The user instructions for your device are in part provided in electronic format and in part in printed format. However, you can request a free printed copy of the user instructions via the website, by telephone or in writing.

The electronic user instructions are available in PDF format (Portable Document Format) and you will need to have a PDF file reader software installed to read the instructions.

The device user instructions can be consulted at the following addresses: www.ultradent.com and www.satelec.com.

It is important for you to have read and understood the content of the user instructions relating to the use of your device and its accessories prior to use.

We recommend that you visit the website regularly to consult and/or to download the latest version of your device's user instructions.

J12109 • V1 • (14) • 03/2014 • NB22US010B

Electronic user
informations

2 Required information

2.1 Indications for use

The treatments performed with this medical device are those described in the User Manual for your ultrasonic generator, in the chapter *Associated documentation page 4*.

2.2 Operating principle

An electrical signal emitted by the medical device is supplied to the dental ultrasonic handpiece. This is connected to the medical device via a cord. The handpiece comprises a piezoelectric ceramic transducer, which transforms the electrical signal into ultrasonic vibrations. Mechanical vibrations are transmitted to a tip or a dental file attached to the end of the ultrasonic handpiece.

2.3 Date of inclusion of EC marking

20(14)

2.4 Latest document update

03/2014

2.5 Repairing or modifying the device

Contact the supplier of your device. Using the services of an unapproved repairer could render your device dangerous for you and your patients. Do not repair or modify the medical device without seeking the prior permission of Ultradent Products Inc..
If the device is modified or repaired, specific checks and tests must be carried out to ensure that the device is still safe to use.

In the event of doubt, contact an approved dealer or the customer service team:
customer.support@ultradent.com
www.ultradent.com

Ultradent Products Inc. at the request of technical personnel working for the network of dealers approved by Ultradent Products Inc., provides all information required to repair the faulty parts on which they may perform repairs.

2.6 Accessory usage conditions

Accessories and handpiece must be cleaned, disinfected and sterilized prior to use.

3 Removal from packaging, installation, connections

3.1 Unpacking the device

When you receive your device, check for any damage that may have occurred during transportation. If you have received this device by mistake, please contact the supplier to arrange for it to be collected. If you have any questions or requirements, contact your supplier.

4 Description of the medical device

4.1 Attach a tip or a file

A tip or a file vibrates correctly when it is perfectly tightened without being forced beyond its stop point. Tighten it moderately using the wrench provided to ensure optimum ultrasound operation. Over-tightening of the tip or file can result in breakage of the tip, file or handpiece.

To prevent self-locking of the tip or the file, the latter must be removed after each use.

4.2 Connecting and disconnecting accessories during use

- Do not connect/disconnect the cord(s) or the handpiece when the medical device is switched on and your foot is on the pedal.
- Do not tighten or loosen the tips when the handpiece is activated.

4.3 Connecting the handpiece

Check for the absence of signs of moisture at the handpiece connections, and eliminate them if necessary (wipe and blow using a multipurpose syringe).

- Lubricate the irrigation system seal located behind the handpiece with dental instrument lubricant to extend its effectiveness and prevent leaks.

Connect the handpiece to the sleeve, by avoiding rotation movement.
Install the handpiece on the support.

4.4 Handpiece cord

The cord is only compatible with Ultradent Products Inc. handpieces .

The Ultrawave™ cord ensures irrigation circulation and electrical connection between the medical device and the handpiece.

- disconnect the handpiece and refer to the handpiece cleaning, disinfection and sterilization protocols by Ultradent Products Inc. J12949;

5 Cleaning, disinfecting and sterilizing

The instructions relating to the accessory cleaning, disinfection and sterilization protocols provided by Ultradent Products Inc. have been approved for each medical device and accessory. The applicable guides are listed in chapter *Associated documentation page 4*

They can be downloaded at the following addresses www.satelec.com/documents or www.ultradent.com

In all cases, the local regulations in force relating to the accessory cleaning, disinfection and sterilization protocols take precedence over the information provided by Ultradent Products Inc.

6 Monitoring and maintenance of the medical device

Before and after each use, check that the device and its accessories are not faulty in any way. This is necessary to detect any insulation fault or damage. If necessary, replace damaged parts.

7 Maintenance

7.1 Identifying incorrect operation

In the event of incorrect operation, refer to the tables below to quickly identify and repair the non-complex parts of the medical device. If the incorrect operation is not described in the tables below, please contact your supplier or the After-Sales team at Ultradent Products Inc.. Do not use the medical device if it appears to be damaged or faulty. Isolate the medical device and make sure that it cannot be used.

7.1.1 No spray

Symptoms: There is no water spray at the tip

Possible causes	Solutions
Flow configuration button on minimum	Adjust the flow configuration button
Tip or file blocked	Unblock the tip or file

Possible causes	Solutions
Incorrect choice of tip	Check the tip
Inadequate amount of spray	Adjust the spray

Do not ever try to unblock a tip or file using a probe.

7.1.2 Ultrasounds not working

Symptoms: the tip does not vibrate, vibration cannot be heard.

Possible causes	Solutions
Tip loose	Tighten the tip using the wrench
Faulty connector contact	Clean the cord contacts
Handpiece cord wire(s) cut	Return to Ultradent Products Inc. After-Sales Department to replace the cord

7.1.3 Water leakage

Symptoms: Water is leaking from one of the following places:


- between the base of the handpiece and its cord

Possible causes	Solutions
Wear of 1.15 mm x1 mm O'ring handpiece seal	Replace the seal using F12304 kit

7.1.4 Replacing the handpiece seal

Remove the faulty seal from the water connector using a thin pair of pliers.

Install the new seal using kit F12304 as shown below.

8 Technical specifications for the medical device

8.1 Identification

Manufacturer	SATELEC®
Name of the medical device	Ultrawave™ handpiece

8.2 Handpiece

Length (in mm)	111,7
Maximum diameter (in mm)	20,5
Weight (in g)	53

8.3 Environmental characteristics

Operating temperature	+10°C to +30°C
-----------------------	----------------

Storage temperature	-0 °C to +60 °C
Operating humidity	30% to 75%
Storage humidity	10% to 70%, including condensation
Atmospheric pressure	Between 800 hPa and 1060 hPa
Altitude	Less than or equal to 2000 meters

8.4 Environmental restrictions

Usage premises	Can be used at all medical premises. The medical device must not be used in an operating theater, or outside.
Use in gas-filled atmosphere	The medical device is not designed for use in a type AP or APG gas-filled atmosphere or in the presence of anesthetic gases.
Immersion	The handpiece must not be immersed.

8.5 Main performance characteristics

Ultrasonic vibrations of the tip or file fitted to the end of the conventional dental ultrasonic handpiece.

- Tip amplitude $\leq 200 \mu\text{m}$.

9 Regulations and standards

9.1 Official Texts

This medical device complies with the essential requirements of European Directive 93/42/EEC. This equipment is designed and developed in compliance with Electrical Safety standard IEC60601-1 in force. It was designed and manufactured in accordance with an EN ISO 13485-certified quality assurance system.

9.2 Medical class of the device

This medical device is a class IIa device according to European Directive 93/42/EEC.

9.3 Standardized Symbols

Symbols	Meaning
	<p>Refer to the accompanying documentation</p>
	<p>Consult the User Manual</p>
 <p>Electronic user informations</p>	<p>Accompanying documentation in electronic format</p>

Symbols	Meaning
	Sterilization at 132 °C in an autoclave
	Washer disinfectant for thermal disinfection
	EC marking

Symbols	Meaning
	Do not dispose of as household waste
	Year of manufacture
Rx Only	Caution: US Federal law restricts this device to sale by or on the order of a Physician

9.4 Manufacturer identification

SATELEC

A Company of ACTEON Group

17, avenue Gustave Eiffel

BP 30216

33708 MERIGNAC cedex

FRANCE

Tel. +33 (0) 556.34.06.07

Fax. +33 (0) 556.34.92.92

E.mail: satelec@acteongroup.com.

www.acteongroup.com

9.5 Distributor

Ultradent Products, Inc.

505 West 10200 South

South Jordan, UT 84095 USA

E-mail : customer.support@ultradent.com

www.ultradent.com

9.6 Disposal and recycling

As an item of Electrical and Electronic Equipment, the medical device must be disposed of via a specialist collection, removal, recycling or destruction channel. This applies in particular to the European market, in reference to Directive 2002/96/EC dated 27/01/2003.

When your medical device has reached the end of its service life, contact Ultradent Products Inc. to find out how to proceed. Relevant contact details are given in the chapter *Distributor page 18*

9.7 Disposing of accessories

An accessory that has reached the end of its service life must be disposed of in infectious clinical waste

10 Index

A

Altitude 13
Amplitude 14
approved dealers 6-7

C

cleaned 7

D

damage 9
dental file 6
dental ultrasonic generators 5
disinfected 7
disposal 18

E

Electrical Safety 14
electronic user instructions 5
European Directive 14

F

fault 9
first inclusion of EC marking 6

G

gas-filled atmosphere 13

H

Handpieces 5, 8

I

incorrect operation 9

K

kit F12304 11

L

latest update 6

M

Manufacturer 12
Medical class 14
moisture 8

10 Index

P

Pressure 13

R

Regulatory 5, 18
repair 6-7

S

seal 11
spray 9
sterilized 7

T

Temperature 12-13
tip 5, 6, 9-10

U

ultrasonic vibrations 6
User Manual 5

W

water leak 11
Wrenches 4

Ref: J12109 • V1 • (14) • 03/2014 • NB22US010B

Ultradent Products, Inc. • 505 West 10200 South • South Jordan • UT 84095 • USA
www.ultradent.com • customer.support@ultradent.com

 **ULTRADENT**
PRODUCTS, INC.

J12949

Manual Pre-disinfection/Cleaning and Sterilization Instructions for Ultradent Scaler Handpieces

Cautions

Do not use steel wool or abrasive cleaners.

Avoid solutions containing iodine or with a high chlorine content.

The pH of detergents/disinfectants must be between 7 and 11.

The cleaning method for scaler Handpieces recommended by Ultradent is manual or automated.

Scaler Handpieces and their accessories (nosepiece, LED ring, optical guide, cord pass-through depending on model) must be cleaned with extreme care.

Do not immerse or place scaler handpieces in an ultrasonic cleaner.

Only the Nosepiece disassembled from Ultrawave™ or Ultrawave™ XS LED Handpieces can be immersed and treated using the ultrasonic procedure.

If the surface of the devices is damaged, ultrasonic cleaning may cause further damage.

It is the responsibility of the end user to ensure that all equipment used to treat Ultradent products is properly installed, validated, maintained and calibrated.

Whenever possible, a washer/disinfector should be used for Handpieces and their accessories. Prevent the overloading of wash baskets for ultrasonic cleaning or cleaning in a washer/disinfector.

Cleaning/sterilization cycle limits

Repeated treatment cycles including manual washing and sterilization have minimal effects on Ultrawave™ XS LED Handpiece scaler Handpieces, LED ring and Optical guide.

Repeated treatment cycles that include ultrasonic cleaning, manual washing and sterilization have minimal effects on scaler Handpiece Nosepieces.

End of service life is normally determined by wear and damage due to use.

Important

- Throughout the procedure, wipe blood and/or debris from the Handpiece to prevent it drying on the surfaces.
- In case of external irrigation tubing or configuration of the device with tank, after each use, rinse the irrigation tubing, Handpiece and tip for 1 minute with an enzymatic or alkaline cleaning solution to prevent debris and contamination from drying on the device. Rinse with distilled or purified water for 1 minute to remove the enzymatic or alkaline cleaner. Flush the irrigation tubing to remove liquid from the inside.
- In the case of irrigation coming from the dental practice water network, clean and rinse the cord, the Handpiece and the tip for 1 minute after each use.
- Soiled devices should be separated from non-contaminated devices to avoid contamination of personnel or surroundings.
- The Handpiece must be covered with a cloth dampened with purified water to prevent blood and/or debris from drying on the surface.

Containment and transportation

Soiled devices must be transported separately from non-contaminated devices to avoid contamination.

Preparation for pre-disinfection/cleaning

It is advisable to treat devices immediately or as soon as possible after use (Ultradent devices must be treated within two hours of use). After each use and prior to cleaning the device:

- Unscrew the tip at the front of the Handpiece-cord assembly and place it in the appropriate container. The cleaning of the tip is subject to a separate protocol.
- Remove the single-use or autoclavable irrigation tubing and autoclavable clips from the handpiece-cord assembly, if applicable. Their cleaning is subject to a separate protocol. Always discard the single-use perforator and/or single-use irrigation tubing.
- Disconnect the Handpiece from the cord.
 - With regards to Ultrawave™ and Ultrawave™ XS LED Handpieces: Unscrew the Nosepiece.
 - With regards to the Ultrawave™ XS LED handpiece: remove the optical guide
 - With regards to the Ultrawave™ XS LED handpiece: disconnect the LED ring from the Handpiece, pulling it gently.

Pre-disinfection/Cleaning - Manual method

Equipment: soft-bristled brush, soft and lint-free swab, lint-free cloth, syringe, pipette and/or water jet, enzymatic or alkaline cleaning solution.

Step duration (minimum)	Cleaning instructions
1 minute	Rinse soiled device under cold running water for at least 1 minute. Use a soft-bristled brush and a lint-free swab to remove most of the contamination and debris. Use a syringe, a pipette or a water jet with an enzymatic or alkaline cleaning solution to rinse the Handpiece cannulation.

Step duration (minimum)	Cleaning instructions
15 minutes	<p>Only the Nosepiece disassembled from Ultrawave™ or Ultrawave™ XS LED handpieces can be treated using the ultrasonic procedure.</p> <p>Wash the Nosepiece using an ultrasound process for 15 minutes in an enzymatic or alkaline cleaning solution.</p>
4 minutes	<p>Clean and disinfect each element separately: the Handpiece body, in particular the screw threading, the nosepiece for Ultrawave™ handpieces, the LED ring and optical guide for the Ultrawave™ XS LED handpiece, with an enzymatic or alkaline cleaning solution for a minimum of 4 minutes.</p>
1 minute	<p>Rinse each component separately with deionized or purified water. Use a syringe, a pipette or a water jet to flush the handpiece cannulation with deionized or purified water.</p>
	<p>Visually inspect device. Repeat the cleaning procedure until no visible contamination remains on the device.</p>
	<p>Dry components using a soft, lint-free cloth or clean compressed air.</p>

Inspection

- Handpieces and Nosepieces must be inspected to check that no contamination remains, that they are not corroded, dulled, discoloured or damaged.
- Before packaging and sterilizing the cleaned products, check they are clean, undamaged and function properly. Repeat cleaning steps until no visible contamination remains on the device.
- Damaged devices must be discarded. Handpieces must not be lubricated.

Packaging

Use suitable packaging for sterilization, such as a Sterile Barrier System in compliance with ISO 11607. Care should be taken to protect Handpieces from contact with other objects that may damage their surface or the Sterile Barrier System.

Sterilization

Unless otherwise specified, non-sterile products can be re-sterilized using validated steam sterilization methods (ISO 17665 or national standards).

Ultradent recommendations for packed Handpieces and their accessories are as follows:

Cycle type	Sterilization exposure time	Sterilization exposure time	Sterilization exposure time
Saturated steam - forced air removal (pre-vacuum)	4 minutes	132 °C	Minimum 20 minutes

Dry times generally range from 20 to 60 minutes due to the difference in packaging materials, steam quality, device materials, total mass, sterilizer performance, and varying cool-down time.

The distributor and manufacturer accept no responsibility for sterilization procedures performed by the customer that are not performed according to these Ultradent recommendations.

Storage

Storage conditions for products labelled “STERILE” are printed on the packaging label. Packaged products should be stored in a dry, clean environment, protected from direct sunlight, pests, and extremes of temperature and humidity. Use products in the order in which they are received (“first in, first out” principle), taking into account the expiry date indicated on the label.

Automated Pre-disinfection/Cleaning and Sterilization Instructions for Ultradent Scaler Handpieces

Cautions

Do not use steel wool or abrasive cleaners.

Avoid solutions containing iodine or with a high chlorine content.

The pH of detergents/disinfectants must be between 7 and 11.

The cleaning method for scaler handpieces recommended by Ultradent is manual or automated.

Scaler handpieces and their accessories (nosepiece, LED ring, optical guide, cord pass-through depending on model) must be cleaned with extreme care.

Do not immerse or place scaler handpieces in an ultrasonic cleaner.

Only the Nosepiece disassembled from Ultrawave™ or Ultrawave™ XS LED handpieces can be immersed and treated using the ultrasonic procedure.

If the surface of the devices is damaged, ultrasonic cleaning may cause further damage.

It is the responsibility of the end user to ensure that all equipment used to treat Ultradent products is properly installed, validated, maintained and calibrated.

Whenever possible, a washer/disinfector should be used for Handpieces and their accessories. Prevent the overloading of wash baskets for ultrasonic cleaning or cleaning in a washer/disinfector.

Cleaning/sterilization cycle limits

Repeated treatment cycles including manual washing and sterilization have minimal effects on Ultrawave™ LED handpiece scaler handpieces, LED ring and Optical guide.

Repeated treatment cycles that include ultrasonic cleaning, manual washing and sterilization have minimal effects on scaler Handpiece Nosepieces.

End of service life is normally determined by wear and damage due to use.

Important

- Throughout the procedure, wipe blood and/or debris from the Handpiece to prevent it drying on the surfaces.
- In case of external irrigation tubing or configuration of the device with tank, after each use, rinse the irrigation tubing, Handpiece and tip for 1 minute with an enzymatic or alkaline cleaning solution to prevent debris and contamination from drying on the device. Rinse with distilled or purified water for 1 minute to remove the enzymatic or alkaline cleaner. Flush the irrigation tubing to remove liquid from the inside.

- In the case of irrigation coming from the dental practice water network, clean and rinse the cord, the Handpiece and the tip for 1 minute after each use.
- Soiled devices should be separated from non-contaminated devices to avoid contamination of personnel or surroundings.
- The Handpiece must be covered with a cloth dampened with purified water to prevent blood and/or debris from drying on the surface.

Containment and transportation

Soiled devices must be transported separately from non-contaminated devices to avoid contamination.

Preparation for pre-disinfection/cleaning

It is advisable to treat devices immediately or as soon as possible after use. Ultradent devices must be treated within two hours of use.

After each use and prior to cleaning the device:

- Unscrew the tip at the front of the Handpiece-cord assembly and place it in the appropriate container. The cleaning of the tip is subject to a separate protocol.
- Remove the single-use or autoclavable irrigation tubing and autoclavable clips from the handpiece-cord assembly, if applicable. Their cleaning is subject to a separate protocol. Always discard the single-use perforator and/or single-use irrigation tubing.
- Disconnect the Handpiece from the cord.
 - With regards to Ultrawave™ and Ultrawave™ XS LED Handpieces: unscrew the Nosepiece.
 - With regards to the Ultrawave™ XS LED handpiece: remove the optical guide
 - With regards to the Ultrawave™ XS LED handpiece: disconnect the LED ring from the Handpiece, pulling it gently.

Pre-disinfection/Cleaning - Automated method

Pre-disinfection/pre-cleaning method

Note: The pre-disinfection/pre-cleaning method must be performed on the Handpiece and accessories prior to automated cleaning.

Equipment: soft-bristled brush, soft and lint-free swab, lint-free cloth, syringe, pipette and/or water jet, enzymatic or alkaline cleaning solution.

Step duration (minimum)	Cleaning instructions
1 minute	Rinse the soiled device under cold running tap water for at least 1 minute. Remove most of the contamination and debris using a soft brush, a lint-free swab or a lint-free cloth. Use a syringe, a pipette or a water jet with an enzymatic or alkaline cleaning solution to rinse the handpiece cannulation.
2 minutes	Only the Nosepiece disassembled from Ultrawave™ or Ultrawave™ XS LED handpieces can be treated using the ultrasonic procedure. Wash the Nosepiece using an ultrasound process for 2 minutes in an enzymatic or alkaline cleaning solution.
2 minutes	Clean and disinfect each element separately: the handpiece body, in particular the screw threading, the nosepiece for Ultrawave™ handpieces, the LED ring and optical guide for the Ultrawave™ XS LED handpiece, with an enzymatic or alkaline cleaning solution for a minimum of 2 minutes.

Step duration (minimum)	Cleaning instructions
1 minute	Rinse the components separately under cold running tap water for at least 1 minute. Use a syringe, a pipette or a water jet with an enzymatic or alkaline cleaning solution to rinse the handpiece cannulation.
	Visually inspect device. Repeat the cleaning procedure until no visible contamination remains on the device.

Pre-disinfection/Cleaning - Automated method

Step	Duration (minimum)	Cleaning/decontamination instructions
Pre-washing	2 minutes	Cold tap water
Washing	10 minutes	Warm tap water (> 40°C); use alkaline or enzymatic cleaning solution
Neutralization	2 minutes	Warm tap water with neutralizer, if necessary

Step	Duration (minimum)	Cleaning/decontamination instructions
Rinsing	2 minutes	Rinse with warm distilled or purified water (>40° C)
Drying	40 minutes	90° C

Thermal disinfection

For automated cleaning: thermal disinfection at 90° C for a minimum of 5 minutes.

Inspection

- Handpieces and Nosepieces must be inspected to check that no contamination remains, that they are not corroded, dulled, discoloured or damaged.
- Before packaging and sterilizing the cleaned products, check they are clean, undamaged and function properly. Repeat cleaning steps until no visible contamination remains on the device.
- Damaged devices must be discarded. Handpieces must not be lubricated.

Packaging

Use suitable packaging for sterilization, such as a Sterile Barrier System in compliance with ISO 11607. Care should be taken to protect handpieces from contact with other objects that may damage their surface or the Sterile Barrier System.

Sterilization

Unless otherwise specified, non-sterile products can be re-sterilized using validated steam sterilization methods (ISO 17665 or national standards).

Ultradent recommendations for packed handpieces and their accessories are as follows:

Cycle type	Sterilization exposure time	Sterilization exposure temperature	Drying time
Saturated steam - forced air removal (pre-vacuum)	4 minutes	132° C	Minimum 20 minutes

Dry times generally range from 20 to 60 minutes due to the difference in packaging materials, steam quality, device materials, total mass, sterilizer performance, and varying cool-down time.

The distributor and manufacturer accept no responsibility for sterilization procedures performed by the customer that are not performed according to these Ultradent recommendations.

Storage

Storage conditions for products labelled "STERILE" are printed on the packaging label. Packaged products should be stored in a dry, clean environment, protected from direct sunlight, pests, and extremes of temperature and humidity. Use products in the order in which they are received ("first in, first out" principle), taking into account the expiry date indicated on the label.

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