
Panasil®

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Product Description

Panasil® impression materials are addition-curing, elastomeric materials with hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation. Panasil® initial contact displays a high initial hydrophilicity.

The Panasil® family consists of five different viscosities (putty, heavy bodied, medium bodied, light bodied, x-light bodied), available in an assortment of delivery systems: standard 1:1 (50 ml automix cartridges) and 5:1 (380 ml cartridges) for use in most commercially available other automatic dispensing and mixing machines (e.g., Sympress), and traditional 1:1 putty jars.

Indications / Techniques

Panasil® Putty (Putty, Putty Fast Set, Putty Soft) and Panasil® binetics Putty (Putty Fast, Putty Soft) are to be used as preliminary materials for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique

- Two-step putty-wash impression technique using a foil (plastic putty spacer)
- One-step putty impression technique for forming functional peripheries

Panasil® tray Fast Heavy, Panasil® tray Soft Heavy and Panasil® tray Soft Heavy Fast are to be used as heavy bodied materials for:

- One-step impression technique (simultaneous technique) using single or dual viscosities
- Two-step impression technique using dual viscosities (Panasil® tray Fast Heavy)
- Functional impressions (one-step: Panasil® tray Soft Heavy)

Panasil® monophasic Medium is to be used as a medium bodied tray or syringeable impression material for:

- Taking impressions over fixed/removable restorations and implants (i.e., transferring impression posts and bridge components)
- Functional impressions
- Fabricating crown and bridgework or inlays

- Fabricating full or partial dentures
- Reline impressions
- Use in the simultaneous mixing technique as well as the putty-wash and Multi tray techniques
- Transferring root posts when fabricating posts and cores indirectly

Panasil® initial contact, Panasil® contact plus X-Light and Panasil® contact two in one Light are to be used as syringeable impression materials for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique
- Two-step impression technique using a spacer
- One-step impression technique (simultaneous technique) using dual viscosities
- Reline impressions
- Fabricating full or partial dentures

Warnings

Do not use Panasil® impression materials as a temporary reliner.

Do not use Panasil® impression materials with condensation-curing silicones, polyether, Vinylsiloxanether® or polysulfide materials.

Panasil® Putty impression materials are not suitable for detailed impressions when used alone.

Cautions

Do not use after expiration date.

Do not leave any residual impression material in the sulcus or oral cavity.

Do not swallow impression material! If swallowed: consult a medical doctor if problems arise or persist.

Avoid contact with eyes. If accidental contact with the eyes occurs, rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

These products do not normally cause allergic reactions. However, for sensitive persons, an allergic reaction to the products cannot be ruled out.

Use of products containing active sulfur, aluminium chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide impression materials, etc.) in conjunction with this product will interfere with the setting reaction of the vinyl polysiloxane material. Use of the materials requires the area to be rinsed thoroughly to remove all residue, before an impression is taken. Do not use latex gloves.

Do not interchange the base paste and catalyst paste lids or scoops for hand-mixed jar putty.

To ensure an optimal connection between the two materials, after the trial positioning, the preliminary impression should also be rinsed with water and dried before the correction material is applied. Both steps of the putty-wash impression technique (preliminary impression and wash impression) should also be carried out in quick succession.

When taking impressions of areas with severe undercuts and wide interdental spaces, use standard block-out techniques.

When using a custom impression tray, ensure that sufficient space remains between the side of the tray and the teeth/jaw.

Do not allow the material to enter the sewer or water system to avoid environmental contamination.

Avoid contact with clothing, since the material cannot be removed by dry cleaning.

Note

To ensure optimal impressions, the temperature of the material should not deviate from 23 °C (74 °F) before applying. Otherwise, working and setting times will be affected.

Impression Tray: Preparation and Adhesives

In general, all conventional impression trays can be used if the necessary pressure built up can be generated. If it is not possible to achieve sufficient retention with the impression material, brush the impression tray with Panasil® adhesive.

The surface of the impression tray must be repeatedly and fully brushed so that the entire surface of the tray is visibly blue. Allow the adhesive to dry in accordance with the information in the direction for use and cover with sufficient impression material in the impression tray.

Instructions for Use: Jar Material for Hand-Mix Putty Preparation

Only use materials with the same lot number.

Do not interchange the base paste and catalyst paste lids or scoops for hand-mixed jar putty.

Use the scoops to dispense equal amounts of base paste and catalyst paste.

The different color scoops must only be used for the materials with corresponding colors.

Close the containers carefully after use, and ensure that the lids are not interchanged.

Contamination of base paste and catalyst paste in the container renders the material unusable.

Knead the base paste and catalyst paste for 45 seconds until the material is homogeneously colored.

Should gloves be worn, test them for compatibility prior to mixing with a sample of the material to be mixed. Certain types of gloves can prevent the polymerization (e.g., latex gloves)

or impair the adhesive bond (e.g., nitrile gloves). Using vinyl gloves is recommended.

Prior to taking the wash impression, the preliminary impression must be cleaned with water, dried and trimmed using standard methods.

Instructions for the Safe Handling of the Kettenbach Plug & Press® System (380 ml Cartridge System; 5:1)

Prior to first use please unscrew the blue plug counterclockwise and pull it out (figure 1, page 146). Check the outlet opening. Then screw on the dynamic mixer (figure 2, page 146). Ensure that the positioning aid is on the top.

Please use only the dispensing and mixing components compatible with the Kettenbach 380 ml cartridge system for dispensing and mixing impression material.

The Kettenbach 380 ml cartridge is heavy and could be damaged if dropped. Damaged cartridges must not be used under any circumstances.

When using the material in a Sympress II dispensing device, please activate the speed setting indicated on the product label. (figure 6, page 146) Failure to comply with this instruction may modify the curing characteristics of the impression materials.

If you are using a unit other than the Sympress mixing device, please follow the instructions for the automatic dispensing and mixing system you are using.

When changing to a different mixing device, always discard at least 3 cm of mixed material.

Using the Kettenbach Plug & Press® System (380 ml Cartridge System; 5:1) in Dispensing Device Sympress / Sympress II

By pressing the plunger return knob, allow the plunger to return to its original position automatically (while the cover is closed; figure 3, page 146). Next, open the cover (figure 4, page 146).

Insert the Kettenbach 380 ml cartridge with the preassembled mixer in the Sympress / Sympress II dispensing and mixing unit (figure 5, page 146). If you are working with a previously used cartridge, the filled mixer will already be located on the cartridge as a cap.

If you are using a Sympress II dispensing device, please activate the speed setting indicated on the product label (figure 6, page 146). A light signal will confirm your selection.

Close the device cover (figure 7, page 147).

Using the Cartridge for the First Time

Dispense the impression material (ensure that the dynamic mixer is wound into place fully) until the base material and catalyst are mixed homogeneously (see coloring of the mixed impression material); discard at least the first 3 cm when dispensing (figure 8, page 147).

At least 3 cm of the mixed material must be discarded each time the plunger is restarted or when using the cartridge in different mixing units.

Load the impression tray (or application syringe) with the required amount of material (figure 9, page 147) by pressing the appropriate feed key. Hold the impression at an angle when loading and press against it lightly. Leave the mixer in the material when dispensing to avoid air bubbles.

Release the feed key as soon as the required amount of material has been dispensed. The filled dynamic mixer should remain on the cartridge as a seal.

Remove the used mixer before the next application by turning it to the left (figure 10, page 147), visually check the outlet opening for blockage and attach a new dynamic mixer (figure 10, page 147).

Remove the completely empty cartridge from the dispensing and mixing unit and dispose of it according to regional regulations.

Further Use of a Cartridge Already in Use

Wind a new dynamic mixer onto the same mixing unit and begin the new application (material no longer has to be discarded). Proceed with the application as described above.

Instructions for Use: Cartridge Material for Use in Applyfix® 4 Dispensing Gun DS-50 1:1/2:1

Insert the cartridge into the Applyfix® 4 dispensing gun DS-50 1:1/2:1. Ensure that the notches on the cartridge base are pointing down. The clasp will not close if the cartridge is not properly inserted into the dispensing gun.

Remove the cartridge cap. The cap can be replaced after initial use.

Advance plunger into cartridge and dispense a small amount of impression material until equal amounts are extruded at the same rate.

Install a mixing tip on the cartridge and turn cap 90 degrees in a clockwise direction to lock in place.

Inject the required amount of material, e.g., directly into the tray or the preliminary impression, or into an application syringe. To apply the material around the prepared teeth use the application syringe or the cartridge with an intraoral tip on the mixer. Check to ensure materials are set prior to mouth removal.

Replace the plug after use. Before using the cartridge again, remove the cartridge cap or the filled mixing tip and check the cartridge outlet openings to ensure that no polymerized materials have sealed the openings. Clear any clogged openings and dispense a small amount of impression material until equal amounts are extruded at the same rate.

Install a new mixing tip and continue as described above.

Cleaning and Disinfection of the Applyfix® 4

Remove the cartridge.

Remove any uncured paste using a cloth soaked in alcohol.

The Applyfix® 4 dispensing gun can be used several times if properly disinfected. Kettenbach recommends exchanging the gun if there are visible traces of use. Disinfection can be carried out using a commercially available dipping disinfectant. Kettenbach recommends the use of glutaraldehyde-based solutions.

Storage Conditions

Store in a dry place at room temperature; do not expose to sunlight and avoid exposure to extreme temperature fluctuations.

Disinfection

Impressions can be disinfected using, for example, a 2 % acidic glutaraldehyde solution. Use disinfection solutions specific for impression materials. See manufacturer instructions.

Electroplating

Impressions may be silver-plated with an alkaline-plating solution only.

Model Preparation

Impressions with Panasil® may be poured immediately after disinfection of the impression up to several weeks with high quality dental stone (class III - V). Impressions can also be poured several times.

Trademarks

Applyfix® 4 and Panasil® are registered trademarks of Kettenbach GmbH & Co. KG.

Note

Selected Kettenbach materials are available in certain markets only.

For professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

Limitation of Liability

Except where prohibited by law, Kettenbach GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development.

However, users must comply with and consider all recommendations and information in connection with any use.

Updated: October 6, 2020

Product Overview

Product name	ISO 4823	Consistency (approx.) mm	Mixing ratio and total content volume per unit	Mixing technique	Mixing element
Panasil® Putty	Type 0, Putty	24	1:1, 900 ml in jars	Manual, 1 scoop (12.3 ml) per component (base paste and catalyst paste)	Hand mix
Panasil® Putty Fast Set	Type 0, Putty	24	1:1, 900 ml in jars	Manual, 1 scoop (12.3 ml) per component (base paste and catalyst paste)	Hand mix
Panasil® Putty Soft	Type 0, Putty	24	1:1, 900 ml in jars	Manual, 1 scoop (12.3 ml) per component (base paste and catalyst paste)	Hand mix
Panasil® binetics Putty Fast	Type 0, Putty	23	5:1, 380 ml cartridge	Sympress/Sympress II or other automatic dispensing and mixing unit	Dynamic mixer, blue
Panasil® binetics Putty Soft	Type 0, Putty	23	5:1, 380 ml cartridge	Sympress/Sympress II or other automatic dispensing and mixing unit	Dynamic mixer, blue
Panasil® tray Fast Heavy	Type 1, Heavy bodied	32	5:1, 380 ml cartridge	Sympress/Sympress II or other automatic dispensing and mixing unit	Dynamic mixer, blue
Panasil® tray Fast Heavy	Type 1, Heavy bodied	34	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm
Panasil® tray Soft Heavy	Type 1, Heavy bodied	32	5:1, 380 ml cartridge	Sympress/Sympress II or other automatic dispensing and mixing unit	Dynamic mixer, blue
Panasil® tray Soft Heavy	Type 1, Heavy bodied	34	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm
Panasil® tray Soft Heavy Fast	Type 1, Heavy bodied	34	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm
Panasil® monophase Medium	Type 2, Medium bodied	35	5:1, 380 ml cartridge	Sympress/Sympress II or other automatic dispensing and mixing unit	Dynamic mixer, blue
Panasil® monophase Medium	Type 2, Medium bodied	35	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm

Product Overview

Product name	ISO 4823	Consistency (approx.) mm	Mixing ratio and total content volume per unit	Mixing technique	Mixing element
Panasil® initial contact Regular	Type 2, Medium bodied	37	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm
Panasil® initial contact Regular Fast	Type 2, Medium bodied	37	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm
Panasil® initial contact Light	Type 3, Light bodied	41	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, yellow, Ø 4.2 mm
Panasil® initial contact Light Fast	Type 3, Light bodied	41	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, yellow, Ø 4.2 mm
Panasil® initial contact X-Light	Type 3, Light bodied	44	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, yellow, Ø 4.2 mm
Panasil® initial contact X-Light Fast	Type 3, Light bodied	44	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, yellow, Ø 4.2 mm
Panasil® contact two in one Light	Type 3, Light bodied	38	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, green, Ø 6.5 mm
Panasil® contact plus X-Light	Type 3, Light bodied	42	1:1, 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	Mixing tip, yellow, Ø 4.2 mm

Technical Data

Product name	Mixing ratio	Working time at 23 °C/ 74 °F ≤	Working time at 35 °C/ 95 °F ≤	Intraoral setting time at 35 °C/ 95 °F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (approx.) %	Elastic Recovery test (approx.) %	Strain in compression (approx.) %
Panasil® Putty	1:1	2 minutes	not applicable	2 minutes	4 minutes	A 66	– 0.20	99.0	2.7
Panasil® Putty Fast Set	1:1	1 minute 30 seconds	not applicable	2 minutes	3 minutes 30 seconds	A 66	– 0.20	99.0	2.7
Panasil® Putty Soft	1:1	2 minutes	not applicable	2 minutes	4 minutes	A 60	– 0.20	99.0	2.7
Panasil® binetics Putty Fast	5:1	1 minute 30 seconds	not applicable	2 minutes 30 seconds	4 minutes	A 63	– 0.20	99.5	2.5
Panasil® binetics Putty Soft	5:1	2 minutes	not applicable	3 minutes	5 minutes	A 56	– 0.20	99.5	3.5
Panasil® tray Fast Heavy	5:1	1 minute 20 seconds	not applicable	2 minutes	3 minutes 20 seconds	A 62	– 0.20	99.7	2.5
Panasil® tray Fast Heavy	1:1	1 minute	not applicable	2 minutes	3 minutes	A 62	– 0.20	99.5	3.0
Panasil® tray Soft Heavy	5:1	2 minutes	not applicable	2 minutes	4 minutes	A 55	– 0.20	99.7	3.0
Panasil® tray Soft Heavy	1:1	2 minutes	not applicable	2 minutes	4 minutes	A 55	– 0.20	99.5	3.0
Panasil® tray Soft Heavy Fast	1:1	1 minute	not applicable	2 minutes	3 minutes	A 55	– 0.20	99.5	3.0
Panasil® monophasic Medium	5:1	2 minutes	1 minute	2 minutes	4 minutes	A 60	– 0.20	99.7	3.5

Technical Data

Product name	Mixing ratio	Working time at 23 °C/ 74 °F ≤	Working time at 35 °C/ 95 °F ≤	Intraoral setting time at 35 °C/ 95 °F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (approx.) %	Elastic Recovery test (approx.) %	Strain in compression (approx.) %
Panasil® monophase Medium	1:1	2 minutes	1 minute	2 minutes	4 minutes	A 60	- 0.20	99.7	3.5
Panasil® initial contact Regular	1:1	1 minute 30 seconds	1 minute	2 minutes 30 seconds	4 minutes	A 46	- 0.20	99.7	3.0
Panasil® initial contact Regular Fast	1:1	1 minute	30 seconds	2 minutes	3 minutes	A 46	- 0.20	99.7	3.0
Panasil® initial contact Light	1:1	1 minute 30 seconds	1 minute	2 minutes 30 seconds	4 minutes	A 46	- 0.20	99.3	3.5
Panasil® initial contact Light Fast	1:1	1 minute	30 seconds	2 minutes	3 minutes	A 46	- 0.20	99.3	3.5
Panasil® initial contact X-Light	1:1	1 minute 30 seconds	1 minute	2 minutes 30 seconds	4 minutes	A 46	- 0.20	99.3	3.5
Panasil® initial contact X-Light Fast	1:1	1 minute	30 seconds	2 minutes	3 minutes	A 46	- 0.20	99.3	3.5
Panasil® contact two in one Light	1:1	2 minutes	1 minute	2 minutes	4 minutes	A 46	- 0.20	99.7	3.5
Panasil® contact plus X-Light	1:1	2 minutes	1 minute	2 minutes	4 minutes	A 46	- 0.20	99.7	5.0

*Total setting time (removal time from mouth) from start of mix.

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