Ultradent®  
Clone™ Bite

**Description/Indications:**

Ultradent® Clone™ Bite is an orange colored addition vulcanizing silicone for all bite registration requirements with a very high final hardness and is non-slumping. It is also a fast setting material perfect for preliminary impressions for crown and bridge procedures and temporaries.

**Technical Data:**

- Working time is 20 seconds.
- Material is completely set in 90 seconds from mix.
- Final hardness is 95 Shore A.

**Procedure:**

1. Load cartridge in dispensing gun.
2. Immediately prior to taking an impression, remove the cartridge cap and discard.
3. Gently squeeze dispenser to bleed enough material from the cartridge to ensure an even flow of the base and catalyst, (about 1/4 of an inch).
4. Attach the mixing tip with or without intra-oral tip. Intra-oral tip is available if the operator prefers a change in angulation of the material being dispensed.
5. Squeeze the dispenser to extrude the product.
6. If taking a bite registration, ensure the patient is brought into maximum intercuspation during the initial working time of 20 seconds.
7. Once maximum intercuspation has been achieved, have the patient remain in this position for 70 seconds to ensure a total set of the material.
8. Remove the material from the patient's mouth. Due to a linear dimensional change of <0.1% the practitioner can rest assured the material is highly accurate.
9. After use, store the cartridges with used mixing tips attached to protect the remaining materials from contamination. The material can be disinfected immediately.

**Precautions:**

1. Due to rapid polymerization, the patient may experience a physiologically harmless temperature increase.
2. Never attempt to extrude the material through a cured mixing canula.
3. Always disinfect cartridges between patients.
4. Store at room temperature.
5. Latex or Sulfur may affect the set of material.
6. Clone Bite may stain clothing.

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Ultradent cartridges have an expiration date and lot number printed on the label.

The lot number is a five digit number used for manufacturing purposes.


The expiration date will read YYYY-MM.

Keep out of reach of children.

For immediate reorder and/or complete descriptions of Ultradent's product line, refer to Ultradent's catalog or call Toll Free **1-800-552-5512**.

Outside U.S. call (801) 572-4200 or visit our website at [www.ultradent.com](http://www.ultradent.com).


### Ultradent® Clone™ Bite

### Hazard Rating

4 = Severe  
3 = Serious  
2 = Moderate  
1 = Slight  
0 = Minimal

For product MSDS please see our website: [www.ultradent.com](http://www.ultradent.com)

EN - Recommended Storage Temperature  
DE - Bei °C lagern  
FR - Température de stockage recommandée  
NL - Aanbevolen bewaartemperatuur  
IT - Temperatura di conservazione suggerita  
ES - Temperatura de Almacenamiento Recomendada  
PT - Temperatura de armazenamento recomendada  
SE - Rekommenderad förvarningstemperatur  
DA - Anbefalet opbevaringstemperatur  
FI - Suositeltava säilytyslämpötila  
EL - Συνιστώμενη θερμοκρασία αποθήκευσης

R<sub>x</sub>

EN-For Professional use only.  
DE-Nur für professionelle Anwendungen.  
FR-Destiné à une utilisation professionnelle uniquement.  
NL-Uitsluitend voor professioneel gebruik.  
IT-Solo per uso professionale.  
ES-Sólo para uso profesional.  
PT-Destina-se apenas à utilização profissional.  
SE-Endast för yrkesmässigt bruk.  
DA-Må kun anvendes af fagligespersonale.  
FI-Vain ammattikäyttöön.  
EL-Μόνο για επαγγελματική χρήση.

Do not reuse to avoid cross contamination.  
DE/ Nicht wiederholt verwenden, um Kreuzkontamination zu verhindern.  
FR/ Ne pas réutiliser pour éviter la contamination croisée.  
NL/ Niet opnieuw gebruiken om kruisbesmetting te vermijden.  
IT/ Per evitare la contaminazione crociata non riutilizzare.  
ES/ No lo reutilice, así evitará la contaminación cruzada.  
PT/ Não reutilizar para evitar a contaminação cruzada.  
SV/ För undvikande av korskontamination får produkten inte återanvändas.  
DA/ Må ikke genbruges for at undgå krydskontaminering.  
FI/ Älä käyttää uudelleen riskikontaminaation välttämiseksi.  
EL/ Μην το επαναχρησιμοποιείτε για να αποφύγετε τη διασταυρωτική μόλυνση.

EC REP

UP Dental GmbH  
Am Westhover Berg 30, 51149  
Cologne Germany

© 2013 Ultradent Products, Inc. All Rights Reserved. Made in Italy.

Manufactured for Ultradent Products, Inc. 505 West 10200 South, South Jordan, Utah 84095 68426.6 112113

**ULTRADENT**  
PRODUCTS, INC.