<b>Endo-Eze™ MTAFlow™ Cement</b> (Root and Pulp Treatment Material)
---

**Description:**
**The Endo-Eze™ MTAFlow™** powder and liquid/gel system consists of an extremely fine, radiopaque, inorganic powder of tricalcium and dicalcium silicate, which sets with a water-based gel for improved placement.

**Indications for Use:**

The **Endo-Eze MTAFlow** cement is intended for use in dental procedures that contact pulp and periarthral tissues. Dental procedures contacting vital pulp tissue are:

- Pulpotomy
- Root capping
- Dental procedures possibly contacting the periodaricular tissues are:
  - Root-end filling
  - Apexification
  - Perforation repair
  - Root resorption

**Contra-Indications**
Hypersensitivity against caustic (high pH) solutions, not for root canal disinfection.

**Direction for Use:**
**Powder Dispensing**
1. Shake the bottle to fluff the powder.
2. Using the enclosed spoon, place 1 level spoonful of powder (see proportions recommended in Table 1) onto the mixing paper pad (a glass slab can be used). Use a spatula educe to level the powder in the spoon.

**Gel Dispensing**
3. Shake the bottle 3 times. Make sure that all gel is concentrated in the spout of the bottle before dropping. Open the cap and dispense drops of gel (see the proportions recommended in Table 1) from the gel bottle.

**NOTE:** The plastic bottle is clear enough to see the gel inside. Avoid air bubbles inside the gel doing the movements above mentioned and verifying the gel at the tip of the bottle before dispensing the drops.

4. For some procedures, a thin consistency may be desired. Table 1 below may be used as a reference to achieve various consistencies of the **Endo-Eze MTAFlow** material for placement.

Table 1. Suggested approximated proportions of powder and gel to be used with different consistencies to be achieved

Applications	Pulp Capping, Pulpa Perforation and Pulpotomy	Resorption, Apexification, Apical Plug	Root End Filling
Powder (measuring spoon)	2 big ends (0,26g)	1 big end plus 1 small end (0,19g)	1 big end plus 1 small end (0,19g)
Gel drops	2 drops	3 drops	1 drop**
Consistency	Thick	Thin	Putty
Deliver tip	Micro Tip (20ga)	NavITip™ (29ga)*	Non-syringe delivery

\* NavITip™ tip not included in the kit
\*\* Starting with 1 drop of gel, but adding more until the desired consistency is achieved.

**NOTE:** More powder or gel can be added during or soon after the mixing to achieve the desired consistency.
5. Using a metal cement spatula, gradually mix the gel into powder until the desired consistency is obtained. Thoroughly mix with the spatula to ensure all the powder particles are hydrated. Unlike other MTA cements, that require mixing in water, due to its consistency, **Endo-Eze MTAFlow** cement will be creamy and homogeneous after mixing.

**NOTE:** **Endo-Eze MTAFlow** cement immediately after mixing, or insert into the clear Skini deliver syringe. Do not leave the mixed cement on the mixing pad for a prolonged period in contact with air, due to evaporating.
6. Remove the cap and plunger of the clear Skini syringe. Take small portions of the mixed **Endo-Eze MTAFlow** cement with the mixing spatula and insert the cement in the back part of the barrel. Replace the plunger. Attach the selected tip securely onto Skini syringe. Gently move the plunger to force the air inside the syringe. Verify flow of cement prior to applying intracanal. If resistance is met, replace the tip and re-check flow. **Use only recommended tips and do not reuse.**

**NOTE:** **Endo-Eze MTAFlow** cement mixed and placed inside the clear Skini syringe can be used for up to **15 minutes** after mixing. **Be sure no air is present inside the syringe.**
7. Place the cement, using the designated tip, in the treatement site.

**Applications:**
**Pulp Capping and Pulpotomy**
1. Complete a cavity preparation outline under rubber dam isolation.
2. Excavate all carious tooth structure using a round bur at low speed, or use hand instruments.

**For pulp capping:**
a. Control rise the exposed pulp with a sterile saline solution.
b. Control hemorrhage with pressure on the exposed pulp using a cotton pellet moistened with saline.

**For pulpotomy:**
a. Remove all remnants of coronal pulp tissue to the level of the orifice of each root canal in multi-rooted teeth. In single-rooted teeth, remove the pulp to 2mm below the level of the cemento-enamel junction.
b. If hemostasis is not achieved with cotton compression, **Astringent®** or **ViscoStat®** hemostatic solutions are recommended.
c. Dry the area.

3. Apply **Consepsis®** antibacterial solution for 60 seconds. Do not rinse. Gently dry air.

4. Insert **Endo-Eze MTAFlow** cement on the pulp exposure, or over the floor of the cavity preparation.

5. Remove excess material at the site with a dry cotton pellet.
6. After application, dry with a small blast of air. Place a flowable composite material (PermaFlow®) or a resin-reinforced glass ionomer cement (UltraCem® cement) over the **Endo-Eze MTAFlow** cement.

**NO:** A composite can be occlusal. When etching for composite placement, do not etch the **Endo-Eze MTAFlow** cement; etch only the tooth.

7. Assess the pulp vitality at three-month intervals or as needed. Periapical status should be assessed by radiographic exams.

**Perforation repair (Pulp Chamber)**
1. Apply **Consepsis®** antibacterial solution for 60 seconds. Do not rinse. Gently dry air.
2. Place Endo-Eze MTAFlow cement into the defect. Depending on the access to the defect, a NavITip tip may be used.
3. Verify the placement with a radiograph.

4. Remove excess material at the site with a sterile dry cotton pellet and rinse gently.
5. After application, dry with a gentle blast of air. Place a composite material or a glass ionomer restorative material over the **Endo-Eze MTAFlow** cement.

**Root Resorption** (External resorption, initiated in the periodontium and invading the dental pulp space)
1. Obtainate the canal space apical to the defect with **Endo-Eze MTAFlow** cement in a thin consistency (see Table 1) and a guta percha cone.

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**For pulp capping:**
a. Control the exposed pulp with a sterile saline solution.
b. Control hemorrhage with pressure on the exposed pulp using a cotton pellet moistened with saline.

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3. Apply **Consepsis®** antibacterial solution for 60 seconds. Do not rinse. Gently dry air.
4. Insert **Endo-Eze MTAFlow** cement on the pulp exposure, or over the floor of the cavity preparation.
5. Remove excess material at the site with a dry cotton pellet.
6. After application, dry with a gentle blast of air. Place a flowable composite material (PermaFlow®) or a resin-reinforced glass ionomer cement (UltraCem® cement) over the **Endo-Eze MTAFlow** cement.

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**Perforation repair (Pulp Chamber)**
1. Apply **Consepsis®** antibacterial solution for 60 seconds. Do not rinse. Gently dry air.
2. Place Endo-Eze MTAFlow cement in areas of the tooth in contact with the gingival sulcus.
3. Verify the placement with a radiograph.

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**Root Resorption** (External resorption, initiated in the periodontium and invading the space of the pulp dentare)
1. Obtainate the canal space apical to the defect with the ciment **Endo-Eze MTAFlow** in a consistance fluide (voir le tableau 1) et un cône de guta percha.

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