





SUGGESTED CURING TIMES

STANDARD

XTRA POWER

XTRA POWER QUADRANT

POWER (mW/cm ²)	1200* ●	3200* ☀	3200* ☀
 <p>POWER BUTTON</p> <p>TIMING INTERVAL LIGHTS</p> <p>MODE/STATUS LIGHT</p> <p>TIME/MODE CHANGE BUTTON</p>	 <p>5 sec. 10 sec. 15 sec. 20 sec.</p>	 <p>1 sec. 2 sec. 3 sec.</p>	 <p>5 consecutive bursts at 3 seconds each</p> <p>3 sec.</p>
METAL & CERAMIC BRACKETS (CURE TIME PER BRACKET)	1 X 10 SECONDS	2 X 3 SECONDS	2 X 3 SECONDS
RESTORATIVE CURING (2MM INCREMENTS)			
PER LAYER	1 X 10 SECONDS	1 X 3 SECONDS	—
FINAL CURE	2 X 10 SECONDS	2 X 3 SECONDS	—

VALO on Xtra Power setting demonstrates that one 3-second exposure directly over the labial face of the bracket, followed later by a second exposure, for a total of 6 seconds, is all that is required.

Applicable for all current and future Opal Orthodontics and Ultradent light-cured materials.

*As measured by a properly calibrated Demetron® LED Radiometer

†As measured by an integrating sphere



888.863.5883
opalorthodontics.com



Isopropyl alcohol will not damage the curing light and is the best option to use when cleaning and wiping down your VALO.