Citric Acid 20% Solution

Indications for use: Citric Acid 20%, a root canal cleanser, is intended to etch root canal walls just prior to obturation to allow an optimum seal.

Description: Citric Acid 20% is a slightly thickened, water-soluble solution. It is tinted slightly orange for easy identification and is provided in a 30ml IndiSpense® Syringe.

Preparation: 1. An empty 1,2ml or 5ml syringe is attached to the end of the large, no-waste IndiSpense syringe by turning the luer lock of the unit dose syringe snugly onto the male thread of the IndiSpense as show in Fig. 1.

2. With palm grasp (Fig. 2) the plunger of the IndiSpense syringe is depressed while the unit dose syringe is stabilized at the plunger with the nondominant hand. Slightly pull on plunger of unit dose syringe to facilitate loading. Do not express contents of unit dose syringe back into IndiSpense syringe. To prevent cross-contamination, a used syringe should never be attached to the IndiSpense syringe.

3. Remove luer lock cap from syringe. Twist a NavixF® Tip or Navitip® firmly onto syringe. Gauge NavixF Tip with the rubber stop at least 2mm short of the root terminus as show in Fig. 3.

Procedure:

Warning: Do not follow Citric Acid rinse with sodium hypochlorite or vice versa, as this can lead to the production of chlorine gas. Always rinse with clean, bacteria-free water between Citric Acid and sodium hypochlorite irrigation.

1. Make sure NavixF Tip is loose in canal and at least 2mm short of terminus. Gently deliver the Citric Acid with NavixF Tip. Move brush tip up and down to agitate and brush solution, cleaning walls to dissolve and remove calcium hydroxide with Citric Acid as shown in Fig. 4. Do not express Citric Acid when canula tip is at full length or if it fits tightly against the canal walls.

2. Flush out canal with clean, bacteria-free water, followed by copious sodium hypochlorite (e.g. ChloroCid®). When Citric Acid is used in lieu of EDTA liquid to remove smear layer, let Citric Acid dwell in canal for 3 minutes. Perform a final rinse with clean, bacteria-free water, plain EDTA or local anesthetic. Quick and effective removal of most of the Citric Acid solution or other irrigants is facilitated by using a Luer Vacuum Adapter with a Capillary Tip or 10mm Micro Capillary Tip.

3. One or two paper points are usually all that are required for final drying of canals if canals are first vacuumed out as described in step 2 above.

4. Fit master cone, then obturate. Ultratrend recommends the EndoREZ® Obturation System including EndoREZ root canal sealer, EndoREZ resin-coated gutta percha points, and specially designed delivery systems.

Warnings: Do not follow Citric Acid rinse with sodium hypochlorite or vice versa, as this can lead to the production of chlorine gas. Always rinse with clean, bacteria-free water between Citric Acid and sodium hypochlorite irrigation.

Precautions: 1. Always use a rubber dam for endodontic therapy. Prevent ingestion of Citric Acid solution. 2. Prevent getting solution in the eyes or leaving it on the soft tissues for an extended period of time. 3. Do not extend tip past the root apex or express Citric Acid past the apex.

4. Do not express Citric Acid from the canula tip at the root terminus or if the tip fits tightly against the canal walls. This may force solution past the apex.

5. Do not use this product on patients with a known sensitivity to Citric Acid.

6. The tips and delivery syringes are disposable and intended for single use.

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EN


Fig. 1. Empty 1.5 ml or 5 ml syringe is inserted into the end of the IndiSpense syringe. Push up and down. Rinse barrel thoroughly with plain EDTA or clean bacteria the water or local anesthetic.


Fig. 2. Pull the plunger back to the stop. Push the plunger back to the stop. Push the plunger back to the stop. Push the plunger back to the stop.

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Fig. 3. Measure NavifX tip 2mm short of root apex and loosely fit into canal.

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