ENI Description and Indications for Use: FORMA™ restorative composite contains bis-GMA, TEGDMA, Bis-EMA, and UDMA. FORMA composite is designed for restorations in anterior and posterior teeth and is suitable for occlusal surfaces. Its formula contains a zirconia/silica inorganic filler combination and barium glass. All shades of FORMA composite are based on the Vita® shade guide and are radiopaque. FORMA composite has a radiopacity value greater than 2 mm of aluminum. Aluminum has a radiopacity equivalent to that of dentine. Thus 1 mm of material having a radiopacity equivalent to 1 mm of aluminum has a radiopacity equivalent to that of dentine and 2 mm of aluminum is equivalent to enamel. Filler load is 64.5% by volume for the body shades, 64.8% for the enamel shades, and 63.5% for the dentin shades. FORMA composite is a nano-hybrid composite with nanoparticles in the 5–50 nm range. For all products described, carefully read and understand all instructions and SDS information prior to use.

List of supporting materials accompanying the product:

The respective Instructions for Use manual accompanies the product.

FORMA Syringe

FORMA Capsule

Directions for use:

Preparation:

1. Apply dentin bonding agent to preparation according to manufacturer's instructions.

Direct Restoration - Syringes:

1. Open cap and twist handle (stem) clockwise to express desired amount of composite. Resistance will lessen with each twist.
2. Place FORMA composite in incremental layers applying the first layer approximately 0.5 – 1.5 mm thick and cure. Follow with up to 2 mm increments; curing after each increment. For curing lights with an intensity of ≥ 800 mW/cm² in the wavelength region of 385 – 515 nm, cure each layer for 10 seconds. A final cure of 20 seconds is recommended.

3. After final increment is placed, sculpt and light cure.
4. Shape, finish, and adjust occlusion using preferred method.
5. Immediately replace cap after dispensing composite to prevent polymerization from ambient light.

Direct Restoration - Singles (capsule):

1. Keep capsule caps on until use.
2. Snap capsule firmly into composite gun before trying to express.
3. Unit dose capsules are designed for single-patient use only. Do not reuse to avoid the risk of cross-contamination. Discard capsule after use.
4. Place FORMA composite in incremental layers applying the first layer approximately 0.5 – 1.5mm thick and cure. Follow with up to 2 mm increments; curing after each increment. For curing lights with an intensity of ≥ 800 mW/cm² in the wavelength region of 385 – 515 nm, cure each layer for 10 seconds. A final cure of 20 seconds is recommended.
5. After final increment is placed, sculpt and light cure.
6. Shape, finish, and adjust occlusion using preferred method.

Indirect Restorations:

1. Make impression of prepared tooth and opposing dentition.
2. Prepare master cast and articulate opposing cast using preferred technique.
3. Place FORMA composite in incremental layers applying the first layer approximately 0.5 – 1.5 mm thick and cure. Follow with up to 2 mm increments; curing after each increment. For curing lights with an intensity of ≥ 800 mW/cm² in the wavelength region of 385 – 515 nm, cure each layer for 10 seconds.
4. Use final increments to create proximal contours and ensure marginal adaptation, a final cure of 20 seconds is recommended.
5. Remove restoration. Clean off residual die stone and trim marginal flash. Additions can be made by roughening restoration surface, applying adhesive, and placing FORMA composite as needed.
6. Sandblast internal aspect of restoration with 50 µm alumina abrasive.
7. Seat restoration onto master cast, evaluate fit, occlusion, and contours. Shape, finish, and adjust occlusion using preferred method.
8. Bond restoration intraorally per adhesive manufacturer's instructions.
9. Following cementation, make necessary occlusal adjustment and polish adjusted surfaces.

Precautions and warnings:

1. Read and understand all instructions before using product.
2. Store at temperature within 2–27°C.
3. Use a quality curing light to ensure maximum physical properties of the composite. Uncured material will result in an inferior restoration and may also cause pulpal irritation.
4. Clinician and patient should wear UV and blue light protective eyewear when curing composite.
5. To lessen pre-polymerization of the composite, minimize exposure to bright lights during the procedure.
6. To optimize bond to each succeeding layer, do not remove oxygen inhibition layer between incremental layers of composite.
7. Ensure composite is totally polymerized between layers and interproximally.
8. Syringes are not for intraoral use.
9. Resins can be sensitizing. Avoid repeated contact of uncured dental resin with skin. Do not use on patients with a known sensitivity to acrylates or other resins. If allergic reaction, dermatitis or rash develops, consult a physician.
10. To avoid pre-polymerization of the composite by ambient light, make sure to always cover the tip of the syringe with its cap after composite has been removed.
11. To avoid cross-contamination, use a clean instrument every time you remove the composite from the syringe.
12. Disinfect contaminated syringes with an intermediate level disinfectant between uses.
13. Properly dispose of syringes, according to current legislation.
14. Keep away from extreme heat and direct sunlight.
15. Do not use after expiration date noted on containers.
16. Keep out of reach of children.

NOTE: for product SDS and additional information about FORMA composite, please see website: www.ultradent.com.br.

FOR PROFESSIONAL USE ONLY

USE BY DATE

CATALOG NUMBER

LOT NUMBER

DO NOT RE-USE TO AVOID
CROSS CONTAMINATION

KEEP AWAY FROM
HEAT/SUNLIGHT

FOR PROFESSIONAL
USE ONLY

RECOMMENDED STORAGE
TEMPERATURE (2-27°C)

REFER TO INSTRUCTIONS
FOR INFORMATION ON
HOW TO USE

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