# User Manual

J58349

**ultrawave™**  
PIEZO ULTRASONIC

# Contents

---

## **1 Documentation**

- 1.1 Associated documentation
- 1.2 Electronic documentation

## **2 Required information**

- 2.1 Indications for use
- 2.2 Operating principle
- 2.3 Date of inclusion of EC marking
- 2.4 Latest document update
- 2.5 Repairing or modifying the device
- 2.6 Warranty
- 2.7 Accessory usage conditions

## **3 Removal from packaging, installation, connections**

- 3.1 Unpacking your medical device
- 3.2 Positioning the medical device
- 3.3 Installing cords
- 3.4 Connecting the medical device to the water system
- 3.5 Connecting the medical device to the electrical network
- 3.6 Plugging the medical device to the electrical network
- 3.7 Installing the control pedal

## **4 Description of the medical device**

- 4.1 Control unit
- 4.2 Power configuration dial
- 4.3 Adjusting ultrasound power
- 4.4 Attach a tip or a file
- 4.5 Connecting and disconnecting accessories during use
- 4.6 Connecting the handpiece
- 4.7 Handpiece cord
- 4.8 Light indicator
- 4.9 Irrigation flow configuration button
- 4.10 Adjusting the irrigation
- 4.11 Initiating irrigation
- 4.12 Air inlets
- 4.13 Connection to water system
- 4.14 Control pedal
- 4.15 Activating ultrasounds using the pedal
- 4.16 Cleaning the irrigation system

## **5 Cleaning, disinfecting and sterilizing**

- 5.1 Cleaning and disinfection of the medical device
- 5.2 Cleaning and disinfecting accessories

## **6 Monitoring and maintenance of the medical device**

### **7 Maintenance**

- 7.1 Replacing the water filter
- 7.2 Identifying incorrect operation
  - 7.2.1 No operation
  - 7.2.2 No spray
  - 7.2.3 The power is not as expected
  - 7.2.4 Ultrasounds not working
  - 7.2.5 Water leakage

### **8 Technical specifications for the medical device**

- 8.1 Identification
- 8.2 Control unit

# Contents

---

|           |                                  |
|-----------|----------------------------------|
| 8.3       | Ultrasonic generator             |
| 8.4       | Length of cords                  |
| 8.5       | Irrigation                       |
| 8.6       | Control pedal                    |
| 8.7       | Environmental characteristics    |
| 8.8       | Environmental restrictions       |
| 8.9       | Main performance characteristics |
| <b>9</b>  | <b>Regulations and standards</b> |
| 9.1       | Official Texts                   |
| 9.2       | Medical class of the device      |
| 9.3       | Standardized Symbols             |
| 9.4       | Manufacturer identification      |
| 9.5       | Distributor                      |
| 9.6       | Disposal and recycling           |
| <b>10</b> | <b>Index</b>                     |

## Foreword

The Ultradent Products Inc. medical device that you are about to install and use in your practice is a medical device designed for professional use. It comprises the chosen tool with which you will provide treatment within the context of your work.

To ensure optimum safety for yourself and your patients, comfort in your daily practice and to benefit fully from the technology of your medical device, please read the documentation provided carefully.

If you have received this medical device by mistake, please contact the supplier to arrange for it to be collected.

Please refer to the general instructions relating to the comprehensive range of dental ultrasonic generators by Ultradent Products Inc. for information about the following:

- documentation format;
- the documentation archiving period;
- warnings concerning user population;
- the treatment area;
- the medical device usage interactions, contraindications and prohibitions;
- electromagnetic compatibility;
- disposal and recycling of the medical device;
- manufacturer responsibility.

Please refer to the accessory cleaning, disinfection and sterilization protocols and the handpiece cleaning, disinfection and sterilization protocols for information about the following:

- preparation of parts for sterilization;
- detailed manual and automatic protocols;
- information concerning the sterilization;
- recommendations for the inspection of parts.

## 1 Documentation

This document contains the following information:

- indications for use;
- description of the medical device;
- installation of the medical device;
- use of the medical device;
- preparation for cleaning and disinfection of the medical device;
- monitoring and general maintenance of the medical device;
- maintenance to be performed by the user.

## 1.1 Associated documentation

This document must be used in association with the following documents:

| Document title                                                                                              | References |
|-------------------------------------------------------------------------------------------------------------|------------|
| Cleaning, disinfection and sterilization protocols for Wrenches                                             | J81009     |
| Cleaning, disinfection and sterilization protocols for Tips                                                 | J02009     |
| Cleaning, disinfection and sterilization protocols for Ultradent Products Inc.Handpieces                    | J12949     |
| General instructions relating to the complete range of Ultradent Products Inc. dental ultrasonic generators | J00039     |
| Ultrawave™ User Manual                                                                                      | J58349     |
| Ultrawave™ Quick Start                                                                                      | J58350     |

## 1.2 Electronic documentation

The user instructions for your device are in part provided in electronic format and in part in printed format. However, you can request a free printed copy of the user instructions via the website, by telephone or in writing.

The electronic user instructions are available in PDF format (Portable Document Format) and you will need to have a PDF file reader software installed to read the instructions.

The device user instructions can be consulted at the following addresses: [www.ultradent.com](http://www.ultradent.com) and [www.satelec.com](http://www.satelec.com).

Electronic user  
informations

It is important for you to have read and understood the content of the user instructions relating to the use of your device and its accessories prior to use.  
We recommend that you visit the website regularly to consult and/or to download the latest version of your device's user instructions.

## 2 Required information

### 2.1 Indications for use

This medical device is intended for prophylaxis, including scaling, for periodontics, for endodontics, and for preservation and restoration dentistry, including prosthesis.

### 2.2 Operating principle

An electrical signal emitted by the medical device is supplied to the dental ultrasonic handpiece. This is connected to the medical device via a cord. The handpiece comprises a piezoelectric ceramic transducer, which transforms the electrical signal into ultrasonic vibrations. Mechanical vibrations are transmitted to a tip or a dental file attached to the end of the ultrasonic handpiece.

### 2.3 Date of inclusion of EC marking

20(14)

## 2.4 Latest document update

03/2014

## 2.5 Repairing or modifying the device

Contact the supplier of your device. Using the services of an unapproved repairer could render your device dangerous for you and your patients.

Do not repair or modify the medical device without seeking the prior permission of Ultradent Products Inc..

If the device is modified or repaired, specific checks and tests must be carried out to ensure that the device is still safe to use.


In the event of doubt, contact an approved dealer or the customer service team:

[customer.support@ultradent.com](mailto:customer.support@ultradent.com)

[www.ultradent.com](http://www.ultradent.com)

Ultradent Products Inc. at the request of technical personnel working for the network of dealers approved by Ultradent Products Inc., provides all information required to repair the faulty parts on which they may perform repairs.

## 2.6 Warranty

The screws marked A, B and C must not be unscrewed by the user under any circumstances as this may void the warranty for the medical device.

## 2.7 Accessory usage conditions

Accessories and handpiece must be cleaned, disinfected and sterilized prior to use.

# 3 Removal from packaging, installation, connections

## 3.1 Unpacking your medical device

When you receive your medical device, check for any damage that may have occurred during transportation. If you have received this medical device by mistake, please contact the supplier to arrange for it to be collected.

If you have any questions or requirements, contact your supplier.

The Ultrawave™ includes the following items:

- a Ultrawave™ unit with non-detachable pedal cord, a non-detachable cord and a handpiece support ;
- tips and wrenches depending on selected options;
- a Quick Start Ultrawave™ [J58350].

## 3.2 Positioning the medical device

Place the control unit in the position that is suitable for your activity.

The medical device must be placed on a secure and flat surface or a surface with a maximum slope of 5 degrees.

Adjust the position of your medical device to correspond to your angle of vision and the characteristics of your workstation, e.g. lighting or distance between the user and the medical device.

Ensure that you can access your medical device quickly.

- | Do not install your medical device near or on another device.

## 3.3 Installing cords

Check that the cords do not hinder the movement or free circulation of anyone.

Never rotate the handpiece connector on its cord as this can damage your medical device.

Never wrap the handpiece cord around the medical device.

Make sure that it is not possible to wheel over or walk on the different cords.

The cord attached to its handpiece must be easily accessible. Make sure that the cord is slack during use.

- | Do not put the medical device cords in a cable cover or a cable tray.

## 3.4 Connecting the medical device to the water system

- | The information below only applies to a medical device that needs to be connected to the water system to operate.

Ask an approved dental installation technician to connect your medical device to the water system.

The water supply system pressure may vary throughout the day. The water supply system pressure must be adapted to the values recommended for your medical device. It is very important to make sure that the maximum pressure permitted for the medical device is never reached or exceeded. If in doubt, you are strongly advised to install or arrange for installation of a water pressure limiting system.

- | The water supply system must comply with the quality criteria compatible with the practice of dental treatments.

## 3.5 Connecting the medical device to the electrical network

Set the medical device to OFF position 0 and check that the power voltage is compatible with that indicated on the medical device or its power plug. Next, connect the cord to the wall socket in compliance with the standards in force in the country of use.

A different voltage would cause damage to the medical device and could injure the patient and/or user. Any variation in the electrical network voltage or electromagnetic field that is non-compliant with the limits in force, could interfere with the medical device's operation.

The medical devices' ground must be connected to a supply network equipped with a protective ground.

- | Do not plug the medical device into an extension lead and do not put the power cord in a cable cover or cable tray.

## 3.6 Plugging the medical device to the electrical network

1. Connect the power adapter to the medical device connector;

## 3.7 Installing the control pedal

The control pedal must be positioned near the feet of the operator and must be readily accessible.

# 4 Description of the medical device

## 4.1 Control unit

The control unit incorporates Newtron® technology patented by SATELEC® S.A. N°US6,765,333B1.

The patented Newtron® technology controls the tips by Cruise Control®, an automatic system for setting the frequency and power in real time. This ensures that Satelec® tip vibration is gentle, regular and controlled.

The control unit incorporates an ultrasound generator equipped with a piezoelectric command.

## 4.2 Power configuration dial

To ensure a quality treatment, you will need to use the tips at the power and irrigation flow settings recommended by Ultradent Products Inc..

The ultrasound power configuration dial ensures:

- Switching On/Switching off of the device (OI symbol).
- Configuration of the operating power: 1 to 20.

The configuration dial has 20 graduations.

- Green: 1 to 9 : very low to low power, used mainly for periodontics.
- Yellow: 6 to 14 : medium power, used mainly for endodontics.
- Blue: 12 to 18 : high power, used mainly for scaling.
- Orange: 17 to 20 : very high power, used mainly for preservation and restoration dentistry.

The ultrasound power configuration dial is not designed to be removed.

## 4.3 Adjusting ultrasound power

Adjust the ultrasound power using the ultrasound power configuration dial. The ultrasound power must be adjusted in accordance with the tip used and the required treatment. The operating power of the tips must be selected in compliance with the Satelec tips color coding system (CCS tips). Details of these indications are given in the adjustment table available at the address [www.satelec.com/documents](http://www.satelec.com/documents) and also via the Tips recommendations for use.

## 4.4 Attach a tip or a file

A tip or a file vibrates correctly when it is perfectly tightened without being forced beyond its stop point. Tighten it moderately using the wrench provided to ensure optimum ultrasound operation. Over-tightening of the tip or file can result in breakage of the tip, file or handpiece.

| To prevent self-locking of the tip or the file, the latter must be removed after each use.

## 4.5 Connecting and disconnecting accessories during use

| Do not connect/disconnect the cord(s) or the handpiece when the medical device is switched on and your foot is on the pedal.

| Do not tighten or loosen the tips when the handpiece is activated.

## 4.6 Connecting the handpiece

Check for the absence of signs of moisture at the handpiece connections, and eliminate them if necessary (wipe and blow using a multipurpose syringe).

| Lubricate the irrigation system seal located behind the handpiece with dental instrument lubricant to extend its effectiveness and prevent leaks.

Connect the handpiece to the sleeve, by avoiding rotation movement.

Install the handpiece on the support.

## 4.7 Handpiece cord

The cord is only compatible with Ultradent Products Inc. handpieces .

The cord ensures irrigation circulation and electrical connection between the medical device and the Ultrawave™ handpiece.

## 4.8 Light indicator

The light indicator is designed to provide information about the status of the device.

When the light indicator is illuminated, the medical device is on and ready to use.

## 4.9 Irrigation flow configuration button

The irrigation flow configuration button starts and stops the irrigation function as well as adjusts the irrigation flow.

The irrigation flow configuration button is not designed to be removed.

## 4.10 Adjusting the irrigation

Adjust the irrigation flow using the irrigation flow configuration button. This adjustment depends on the tip and the treatment.

As work habits, feedback and professional training differ from one professional to another, the user must make sure that the irrigation flow is perfectly adapted to the treatment to be carried out to avoid burning the treatment area.

## 4.11 Initiating irrigation

The medical device must be set to minimum power depending on the required irrigation flow rate. Press the pedal until a spray appears.

## 4.12 Air inlets

Air inlets ensure correct ventilation of the control unit. Leave them uncovered to allow air to circulate.

## 4.13 Connection to water system

The supply pipe connector is used to connect the medical device to the domestic water distribution system.

The connector is extended by a pipe to which a filter is attached. The filter needs to be cleaned and/or replaced regularly as specified in the chapter Replacing the water filter page 19.

| The water quality must meet the criteria required to perform dental treatments.

## 4.14 Control pedal

The ON/OFF type pedal is used by the practitioner to operate the medical device.

Pressing the pedal automatically activates the handpiece ultrasounds, and the irrigation function if it is not in 0 position.

The control pedal equipped with its cord cannot be disconnected. Its weight and antislip pad ensure good stability.

## 4.15 Activating ultrasounds using the pedal

To activate the ultrasounds on your medical device, press the control pedal.

## 4.16 Cleaning the irrigation system

After installation and before first use, at the end of the day and following a period of prolonged non-use of the medical device, it is important to clean the irrigation system.

When the irrigation system has been cleaned, perform the following operations:

1. disconnect the handpiece and refer to the handpiece cleaning, disinfection and sterilization protocols by Ultradent Products Inc. J12949;

2. clean and disinfect the medical device as indicated in the chapter *Cleaning and disinfection of the medical device* page 9

3. follow the instructions for accessory cleaning, disinfection and sterilization protocols as outlined in the documents [J81009] and [J02009].

# 5 Cleaning, disinfecting and sterilizing

The instructions relating to the accessory cleaning, disinfection and sterilization protocols provided by Ultradent Products Inc. have been approved for each medical device and accessory. The applicable guides are listed in chapter *Associated documentation page 4*. They can be downloaded at the following addresses [www.satelec.com/documents](http://www.satelec.com/documents) or [www.ultradent.com](http://www.ultradent.com). In all cases, the local regulations in force relating to the accessory cleaning, disinfection and sterilization protocols take precedence over the information provided by Ultradent Products Inc.

## 5.1 Cleaning and disinfection of the medical device

The medical device must be in OFF or O stop position during cleaning and disinfecting procedures.

Refer to the instructions detailed in the chapter *Cleaning the irrigation system page 8*

Avoid using cleaning and disinfection products that contain flammable agents.

Otherwise, ensure that the product has completely evaporated from or that there is nothing flammable left on the medical device and its accessories before switching it on.

- | Do not use abrasive product to clean the medical device.
  - | Never apply sprays directly to the medical device to clean it. Always spray the product onto a wipe, then clean the medical device.
- The medical device control unit, handpiece cord and control pedal must be cleaned and disinfected daily.
- | Use alcohol-impregnated wipes for disinfection.

## 5.2 Cleaning and disinfecting accessories

Refer to the accessory cleaning, disinfection and sterilization protocols listed in the chapter *Associated documentation page 4*.

# 6 Monitoring and maintenance of the medical device

Before and after each use, check that the device and its accessories are not faulty in any way. This is necessary to detect any insulation fault or damage. If necessary, replace damaged parts.

Monitor the cleanliness of the air inlets on the control unit to prevent any heating.

# 7 Maintenance

The only preventive maintenance the medical device requires is:


- checking of accessories;
- everyday cleaning, disinfection and sterilization procedures;
- cleaning.
- replacement of the water filter cartridge.



## 7.1 Replacing the water filter

The water filter must be replaced or at least cleaned every 6 months.

Proceed as follows:

- shut off the water supply;
- stop the medical device (position O);
- disconnect the network plug;
- unscrew the two filter sections;

- check that the spray works correctly and that there is no leakage.  
| A damaged or blocked cartridge must be replaced.

## 7.2 Identifying incorrect operation

In the event of incorrect operation, refer to the tables below to quickly identify and repair the non-complex parts of the medical device.

If the incorrect operation is not described in the tables below, please contact your supplier or the After-Sales team at Ultradent Products Inc..

Do not use the medical device if it appears to be damaged or faulty. Isolate the medical device and make sure that it cannot be used.

### 7.2.1 No operation

Symptoms: the indicator light on the medical device is off and the medical device is not working.

| Possible causes                                                         | Solutions                                            |
|-------------------------------------------------------------------------|------------------------------------------------------|
| No electrical current                                                   | Contact your electrician                             |
| Internal fuse not working                                               | Return to Ultradent Products Inc. After-Sales team   |
| Faulty connection between the power cord and the electrical wall socket | Connect the power cord to the electrical wall socket |

| The medical device also has internal fuses F1, F2 and F3 that cannot be accessed by the user.

### 7.2.2 No spray

Symptoms: There is no water spray at the tip.

| Possible causes                      | Solutions                            |
|--------------------------------------|--------------------------------------|
| Flow configuration button on minimum | Adjust the flow configuration button |
| Faulty water pipe connection         | Check the water inlet                |
| Low water pressure                   | Check the water system pressure      |
| Blocked filter                       | Clean or change the filter           |
| Tip or file blocked                  | Unblock the tip or file              |
| Incorrect choice of tip              | Check the tip                        |
| Inadequate amount of spray           | Adjust the spray                     |

## 7.2.3 The power is not as expected

Symptoms: the tip does not vibrate at the expected frequency, the treatment is not progressing as normal and is taking longer or is at a standstill.

| Possible causes                                                             | Solutions                                                                                                              |
|-----------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|
| Worn or distorted tip                                                       | Replace the tip                                                                                                        |
| Incorrect use: incorrect approach angle or inadequate pressure on the tooth | Refer to the configuration table available at <a href="http://www.satelec.com/documents">www.satelec.com/documents</a> |
| Presence of liquid or moisture between the handpiece and cord               | Thoroughly dry the electrical contacts                                                                                 |
| Wear of the O'ring handpiece seal                                           | Replace the O'ring handpiece seal                                                                                      |

## 7.2.4 Ultrasounds not working

Symptoms: the tip does not vibrate, vibration cannot be heard.

| Possible causes            | Solutions                                                                    |
|----------------------------|------------------------------------------------------------------------------|
| Tip loose                  | Tighten the tip using the wrench                                             |
| Faulty connector contact   | Clean the cord contacts                                                      |
| Handpiece cord wire(s) cut | Return to Ultradent Products Inc. After-Sales Department to replace the cord |

## 7.2.5 Water leakage

Symptoms: Water is leaking from one of the following places:

- between the base of the handpiece and its cord

| Possible causes                             | Solutions                         |
|---------------------------------------------|-----------------------------------|
| Wear of 1.15 mm x1 mm O'ring handpiece seal | Replace the seal using F12304 kit |

# 8 Technical specifications for the medical device

## 8.1 Identification

|                            |            |
|----------------------------|------------|
| Manufacturer               | SATELEC®   |
| Name of the medical device | Ultrawave™ |

## 8.2 Control unit

|                           |      |
|---------------------------|------|
| Width (in mm)             | 150  |
| Height (in mm)            | 80   |
| Depth (in mm)             | 170  |
| Weight (in g)             | 1380 |
| Ingress protection rating | IPX0 |

## 8.3 Ultrasonic generator

|                                           |                                                                                                                                                                                   |
|-------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Supply voltage                            | 110V - 120V                                                                                                                                                                       |
| Power consumption                         | 30 VA                                                                                                                                                                             |
| Output frequency                          | Minimum 28 kHz                                                                                                                                                                    |
| Operating mode                            | Intermittent: 10 minutes ON / 5 minutes OFF                                                                                                                                       |
| Type of leakage currents                  | BF                                                                                                                                                                                |
| Electrical rating                         | 2                                                                                                                                                                                 |
| Internal fuses not accessible to the user | Ref FU1 et FU2: 315 mA <sub>T</sub> à 230 VAC<br>Ref FU1 et FU2: 400 mA <sub>T</sub> à 115 VAC<br>Ref FU1 et FU2: 400 mA <sub>T</sub> à 100 VAC<br>Ref FU3: CMA / 1,5AT / 125 VAC |

## 8.4 Length of cords

|                               |      |
|-------------------------------|------|
| Scaler handpiece cord (in mm) | 2000 |
| Control pedal cord (in mm)    | 2000 |

## 8.5 Irrigation

|                                  |                                                 |
|----------------------------------|-------------------------------------------------|
| Water pressure at inlet          | 1 to 5 bar                                      |
| Water output adjustable with tap | Min: 0 ml/min<br>Max: 90 ml/min<br>on half turn |

## 8.6 Control pedal

|                           |      |
|---------------------------|------|
| Width (in mm)             | 70   |
| Height (in mm)            | 30   |
| Depth (in mm)             | 95   |
| Weight (in g)             | 150  |
| Ingress protection rating | IPX1 |

## 8.7 Environmental characteristics

|                       |                                    |
|-----------------------|------------------------------------|
| Operating temperature | +10 °C to +30 °C                   |
| Storage temperature   | -0 °C to +60 °C                    |
| Operating humidity    | 30% to 75%                         |
| Storage humidity      | 10% to 70%, including condensation |
| Atmospheric pressure  | Between 800 hPa and 1060 hPa       |
| Altitude              | Less than or equal to 2000 meters  |

## 8.8 Environmental restrictions

|                              |                                                                                                                              |
|------------------------------|------------------------------------------------------------------------------------------------------------------------------|
| Usage premises               | Can be used at all medical premises. The medical device must not be used in an operating theater, or outside.                |
| Use in gas-filled atmosphere | The medical device is not designed for use in a type AP or APG gas-filled atmosphere or in the presence of anesthetic gases. |
| Immersion                    | The handpiece must not be immersed.                                                                                          |

## 8.9 Main performance characteristics

Ultrasonic vibrations of the tip or file fitted to the end of the conventional dental ultrasonic handpiece.

- Vibration frequency  $\geq 28$  kHz.
- Tip amplitude  $\leq 200$   $\mu\text{m}$ .

# 9 Regulations and standards

## 9.1 Official Texts

This medical device complies with the essential requirements of European Directive 93/42/EEC. This equipment is designed and developed in compliance with Electrical Safety standard IEC60601-1 in force. It was designed and manufactured in accordance with an EN ISO 13485-certified quality assurance system.

## 9.2 Medical class of the device

This medical device is a class IIa device according to European Directive 93/42/EEC.

## 9.3 Standardized Symbols

| Symbols                                                                                                           | Meaning                                         |
|-------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|
|                                  | Refer to the accompanying documentation         |
|                                  | Consult the User Manual                         |
| <br>Electronic user informations | Accompanying documentation in electronic format |
|                                  | BF type                                         |
|                                  | Class 2                                         |
|                                  | Sterilization at 132 °C in an autoclave         |
|                                  | Washer disinfector for thermal disinfection     |
|                                  | EC marking                                      |

| Symbols                                                                                | Meaning                                                                                                                                                                                   |
|----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|       | Do not dispose of as household waste                                                                                                                                                      |
| YYYY  | Year of manufacture                                                                                                                                                                       |
|       | Control pedal                                                                                                                                                                             |
| IPX1                                                                                   | IP: ingress protection ratings procured by a range<br>X: no ingress protection rating claim against the penetration of solids<br>1: protects against the vertical falls of drops of water |
| Rx Only                                                                                | Caution: US Federal law restricts this device to sale by or on the order of a Physician                                                                                                   |

## 9.4 Manufacturer identification

SATELEC

A Company of ACTEON Group

17, avenue Gustave Eiffel

BP 30216

33708 MERIGNAC cedex

FRANCE

Tel. +33 (0) 556.34.06.07

Fax. +33 (0) 556.34.92.92

E.mail: satelec@acteongroup.com.

www.acteongroup.com

## 9.5 Distributor

Ultradent Products, Inc.

505 West 10200 South

South Jordan, UT 84095 USA

E-mail : customer.support@ultradent.com

www.ultradent.com

## 9.6 Disposal and recycling

As an item of Electrical and Electronic Equipment, the medical device must be disposed of via a specialist collection, removal, recycling or destruction channel. This applies in particular to the European market, in reference to Directive 2002/96/EC dated 27/01/2003.

When your medical device has reached the end of its service life, contact Ultradent Products Inc. to find out how to proceed. Relevant contact details are given in the chapter *Distributor page 16*

# 10 Index

## A

After installation 8  
air circulation 8  
air inlets 8, 9  
Altitude 13  
Amplitude 13  
approved dealers 5

## C

clean and disinfect the device 8-9  
clean irrigation system 8-9  
cleaned 5  
color coding system 7  
control pedal 7-8  
control unit 7  
Cruise Control® 7

## D

damage 9  
dental file 4  
dental ultrasonic generators 4  
disinfected 5  
disposal 16

## E

Electrical Safety 13  
electronic user instructions 4

end of the day 8  
endodontics 4  
European Directive 13

## F

fault 9  
filter 10  
first inclusion of EC marking 4  
first use 8  
fuse 10

## G

gas-filled atmosphere 13  
general instructions relating to the comprehensive  
range of dental ultrasonic generators 3

## H

Handpieces 4, 7

## I

incorrect operation 10  
indicator light 10  
irrigation flow 7-8

## K

kit F12304 11

## L

light indicator 8  
Liquids 8

## M

Manufacturer 11  
Medical class 13  
moisture 7

## N

non-use 8

## P

pedal 8  
periodontics 4  
power settings 7  
preservation and restoration dentistry 4  
Pressure 13  
prophylaxis 4  
prosthesis 4

## R

Regulatory 4-5, 16  
repair 5

## S

scaling 4  
seal 11  
spray 8, 10  
sterilized 5  
supply pipe 8

## T

Temperature 13  
tip 4, 10-11  
treatment 7-8

## U

ultrasonic vibrations 4  
ultrasound power 7  
User Manual 4

## V

Vibration frequency 13

## W

water filter 9  
water leak 11  
water system 6  
Wrenches 4

# General Instructions

Range of dental ultrasonic generators

J00039

# Contents

---

- 1 Documentation**
  - 1.1 Electronic documentation
  - 1.2 Archiving duration
- 2 Required information**
  - 2.1 Latest document update
  - 2.2 Official Texts
  - 2.3 Standardized Symbols
  - 2.4 Manufacturer identification
  - 2.5 Manufacturer responsibility
  - 2.6 Distributor
- 3 Warnings**
  - 3.1 Federal Law
  - 3.2 Warning applicable to all countries in which the device is sold
  - 3.3 User population
  - 3.4 Specific user training
  - 3.5 Patient population restriction
  - 3.6 Parts of the body or types of tissues treated
  - 3.7 Elements applied
  - 3.8 Basic safety in normal use
  - 3.9 Normal usage conditions
  - 3.10 Irrigation spray
  - 3.11 Tip service life
  - 3.12 Broken tips
  - 3.13 Stopping the medical device
- 4 Interactions, contraindications, prohibitions**
  - 4.1 Interferences with other medical devices
  - 4.2 Using accessories not supplied by
  - 4.3 Prohibited uses
  - 4.4 Irrigation systems of the medical devices
  - 4.5 Assembly and disassembly
- 5 Electromagnetic compatibility**
  - 5.1 Cable length
  - 5.2 Recommended separation distances
  - 5.3 Electromagnetic emissions
  - 5.4 Electromagnetic immunity
  - 5.5 Electromagnetic immunity, handheld radiofrequency equipment
- 6 Cleaning, disinfecting and sterilizing**
  - 6.1 Disposal and recycling
  - 6.2 Disposing of accessories
- 7 Index**

## Foreword

The Ultradent Products Inc. medical device that you are about to install and use in your practice is a medical device designed for professional use. It comprises the chosen tool with which you will provide treatment within the context of your work.

To ensure optimum safety for yourself and your patients, comfort in your daily practice and to benefit fully from the technology of your medical device, please read the documentation provided carefully.

If you have received this medical device by mistake, please contact the supplier to arrange for it to be collected.

Please refer to the accessory cleaning, disinfection and sterilization protocols and the handpiece cleaning, disinfection and sterilization protocols for information about the following:

- preparation of parts for sterilization;
- detailed manual and automatic protocols;
- information concerning the sterilization;
- recommendations for the inspection of parts.

Please refer to the User Manuals, Quick Start Guides and Quick Clean Guides for each medical device for information about the following:

- unpacking and installing the medical device;
- using the medical device;
- monitoring and maintaining the medical device;
- technical specifications for the medical device.

# 1 Documentation

This document contains the following information:

- relating to patient, practitioner and environmental safety;
- required to install your medical device in optimum conditions;
- required to contact the manufacturer or the latter's representatives if necessary;

## 1.1 Electronic documentation

The user instructions for your device are in part provided in electronic format and in part in printed format. However, you can request a free printed copy of the user instructions via the website, by telephone or in writing.

The electronic user instructions are available in PDF format (Portable Document Format) and you will need to have a PDF file reader software installed to read the instructions.

The device user instructions can be consulted at the following addresses: [www.ultradent.com](http://www.ultradent.com) and [www.satelec.com](http://www.satelec.com).

Electronic user  
informations

It is important for you to have read and understood the content of the user instructions relating to the use of your device and its accessories prior to use. We recommend that you visit the website regularly to consult and/or to download the latest version of your device's user instructions.

## 1.2 Archiving duration

Users are asked to keep documentation on hand so that it can be consulted when necessary. You are asked to print and/or to download all documents or sections of documents that you may need to consult in the event of an emergency, if you are unable to connect to the internet or if your electronic display tool is not working (computer, tablet, etc.). All printed and electronic format documentation relating to your medical device must be kept for the device's entire service life. When loaning out or selling the medical device, the documentation must be provided with it.

# 2 Required information

## 2.1 Latest document update

12/2013

## 2.2 Official Texts

This medical device complies with the essential requirements of European Directive 93/42/EEC. This equipment is designed and developed in compliance with Electrical Safety standard IEC60601-1 in force. It was designed and manufactured in accordance with an EN ISO 13485-certified quality assurance system.

## 2.3 Standardized Symbols

| Symbols                                                                                                           | Meaning                                         |
|-------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|
|                                  | Refer to the accompanying documentation         |
|                                  | Consult the User Manual                         |
| <br>Electronic user informations | Accompanying documentation in electronic format |
|                                  | BF type                                         |
|                                  | Alternating voltage                             |
|                                  | Sterilization at 132 °C in an autoclave         |
|                                  | Washer disinfector for thermal disinfection     |
|                                  | EC marking                                      |
|                                  | Do not dispose of as household waste            |

| Symbols                                                                               | Meaning                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| YYYY  | Year of manufacture                                                                                                                                                                       |
|      | Control pedal                                                                                                                                                                             |
| O                                                                                     | Device OFF                                                                                                                                                                                |
| I                                                                                     | Device ON                                                                                                                                                                                 |
| IPX1                                                                                  | IP: ingress protection ratings procured by a range<br>X: no ingress protection rating claim against the penetration of solids<br>1: protects against the vertical falls of drops of water |
| Rx Only                                                                               | Caution: US Federal law restricts this device to sale by or on the order of a Physician                                                                                                   |

## 2.4 Manufacturer identification

SATELEC

A Company of ACTEON Group

17, avenue Gustave Eiffel

BP 30216

33708 MERIGNAC cedex

FRANCE

Tel. +33 (0) 556.34.06.07

Fax. +33 (0) 556.34.92.92

E.mail: [satelec@acteongroup.com](mailto:satelec@acteongroup.com).

[www.acteongroup.com](http://www.acteongroup.com)

## 2.5 Manufacturer responsibility

The manufacturer shall under no circumstances be liable for:

- non-compliance with manufacturer recommendations during installation, whether this is the network voltage or the electromagnetic environment,
- maintenance or repair procedures performed by people who are unauthorized by the manufacturer,
- use on an electrical fixture that is not compliant with regulations in force,
- uses other than those specified in this manual,
- use of accessories, tips or handpieces other than those supplied by SATELEC® or Ultradent Products Inc.,
- non-compliance with the instructions contained in this document.

Note: the manufacturer reserves the right to modify the medical device and/or any documentation without notice.

## 2.6 Distributor

Ultradent Products, Inc.

505 West 10200 South

South Jordan, UT 84095 USA

E-mail : [customer.support@ultradent.com](mailto:customer.support@ultradent.com)

[www.ultradent.com](http://www.ultradent.com)

# 3 Warnings

## 3.1 Federal Law

• The indication below applies to the United States of America only.  
The United States Federal Law restricts the use of this medical device in its territory to qualified dental health professionals, fit and certified to perform and manage their professional duties.

## 3.2 Warning applicable to all countries in which the device is sold

The information below is based on the requirements of standards to which the manufacturers of medical devices must adhere (as stated in standard IEC62366).

## 3.3 User population

This medical device must only be used by qualified dental health practitioners, fit and certified to perform their professional duties.

Users must know and comply with the rules of dental practice in compliance with knowledge acquired in the field and the key medical hygiene principles including cleaning, disinfection and sterilization of medical devices.

The medical device can be used by any adult dental practitioner of any weight, age, height, gender and nationality.

The user must wear gloves.

The user is not the patient.

The user must not be prone to any of the following:

- visual impairments: any vision problems must be corrected by glasses or lenses.
- disability of arms that may prevent the user from holding a handpiece, or legs that may prevent use of a control pedal;
- hearing difficulties that could prevent the user hearing audible alarms depending on medical devices;
- difficulty memorizing or concentrating that could affect the setting of sequences or the performance of treatment protocols.

## 3.4 Specific user training

No specific training other than initial professional training is required to use this medical device.

The practitioner is responsible for performing clinical treatments and for dangers that may arise due to a lack of skill and/or training.

## 3.5 Patient population restriction

This medical device must not be used on the following patient population:

- pregnant or breastfeeding women due to restrictions associated with the possible use of medical solutions such as anesthetics;
- patients with medical issues;
- patients with allergies;
- patients with a clinical site not suitable for treatment.

The patient must be calm, relaxed, still, ideally lying flat on a dental chair.

The user is the only person who can decide whether or not to treat his/her patients.

## 3.6 Parts of the body or types of tissues treated

Treatments must only be performed on the patient's oral environment.

## 3.7 Elements applied

|                                               |                             |
|-----------------------------------------------|-----------------------------|
| Elements in direct contact with the patient   | Tip                         |
| Elements in indirect contact with the patient | Handpiece nose<br>Handpiece |

## 3.8 Basic safety in normal use

The active part, the handpiece is held by the practitioner throughout the treatment. Being medically qualified, the practitioner is qualified to immediately detect any problem at the treatment site and to react accordingly. The force applied to the handpiece fitted with its tip must be managed by the practitioner in compliance with good dental practices. Basic safety is ensured by the practitioner. It is advisable to have a spare medical device or an alternative means with which to perform the medical treatment in the event of device failure.

## 3.9 Normal usage conditions

The normal usage conditions are as follows:

- Storage;
- Installation;
- Use;
- Maintenance;
- Disposal.

## 3.10 Irrigation spray

The irrigation spray is required to cool down or rinse the treatment area. Tips must be used with the irrigation function in accordance with manufacturer recommendations.

## 3.11 Tip service life

The shape and weight of the tips are determined to obtain maximum performance of the ultrasonic generator. The medical device will perform best if the user pays attention to these two characteristics. Therefore, we strongly advise against the modification of the structure of the tips by filing, twisting or by performing other types of modification. In addition, tip ageing leading to normal wear can modify the tip characteristics. Systematically replace a tip that is damaged due to wear or accidental impact (drop, distortion, etc.). As it is not possible to establish a maximum number of uses (that may be determined by many parameters such as duration of use, hardness of enamel, the force applied, wear, etc.), we recommend that tips being used are replaced at least once a year.

## 3.12 Broken tips

A tip is a medical device to which a mechanical force is applied to perform dental treatments. The mechanical force of the user is combined with the ultrasonic force generated by the dental handpiece. The tips have been developed to ensure safe use in association with handpieces in accordance with the power levels defined. However, the tips may break depending on frequency of use, force exerted or by being dropped. To reduce all risk, however minimal, we recommend the use of a suction device such as a saliva suction cannula. You should also encourage your patient to breathe through their nose.

## 3.13 Stopping the medical device

For devices connected to the water network, shut off the dental practice's water inlet. In the event of non-use, storage or prolonged absence, unplug the power plug from the electrical supply system. Before unplugging the power cord, switch the medical device to stop position O. To unplug the power cord, take hold of it and hold the wall socket.

# 4 Interactions, contraindications, prohibitions

This includes information relating to the interactions, contraindications and prohibitions known by on the date on which this document was written.

## 4.1 Interferences with other medical devices

Interferences may occur when the system is used on patients fitted with a pacemaker.

The medical device presents potential risks due to the emission of electromagnetic fields.

It may in particular cause malfunction of implanted devices such as a pacemaker or implantable defibrillator and, generally, of any type of active implant:

- before using this product, check whether patients and practitioners are fitted with a device of this type;
- explain the situation;
- weigh up the benefits versus the risks and contact your patient's cardiologist or another qualified health professional prior to starting treatment;
- keep this system away from implantable devices;
- apply suitable emergency measures and act fast if the patient shows signs of being unwell.
  - | Symptoms such as an increased heart beat, irregular pulse or dizziness may indicate a malfunction of a pacemaker or an implantable defibrillator.

The medical device is not designed to withstand electrical defibrillation shocks.

## 4.2 Using accessories not supplied by

The medical device was designed and developed with its accessories to guarantee maximum safety and performance. The use of accessories from another source could put you and your patients at risk and could damage your medical device. Do not try to connect accessories not provided by to your medical device connector(s) or to the handpiece. Even if the manufacturer or dealer of your accessory claims full compatibility with equipment, it is advisable to exercise caution with regards to the origin and safety of the product offered. Look out in particular for lack of information, information in a foreign language, very attractive prices, suspect appearance, mediocre quality or premature wear. If necessary, contact an approved dealer or the customer services department.

## 4.3 Prohibited uses

- Never cover the medical device and/or obstruct the air inlets.
- Do not immerse or use outside.
- Do not place the medical device next to a source of heat or in direct sunlight.
- Do not expose the medical device to water spray or mist.

The medical device is not designed to operate near a source of ionizing radiation.

A hot/cold temperature contrast can cause condensation to form in the medical device, which may be dangerous. If the medical device needs to be moved from a cold place to a warm place, do not use the device immediately but wait until it reaches room temperature.

The medical device may not be stored or used outside the atmospheric pressure and temperature ranges recommended in the User Manual supplied with your medical device.

Do not touch accessible electrical connections.

## 4.4 Irrigation systems of the medical devices

| This concerns medical devices that use tanks, bottles or bags.

- Never fill an irrigation solution tank or bottle when they are installed on the medical device.
- Always remove the irrigation tanks or bottles from the medical device before filling them.
- Never fill an irrigation solution tank or bottle over the maximum level mark.
- Irrigation tank covers must always be closed when the medical device is in use.
- Use the medical device with bottles or bags that do not exceed the maximum levels recommended for the medical device (depending on option).
- To prevent any interaction such as crystallization or precipitate between the different irrigation solutions used, rinse through the irrigation circuit and clean the tank in accordance with the instructions provided in the medical device User Manual.

## 4.5 Assembly and disassembly

Unless otherwise indicated in the instructions specific to your medical device:

- Control devices are not designed to be removed or disassembled.
- Access doors and/or flaps are not designed to be removed or disassembled.

# 5 Electromagnetic compatibility

| All the information below is based on the requirements of standards to which the manufacturers of electrical medical devices must adhere (as stated in standard IEC60601-1-2).

The medical device complies with the electromagnetic compatibility standards in force. However, the user will make sure that any electromagnetic interference does not create an additional risk, such as radiofrequency transmitters, or other electronic devices.

This chapter contains the information required for you to install and use your medical device in optimum conditions in terms of electromagnetic compatibility.

The different medical device cords must be kept away from each other.

Some types of mobile telecommunication devices such as mobile phones may interfere with the medical device. The separation distances recommended in this chapter **MUST** be complied with.

The medical device must not be used near another device or placed on top of it. If this cannot be avoided, correct operation of the device in operating conditions must be checked prior to use.

The use of accessories other than those specified or sold by SATELEC® or Ultradent Products Inc. as replacement parts, may increase the transmission or reduce the immunity of the medical device.

## 5.1 Cable length

| Cables and accessories | Maximum length | Test type                                                             | In compliance with: |
|------------------------|----------------|-----------------------------------------------------------------------|---------------------|
| Cables/Cords           | < 3m           | RF emission                                                           | CISPR 11, Class B   |
|                        |                | Harmonic current emission                                             | IEC61000-3-2        |
|                        |                | Voltage fluctuation and flickers                                      | IEC61000-3-3        |
|                        |                | Electrostatic discharge immunity                                      | IEC61000-4-2        |
|                        |                | Radiated immunity - Electromagnetic fields                            | IEC61000-4-3        |
|                        |                | Electrical fast transient/burst immunity                              | IEC61000-4-4        |
|                        |                | Surge immunity                                                        | IEC61000-4-5        |
|                        |                | Immunity to conducted disturbances, induced by radio-frequency fields | IEC61000-4-6        |
|                        |                | Radiated immunity - Magnetic fields                                   | IEC61000-4-8        |
|                        |                | Voltage dips, short interruptions and voltage variation immunity      | IEC 61000-4-11      |

## 5.2 Recommended separation distances

The medical device is designed to be used in an electromagnetic environment in which interferences caused by RF radiation are controlled.

The user or installer of the medical device may help to prevent electromagnetic interference by maintaining a minimum distance, depending on the maximum power of the handheld and mobile radiofrequency transmission equipment (transmitters), between the medical device and the equipment as recommended in the table below.

|                                                                                                                                                                                                                                                                                                                  | Separation distance in metres (m) determined by transmitter frequency |                                                    |                                                     |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|----------------------------------------------------|-----------------------------------------------------|
| Max. nominal power of the transmitter in Watts                                                                                                                                                                                                                                                                   | From 150 kHz to 80 Mhz<br>$d = \sqrt{P} \cdot 1.2$                    | From 80 MHz to 800 MHz<br>$d = \sqrt{P} \cdot 1.2$ | From 800 MHz to 2.5 GHz<br>$d = \sqrt{P} \cdot 2.3$ |
| 0.01                                                                                                                                                                                                                                                                                                             | 0.12 m                                                                | 0.12 m                                             | 0.23 m                                              |
| 0.1                                                                                                                                                                                                                                                                                                              | 0.38 m                                                                | 0.38 m                                             | 0.73 m                                              |
| 1                                                                                                                                                                                                                                                                                                                | 1.2 m                                                                 | 1.2 m                                              | 2.3 m                                               |
| 10                                                                                                                                                                                                                                                                                                               | 3.8 m                                                                 | 3.8 m                                              | 7.3 m                                               |
| 100                                                                                                                                                                                                                                                                                                              | 12 m                                                                  | 12 m                                               | 23 m                                                |
| With regards to transmitters for which the maximum power is not listed above, the recommended separation distance (d) in metres (m) can be estimated by using the equation applicable to the transmitter frequency where (P) is the maximum power of the transmitter in watts (W) according to the manufacturer. |                                                                       |                                                    |                                                     |

Note 1: At 80 MHz and 800 MHz, the highest frequency range applies.

Note 2: These specifications cannot be applied to all situations.

The electromagnetic propagation is reduced by the absorption and reflection of structures, objects and people.

## 5.3 Electromagnetic emissions

The medical device is designed for use in the electromagnetic environment described in the table below. The user and/or installer must ensure that the medical device is used in the environment described below.

| Emission test                                   | Conformity | Electromagnetic environment - comments                                                                                                                                                                        |
|-------------------------------------------------|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emission (CISPR 11)                          | Group 1    | The medical device uses RF energy for its internal operation. Consequently, its radiofrequency emissions are very low and are not likely to create any interference with other nearby equipment.              |
| RF emission (CISPR 11)                          | Class B    | The medical device is suitable for use in all establishments, including domestic and those directly connected to the low voltage energy supply public network supplying buildings used for domestic purposes. |
| Harmonic current emission (IEC61000-3-2)        | Class A    |                                                                                                                                                                                                               |
| Voltage fluctuation and flickers (IEC61000-3-3) | Conforming |                                                                                                                                                                                                               |

## 5.4 Electromagnetic immunity

The medical device is designed for use in the magnetic and electromagnetic environment described in the table below. The user and/or installer must ensure conformity of the electromagnetic environment.

| Immunity test                                                            | Test level in accordance with IEC60601                                                                                                                                                              | Conformity level                                                                                                                                                                                    | Electromagnetic environment / comments                                                                                                                                                                                                                                                 |
|--------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD) (IEC61000-4-2)                             | ± 6 kV in contact<br>± 8 kV in the air                                                                                                                                                              | ± 6 kV in contact<br>± 8 kV in the air                                                                                                                                                              | Floors must be wood, concrete, cement or tiled. If floors are covered with synthetic materials (carpet, etc.), the relative humidity must be 30% minimum.                                                                                                                              |
| Electrical fast transient (IEC61000-44)                                  | ± 2 kV for electricity supply lines                                                                                                                                                                 | ± 2 kV for electricity supply lines                                                                                                                                                                 | The quality of the electricity supply must be equivalent to that of a typical commercial environment or hospital establishment (hospital, clinic).                                                                                                                                     |
| Surge (IEC61000-4-5)                                                     | ± 1 kV in differential mode<br>± 2 kV in common mode                                                                                                                                                | ± 1 kV in differential mode<br>± 2 kV in common mode                                                                                                                                                | The quality of the electricity supply must be equivalent to that of a typical commercial environment or a hospital.                                                                                                                                                                    |
| Magnetic field at 50 Hz/60 Hz (IEC61000-4-8)                             | 3A/m                                                                                                                                                                                                | 3A/m                                                                                                                                                                                                | The magnetic field intensity must be equal to the level found in a typical commercial or hospital environment.                                                                                                                                                                         |
| Voltage dips, short interruptions and voltage variation (IEC 61000-4-11) | <5% UT<br>(>95% dip in UT)<br>for 0.5 cycles<br><br>40% UT<br>(60% dip in UT)<br>for 5 cycles<br><br>70% UT<br>(30% dip in UT)<br>for 25 cycles<br><br><5% UT<br>(>95% dip in UT)<br>for 250 cycles | <5% UT<br>(>95% dip in UT)<br>for 0.5 cycles<br><br>40% UT<br>(60% dip in UT)<br>for 5 cycles<br><br>70% UT<br>(30% dip in UT)<br>for 25 cycles<br><br><5% UT<br>(>95% dip in UT)<br>for 250 cycles | The quality of the network supply must be equal to that of a typical commercial or hospital environment. If the use of the system requires continuous operation during mains power outages, it is advisable to supply the medical device using an separate current source (UPS, etc.). |

## 5.5 Electromagnetic immunity, handheld radiofrequency equipment

The medical device is designed for use in the magnetic and electromagnetic environment described in the table below. The user and/or installer must ensure conformity of the electromagnetic environment.

| Immunity test                                                                                                                                                                                                                                                                                                                                                                      | Test level                 | Conformity level | Electromagnetic environment - comments                                                                                                                                                                                                                                         |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Handheld and mobile radiofrequency communication devices must not be used near the device (including cables) at a distance below that recommended and calculated according to the frequency and power of the transmitter.                                                                                                                                                          |                            |                  |                                                                                                                                                                                                                                                                                |
| Radiated radiofrequency electromagnetic field (IEC61000-4-3)                                                                                                                                                                                                                                                                                                                       | 3 V/m<br>80 MHz to 2.5 GHz | 3 V/m            | $d = 1.2 / \sqrt{P}$ 80 MHz to 800 MHz<br>$d = 2.3 / \sqrt{P}$ 800 MHz to 2.5 GHz<br>Where (P) is the maximum nominal power of the transmitter in Watts (W) according to the manufacturer specifications and (d) is the minimum recommended separation distance in meters (m). |
| Radiofrequency conducted disturbance (IEC61000-4-6)                                                                                                                                                                                                                                                                                                                                | 3 V/m<br>150KHz to 80MHz   | 3 V/m            | Recommended separation distance: $d = 1.2 / \sqrt{P}$                                                                                                                                                                                                                          |
| The electromagnetic field intensity of fixed radiofrequency transmitters, as determined by an electromagnetic environment measurement (a), must be less than the conformity level for each frequency range (b). Interferences may occur near the equipment identified by the following symbol:<br> |                            |                  |                                                                                                                                                                                                                                                                                |

Note 1: At 80 MHz and 800 MHz, the highest frequency range applies.

Note 2: These specifications cannot be applied to all situations. The electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.

(a) The electromagnetic field intensity of fixed radiofrequency transmitters, such as base stations for portable phones (mobiles / wireless), mobile radios, radio amateurs, AM/FM radio transmissions and TV transmissions cannot be determined accurately by the theory.

To assess the electromagnetic environment caused by fixed radiofrequency transmitters, an electromagnetic environment measurement must be taken. If the measured intensity of the radiofrequency field in the immediate product use environment exceeds the radiofrequency conformity level specified above, it is necessary to test product performance to check this complies with specifications. If abnormal performance is observed, additional measures may be necessary, such as changing the direction of or moving the product. (b) In the 150 kHz to 80 Mhz frequency range, the electromagnetic fields must be less than 3 V/m.

## 6 Cleaning, disinfecting and sterilizing

The instructions relating to the accessory cleaning, disinfection and sterilization protocols provided by Ultradent Products Inc. have been approved for each medical device and accessory. They can be downloaded at the following addresses [www.satelec.com/documents](http://www.satelec.com/documents) or [www.ultradent.com](http://www.ultradent.com). In all cases, the local regulations in force relating to the accessory cleaning, disinfection and sterilization protocols take precedence over the information provided by Ultradent Products Inc.

## 6.1 Disposal and recycling

As an item of Electrical and Electronic Equipment, the medical device must be disposed of via a specialist collection, removal, recycling or destruction channel. This applies in particular to the European market, in reference to Directive 2002/96/EC dated 27/01/2003.

When your medical device has reached the end of its service life, contact Ultradent Products Inc. to find out how to proceed. Relevant contact details are given in the chapter *Distributor page 5*

| The indication below applies to France only.

In compliance with the provisions of the French Environment Code relating to the disposal of electronic and electrical equipment waste or WEEE (Decree no. 2012-617 dated 2 May 2012), our Company fulfils its obligations to reclaim and dispose of its electrical and electronic equipment through the means established by the approved organisation Recydent (NOR approval: DEVP1229534A).

As a producer, our Company is listed in the National Register of Producers kept by the ADEME (French Environment and Energy Management Agency). Professionals buying our products directly from the distribution chain are responsible for passing on this information about our established recycling methods to the end user.

In addition, the buyer agrees to take back our brand's devices at the end of their service life and to transfer them to one of the collection centres set up by Recydent for recycling (see list of collection centres on the site [Recydent.fr](http://Recydent.fr)). If necessary, Recydent can come and collect these devices from you free of charge once the quantity of devices has reached a certain level in the pallets-containers with which you are provided to store this waste.

## 6.2 Disposing of accessories

An accessory that has reached the end of its service life must be disposed of in infectious clinical waste containers.

# 7 Index

## D

disposal 13

## E

Electrical Safety 3

electronic user instructions 3

European Directive 3

## R

Regulatory 3, 13

# J81009

## Manual Pre-disinfection/Cleaning and Sterilization Instructions for SATELEC wrenches

### Caution

Do not use steel wool or abrasive cleaners.  
Avoid solutions containing iodine or with a high chlorine content.  
The pH of detergents/disinfectants must be between 7 and 11.

The cleaning method for wrenches recommended by Satelec is manual or automated.  
All devices must be carefully cleaned and must be sterilized after cleaning prior to use.  
The sterilization parameters are only valid for correctly cleaned devices.  
Wrenches must be cleaned with extreme care.  
It is the responsibility of the end user to ensure that all equipment used to treat Satelec devices is properly installed, validated, maintained and calibrated.  
Whenever possible, a washer/disinfector should be used for wrenches. Prevent the overloading of wash baskets for ultrasonic cleaning or cleaning in a washer/disinfector.

### Cleaning/sterilization cycle limits

Repeated treatment cycles that include ultrasonic cleaning, manual washing and sterilization have minimal effects on wrenches.  
End of service life is normally determined by wear and damage due to use.

### Important

- Soiled devices should be separated from non-contaminated devices to avoid contamination of personnel or surroundings.
- Wipe blood and/or debris from the wrenches to prevent it from drying onto the surfaces.
- Cover the wrench with a soft lint-free cloth dampened with purified water to prevent blood and/or debris from drying.

### Containment and transportation

Soiled devices must be transported separately from non-contaminated devices to avoid contamination.

### Preparation for pre-disinfection/cleaning

It is advisable to treat devices immediately or as soon as possible after use. Satelec devices must be treated within two hours of use.

### Pre-disinfection/Cleaning - Manual method

**Equipment:** soft-bristled brush, soft lint-free swab, lint-free cloth, enzymatic or alkaline cleaner, ultrasonic cleaner.

| I | Step duration (minimum) | Cleaning instructions                                                                                                                                                                                                   |
|---|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   | 1 minute                | Rinse soiled device under cold running water. Use a soft-bristled brush or a soft and lint-free clean cloth to remove most of the contamination and debris.                                                             |
|   | 10 minutes              | In an ultrasonic cleaner, immerse the device in a freshly prepared enzymatic or alkaline cleaning solution.                                                                                                             |
|   | 1 minute                | Rinse device under cold running water.                                                                                                                                                                                  |
|   | 2 minutes 30 seconds    | Manually wash device in a freshly prepared enzymatic or alkaline cleaning solution. Use a soft-bristled brush to remove soil and debris. Wash the device under running water to prevent aerosolization of contaminants. |
|   | 1 minute                | Thoroughly rinse device with distilled or purified water.                                                                                                                                                               |
|   |                         | Repeat the manual cleaning procedure until no visible contamination remains on device.                                                                                                                                  |
|   |                         | Perform a final rinse of device using distilled or purified water.                                                                                                                                                      |
|   |                         | Dry device using a soft, lint-free cloth or clean compressed air.                                                                                                                                                       |

## Inspection

- Wrenches must be inspected to check that no contamination remains, that they are not corroded, dulled, discoloured or damaged.
- Before packaging and sterilizing the cleaned products, check they are clean, undamaged and function properly.
- Damaged wrenches must be discarded. Wrenches must not be lubricated.

## Packaging

Use suitable packaging for sterilization, such as a Sterile Barrier System in compliance with ISO 11607. Care should be taken to protect wrenches from contact with other objects that may damage their surface or the Sterile Barrier System.

## Sterilization

Unless otherwise specified, non-sterile products can be re-sterilized using validated steam sterilization methods (ISO 17665 or national standards).

Satelec recommendations for packed wrenches are as follows:

| Cycle type                                        | Sterilization exposure time | Sterilization exposure temperature | Dry time           |
|---------------------------------------------------|-----------------------------|------------------------------------|--------------------|
| Saturated steam - forced air removal (pre-vacuum) | 4 minutes                   | 132 °C                             | Minimum 20 minutes |

Dry times generally range from 20 to 60 minutes due to the difference in packaging materials, steam quality, device materials, total mass, sterilizer performance, and varying cool-down time.

The distributor and manufacturer accept no responsibility for sterilization procedures performed by the customer that are not performed according to these Satelec recommendations.

## Storage

Storage conditions for products labelled “STERILE” are printed on the packaging label. Packaged products should be stored in a dry, clean environment, protected from direct sunlight, pests, and extremes of temperature and humidity. Use products in the order in which they are received (“first in, first out” principle), taking into account the expiry date indicated on the label.

# Automated Pre-disinfection/Cleaning and Sterilization Instructions for all SATELEC Wrenches

## Cautions

Do not use steel wool or abrasive cleaners.  
Avoid solutions containing iodine or with a high chlorine content.  
The pH of detergents/disinfectants must be between 7 and 11.

The cleaning method for wrenches recommended by Satelec is manual or automated.  
All devices must be carefully cleaned and must be sterilized after cleaning prior to use.  
The sterilization parameters are only valid for correctly cleaned devices.  
Wrenches must be cleaned with extreme care.

It is the responsibility of the end user to ensure that all equipment used to treat Satelec devices is properly installed, validated, maintained and calibrated.  
Whenever possible, a washer/disinfector should be used for wrenches. Prevent the overloading of wash baskets for ultrasonic cleaning or cleaning in a washer/disinfector.

## Cleaning/sterilization cycle limits

Repeated treatment cycles that include ultrasonic cleaning, manual washing and sterilization have minimal effects on wrenches.  
End of service life is normally determined by wear and damage due to use.

## Important

- Soiled devices should be separated from non-contaminated devices to avoid contamination of personnel or surroundings.
- Wipe blood and/or debris from the wrenches to prevent it from drying onto the surfaces.
- Cover the wrench with a soft lint-free cloth dampened with purified water to prevent blood and/or debris from drying.

## Containment and transportation

Soiled devices must be transported separately from non-contaminated devices to avoid contamination.

## Preparation for pre-disinfection/cleaning

It is advisable to treat devices immediately or as soon as possible after use. Satelec devices must be treated within two hours of use.

## Pre-disinfection/pre-cleaning method

Note: The pre-disinfection/pre-cleaning method has to be performed for the wrench prior to automated cleaning.  
Equipment: soft-bristled brush, soft lint-free swab, lint-free cloth, ultrasonic cleaner, washer/disinfector, enzymatic or alkaline cleaner.

| Step duration (minimum) | Cleaning instructions                                                                                                                                                                                                   |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 30 seconds              | Rinse soiled device under cold running water. Use a soft-bristled brush or a soft and lint-free clean cloth to remove most of the contamination and debris.                                                             |
| 1 minute                | Manually wash device in a freshly prepared enzymatic or alkaline cleaning solution. Use a soft-bristled brush to remove soil and debris. Wash the device under running water to prevent aerosolization of contaminants. |
| 30 seconds              | Rinse device under cold running water.                                                                                                                                                                                  |
| 8 minutes               | In an ultrasonic cleaner, immerse the device in a freshly prepared enzymatic or alkaline cleaning solution.                                                                                                             |
| 1 minute                | Thoroughly rinse device with distilled or purified water.                                                                                                                                                               |
|                         | Repeat the pre-cleaning procedure until no visible contamination remains on device.                                                                                                                                     |

## Automated cleaning method

| Step           | Duration (minimum) | Cleaning/decontamination instructions                                   |
|----------------|--------------------|-------------------------------------------------------------------------|
| Pre-washing    | 2 minutes          | Cold tap water                                                          |
| Washing        | 2 minutes          | Warm tap water (> 40°C); use an alkaline or enzymatic cleaning solution |
| Neutralization | 5 minutes          | Warm tap water with neutralizer, if necessary                           |
| Rinsing        | 2 minutes          | Rinse with warm deionised or purified water (> 40°C)                    |
| Drying         | 40 minutes         | 90°C                                                                    |

## Thermal disinfection

For automated cleaning: thermal disinfection at 90°C for a minimum of 5 minutes.

## Inspection

- Wrenches must be inspected to check that no contamination remains, that they are not corroded, dulled, discoloured or damaged.
- Before packaging and sterilizing the cleaned products, check they are clean, undamaged and function properly.
- Damaged wrenches must be discarded. Wrenches must not be lubricated.

## Packaging

Use suitable packaging for sterilization, such as a Sterile Barrier System in compliance with ISO 11607. Care should be taken to protect wrenches from contact with other objects that may damage their surface or the Sterile Barrier System.

## Sterilization

Unless otherwise specified, non-sterile products can be re-sterilized using validated steam sterilization methods (ISO 17665 or national standards).

Satelec recommendations for packed wrenches are as follows:

| Cycle type                                        | Sterilization exposure time | Sterilization exposure temperature | Drying time        |
|---------------------------------------------------|-----------------------------|------------------------------------|--------------------|
| Saturated steam - forced air removal (pre-vacuum) | 4 minutes                   | 132°C                              | Minimum 20 minutes |

Dry times generally range from 20 to 60 minutes due to the difference in packaging materials, steam quality, device materials, total mass, sterilizer performance, and varying cool-down time.

The distributor and manufacturer accept no responsibility for sterilization procedures performed by the customer that are not performed according to these Satelec recommendations.

## Storage

Storage conditions for products labelled “STERILE” are printed on the packaging label. Packaged products should be stored in a dry, clean environment, protected from direct sunlight, pests, and extremes of temperature and humidity. Use products in the order in which they are received (“first in, first out” principle), taking into account the expiry date indicated on the label.

# J02009

# Manual Pre-disinfection/Cleaning and Sterilization Instructions for SATELEC Tips and Files

## Caution

Do not use steel wool or abrasive cleaners.  
Avoid solutions containing iodine or with a high chlorine content.  
The pH of detergents/disinfectants must be between 7 and 11.

The cleaning method for tips and files recommended by Satelec is manual or automated.  
All devices must be carefully cleaned and must be sterilized after cleaning prior to use.  
The sterilization parameters are only valid for correctly cleaned devices.

Tips and files must be cleaned with extreme care.

During automated cleaning, tips or files must be placed on suitable instrument holders or in small baskets to prevent them from damage during washing.

It is the responsibility of the end user to ensure that all equipment used to treat Satelec devices is properly installed, validated, maintained and calibrated.

Whenever possible, a washer/disinfector should be used for tips. Prevent the overloading of wash baskets for ultrasonic cleaning or cleaning in a washer/disinfector.

## Cleaning/sterilization cycle limits

Repeated treatment cycles that include ultrasonic cleaning, manual washing and sterilization have minimal effects on tips and files.

End of service life is normally determined by wear and damage due to use.

## Important

If the irrigation tubing is external or if configuring the device with tank, press the cleaning button on the device after each use and

- o Remove the perforator from the irrigation solution, and then immerse it in an enzymatic or alka line cleaning solution. By pressing the flush button on the device, rinse the irrigation tubing, the handpiece and the tip (or the file) for 1 minute to prevent the drying of soil and debris inside the device.
- o Immerse the perforator in a distilled or purified water solution, rinse for 1 minute in order to eliminate the enzymatic or alkaline cleaning solution.
- o Remove the perforator from the distilled or purified water solution and purge the irrigation tubing to remove the liquid from the inside.
- If the case of irrigation coming from the dental practice water network, clean and rinse the cord, the handpiece and the tip (or file) for 1 minute after each use.
- Soiled devices should be separated from non-contaminated devices to avoid contamination of personnel or surroundings.
- Wipe blood and/or debris from the tips or files to prevent it from drying onto the surfaces. Cover the tip or the file with a soft lint-free wipe dampened with purified water to prevent blood and/or debris from drying.

## Containment and transportation

Soiled devices must be transported separately from non-contaminated devices to avoid contamination.

## Preparation for pre-disinfection/cleaning

It is advisable to treat devices immediately or as soon as possible after use. Satelec devices must be treated within two hours of use.

Before cleaning, unscrew the tip or the file at the front of the ultrasonic handpiece after each use.

## Pre-disinfection/Cleaning - Manual method

**Equipment:** soft-bristled brush, soft and lint-free swab, lint-free cloth, syringe, pipette and/or water jet, enzymatic or alkaline cleaning solution, ultrasonic cleaner.

| Step duration (minimum)                                                   | Cleaning instructions                                                                                                                                                                                                                                                             |
|---------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 minute                                                                  | Rinse soiled device under cold running water. Use a soft-bristled brush or a soft and lint-free clean cloth to remove most of the contamination and debris. Use a syringe, a pipette or a water jet with an enzymatic or alkaline cleaning solution to rinse the tip cannulation. |
| 10 minutes                                                                | In an ultrasonic cleaner, immerse the device in a freshly prepared enzymatic or alkaline cleaning solution.                                                                                                                                                                       |
| 1 minute                                                                  | Rinse device under cold running water. Use a syringe, a pipette or a water jet to rinse the tip cannulation.                                                                                                                                                                      |
| 4 minutes                                                                 | Manually wash device in a freshly prepared enzymatic or alkaline cleaning solution. Use a soft-bristled brush to remove soil and debris. Wash device under water to prevent aerosolization of contaminants.                                                                       |
| 1 minute                                                                  | Thoroughly rinse device with distilled or purified water. Use a syringe, a pipette or a water jet to rinse the tip cannulation. Repeat rinsing.                                                                                                                                   |
|                                                                           | Repeat the pre-cleaning procedure until no visible contamination remains on device.                                                                                                                                                                                               |
| Perform a final rinse of device using distilled or purified water.        |                                                                                                                                                                                                                                                                                   |
| Dry the device with a soft, lint-free cloth or with clean compressed air. |                                                                                                                                                                                                                                                                                   |

## Inspection

- Devices must be inspected to check that no contamination remains, that they are not corroded, dulled, discoloured or damaged.
- Before packaging and sterilizing the cleaned devices, check they are clean, undamaged and function properly.
- Damaged devices must be discarded. Tips must not be lubricated.

## Packaging

Use suitable packaging for sterilization, such as a Sterile Barrier System in compliance with ISO 11607. Care should be taken to protect tips (or files) from contact with other objects that may damage their surface or the Sterile Barrier System.

## Sterilization

Unless otherwise specified, non-sterile products can be re-sterilized using validated steam sterilization methods (ISO 17665 or national standards).

Satelec recommendations for packed tips are as follows:

| Cycle type                                        | Sterilization exposure time | Sterilization exposure temperature | Drying time        |
|---------------------------------------------------|-----------------------------|------------------------------------|--------------------|
| Saturated steam - forced air removal (pre-vacuum) | 4 minutes                   | 132 °C                             | Minimum 20 minutes |

Dry times generally range from 20 to 60 minutes due to the difference in packaging materials, steam quality, device materials, total mass, sterilizer performance, and varying cool-down time.  
The distributor and manufacturer accept no responsibility for sterilization procedures performed by the customer that are not performed according to these Satelec recommendations.

## Storage

Storage conditions for products labelled “STERILE” are printed on the packaging label. Packaged products should be stored in a dry, clean environment, protected from direct sunlight, pests, and extremes of temperature and humidity. Use products in the order in which they are received (“first in, first out” principle), taking into account the expiry date indicated on the label.

# Automated Pre-disinfection/Cleaning and Sterilization Instructions for SATELEC Tips and Files

## Caution

Do not use steel wool or abrasive cleaners.  
Avoid solutions containing iodine or with a high chlorine content.  
The pH of detergents/disinfectants must be between 7 and 11.

The cleaning method for tips and files recommended by Satelec is manual or automated.

All devices must be carefully cleaned and must be sterilized after cleaning prior to use.

The sterilization parameters are only valid for correctly cleaned devices.

Tips and files must be cleaned with extreme care.

During automated cleaning, tips or files must be placed on suitable instrument holders or in small baskets to prevent them from damage during washing.

It is the responsibility of the end user to ensure that all equipment used to treat Satelec devices is properly installed, validated, maintained and calibrated.

Whenever possible, a washer/disinfector should be used for tips. Prevent the overloading of wash baskets for ultrasonic cleaning or cleaning in a washer/disinfector.

## Cleaning/sterilization cycle limits

Repeated treatment cycles that include ultrasonic cleaning, manual washing and sterilization have minimal effects on tips and files.

End of service life is normally determined by wear and damage due to use.

## Important

- If the irrigation tubing is external or if configuring the device with tank, press the cleaning button on the device after each use and
  - o Remove the perforator from the irrigation solution, and then immerse it in an enzymatic or alkaline cleaning solution. By pressing the flush button on the device, rinse the irrigation tubing, the handpiece and the tip (or the file) for 1 minute to prevent the drying of soil and debris inside the device.
  - o Immerse the perforator in a distilled or purified water solution, rinse for 1 minute in order to eliminate the enzymatic or alkaline cleaning solution.
  - o Remove the perforator from the distilled or purified water solution and purge the irrigation tubing to remove the liquid from the inside.

- If the case of irrigation coming from the dental practice water network, clean and rinse the cord, the handpiece and the tip (or file) for 1 minute after each use.
- Soiled devices should be separated from non-contaminated devices to avoid contamination of personnel or surroundings.
- Wipe blood and/or debris from the tips or files to prevent it from drying onto the surfaces.
- Cover the tip or the file with a soft lint-free wipe dampened with purified water to prevent blood and/or debris from drying.

## Containment and transportation

Soiled devices must be transported separately from non-contaminated devices to avoid contamination.

## Preparation for pre-disinfection/cleaning

It is advisable to treat devices immediately or as soon as possible after use. Satelec devices must be treated within two hours of use.

Before cleaning, unscrew the tip or the file at the front of the ultrasonic handpiece after each use.

## Pre-disinfection/Cleaning - Automated method

### Pre-disinfection/pre-cleaning method

Note: The manual pre-disinfection/pre-cleaning method has to be performed for the tips prior to automated cleaning.

**Equipment:** soft-bristled brush, soft and lint-free swab, lint-free cloth, syringe, pipette and/or water jet, ultrasonic cleaner, washer/disinfectant, enzymatic or alkaline cleaning solution.

| Step duration (minimum) | Cleaning instructions                                                                                                                                                                                                                                                             |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 minute                | Rinse soiled device under cold running water. Use a soft-bristled brush or a soft and lint-free clean cloth to remove most of the contamination and debris. Use a syringe, a pipette or a water jet with an enzymatic or alkaline cleaning solution to rinse the tip cannulation. |
| 1 minute 30 seconds     | Manually wash device in a freshly prepared enzymatic or alkaline cleaning solution. Use a soft-bristled brush to remove soil and debris. Wash device under water to prevent aerosolization of contaminants.                                                                       |
| 1 minute                | Rinse device under cold running water. Use a syringe, a pipette or a water jet to rinse the tip cannulation.                                                                                                                                                                      |
| 10 minutes              | In an ultrasonic cleaner, immerse the device in a freshly prepared enzymatic or alkaline cleaning solution.                                                                                                                                                                       |
| 1 minute                | Thoroughly rinse device with distilled or purified water. Use a syringe, a pipette or a water jet to rinse the tip cannulation. Repeat rinsing.                                                                                                                                   |
|                         | Repeat the pre-cleaning procedure until no visible contamination remains on device.                                                                                                                                                                                               |

## Automated cleaning method

| Step           | Duration   | Cleaning instructions                                                |
|----------------|------------|----------------------------------------------------------------------|
| Pre-washing    | 2 minutes  | Cold tap water                                                       |
| Washing        | 10 minutes | Warm tap water (> 40°C); use alkaline or enzymatic cleaning solution |
| Neutralization | 2 minutes  | Warm tap water with neutralizer, if necessary                        |
| Rinsing        | 2 minutes  | Rinse with warm distilled or purified water (>40° C)                 |
| Drying         | 40 minutes | 90° C                                                                |

## Thermal disinfection

For automated cleaning: thermal disinfection at 90° C for a minimum of 5 minutes.

## Inspection

- Tips must be inspected to check that no contamination remains, that they are not corroded, dulled, discoloured or damaged.
- Before packaging and sterilizing the cleaned products, check they are clean, undamaged and function properly.
- Damaged tips must be discarded. Tips must not be lubricated.

## Packaging

Use suitable packaging for sterilization, such as a Sterile Barrier System in compliance with ISO 11607. Care should be taken to protect tips (or files) from contact with other objects that may damage their surface or the Sterile Barrier System.

## Sterilization

Unless otherwise specified, non-sterile products can be re-sterilized using validated steam sterilization methods (ISO 17665 or national standards).

Satelec recommendations for packed tips are as follows:

| Cycle type                                        | Sterilization exposure time | Sterilization exposure temperature | Drying time        |
|---------------------------------------------------|-----------------------------|------------------------------------|--------------------|
| Saturated steam - forced air removal (pre-vacuum) | 4 minutes                   | 132° C                             | Minimum 20 minutes |

Dry times generally range from 20 to 60 minutes due to the difference in packaging materials, steam quality, device materials, total mass, sterilizer performance, and varying cool-down time.  
The distributor and manufacturer accept no responsibility for sterilization procedures performed by the customer that are not performed according to these Satelec recommendations.

## Storage

Storage Conditions for products labelled “STERILE” are printed on the packaging label. Packaged products should be stored in a dry, clean environment, protected from direct sunlight, pests, and extremes of temperature and humidity. Use products in the order in which they are received (“first in, first out” principle), taking into account the expiry date indicated on the label.

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