Opal™ Seal

Fluoride Releasing Primer and Sealant

Description and Indications for Use:

Opal™ Seal primer and sealant is a light cure primer that is used when bonding orthodontic appliances to etched enamel. Opal Seal primer and sealant is 38% filled with proprietary GI particles and nano fillers that release fluoride during orthodontic treatment. When Opal Seal primer and sealant has been applied to the tooth surface the use of a UV "black" light will make the sealant illuminate showing that the resin containing fluoride is present.

Direct Bonding Procedure

NOTE: Working one quadrant at a time is recommended.

- Prepare enamel surface
- Use oil-free, water-based pumice to clean enamel surface
- Rinse, dry and isolate the teeth to be etched.
- Securely attach functional tip to 35% phosphoric acid etchant gel.
- Verify etch flow prior to applying intra-orally.
- Apply etchant to enamel and leave in place 15–30 seconds.

CAUTION: Do not allow etch gel to contact soft tissue.

- Rinse and dry thoroughly and re-isolate. Enamel will appear chalky white.
- Apply Opal Seal primer and sealant
- Remove Luer cap and securely attach brush tip.
- Verify flow prior to applying intra-orally by slightly twisting pink stem and moistening the tip.

NOTE: If too much is expressed, dab excess on a pad.

- Apply a thin layer of Opal Seal primer and sealant to the etched enamel, moving quickly from tooth to tooth.
- Gently air dry for 2 seconds per tooth to thin. Sealant should appear shiny, not pooled on the tooth surface.

OPTIONAL: Before placing bracket, tack-cure Opal Seal primer and sealant for 5 seconds per tooth to avoid "skating" upon placement.

NOTE: Opal Seal primer and sealant is designed to be re-applied periodically during orthodontic treatment. Shining a UV

"black" light on the tooth will indicate the presence of the sealant.

Re-Application Procedure

- With orthodontic wire removed, ensure enamel surface is free of plaque and debris.
NOTE: For optimal treatment: clean enamel surface with an oil-free, water-based pumice.
- Follow Steps 1b-1 applying the 35% phosphoric acid etch to the exposed enamel and up to the base of the bracket for 15–30 seconds.

NOTE: Photocuring acid etchants will not cause de-bonding of any appliances.

- Follow Steps 2-3.

- Light cure 10 sec. per tooth (20 sec. for lights with output<600mW/cm²).

Removal Procedure

- At the conclusion of treatment the remaining sealant should be removed with a fluted bur or an abrasive eraser using caution to avoid overheating the tooth.
- The use of a UV "black" light will illuminate the sealant for ease of removal.

Precautions:

- Use as a personal use only.
- Carefully read and understand all instructions, precautions, and MSDS before using.
- Verify flow prior to applying intra-orally. If resistance is met, replace tip and re-check. Use only recommended tips.
- Dispose of tips and empty syringes properly. Tips are single use only.
- Apply a syringe with Luer Lock cap and disinfect with an intermediate level disinfectant between uses.
- Use quality curing lights to cure product. Uncured resin is deficient in its physical properties.
- Use quality curing lights to cure product. Uncured resin is deficient in its physical properties.
- Do not use if the sealant is past its expiration date.
- Not all bonding agents and adhesives are compatible. Test before using products outside the Ultradent or any specific manufacturer's product family.
5. Skin or soft tissue is exposed to product, thoroughly wash with soap and water. If allergic reactions or rash develops, consult a physician.

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