Consepsis™ & Consepsis™ V Chlorhexidine Antibacterial Solution

Description:

Consepsis™ is a lightly flavored 2.0% chlorhexidine gluconate disinfecting solution.

Note: Chlorhexidine hand soaps can adversely affect bond strength of dentin bonding agents. Consepsis contains no surfactants or emollients that interfere with bond strength.

Consepsis™ V is a lightly flavored 2.0% chlorhexidine gluconate disinfecting solution that is more viscous than regular Consepsis. Consepsis V is specifically designed for endodontic disinfecting purposes only.

Indications:

Consepsis liquid is used before crown cementation (temporary and/or permanent) and for restorative preparations of crowns, inlays, and composites. Consepsis is a demonstrated, quality wetting agent for bonding, and increasing bond strengths.^{1,2} Consepsis is recommended for procedural endodontic disinfection, as a final endodontic rinse prior to canal obturation (see EndoREZ™ instruction) and as an antimicrobial prior to pulp capping.³ Consepsis helps prevent the influx of micro-organisms into dentinal tubules.⁴ A clinician can reduce potential post-op pulpitis and sensitivity by thoroughly cleaning and disinfecting preparations with Consepsis before sealing and restoring.

Consepsis V's viscosity enables the material to be placed

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EN and retained in the canal without quickly running out. Note: Do not use Consepsis or Consepsis V in a root canal which is wet with sodium hypochlorite, as a brown precipitate will form.

Pre-Treatment Procedure:

- Using 2ml prefilled syringes:
 - Remove Luer Lock cap from syringe.
 - Twist disposable tip firmly onto syringe: Consepsis: Use a Black Mini™ Brush tip or a Blue Mini™ Dento-Infusor™ tip.
 - For optimum control, hold syringe in palm of hand rather than in conventional manner with thumb on plunger (Fig. 1).
- Using 30ml IndiSpense™ System:
 - A pre-labeled unit dose syringe is attached to the end of the large, no-waste IndiSpense syringe by turning the Luer Lock of the unit dose syringe snugly onto the male thread of the IndiSpense (Fig. 2).
 - Using the palm grasp (Fig. 3) the plunger of the IndiSpense syringe is depressed while the unit dose syringe is stabilized at the plunger with the non-dominant hand. To prevent cross-contamination, a used syringe should not be reattached to the IndiSpense syringe. Do not express contents of unit dose syringe back into the IndiSpense syringe.
 - If not using the pre-loaded syringe immediately, replace Luer Lock cap until ready for use.
 - Before use, twist recommended tip firmly onto the unit dose syringe. (Tip recommendation: Consepsis: See A.2. Above. Consepsis V: Black Micro™ FX or

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Inspiral™ Brush Tip.)

Procedure:

- Consepsis
 - When placing direct-placed or resin-luted restorations, etch the surfaces to be bonded before applying Consepsis. Dispense a small amount of Consepsis and thoroughly rub preparation (Figs. 4 and 5). Air dry. Do not rinse. Note: Some clinicians instruct to etch, disinfect, and then apply bonding primers. Our tests have demonstrated no reduction in bond strength whether Consepsis is placed before or after etching. Our testing has also shown a slight gain in bond strength by moistening dentin with Consepsis followed by removal of excess moisture on surface only and placement of A/B Primers.
 - Remove excess fluid to obtain a moist dentin surface (for recommended "wet" bonding techniques see PQ1™ or PermaQuick™ bonding systems instructions).
 - Continue with appropriate bonding technique according to instructions.

Consepsis Treatment Notes:

- For composite and amalgam preparations, etc., treat surfaces the same as temporary crown cementation. Place appropriate dentin liner and/or base as needed (Ultra-Blend™ plus™ is recommended).
- For composites, resin luting, amalgams, and glass ionomers, disinfect before or after etching and bonding with a quality bonding agent such as PQ1 or PermaQuick.

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- For disinfecting close to pulp, gently scrub with Consepsis liquid for 60 seconds. Gently air dry (don't rinse). Follow with application of Ultra-Blend plus, Ultra-Etch™ 35% phosphoric acid and a bonding agent to seal the dentin.⁵ Use all materials according to their instructions.

B. Consepsis V

- Upon completion of cleaning and shaping canals, rinse with Consepsis V. Leave Consepsis V in canals for 60 seconds.
- Rinse thoroughly with sterile water, or EDTA.
- Dry completely with Ultradent Luer Vacuum Adapter and paper points.
- Obturate with materials of choice.

Consepsis/Consepsis V Precautions:

- Caution: Chemical may set in tip after a few hours. To avoid clogging, remove tip and re-cap with Luer Lock Caps between each use. Always verify material flow prior to applying intraorally. If resistance is met, replace tip and re-check. This does NOT apply to dual-barrel syringes and mixing tips.
- tips are disposable; do not reuse.
- DO NOT allow volumes of material to be ingested. Collect extraneous material with proper vacuuming techniques.
- DO NOT use on patients with known sensitivity to chlorhexidine.
- Use Ultradent™ Syringe Covers for a protective barrier to prevent cross-contamination.
- Consepsis V is specifically designed for endodontic disinfecting purposes only.

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Consepsis™ & Consepsis™ V

Chlorhexidine Antibacterial Solution

Manufactured by:
Ultradent Products, Inc.
505 West Ultradent Drive (10200 South)
South Jordan, Utah 84095 USA

Manufactured in the USA from globally sourced materials

Australian Sponsor:
Ultradent Australia Pty Ltd
16/1 Market St. Sydney, NSW 2000
1800 29 09 29 | www.ultradent.com.au

ULTRADENT
PRODUCTS, INC.

8 0 0 . 5 5 2 . 5 5 1 2 | | www.ultradent.com
© 2022 Ultradent Products, Inc. All Rights Reserved. 16024AR13 021122

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U.S. Federal law restricts this device to sale by or on the order of a State licensed dental professional.

For immediate reorder and/or complete descriptions of Ultradent's product line, refer to Ultradent's catalog or call Toll Free 1-800-552-5512. Outside U.S. call (801) 572-4200, or visit www.ultradent.com.

- 1 Cao, D., Hollis, R., Christensen, G.J., Christensen, R.P., "Effect of Tooth Disinfecting Procedures on Dentin Shear Bond Strength," Journal of Dental Research, Vol. 74 (AADR Abstract #493) CRA reports to IADR (AADR) San Antonio Meeting, March 1995.
- 2 Dickinson, G., Smith, C., Morris, C., Murphy, S., "New Universal Bonding Systems for Aesthetic Restorations," Practical Periodontics & Aesthetic Dentistry, June/July 1995.
- 3 "Clinically Significant Research," Clinical Research Associates Newsletter, Volume 18, Issue 4, April 1994.
- 4 Brännström, M. "Infection beneath composite resin restoration: Can it be avoided?" Operative Dentistry, Autumn, 1987:12(4); 158-163.
- 5 Pameijer, C., D.M.D., PhD., Histopathology Study with Ultra-Blend, 1998 (unpublished).

- 1 Cao, D., Hollis, R., Christensen, G.J., Christensen, R.P., "Effect of Tooth Disinfecting Procedures on Dentin Shear Bond Strength," Journal of Dental Research, Vol. 74 (AADR Abstract #493) CRA reports to IADR (AADR) San Antonio Meeting, March 1995.
- 2 Dickinson, G., Smith, C., Morris, C., Murphy, S., "New Universal Bonding Systems for Aesthetic Restorations," Practical Periodontics & Aesthetic Dentistry, June/July 1995.
- 3 "Clinically Significant Research," Clinical Research Associates Newsletter, Volume 18, Issue 4, April 1994.
- 4 Brännström, M. "Infection beneath composite resin restoration: Can it be avoided?" Operative Dentistry, Autumn, 1987:12(4); 158-163.
- 5 Pameijer, C., D.M.D., PhD., Histopathology Study with Ultra-Blend, 1998 (unpublished).

- 1 Cao, D., Hollis, R., Christensen, G.J., Christensen, R.P., "Effect of Tooth Disinfecting Procedures on Dentin Shear Bond Strength," Journal of Dental Research, Vol. 74 (AADR Abstract #493) CRA reports to IADR (AADR) San Antonio Meeting, Marzo 1995.
- 2 Dickinson, G., Smith, C., Morris, C., Murphy, S., "New Universal Bonding Systems for Aesthetic Restorations," Practical Periodontics & Aesthetic Dentistry, Junio/Julio 1995.
- 3 "Clinically Significant Research," Clinical Research Associates Newsletter, Volume 18, Issue 4, Abril 1994.
- 4 Brännström, M. "Infection beneath composite resin restoration: Can it be avoided?" Operative Dentistry, Autumn, 1987:12(4); 158-163.
- 5 Pameijer, C., D.M.D., PhD., Histopathology Study with Ultra-Blend, 1998 (no publicado).